

*2015-2016* | PRODUCT CATALOG

## Preface

TransGen Biotech, Inc. is a researcher, developer, manufacturer and distributor of more than 200 molecular and cellular biology products and kits for life science research and molecular diagnostics. In 2001, the company was founded by three scientists with a mission to produce innovative and cost-effective products for life science research.

In March 21, 2006, TransGen Biotech was incorporated in Beijing, China. The company's headquarters, R&D, and manufacturing facility are located in Beijing. To date, the company has more than 200 scientists in Beijing, and has more than 30 distribution centers covering all major cities in China. Our extensive R&D experience and state-of-the-art facilities enable us to keep generating the most innovative and the highest quality products. Since 2006, the company was consecutively awarded as one of the "High Tech Corporation in Beijing" by Beijing local government.

Currently, our products cover: plasmid based DNA markers, high efficiency chemically competent cells, 5 minutes PCR product cloning and expression vectors, a variety of PCR enzymes and SuperMixes, RNase H deficient and high temperature RT enzymes, qPCR and qRT-PCR SuperMixes, 5 minutes fast restriction enzymes the highest efficiency mutagenesis kits, high quality nucleic acid extraction and purification kits, unstained and prestained protein markers, Western blot markers, and protein purification resins, cell culture and transfection reagents, antibodies.

As the leading bioreagent company in China, we are looking forward to partner with you in your quest for ground-breaking life science discoveries.

CEO: 


Certain products will not be sold in some countries. Please contact TransGen for detailed information.




# TransStart® Hot Start ( Double Blocking )



## PCR Preparation without Hot Start



## PCR Preparation with *TransStart*® Method



- Blocking efficiency up to 100%.
- Different from *Taq* antibody blocking, risks of contamination from mammals DNA are avoided.
- Different from chemical modified blocking, long denaturing step is not needed.

- *TransStart*® *Taq* DNA Polymerase
- *TransStart*® *TopTaq* DNA Polymerase
- *TransStart*® *FastPfu* DNA Polymerase
- 2x *TransStart*® *FastPfu* PCR SuperMix
- *TransStart*® *FastPfu* Fly DNA Polymerase
- *TransStart*® *KD Plus* DNA Polymerase
- *TransStart*® Green qPCR SuperMix
- *TransStart*® Green qPCR SuperMix UDG
- *TransStart*® Top Green qPCR SuperMix
- *TransStart*® Tip Green qPCR SuperMix
- *TransStart*® Probe qPCR SuperMix

# TransStart® FastPfu DNA Polymerase

## TransStart® FastPfu Fly DNA Polymerase



### Fast, high fidelity, hot start DNA polymerase

#### Fast extension rate

TransStart® FastPfu DNA polymerase has an extension rate of up to 4 kb/min.

TransStart® FastPfu Fly DNA polymerase has an extension rate of up to 6 kb/min.

- ◆ Fast
- ◆ Highest Fidelity
- ◆ High Sensitivity

#### High fidelity


TransStart® FastPfu DNA Polymerase offers 54-fold fidelity as compared to EasyTaq® DNA Polymerase.

TransStart® FastPfu Fly DNA Polymerase offers 108-fold fidelity as compared to EasyTaq® DNA Polymerase.

#### Better amplification efficiency

Suitable for long fragment or low copy gene amplification

Amplification using TransStart® FastPfu DNA Polymerase



M: 1 Kb Plus DNA Ladder

|             |        |            |
|-------------|--------|------------|
| 1: NCBP     | 2.5 kb | 2 h 20 min |
| 2: ACTR     | 3 kb   | 2 h 20 min |
| 3: HDP      | 3.5 kb | 2 h 20 min |
| 4: β-globin | 3 kb   | 1 h 27 min |
| 5: Rhod     | 4.1 kb | 1 h 27 min |
| 6: β-globin | 4.1 kb | 1 h 27 min |
| 7: UDG      | 7 kb   | 1 h 36 min |
| 8: LN       | 10 kb  | 1 h 55 min |





4 kb: Genomic DNA;  
7 kb and 10 kb: Plasmid DNA



# ***TransDirect***<sup>®</sup> **Blood PCR Kit**





 High resistance to inhibitors and impurities.

 Direct PCR amplification using the whole blood or cell culture as template without DNA extraction.




# ***TransScript***<sup>®</sup> **Fly** **First-Strand cDNA** **Synthesis SuperMix**

-  Fast: 5 minutes reverse transcription.
-  cDNA up to 12 kb.





# One-Step gDNA Removal and cDNA Synthesis SuperMix



- Simultaneous genomic DNA removal and cDNA synthesis.
- Easy to use SuperMix.

• *TransScript*<sup>®</sup> -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix (42°C-65°C)

• *EasyScript*<sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix (42°C)


• *TransScript*<sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix (42°C)

• *TransScript*<sup>®</sup> II One-Step gDNA Removal and cDNA Synthesis SuperMix (42°C-55°C)



# **TransScript®-Uni Cell to cDNA Synthesis SuperMix for qPCR**

- Resulting lysate without purification can be directly used for reverse transcription.
- Simultaneous genomic DNA removal and cDNA synthesis.
- Suitable for qPCR directly from cells.





# RT “All-in-One” SuperMix

Primers

dNTPs

RI

All-Mix

RTase

Buffer



- Easy: All components (except RNA template) are premixed.
- Fast: 30 minutes RT reaction for PCR template; 15 minutes one step gDNA removal and RT reaction for qPCR template.

- *TransScript*® All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)
- *TransScript*® II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)
- *TransScript*® All-in-One First-Strand cDNA Synthesis SuperMix for PCR
- *TransScript*® II All-in-One First-Strand cDNA Synthesis SuperMix for PCR

# FlyCut™ Restriction Enzymes



- Fast, 5 Minutes Digestion
- High Efficiency
- No Star Activity
- Universal Buffer

- *FlyCut™* Avr II-HF
- *FlyCut™* Bam HI-HF
- *FlyCut™* Bgl II-HF
- *FlyCut™* Bsg I-HF
- *FlyCut™* Eag I-HF

- *FlyCut™* EcoR I-HF
- *FlyCut™* EcoR V-HF
- *FlyCut™* Hind III-HF
- *FlyCut™* Kpn I-HF
- *FlyCut™* Nco I-HF
- *FlyCut™* Nde I-HF
- *FlyCut™* Nhe I-HF
- *FlyCut™* Not I-HF
- *FlyCut™* Pst I-HF
- *FlyCut™* Pvu I-HF

- *FlyCut™* Sac I-HF
- *FlyCut™* Sac II-HF
- *FlyCut™* Sal I-HF
- *FlyCut™* Sca I-HF
- *FlyCut™* Sma I-HF
- *FlyCut™* Spe I-HF
- *FlyCut™* Sph I-HF
- *FlyCut™* Xba I-HF
- *FlyCut™* Xho I-HF
- *FlyCut™* Xma I-HF






# TransStart® Tip Green qPCR SuperMix




- ◆ *High specificity*
- ◆ *High sensitivity*

- 🔬 A combination of chemical blocking technique and *TransStart*® hot start technique to achieve complete blocking. Compared with double blocking *TransStart*® *TopTaq*, this method provides higher sensitivity, and better amplification.
- 🔬 Double cation ( $K^+$ ,  $NH_4^+$ ) buffer enhances specificity and reduces primer-dimers formation.
- 🔬 Passive reference dyes for different qPCR instruments.




High reproducibility (Roche LightCycler480)



High consistency (Bio-Rad CFX96)



# **pEASY<sup>®</sup> Cloning** **Room Temperature** **5 Minutes Fast Cloning**



- pEASY<sup>®</sup>-T1 Cloning Kit
- pEASY<sup>®</sup>-Blunt Cloning Kit
- pEASY<sup>®</sup>-T1 Simple Cloning Kit
- pEASY<sup>®</sup>-Blunt Simple Cloning Kit
- pEASY<sup>®</sup>-T3 Cloning Kit
- pEASY<sup>®</sup>-Blunt3 Cloning Kit
- pEASY<sup>®</sup>-T5 Zero Cloning Kit
- pEASY<sup>®</sup>-Blunt Zero Cloning Kit


- pEASY<sup>®</sup>-Blunt E1 Expression Kit
- pEASY<sup>®</sup>-Blunt E2 Expression Kit
- pEASY<sup>®</sup>-Blunt M2 Expression Kit
- pEASY<sup>®</sup>-Blunt M3 Expression Kit



# pEASY<sup>®</sup>-Uni Seamless Cloning and Assembly Kit

- Fast:** 15 minutes.
- Broad:** no restriction enzyme digestions. Can be cloned into any sites.
- High efficiency:** up to 95% cloning efficiency.
- Seamless:** no extra sequences introduced; up to 5 fragments assembly.

1. Prepare linearized vector by PCR/Enzyme digestion
2. PCR amplify inserts with 15-25 bp overlapping sequences
3. Mix vector, DNA fragments and Assembly Mix together, incubate at 50°C for 15 minutes
4. Transformation



# Fast Mutagenesis System

## High fidelity and fast amplification

2xTransStart® FastPfu PCR SuperMix improves the fidelity and shorts the amplification time.


## Visible

Amplification products can be visualized on agarose gel.

## High efficiency

Both primers have the desired mutation sites providing higher mutation efficiency. DMT enzyme digests parental plasmids *in vitro* and DMT competent cell digests parental plasmids *in vivo* providing much lower background.

- ◆ Fast
- ◆ Convenient
- ◆ High Performance



### PCR amplification

Mutagenesis by PCR amplification with two overlapping primers. Both primers contain the target mutations.

### Digestion of parental plasmid


DMT enzyme digests and DMT competent cell further digests parental plasmids.

X= mutation





# ProteinIso<sup>®</sup> Ni-NTA Resin



- High selectivity for high purity.
- Binding under denaturing and non-denaturing conditions.
- Easy to regenerate.


- ◆ *High adsorption capacity*
- ◆ *Good selectivity*
- ◆ *Strong permeability*
- ◆ *Easy regeneration*



# Blue Plus® IV Protein Marker

## Visible estimation of protein electrophoresis and membrane transfer efficiency

It is composed of prestained proteins from 10 kDa to 180 kDa. Different color bands are favorable to monitor electrophoresis and estimate membrane transfer efficiency.



## Convenience

Ready-to-use format.





# EasySee® Western Marker





## Visible estimation of protein electrophoresis and membrane transfer efficiency

It is composed of prestained and unstained proteins from 25 kDa to 150 kDa. Different color bands are favorable to monitor electrophoresis, estimate membrane transfer efficiency and determine direction of membrane transfer.

## Real visualization and accuracy

Bands from unstained proteins are visible by alkaline phosphatase and horseradish peroxidase chemiluminescence detection, providing more accurate molecular weight estimation than dye-detached protein markers.

- ◆ *Visible Western Blot*
- ◆ *High sensitivity*




## Convenience

Ready-to-use format.

# TransSerum® HQ Fetal Bovine Serum



- Low toxicity
- Better cell growth
- Suitable for a broad range of cells



- ◆ *Fast growth*
- ◆ *Short doubling time*




# ***TransLipid***<sup>®</sup> HL Transfection Reagent



## **High efficiency and Low cytotoxicity**


- Transfect DNA, RNA, siRNA.
- Adherent or suspension cells.
- Can be used in the presence of serum and antibiotics.



# ***TransIn™* EL Transfection Reagent**


**High efficiency and Low cytotoxicity**

- Non-liposomal Transfection Reagent.
- Adherent or suspension cells.
- Can be used in the presence of serum and antibiotics.






# ***TransDetect***<sup>®</sup> **Double-Luciferase Reporter Assay Kit**



- ◆ *Fast detection*
- ◆ *High sensitivity*
- ◆ *Broad detection range*
- ◆ *No endogenous activity*




# TransDetect® Cell Counting Kit (CCK)



- ◆ *Fast and sensitive*
- ◆ *Minimal cytotoxicity*
- ◆ *Broad linear range*
- ◆ *High reproducibility*

## Detection of CCK



Low cytotoxicity





# TransDetect® Annexin V -FITC/PI Cell Apoptosis Detection Kit

# TransDetect® Annexin V -EGFP/PI Cell Apoptosis Detection Kit


- ◆ High sensitivity
- ◆ High specificity



- 🔬 Rapid fluorescent detection of annexin V bound cells.
- 🔬 No cell fixation, the cells can be used for further study after this assay.
- 🔬 Propidium iodide provided to differentiate apoptotic cells from viable and necrotic cells.



# Primary Antibodies Fluorescent-labeled Secondary Antibodies



Fusion protein



Primary antibodies



XX Conjugated  
secondary antibodies

Protein detection



- ◆ *High sensitivity*
- ◆ *High specificity*



# Index

## Chapter 1 PCR, RT-PCR, qPCR and qRT-PCR

### PCR Enzyme

|  |     |
|--|-----|
| <i>TransFast</i> <sup>®</sup> Taq DNA Polymerase .....                 | 006 |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase .....                       | 007 |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase for PAGE .....              | 009 |
| <i>TransTaq</i> <sup>®</sup> -T DNA Polymerase .....                   | 010 |
| <i>TransTaq</i> <sup>®</sup> DNA Polymerase High Fidelity (HiFi) ..... | 011 |
| <i>TransStart</i> <sup>®</sup> Taq DNA Polymerase .....                | 014 |
| <i>TransStart</i> <sup>®</sup> TopTaq DNA Polymerase .....             | 016 |
| <i>EasyPfu</i> DNA Polymerase .....                                    | 018 |
| <i>TransStart</i> <sup>®</sup> FastPfu DNA Polymerase .....            | 019 |
| <i>TransStart</i> <sup>®</sup> FastPfu Fly DNA Polymerase .....        | 021 |
| <i>TransStart</i> <sup>®</sup> KD Plus DNA Polymerase .....            | 023 |
| GC Enhancer .....  | 025 |
| PCR Stimulant .....  | 025 |

### PCR SuperMix

|   |     |
|---|-----|
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix .....                       | 027 |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix for PAGE .....              | 028 |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix .....                   | 029 |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix ..... | 030 |
| 2× <i>EasyPfu</i> PCR SuperMix .....                                    | 032 |
| 2× <i>TransStart</i> <sup>®</sup> FastPfu PCR SuperMix .....            | 033 |

### Direct PCR

|   |     |
|---|-----|
| <i>TransDirect</i> <sup>®</sup> Animal Tissue PCR Kit ..... | 034 |
| <i>TransDirect</i> <sup>®</sup> Plant Tissue PCR Kit .....  | 036 |
| <i>TransDirect</i> <sup>®</sup> Blood PCR Kit .....         | 037 |

### RT-PCR

|   |     |
|---|-----|
| <i>EasyScript</i> <sup>®</sup> Reverse Transcriptase .....  | 041 |
| <i>TransScript</i> <sup>®</sup> Reverse Transcriptase .....   | 043 |
| <i>TransScript</i> <sup>®</sup> II Reverse Transcriptase .....  | 044 |
| <i>EasyScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix .....   | 045 |
| <i>EasyScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix .....                                    | 047 |
| <i>TransScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix .....  | 048 |
| <i>TransScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix .....                                   | 049 |
| <i>TransScript</i> <sup>®</sup> Fly First-Strand cDNA Synthesis SuperMix .....  | 050 |
| <i>TransScript</i> <sup>®</sup> -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix .....                              | 051 |
| <i>TransScript</i> <sup>®</sup> -Uni Cell to cDNA Synthesis SuperMix for qPCR .....                                       | 053 |
| <i>TransScript</i> <sup>®</sup> miRNA First-Strand cDNA Synthesis SuperMix .....  | 055 |
| <i>TransScript</i> <sup>®</sup> II First-Strand cDNA Synthesis SuperMix .....   | 057 |
| <i>TransScript</i> <sup>®</sup> II One-Step gDNA Removal and cDNA Synthesis SuperMix .....                                | 058 |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for PCR .....                             | 059 |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) .....    | 060 |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for PCR .....                          | 062 |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) ..... | 063 |



|   |     |
|---|-----|
| <i>TransScript</i> <sup>®</sup> Two-Step RT-PCR SuperMix .....              | 064 |
| <i>TransScript</i> <sup>®</sup> II Two-Step RT-PCR SuperMix .....           | 065 |
| <i>EasyScript</i> <sup>®</sup> One-Step RT-PCR SuperMix .....               | 066 |
| <i>TransScript</i> <sup>®</sup> One-Step RT-PCR SuperMix .....              | 067 |
| <i>TransScript</i> <sup>®</sup> II One-Step RT-PCR SuperMix .....           | 068 |
| Ribonuclease Inhibitor .....  | 070 |
| <br>  |     |
| qPCR and qRT-PCR SuperMix   |     |
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix .....                    | 071 |
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix UDG .....                | 073 |
| <i>TransStart</i> <sup>®</sup> Top Green qPCR SuperMix .....                | 074 |
| <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix .....                | 075 |
| <i>TransScript</i> <sup>®</sup> Green Two-Step qRT-PCR SuperMix .....       | 076 |
| <i>TransScript</i> <sup>®</sup> Green miRNA Two-Step qRT-PCR SuperMix ..... | 078 |
| <i>TransScript</i> <sup>®</sup> II Green Two-Step qRT-PCR SuperMix .....    | 079 |
| <i>TransScript</i> <sup>®</sup> Green One-Step qRT-PCR SuperMix .....       | 080 |
| <i>TransScript</i> <sup>®</sup> II Green One-Step qRT-PCR SuperMix .....    | 082 |
| <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix .....                    | 084 |
| <i>TransScript</i> <sup>®</sup> Probe One-Step qRT-PCR SuperMix .....       | 085 |
| <i>TransScript</i> <sup>®</sup> II Probe One-Step qRT-PCR SuperMix .....    | 087 |
| <br>  |     |
| High Pure dNTPs   |     |
| High Pure dNTPs .....   | 089 |

## Chapter 2 DNA Molecular Weight Standards

|  |     |
|--|-----|
| <i>Trans2K</i> <sup>®</sup> DNA Marker .....         | 092 |
| <i>Trans2K</i> <sup>®</sup> Plus DNA Marker .....    | 092 |
| <i>Trans2K</i> <sup>®</sup> Plus II DNA Marker ..... | 092 |
| <i>Trans5K</i> DNA Marker .....                      | 093 |
| <i>Trans15K</i> DNA Marker .....                     | 093 |
| 1Kb DNA Ladder .....                                 | 093 |
| 1Kb Plus DNA Ladder .....                            | 094 |
| 100bp DNA Ladder .....                               | 094 |
| 100bp Plus DNA Ladder .....                          | 094 |
| 100bp Plus II DNA Ladder .....                       | 095 |
| GelStain .....                                       | 095 |
| Agarose .....  | 095 |

## Chapter 3 Cloning and Mutagenesis System

|  |     |
|--|-----|
| Cloning Vectors  |     |
| <i>pEASY</i> <sup>®</sup> -T1 Cloning Kit .....                        | 099 |
| <i>pEASY</i> <sup>®</sup> -Blunt Cloning Kit .....                     | 102 |
| <i>pEASY</i> <sup>®</sup> -T1 Simple Cloning Kit .....                 | 103 |
| <i>pEASY</i> <sup>®</sup> -Blunt Simple Cloning Kit .....              | 104 |
| <i>pEASY</i> <sup>®</sup> -T3 Cloning Kit .....                        | 105 |
| <i>pEASY</i> <sup>®</sup> -Blunt3 Cloning Kit .....                    | 106 |
| <i>pEASY</i> <sup>®</sup> -T5 Zero Cloning Kit .....                   | 107 |
| <i>pEASY</i> <sup>®</sup> -Blunt Zero Cloning Kit .....                | 108 |
| <i>pEASY</i> <sup>®</sup> -Uni Seamless Cloning and Assembly Kit ..... | 109 |



|   |     |
|---|-----|
| Cloning Competent Cells   |     |
| <i>Trans</i> 10 Chemically Competent Cell .....                   | 112 |
| <i>Trans</i> 5 $\alpha$ Chemically Competent Cell .....           | 112 |
| <i>Trans</i> 109 Chemically Competent Cell .....                  | 113 |
| <i>Trans</i> 110 Chemically Competent Cell .....                  | 113 |
| <i>Trans</i> 1-Blue Chemically Competent Cell .....               | 113 |
| <i>Trans</i> 2-Blue Chemically Competent Cell .....               | 114 |
| <i>Trans</i> 1-T1 Phage Resistant Chemically Competent Cell ..... | 114 |
| DMT Chemically Competent Cell .....                               | 114 |
| <i>Trans</i> Stbl3 Chemically Competent Cell .....                | 115 |
| <i>Trans</i> DB3.1 Chemically Competent Cell .....                | 115 |

|  |     |
|--|-----|
| Mutagenesis System                             |     |
| <i>Fast</i> Mutagenesis System .....           | 116 |
| <i>Fast</i> MultiSite Mutagenesis System ..... | 117 |

## Chapter 4 Nucleic Acid Purification

|   |     |
|---|-----|
| Genomic DNA Purification  |     |
| <i>BloodZol</i> .....   | 121 |
| <i>PlantZol</i> .....   | 122 |
| <i>EasyPure</i> <sup>®</sup> Genomic DNA Kit .....                          | 123 |
| <i>EasyPure</i> <sup>®</sup> Plant Genomic DNA Kit .....                    | 125 |
| <i>EasyPure</i> <sup>®</sup> Blood Genomic DNA Kit .....                    | 126 |
| <i>EasyPure</i> <sup>®</sup> Marine Animal Genomic DNA Kit .....            | 127 |
| <i>EasyPure</i> <sup>®</sup> Bacteria Genomic DNA Kit .....                 | 128 |
| <i>EasyPure</i> <sup>®</sup> Food and Fodder Security Genomic DNA Kit ..... | 129 |
| <i>EasyPure</i> <sup>®</sup> Micro Genomic DNA Kit .....                    | 131 |

|  |     |
|--|-----|
| Plasmid DNA Purification and <i>E.coli</i> Medium              |     |
| <i>EasyPure</i> <sup>®</sup> Plasmid MiniPrep Kit .....        | 132 |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MiniPrep Kit ..... | 133 |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MaxiPrep Kit ..... | 134 |
| <i>ArtMedia</i> <sup>®</sup> Plasmid Culture .....             | 135 |

|   |     |
|---|-----|
| DNA Purification  |     |
| <i>EasyPure</i> <sup>®</sup> PCR Purification Kit .....     | 136 |
| <i>EasyPure</i> <sup>®</sup> Quick Gel Extraction Kit ..... | 137 |

|   |     |
|---|-----|
| RNA Purification  |     |
| <i>TransZol</i> .....                                   | 138 |
| <i>TransZol</i> Up .....                                | 139 |
| <i>TransZol</i> Plant .....                             | 140 |
| <i>EasyPure</i> <sup>®</sup> RNA Kit .....              | 141 |
| <i>EasyPure</i> <sup>®</sup> Viral DNA/RNA Kit .....    | 142 |
| <i>EasyPure</i> <sup>®</sup> Plant RNA Kit .....        | 143 |
| <i>EasyPure</i> <sup>®</sup> Blood RNA Kit .....        | 144 |
| <i>TransZol</i> Up Plus RNA Kit .....                   | 145 |
| <i>EasyPure</i> <sup>®</sup> miRNA Kit .....            | 146 |
| <i>EasyPure</i> <sup>®</sup> RNA Purification Kit ..... | 147 |
| <i>RNAhold</i> <sup>®</sup> .....                       | 148 |



## Chapter 5 Gene Expression

|  |     |
|--|-----|
| Prokaryotic Expression Vectors                           |     |
| <i>pEASY</i> <sup>®</sup> -Blunt E1 Expression Kit ..... | 151 |
| <i>pEASY</i> <sup>®</sup> -Blunt E2 Expression Kit ..... | 154 |
| Expression Medium  |     |
| <i>ArtMedia</i> <sup>®</sup> Protein Expression .....    | 155 |
| Expression Competent Cells                               |     |
| BL21(DE3) Chemically Competent Cell .....                | 156 |
| BL21(DE3)pLysS Chemically Competent Cell .....           | 156 |
| <i>Transetta</i> (DE3) Chemically Competent Cell .....   | 157 |
| <i>TransB</i> (DE3) Chemically Competent Cell .....      | 157 |
| BL21 Chemically Competent Cell .....                     | 157 |
| Mammalian Expression Vectors                             |     |
| <i>pEASY</i> <sup>®</sup> -Blunt M2 Expression Kit ..... | 158 |
| <i>pEASY</i> <sup>®</sup> -Blunt M3 Expression Kit ..... | 161 |

## Chapter 6 Protein Extraction, Purification and Detection

|   |     |
|---|-----|
| Protein Extraction  |     |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Total Protein Extraction Kit .....                   | 164 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Nuclear and Cytoplasmic Protein Extraction Kit ..... | 165 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Membrane Protein Extraction Kit .....                | 166 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit for Cultured Cells .....  | 167 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit for Tissue .....          | 168 |
| Protein Purification  |     |
| <i>ProteinIso</i> <sup>®</sup> Ni-NTA Resin .....   | 169 |
| <i>ProteinIso</i> <sup>®</sup> Ni-IDA Resin .....   | 171 |
| <i>ProteinIso</i> <sup>®</sup> GST Resin .....  | 173 |
| <i>ProteinIso</i> <sup>®</sup> Protein A Resin .....  | 175 |
| <i>ProteinIso</i> <sup>®</sup> Protein G Resin .....  | 177 |
| Unstained Protein Marker  |     |
| <i>ProteinRuler</i> <sup>®</sup> I .....  | 180 |
| <i>ProteinRuler</i> <sup>®</sup> II .....   | 180 |
| <i>ProteinRuler</i> <sup>®</sup> IV .....   | 181 |
| Prestained Protein Marker   |     |
| <i>Blue Plus</i> <sup>®</sup> Protein Marker .....  | 182 |
| <i>Blue Plus</i> <sup>®</sup> II Protein Marker .....   | 182 |
| <i>Blue Plus</i> <sup>®</sup> III Protein Marker .....  | 183 |
| <i>Blue Plus</i> <sup>®</sup> IV Protein Marker .....   | 183 |
| Western Protein Marker  |     |
| <i>EasySee</i> <sup>®</sup> Western Marker .....  | 184 |
| <i>EasySee</i> <sup>®</sup> II Western Marker .....   | 185 |



## Related Products

|   |     |
|---|-----|
| <i>EasySee</i> <sup>®</sup> Western Blot Kit .....            | 186 |
| 6xProtein Loading Buffer .....                                | 186 |
| <i>Easy</i> Protein Quantitative Kit (Bradford) .....         | 187 |
| <i>Easy</i> II Protein Quantitative Kit (BCA) .....           | 188 |
| <i>ProteinEle</i> <sup>™</sup> Precast Tris-Glycine Gel ..... | 190 |

## Chapter 7 Cell Culture and Detection

### Cell Culture

|  |     |
|--|-----|
| <i>TransSerum</i> <sup>®</sup> HQ Fetal Bovine Serum .....   | 192 |
| <i>TransLipid</i> <sup>®</sup> HL Transfection Reagent ..... | 193 |
| <i>TransIn</i> <sup>™</sup> EL Transfection Reagent .....    | 195 |
| Penicillin-Streptomycin (100x) .....                         | 197 |
| L-Glutamine (100x) .....                                     | 197 |
| Trypsin .....  | 197 |
| G418 .....   | 197 |
| PBS (1x) .....   | 198 |

### Cell Detection

|   |     |
|---|-----|
| <i>TransDetect</i> <sup>®</sup> Double-Luciferase Reporter Assay Kit .....                          | 198 |
| <i>TransDetect</i> <sup>®</sup> Cell Counting Kit (CCK) .....                                       | 199 |
| <i>TransDetect</i> <sup>®</sup> Annexin V-FITC/PI Cell Apoptosis Detection Kit .....                | 201 |
| <i>TransDetect</i> <sup>®</sup> Annexin V-EGFP/PI Cell Apoptosis Detection Kit .....                | 202 |
| <i>TransDetect</i> <sup>®</sup> <i>In Situ</i> Fluorescein TUNEL Cell Apoptosis Detection Kit ..... | 203 |

## Chapter 8 Antibodies

### Primary Antibodies

|   |     |
|---|-----|
| <i>ProteinFind</i> <sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody .....        | 205 |
| <i>ProteinFind</i> <sup>®</sup> Anti-DYKDDDDK Tag Mouse Monoclonal Antibody ..... | 205 |
| <i>ProteinFind</i> <sup>®</sup> Anti-HA Mouse Monoclonal Antibody .....           | 206 |
| <i>ProteinFind</i> <sup>®</sup> Anti-V5 Mouse Monoclonal Antibody .....           | 207 |
| <i>ProteinFind</i> <sup>®</sup> Anti-His Mouse Monoclonal Antibody .....          | 207 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GST Mouse Monoclonal Antibody .....          | 208 |
| <i>ProteinFind</i> <sup>®</sup> Anti-MBP Mouse Monoclonal Antibody .....          | 208 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GFP Mouse Monoclonal Antibody .....          | 209 |

### Control Antibodies

|  |     |
|--|-----|
| <i>ProteinFind</i> <sup>®</sup> Anti- $\beta$ -Tubulin Mouse Monoclonal Antibody ..... | 209 |
| <i>ProteinFind</i> <sup>®</sup> Anti- $\beta$ -Actin Mouse Monoclonal Antibody .....   | 210 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GAPDH Mouse Monoclonal Antibody .....             | 210 |

### Secondary Antibodies

|  |     |
|--|-----|
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), HRP Conjugate .....   | 211 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), FITC Conjugate .....  | 211 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), PE Conjugate .....    | 212 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), AF488 Conjugate ..... | 212 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), HRP Conjugate .....    | 213 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), FITC Conjugate .....   | 214 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), PE Conjugate .....     | 214 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), AF488 Conjugate .....  | 215 |



---

### Related Products

|                                 |     |
|---------------------------------|-----|
| TMB ELISA Substrate .....       | 216 |
| Super TMB ELISA Substrate ..... | 216 |

## Chapter 9 Other Products

|  |     |
|--|-----|
| T4 DNA Ligase .....                        | 218 |
| DMT Enzyme .....                           | 219 |
| DNase I (RNase-free) .....                 | 219 |
| RNase A .....                              | 220 |
| Proteinase K .....                         | 220 |
| IPTG .....                                 | 220 |
| X-gal .....                                | 221 |
| Ampicillin .....                           | 221 |
| Kanamycin .....                            | 221 |
| Chloramphenicol .....                      | 221 |
| 6×DNA Loading Buffer .....                 | 222 |
| 2×RNA Loading Buffer .....                 | 222 |
| ddH <sub>2</sub> O .....                   | 222 |
| RNase-free Water .....                     | 222 |
| T7 High Efficiency Transcription Kit ..... | 223 |
| shRNA Synthesis Kit .....                  | 225 |

|                |     |
|----------------|-----|
| Services ..... | 227 |
|----------------|-----|



---

## Chapter 1 PCR, RT-PCR, qPCR and qRT-PCR

---

### PCR Enzyme

|  |     |
|--|-----|
| <i>TransFast</i> <sup>®</sup> <i>Taq</i> DNA Polymerase .....          | 006 |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase .....                       | 007 |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase for PAGE .....              | 009 |
| <i>TransTaq</i> <sup>®</sup> -T DNA Polymerase .....                   | 010 |
| <i>TransTaq</i> <sup>®</sup> DNA Polymerase High Fidelity (HiFi) ..... | 011 |
| <i>TransStart</i> <sup>®</sup> <i>Taq</i> DNA Polymerase .....         | 014 |
| <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> DNA Polymerase .....      | 016 |
| <i>EasyPfu</i> DNA Polymerase .....                                    | 018 |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> DNA Polymerase .....     | 019 |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> Fly DNA Polymerase ..... | 021 |
| <i>TransStart</i> <sup>®</sup> <i>KD Plus</i> DNA Polymerase .....     | 023 |
| GC Enhancer .....  | 025 |
| PCR Stimulant .....  | 025 |

### PCR SuperMix

|   |     |
|---|-----|
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix .....                       | 027 |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix for PAGE .....              | 028 |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix .....                   | 029 |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix ..... | 030 |
| 2× <i>EasyPfu</i> PCR SuperMix .....                                    | 032 |
| 2× <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix .....     | 033 |

### Direct PCR

|   |     |
|---|-----|
| <i>TransDirect</i> <sup>®</sup> Animal Tissue PCR Kit ..... | 034 |
| <i>TransDirect</i> <sup>®</sup> Plant Tissue PCR Kit .....  | 036 |
| <i>TransDirect</i> <sup>®</sup> Blood PCR Kit .....         | 037 |

### RT-PCR

|   |     |
|---|-----|
| <i>EasyScript</i> <sup>®</sup> Reverse Transcriptase .....                                | 041 |
| <i>TransScript</i> <sup>®</sup> Reverse Transcriptase .....                               | 043 |
| <i>TransScript</i> <sup>®</sup> II Reverse Transcriptase .....                            | 044 |
| <i>EasyScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix .....                 | 045 |
| <i>EasyScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA<br>Synthesis SuperMix ..... | 047 |

---

## Chapter 1 PCR, RT-PCR, qPCR and qRT-PCR

---

### RT-PCR

|   |     |
|---|-----|
| <i>TransScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix .....  | 048 |
| <i>TransScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix .....                                   | 049 |
| <i>TransScript</i> <sup>®</sup> Fly First-Strand cDNA Synthesis SuperMix .....  | 050 |
| <i>TransScript</i> <sup>®</sup> -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix .....                              | 051 |
| <i>TransScript</i> <sup>®</sup> -Uni Cell to cDNA Synthesis SuperMix for qPCR .....                                       | 053 |
| <i>TransScript</i> <sup>®</sup> miRNA First-Strand cDNA Synthesis SuperMix .....  | 055 |
| <i>TransScript</i> <sup>®</sup> II First-Strand cDNA Synthesis SuperMix .....   | 057 |
| <i>TransScript</i> <sup>®</sup> II One-Step gDNA Removal and cDNA Synthesis SuperMix .....                                | 058 |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for PCR .....                             | 059 |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) .....    | 060 |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for PCR .....                          | 062 |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) ..... | 063 |
| <i>TransScript</i> <sup>®</sup> Two-Step RT-PCR SuperMix .....  | 064 |
| <i>TransScript</i> <sup>®</sup> II Two-Step RT-PCR SuperMix .....   | 065 |
| <i>EasyScript</i> <sup>®</sup> One-Step RT-PCR SuperMix .....   | 066 |
| <i>TransScript</i> <sup>®</sup> One-Step RT-PCR SuperMix .....  | 067 |
| <i>TransScript</i> <sup>®</sup> II One-Step RT-PCR SuperMix .....   | 068 |
| Ribonuclease Inhibitor .....  | 070 |

### qPCR and qRT-PCR SuperMix

|  |     |
|--|-----|
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix .....     | 071 |
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix UDG ..... | 073 |
| <i>TransStart</i> <sup>®</sup> Top Green qPCR SuperMix ..... | 074 |
| <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix ..... | 075 |

---

## Chapter 1 PCR, RT-PCR, qPCR and qRT-PCR

---

### qPCR and qRT-PCR SuperMix

|   |     |
|---|-----|
| <i>TransScript</i> <sup>®</sup> Green Two-Step qRT-PCR SuperMix .....       | 076 |
| <i>TransScript</i> <sup>®</sup> Green miRNA Two-Step qRT-PCR SuperMix ..... | 078 |
| <i>TransScript</i> <sup>®</sup> II Green Two-Step qRT-PCR SuperMix .....    | 079 |
| <i>TransScript</i> <sup>®</sup> Green One-Step qRT-PCR SuperMix .....       | 080 |
| <i>TransScript</i> <sup>®</sup> II Green One-Step qRT-PCR SuperMix .....    | 082 |
| <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix .....                    | 084 |
| <i>TransScript</i> <sup>®</sup> Probe One-Step qRT-PCR SuperMix .....       | 085 |
| <i>TransScript</i> <sup>®</sup> II Probe One-Step qRT-PCR SuperMix .....    | 087 |

### High Pure dNTPs



|                       |     |
|-----------------------|-----|
| High Pure dNTPs ..... | 089 |
|-----------------------|-----|





# PCR Enzymes

PCR enzymes are purified from *E. coli* strains carrying genes for specific DNA polymerase. The choice of DNA polymerase for PCR application highly depends on the characteristics of the system as well as the desired results. The following is PCR Enzyme Selection Guide.



## Feature of PCR Enzymes

| PCR Enzyme                                      | <i>TransFast</i> <sup>®</sup> <i>Taq</i>  | <i>EasyTaq</i> <sup>®</sup> | <i>TransTaq</i> <sup>®</sup> -T | <i>TransTaq</i> <sup>®</sup> HiFi | <i>TransStart</i> <sup>®</sup> <i>Taq</i> | <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> |
|---|---|-----------------------------|---------------------------------|-----------------------------------|---|--|
| Amplification Efficiency                        | <i>TransFast</i> <sup>®</sup> <i>Taq</i> = <i>EasyTaq</i> <sup>®</sup> < <i>TransTaq</i> <sup>®</sup> -T < <i>TransTaq</i> <sup>®</sup> HiFi < <i>TransStart</i> <sup>®</sup> <i>Taq</i> < <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> |                             |                                 |                                   |   |  |
| Specificity                                     | <i>TransFast</i> <sup>®</sup> <i>Taq</i> = <i>EasyTaq</i> <sup>®</sup> < <i>TransTaq</i> <sup>®</sup> -T < <i>TransTaq</i> <sup>®</sup> HiFi < <i>TransStart</i> <sup>®</sup> <i>Taq</i> < <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> |                             |                                 |                                   |   |  |
| Fidelity (vs. <i>EasyTaq</i> <sup>®</sup> )     | 1x  | 1x                          | 18x                             | 18x                               | 18x                                       | 18x  |
| Extension Rate                                  | 6 kb/min  | 1-2 kb/min                  | 1-2 kb/min                      | 1-2 kb/min                        | 1-2 kb/min                                | 1-2 kb/min                                   |
| Hot Start                                       | -   | -                           | +                               | +                                 | +   | +  |
| "A" at 3' end                                   | +   | +                           | +                               | +                                 | +   | +  |
| Product Size<br>(human genomic DNA as template) | ≤4 kb   | ≤4 kb                       | ≤8 kb                           | ≤15 kb                            | ≤15 kb                                    | ≤15 kb                                       |

High quality products



| PCR Enzyme                                      | <i>EasyPfu</i>   | <i>TransStart<sup>®</sup> KD Plus</i> | <i>TransStart<sup>®</sup> FastPfu</i> | <i>TransStart<sup>®</sup> FastPfu Fly</i> |
|---|--|---------------------------------------|---------------------------------------|---|
| Amplification Efficiency                        | <i>EasyPfu</i> < <i>TransStart<sup>®</sup> KD Plus</i> = <i>TransStart<sup>®</sup> FastPfu</i> < <i>TransStart<sup>®</sup> FastPfu Fly</i> |                                       |                                       |   |
| Specificity                                     | <i>EasyPfu</i> < <i>TransStart<sup>®</sup> KD Plus</i> = <i>TransStart<sup>®</sup> FastPfu</i> < <i>TransStart<sup>®</sup> FastPfu Fly</i> |                                       |                                       |   |
| Fidelity (vs. <i>EasyTaq<sup>®</sup></i> )      | 18x  | 108x                                  | 54x                                   | 108x                                      |
| Extension Rate                                  | 0.5 kb/min   | 1 kb/min                              | 2-4 kb/min                            | 2-6 kb/min                                |
| Hot Start                                       | +  | +                                     | +                                     | +   |
| "A" at 3' end                                   | -  | -                                     | -                                     | -   |
| Product Size<br>(human genomic DNA as template) | ≤6 kb  | ≤15 kb                                | ≤15 kb                                | ≤15 kb                                    |
| Product Size<br>(plasmid DNA as template)       | ≤10 kb   | ≤20 kb                                | ≤20 kb                                | ≤20 kb                                    |

## Applications

| PCR Enzyme   | Application  |
|--|--|
| <i>TransFast<sup>®</sup> Taq</i> DNA Polymerase  | routine PCR, fast PCR, colony PCR  |
| <i>EasyTaq<sup>®</sup></i> DNA Polymerase  | routine PCR, colony PCR  |
| <i>EasyTaq<sup>®</sup></i> DNA Polymerase for PAGE   | short fragment PCR   |
| <i>TransTaq<sup>®</sup> -T</i> DNA Polymerase  | complex templates, TA cloning  |
| <i>TransTaq<sup>®</sup> HiFi</i> DNA Polymerase<br><i>TransStart<sup>®</sup> Taq</i> DNA Polymerase<br><i>TransStart<sup>®</sup> TopTaq</i> DNA Polymerase | GC/AT-rich templates, complex templates, long PCR, TA cloning                            |
| <i>TransStart<sup>®</sup> Taq</i> DNA Polymerase<br><i>TransStart<sup>®</sup> TopTaq</i> DNA Polymerase  | GC/AT-rich templates, complex templates, qPCR, multiplex PCR, TA cloning                 |
| <i>EasyPfu</i> DNA Polymerase<br><i>TransStart<sup>®</sup> KD Plus</i> DNA Polymerase  | high fidelity PCR, blunt cloning; site-directed mutagenesis                              |
| <i>TransStart<sup>®</sup> FastPfu</i> DNA Polymerase<br><i>TransStart<sup>®</sup> FastPfu Fly</i> DNA Polymerase   | high fidelity PCR, fast PCR, complex templates, blunt cloning, site-directed mutagenesis |



# TransFast<sup>®</sup> Taq DNA Polymerase

|                   |          |             |
|-------------------|----------|-------------|
| dNTPs-free        | AP101-01 | 500 units   |
|                   | AP101-02 | 6×500 units |
| dNTPs<br>(2.5 mM) | AP101-11 | 500 units   |
|                   | AP101-12 | 6×500 units |

## Concentration

5 units/μl

## Contents

- TransFast<sup>®</sup> Taq DNA Polymerase
- 10×TransFast<sup>®</sup> Taq Buffer  
(200 mM Tris-HCl pH 8.4; 100 mM KCl;  
100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; others)
- 6×DNA Loading Buffer

## Storage

at -20°C for two years

## Description

TransFast<sup>®</sup> Taq DNA Polymerase is an engineered version of Taq DNA Polymerase. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kDa. TransFast<sup>®</sup> Taq DNA Polymerase has 5'-3' DNA polymerase activity and 5'-3' exonuclease activity. It lacks 3'-5' exonuclease activity.

- Extension rate is about 6 kb/min.
- Template-independent "A" can be generated at the 3' end of the PCR product. PCR products can be directly cloned into pEASY<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 4 kb.

## Applications

- Routine PCR
- Fast PCR
- Colony PCR

## Unit Definition

One unit of TransFast<sup>®</sup>Taq DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

## Quality Control

TransFast<sup>®</sup> Taq DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransFast<sup>®</sup> Taq DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

## PROTOCOL

### Reaction Components

| Component                                 | Volume    | Final Concentration |
|---|-----------|---------------------|
| Template                                  | Variable  | as required         |
| Forward Primer (10 μM)                    | 1 μl      | 0.2 μM              |
| Reverse Primer (10 μM)                    | 1 μl      | 0.2 μM              |
| 10×TransFast <sup>®</sup> Taq Buffer      | 5 μl      | 1×                  |
| 2.5 mM dNTPs                              | 4 μl      | 0.2 mM              |
| TransFast <sup>®</sup> Taq DNA Polymerase | 0.5 -1 μl | 2.5-5 units         |
| ddH <sub>2</sub> O                        | Variable  | -                   |
| Total volume                              | 50 μl     | -                   |







**Thermal cycling conditions**

|         |          |                |
|---------|----------|----------------|
| 94°C    | 2-5 min  | } 30-35 cycles |
| 94°C    | 5 sec    |                |
| 50-60°C | 15 sec   |                |
| 72°C    | x sec    |                |
| 72°C    | 5-10 min |                |

**Target                      Extension time**

|        |           |
|--------|-----------|
| 0-2 kb | 10 sec/kb |
| 2-3 kb | 20 sec/kb |
| >3 kb  | 30 sec/kb |



# EasyTaq® DNA Polymerase

|                |          |                |
|----------------|----------|----------------|
| dNTPs-free     | AP111-01 | 500 units      |
|                | AP111-02 | 6×500 units    |
|                | AP111-03 | 4×2,500 units  |
|                | AP111-04 | 10×5,000 units |
| dNTPs (2.5 mM) | AP111-11 | 500 units      |
|                | AP111-12 | 6×500 units    |
|                | AP111-13 | 4×2,500 units  |

**Concentration**

5 units/μl

**Contents**

- EasyTaq® DNA Polymerase
- 10×EasyTaq® Buffer (200 mM Tris-HCl pH 8.3; 200 mM KCl; 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; others)
- 6×DNA Loading Buffer

**Storage**

at -20°C for two years

**Description**

EasyTaq® DNA Polymerase is purified from *E. coli* expressing a cloned DNA polymerase from *Thermus aquaticus*. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kDa. EasyTaq® DNA Polymerase has 5'-3' DNA polymerase activity and 5'-3' exonuclease activity. It lacks 3'-5' exonuclease activity. EasyTaq® DNA Polymerase is suitable for routine amplification. PCR products are not suitable for PAGE.

- Extension rate is about 1-2 kb/min.
- Template-independent "A" can be generated at the 3' end of the PCR product. PCR products can be directly cloned into pEASY®-T vectors.
- Amplification of genomic DNA fragment up to 4 kb.

**Applications**

- Routine PCR
- Colony PCR

**Unit Definition**

One unit of EasyTaq® DNA Polymerase incorporates 10 nmol of



deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

**Quality Control**

*EasyTaq*® DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of *EasyTaq*® DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

**PROTOCOL**


**Reaction Components**

| Component                       | Volume   | Final Concentration |
|---------------------------------|----------|---------------------|
| Template                        | Variable | as required         |
| Forward Primer (10 µM)          | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)          | 1 µl     | 0.2 µM              |
| 10× <i>EasyTaq</i> ® Buffer     | 5 µl     | 1×                  |
| 2.5 mM dNTPs                    | 4 µl     | 0.2 mM              |
| <i>EasyTaq</i> ® DNA Polymerase | 0.5-1 µl | 2.5-5 units         |
| ddH <sub>2</sub> O              | Variable | -                   |
| Total volume                    | 50 µl    | -                   |

**Thermal cycling conditions**

94°C            2-5 min  
 94°C            30 sec  
 50-60°C       30 sec  
 72°C           1-2 kb/min  
 72°C           5-10 min

30-35 cycles



M: 1Kb Plus DNA Ladder  
 1: CCRD 0.5 kb; 2: BDNF 0.8 kb;  
 3: Rhod 1.2 kb; 4: Rhod 2 kb;  
 5: Rhod 4.17 kb.  
 50 ng of Human Genomic DNA as templates



## EasyTaq<sup>®</sup> DNA Polymerase for PAGE

|                   |          |               |
|-------------------|----------|---------------|
| dNTPs-free        | AP112-01 | 2,500 units   |
|                   | AP112-02 | 4x2,500 units |
| dNTPs<br>(2.5 mM) | AP112-11 | 2,500 units   |
|                   | AP112-12 | 4x2,500 units |

### Concentration

5 units/μl

### Contents

- EasyTaq<sup>®</sup> DNA Polymerase for PAGE
- 10xEasyTaq<sup>®</sup> Buffer for PAGE  
(200 mM Tris-HCl pH 8.3; 200 mM KCl; 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; others)
- 6xDNA Loading Buffer

### Storage

at -20°C for two years

### Description

EasyTaq<sup>®</sup> DNA Polymerase for PAGE is purified from *E. coli* expressing a cloned DNA polymerase from *Thermus aquaticus*. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kDa. EasyTaq<sup>®</sup> DNA Polymerase for PAGE has 5'-3' DNA polymerase activity and 5'-3' exonuclease activity. It lacks 3'-5' exonuclease activity.

- Extension rate is about 1-2 kb/min.
- Unique buffer system compatible with PAGE.
- Template-independent "A" can be generated at the 3' end of the PCR product. PCR products can be directly cloned into pEASY<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 3 kb.

### Application

Short fragment PCR

### Unit Definition

One unit of EasyTaq<sup>®</sup> DNA Polymerase for PAGE incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

EasyTaq<sup>®</sup> DNA Polymerase for PAGE has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of EasyTaq<sup>®</sup> DNA Polymerase for PAGE has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

## PROTOCOL

### Reaction Components

| Component                                    | Volume   | Final Concentration |
|--|----------|---------------------|
| Template                                     | Variable | as required         |
| Forward Primer (10 μM)                       | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                       | 1 μl     | 0.2 μM              |
| 10xEasyTaq <sup>®</sup> Buffer for PAGE      | 5 μl     | 1x                  |
| 2.5 mM dNTPs                                 | 4 μl     | 0.2 mM              |
| EasyTaq <sup>®</sup> DNA Polymerase for PAGE | 0.5-1 μl | 2.5-5 units         |
| ddH <sub>2</sub> O                           | Variable | -                   |
| Total volume                                 | 50 μl    | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |





# TransTaq<sup>®</sup>-T DNA Polymerase

|                   |          |             |
|-------------------|----------|-------------|
|                   | AP122-01 | 250 units   |
| dNTPs-free        | AP122-02 | 500 units   |
|                   | AP122-03 | 6×500 units |
|                   | AP122-11 | 250 units   |
| dNTPs<br>(2.5 mM) | AP122-12 | 500 units   |
|                   | AP122-13 | 6×500 units |

## Concentration

5 units/μl

## Contents

- TransTaq<sup>®</sup>-T DNA Polymerase
- 10×TransTaq<sup>®</sup>-T Buffer  
(200 mM Tris-HCl pH 9.0; 100 mM KCl; 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; others)
- 6×DNA Loading Buffer

## Storage

at -20°C for two years

## Description

TransTaq<sup>®</sup>-T DNA Polymerase is a mixture of EasyTaq<sup>®</sup> DNA Polymerase with a proofreading 3'-5' exonuclease. The fidelity is equal to EasyPfu DNA Polymerase. The yield is equal to that from EasyTaq<sup>®</sup> DNA Polymerase. It is more suitable for high fidelity TA cloning.

- TransTaq<sup>®</sup>-T DNA Polymerase offers 18-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent "A" can be generated at the 3' end of the PCR product. PCR products can be directly cloned into pEASY<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 8 kb.

## Applications

- Complex templates
- TA cloning

## Unit Definition

One unit of TransTaq<sup>®</sup>-T DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

## Quality Control

TransTaq<sup>®</sup>-T DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransTaq<sup>®</sup>-T DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

## PROTOCOL


### Reaction Components

| Component                               | Volume   | Final Concentration |
|---|----------|---------------------|
| Template                                | Variable | as required         |
| Forward Primer (10 μM)                  | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                  | 1 μl     | 0.2 μM              |
| 10×TransTaq <sup>®</sup> -T Buffer      | 5 μl     | 1×                  |
| 2.5 mM dNTPs                            | 4 μl     | 0.2 mM              |
| TransTaq <sup>®</sup> -T DNA Polymerase | 0.5-1 μl | 2.5-5 units         |
| ddH <sub>2</sub> O                      | Variable | -                   |
| Total volume                            | 50 μl    | -                   |



### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



M: 1Kb Plus DNA Ladder  
 1: BDNF 0.8 kb; 2: Rhod 1.2 kb;  
 3:  $\beta$ -globin 1.3 kb; 4: Rhod 2.0 kb;  
 5:  $\beta$ -globin 3.0 kb; 6: Rhod 4.17 kb;  
 7: Factor IX 7.5 kb  
 50 ng of Human Genomic DNA as templates

## TransTaq<sup>®</sup> DNA Polymerase High Fidelity (HiFi)

|                   |          |             |
|-------------------|----------|-------------|
|                   | AP131-01 | 250 units   |
| dNTPs-free        | AP131-02 | 500 units   |
|                   | AP131-03 | 6×500 units |
|                   | AP131-11 | 250 units   |
| dNTPs<br>(2.5 mM) | AP131-12 | 500 units   |
|                   | AP131-13 | 6×500 units |

### Concentration

5 units/ $\mu$ l

### Contents

- *TransTaq*<sup>®</sup> HiFi DNA Polymerase
- 10×*TransTaq*<sup>®</sup> HiFi Buffer I, II  
(200 mM Tris-HCl pH 9.0; 100 mM KCl;  
100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>;  
10% glycerol; others)
- GC Enhancer
- 6×DNA Loading Buffer

### Storage

at -20°C for two years

### Description

*TransTaq*<sup>®</sup> DNA Polymerase High Fidelity (*TransTaq*<sup>®</sup> HiFi DNA Polymerase) contains *TransTaq*<sup>®</sup>-T DNA Polymerase and a proofreading 3'-5' exonuclease. *TransTaq*<sup>®</sup> HiFi DNA Polymerase provides higher specificity and higher amplification efficiency than *TransTaq*<sup>®</sup>-T DNA Polymerase. Two different buffers are provided in the kit. *TransTaq*<sup>®</sup> HiFi Buffer I is optimized for the amplification of genomic DNA and *TransTaq*<sup>®</sup> HiFi Buffer II is optimized for the amplification of  $\lambda$ DNA, cDNA or plasmid DNA.

- *TransTaq*<sup>®</sup> HiFi DNA Polymerase offers 18-fold fidelity as compared to *EasyTaq*<sup>®</sup> DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent "A" can be generated at the 3' end of the PCR product. PCR products can be directly cloned into *pEASY*<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 15 kb.

### Applications

- Complex templates
- GC/AT rich templates
- Long PCR
- High yield PCR

### Unit Definition

One unit of *TransTaq*<sup>®</sup> HiFi DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.



### Quality Control

*TransTaq*<sup>®</sup> HiFi DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of *TransTaq*<sup>®</sup> HiFi DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

## PROTOCOL


### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 μM)                            | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                            | 1 μl     | 0.2 μM              |
| 10× <i>TransTaq</i> <sup>®</sup> HiFi Buffer I/II | 5 μl     | 1×                  |
| 2.5 mM dNTPs                                      | 4 μl     | 0.2 mM              |
| <i>TransTaq</i> <sup>®</sup> HiFi DNA Polymerase  | 0.5-1 μl | 2.5-5 units         |
| ddH <sub>2</sub> O                                | Variable | -                   |
| Total volume                                      | 50 μl    | -                   |


### Thermal cycling conditions

94°C            2-5 min  
 94°C            30 sec  
 50-60°C       30 sec  
 72°C            1-2 kb/min  
 72°C            5-10 min

} 30-35 cycles



*TransTaq*<sup>®</sup> HiFi Buffer I  
 M1: 1Kb Plus DNA Ladder  
 M2: *Trans*15K DNA Marker  
 1: BDNF 0.8 kb;  
 2: Rhod 1.2 kb;  
 3: Rhod 2.0 kb;  
 4: β-globin 3.0 kb;  
 5: Rhod 4.17 kb;  
 6: Factor IX 7.5 kb;  
 7: Serum albumin 12.4 kb  
 50 ng of Human Genomic DNA as templates




*TransTaq*<sup>®</sup> HiFi Buffer II  
 M1: 1Kb Plus DNA Ladder  
 M2: *Trans*15K DNA Marker  
 1: REPA 1.8 kb; 2: NCBP 2.5 kb;  
 3: HDP 3.5 kb; 4: VIN 4.6 kb;  
 5: Pol 6.8 kb; 6: APC 8.5 kb;  
 7: Dynein 12.3 kb  
 100 ng of Human total RNA as templates






# TransStart<sup>®</sup> Hot Start (Double Blocking)

## PCR Preparation without Hot Start



## PCR Preparation with TransStart<sup>®</sup> Method



At room temperature, one binding protein binds to double-strand DNA template and another binding protein binds to primer. These unique formulations effectively neutralize the DNA polymerase activity at room temperature. Blocking proteins are released from templates and primers during the initial denaturation. This double blocking method has higher efficiency than antibody based, or chemically modified hot start PCR.



# TransStart<sup>®</sup> Taq DNA Polymerase

|                   |          |             |
|-------------------|----------|-------------|
| dNTPs-free        | AP141-01 | 250 units   |
|                   | AP141-02 | 500 units   |
|                   | AP141-03 | 6x500 units |
| dNTPs<br>(2.5 mM) | AP141-11 | 250 units   |
|                   | AP141-12 | 500 units   |
|                   | AP141-13 | 6x500 units |

## Concentration

2.5 units/μl

## Contents

- TransStart<sup>®</sup> Taq DNA Polymerase
- 10xTransStart<sup>®</sup> Taq Buffer  
(500 mM Tris-HCl pH 9.0; 200 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; 10% glycerol; others)
- GC Enhancer
- 6xDNA Loading Buffer

## Storage

at -20°C for two years

## Description

TransStart<sup>®</sup> Taq DNA Polymerase is a hot start Taq DNA polymerase containing Taq DNA polymerase and two proprietary DNA binding proteins. At room temperature, one binding protein binds to double-strand DNA template and another binding protein binds to primer. These unique formulations effectively neutralize the DNA polymerase activity at room temperature. Blocking proteins are released from templates and primers during the initial denaturation. This double blocking method has higher efficiency than antibody based, or chemically modified hot start PCR.

- TransStart<sup>®</sup> Taq DNA Polymerase offers 18-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be directly cloned into pEASY<sup>®</sup>-T vectors.
- Reduced nonspecific amplification and primer dimer formation.
- Different from Taq antibody, no risk of contamination from mammalian DNA.
- Different from chemical modification, long denaturing step is not needed.
- Amplification of genomic DNA fragment up to 15 kb.

## Applications

- Complex templates
- GC/AT-rich templates
- Multiplex PCR
- High yield PCR

## Unit Definition

One unit of TransStart<sup>®</sup> Taq DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

## Quality Control

TransStart<sup>®</sup> Taq DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransStart<sup>®</sup> Taq DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.



## PROTOCOL


### Reaction Components

| Component  | Volume        | Final Concentration |
|--|---------------|---------------------|
| Template   | Variable      | as required         |
| Forward Primer (10 $\mu$ M)                              | 1 $\mu$ l     | 0.2 $\mu$ M         |
| Reverse Primer (10 $\mu$ M)                              | 1 $\mu$ l     | 0.2 $\mu$ M         |
| 10x <i>TransStart</i> <sup>®</sup> <i>Taq</i> Buffer     | 5 $\mu$ l     | 1x                  |
| 2.5 mM dNTPs   | 4 $\mu$ l     | 0.2 mM              |
| <i>TransStart</i> <sup>®</sup> <i>Taq</i> DNA Polymerase | 0.5-1 $\mu$ l | 1.25-2.5 units      |
| ddH <sub>2</sub> O                                       | Variable      | -                   |
| Total volume   | 50 $\mu$ l    | -                   |

### Thermal cycling conditions

94°C 2-5 min  
 94°C 30 sec  
 50-60°C 30 sec  
 72°C 1-2 kb/min  
 72°C 5-10 min

30-35 cycles



M1: 1Kb Plus DNA Ladder  
 M2: *Trans*15K DNA Marker  
 1: Numb 0.3 kb; 2: CCRD 0.5 kb;  
 3: BDNF 0.8 kb; 4: Rhod 1.2 kb;  
 5: Rhod 2.0 kb; 6:  $\beta$ -globin 3.0 kb;  
 7: Rhod 4.17 kb; 8: Factor IX 7.5 kb;  
 9: Serum albumin 12.4 kb  
 50 ng of Human Genomic DNA as templates





# TransStart<sup>®</sup> TopTaq DNA Polymerase

|                   |          |             |
|-------------------|----------|-------------|
| dNTPs-free        | AP151-01 | 250 units   |
|                   | AP151-02 | 500 units   |
|                   | AP151-03 | 6×500 units |
| dNTPs<br>(2.5 mM) | AP151-11 | 250 units   |
|                   | AP151-12 | 500 units   |
|                   | AP151-13 | 6×500 units |

## Concentration

2.5 units/μl

## Contents

- TransStart<sup>®</sup> TopTaq DNA Polymerase
- 10×TransStart<sup>®</sup> TopTaq Buffer  
(500 mM Tris-HCl (pH 9.0); 200 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 20 mM MgSO<sub>4</sub>; 10% glycerol others)
- GC Enhancer
- 6×DNA Loading Buffer

## Storage

at -20°C for two years

## Description

TransStart<sup>®</sup> TopTaq DNA Polymerase is an engineered version of Taq DNA Polymerase combined with TransStart<sup>®</sup> technique. One binding protein binds to double-strand DNA template, preventing polymerase activity at room temperature. Other two binding proteins bind primers, preventing primer-dimer formation. Blocking proteins are released from primers and templates during the initial denaturation. This double blocking method has higher efficiency than antibody based, or chemically modified hot start PCR.

- Compared with TransStart<sup>®</sup> Taq DNA Polymerase, TransStart<sup>®</sup> TopTaq DNA Polymerase has higher amplification efficiency, specificity and sensitivity.
- TransStart<sup>®</sup> TopTaq DNA Polymerase offers 18-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.
- The specificity is higher than antibody based or chemically modified hot start DNA polymerases.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be directly cloned into pEASY<sup>®</sup>-T vectors.
- Reduced nonspecific amplification and primer dimer formation.
- Different from Taq antibody, no risk of contamination from mammalian DNA.
- Different from chemical modification, long denaturing step is not needed.
- Amplification of genomic DNA fragment up to 15 kb.

## Applications

- Complex templates
- GC/AT-rich templates
- Multiplex PCR
- High yield PCR

## Unit Definition

One unit of TransStart<sup>®</sup> TopTaq DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

## Quality Control

TransStart<sup>®</sup> TopTaq DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransStart<sup>®</sup> TopTaq DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.



## PROTOCOL


### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 µM)                                      | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                                      | 1 µl     | 0.2 µM              |
| 10x <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> Buffer     | 5 µl     | 1x                  |
| 2.5 mM dNTPs  | 4 µl     | 0.2 mM              |
| <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> DNA Polymerase | 0.5-1 µl | 1.25-2.5 units      |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total volume  | 50 µl    | -                   |


### Thermal cycling conditions

94°C            2-5 min  
 94°C            30 sec  
 50-60°C       30 sec  
 72°C           1-2 kb/min  
 72°C           5-10 min

} 30-35 cycles



M: *Trans2K*<sup>®</sup> Plus II DNA Marker  
 1: Rhod 1.2 kb    2: Rhod 2.0 kb    3: β-globin 3.0 kb  
 4: β-globin 4.1 kb    5: Rhod 4.17 kb    6: β-globin 6.1 kb  
 7: Factor IX 7.5 kb    8: IGF2R 8.9 kb    9: Serum albumin 12.4 kb  
 50 ng Human Genomic DNA as templates



M: *Trans2K*<sup>®</sup> Plus DNA Marker  
 1: DMD1 0.3 kb    2: Numb 0.3 kb  
 3: P53 0.5 kb    4: P1P2 1.2 kb  
 A: Competitor A Hot Start DNA Polymerase  
 T: *TransStart*<sup>®</sup> *TopTaq* DNA Polymerase



## EasyPfu DNA Polymerase

|            |          |             |                   |          |             |
|------------|----------|-------------|-------------------|----------|-------------|
| dNTPs-free | AP211-01 | 250 units   | dNTPs<br>(2.5 mM) | AP211-11 | 250 units   |
|            | AP211-02 | 500 units   |                   | AP211-12 | 500 units   |
|            | AP211-03 | 6x500 units |                   | AP211-13 | 6x500 units |

### Concentration

2.5 units/μl

### Contents

- EasyPfu DNA Polymerase
- 10xEasyPfu Buffer  
(200 mM Tris-HCl pH 8.8; 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 100 mM KCl; 20 mM MgSO<sub>4</sub>; others)
- 50 mM MgSO<sub>4</sub>
- 6xDNA Loading Buffer

### Storage

at -20°C for two years

### Description

EasyPfu DNA Polymerase is an engineered version of *pfu* DNA Polymerase with enhanced yield and higher fidelity. EasyPfu DNA Polymerase possesses a proofreading 3'-5' exonuclease activity.

- EasyPfu DNA Polymerase offers 18-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.
- Extension rate is about 0.5 kb/min.
- PCR products can be directly cloned into pEASY<sup>®</sup>-Blunt vectors.
- Amplification of genomic DNA fragment up to 6 kb.
- Amplification of plasmid DNA fragment up to 10 kb.

### Applications

- High fidelity PCR
- Blunt-end cloning
- Site-directed mutagenesis

### Unit Definition

One unit of EasyPfu DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

EasyPfu DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity, >99% homogeneous measured by SDS-PAGE. Each batch of EasyPfu DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

## PROTOCOL


### Reaction Components

| Component              | Volume   | Final Concentration |
|------------------------|----------|---------------------|
| Template               | Variable | as required         |
| Forward Primer (10 μM) | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM) | 1 μl     | 0.2 μM              |
| 10xEasyPfu Buffer      | 5 μl     | 1x                  |
| 2.5 mM dNTPs           | 4 μl     | 0.2 mM              |
| EasyPfu DNA Polymerase | 1 μl     | 2.5 units           |
| ddH <sub>2</sub> O     | Variable | -                   |
| Total volume           | 50 μl    | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 0.5 kb/min |                |
| 72°C    | 5-10 min   |                |





## TransStart<sup>®</sup> FastPfu DNA Polymerase

|                   |          |             |
|-------------------|----------|-------------|
|                   | AP221-01 | 250 units   |
| dNTPs-free        | AP221-02 | 500 units   |
|                   | AP221-03 | 6×500 units |
|                   | AP221-11 | 250 units   |
| dNTPs<br>(2.5 mM) | AP221-12 | 500 units   |
|                   | AP221-13 | 6×500 units |

### Concentration

2.5 units/ $\mu$ l

### Contents

- TransStart<sup>®</sup> FastPfu DNA Polymerase
- 5×TransStart<sup>®</sup> FastPfu Buffer  
(100 mM Tris-SO<sub>4</sub> pH 9.2; 50 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 200 mM KCl; 10 mM MgSO<sub>4</sub>; 10% glycerol; others)
- 50 mM MgSO<sub>4</sub>
- PCR Stimulant
- 6×DNA Loading Buffer

### Storage

at -20°C for two years

### Description

TransStart<sup>®</sup> FastPfu DNA Polymerase is a fast, high fidelity and high processivity hot start DNA polymerase.

- Extension rate is about 2-4 kb/min.
- TransStart<sup>®</sup> FastPfu DNA Polymerase offers 54-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.
- PCR products can be directly cloned into pEASY<sup>®</sup>-Blunt vectors.
- Amplification of genomic DNA fragment up to 15 kb.
- Amplification of plasmid DNA fragment up to 20 kb.

### Applications

- High fidelity PCR
- High yield and fast PCR
- Blunt end cloning
- Site-directed mutagenesis
- Complex templates

### Unit Definition

One unit of TransStart<sup>®</sup> FastPfu DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

TransStart<sup>®</sup> FastPfu DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransStart<sup>®</sup> FastPfu DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.



## PROTOCOL

### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| Template   | Variable | as required         |
| Forward Primer (10 µM)                                       | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                                       | 1 µl     | 0.2 µM              |
| 5× <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> Buffer      | 10 µl    | 1×                  |
| 2.5 mM dNTPs   | 4 µl     | 0.2 mM              |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> DNA Polymerase | 1 µl     | 2.5 units           |
| ddH <sub>2</sub> O   | Variable | -                   |
| Total volume   | 50 µl    | -                   |

### Suggested conditions


(50 µl reaction volume)

| Parameter         | Targets ≤10 kb  | Targets ≥10 kb                                | cDNA  |
|-------------------|---|---|---|
| Template          | 100 ng Genomic DNA<br>5-30 ng Plasmid DNA   | 200-500 ng Genomic DNA<br>5-30 ng Plasmid DNA | 1-2 µl cDNA from RT reaction<br>(50-500 ng RNA for RT reaction) |
| MgSO <sub>4</sub> | Add extra 1-2 µl of 50 mM MgSO <sub>4</sub> to a final concentration of 3-4 mM if the amplified product is larger than 5 kb |   |   |

### Thermal cycling conditions

| Number of cycles   | Temperature | Plasmid or Genomic DNA                                     | cDNA     |
|--|-------------|--|----------|
| 1 cycle  | 95°C        | 2 min  | 1 min    |
| Plasmid or Genomic DNA: 30-35 cycles<br>cDNA: 35-40 cycles | 95°C        | 20 sec   | 20 sec   |
|  | Tm-5°C      | 20 sec   | 20 sec   |
|  | 72°C        | 4 kb/min for targets ≤1 kb<br>2-4 kb/min for targets >1 kb | 2 kb/min |
| 1 cycle  | 72°C        | 5 min  | 5 min    |

### Amplification from cDNA templates with *TransStart*<sup>®</sup> *FastPfu* DNA Polymerase




M1: 1Kb Plus DNA Ladder  
M2: *Trans*15K DNA Marker  
1: ACTR 3.5 kb; 2 hrs 14 min  
2: VIN 4.6 kb; 2 hrs 34 min  
3: Pol 6.8 kb; 3 hrs 09 min  
4: APC 8.5 kb; 3 hrs 41 min  
5: Dynein 12.3 kb; 4 hrs 48 min

High quality products




**Amplification from genomic DNA templates with *TransStart® FastPfu* DNA Polymerase**



M1: 1Kb Plus DNA Ladder  
 M2: *Trans*15K DNA Marker  
 1: Rhod 2.0 kb; 1 hrs 19 min  
 2:  $\beta$ -globin 3.0 kb; 1 hrs 27 min  
 3: Rhod 4.17 kb; 1 hrs 29 min  
 4: Factor IX 7.5 kb; 3 hrs 25 min  
 5: Serum albumin 12.4 kb; 4 hrs 48 min

**Amplification from plasmid DNA templates with *TransStart® FastPfu* DNA Polymerase**



M: *Trans*15K DNA Marker  
 1: UDG 7.0 kb; 1 hrs 36 min  
 2: LN 10.0 kb; 1 hrs 55 min  
 3: Fang 14.7 kb; 2 hrs 26 min

***TransStart® FastPfu* Fly DNA Polymerase**

|                   |          |             |
|-------------------|----------|-------------|
|                   | AP231-01 | 250 units   |
| dNTPs-free        | AP231-02 | 500 units   |
|                   | AP231-03 | 6×500 units |
|                   | AP231-11 | 250 units   |
| dNTPs<br>(2.5 mM) | AP231-12 | 500 units   |
|                   | AP231-13 | 6×500 units |

**Concentration**

2.5 units/ $\mu$ l

**Contents**

- *TransStart® FastPfu* Fly DNA Polymerase
- 5×*TransStart® FastPfu* Fly Buffer (100 mM Tris-SO<sub>4</sub> pH 9.2; 50 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 200 mM KCl; 10 mM MgSO<sub>4</sub>; 10% glycerol; others)
- 50 mM MgSO<sub>4</sub>
- PCR Stimulant
- 6×DNA Loading Buffer

**Storage**

at -20°C for two years


**Description**

*TransStart® FastPfu* Fly DNA Polymerase is a hot start, high fidelity and high processivity DNA Polymerase. *TransStart® FastPfu* Fly DNA Polymerase has an extension rate of up to 6 kb/min. Compared with *TransStart® FastPfu* DNA Polymerase, *TransStart® FastPfu* Fly DNA Polymerase has higher extension rate, higher fidelity, and higher amplification efficiency.

- *TransStart® FastPfu* Fly DNA Polymerase offers 108-fold fidelity as compared to *EasyTaq®* DNA Polymerase.
- Extension rate is about 2-6 kb/min.
- PCR products can be directly cloned into *pEASY®*-Blunt vectors.
- Amplification of genomic DNA fragment up to 15 kb.
- Amplification of plasmid DNA fragment up to 20 kb.

**Applications**

- High fidelity PCR
- High yield and fast PCR
- Blunt end cloning
- Site-directed mutagenesis
- Complex templates



### Unit Definition

One unit of *TransStart<sup>®</sup> FastPfu* Fly DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

*TransStart<sup>®</sup> FastPfu* Fly DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of *TransStart<sup>®</sup> FastPfu* Fly DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| Template   | Variable | as required         |
| Forward Primer (10 µM)                                   | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                                   | 1 µl     | 0.2 µM              |
| 5× <i>TransStart<sup>®</sup> FastPfu</i> Fly Buffer      | 10 µl    | 1×                  |
| 2.5 mM dNTPs   | 4 µl     | 0.2 mM              |
| <i>TransStart<sup>®</sup> FastPfu</i> Fly DNA Polymerase | 1 µl     | 2.5 units           |
| ddH <sub>2</sub> O                                       | Variable | -                   |
| Total volume   | 50 µl    | -                   |

### Suggested conditions

(50 µl reaction volume)

| Parameter         | Targets ≤10 kb  | Targets ≥10 kb                                | cDNA  |
|-------------------|---|---|---|
| Template          | 100 ng Genomic DNA<br>5-30 ng Plasmid DNA   | 200-500 ng Genomic DNA<br>5-30 ng Plasmid DNA | 1-2 µl cDNA from RT reaction<br>(50-500 ng RNA for RT reaction) |
| MgSO <sub>4</sub> | Add extra 1-2 µl of 50 mM MgSO <sub>4</sub> to a final concentration of 3-4 mM if the amplified product is larger than 5 kb |   |   |

### Thermal cycling conditions

| Number of cycles   | Temperature         | cDNA or Genomic DNA  | Plasmid DNA  |
|--|---------------------|--|--|
| 1 cycle  | 95°C                | 2 min  | 2 min  |
| Plasmid or Genomic DNA: 30-35 cycles<br>cDNA: 35-40 cycles | 95°C                | 20 sec   | 20 sec   |
|  | T <sub>m</sub> -5°C | 20 sec   | 20 sec   |
|  | 72°C                | 6 kb/min for targets ≤2 kb<br>2-4 kb/min for targets >2 kb | 6 kb/min for targets ≤6 kb<br>2-4 kb/min for targets >6 kb |
| 1 cycle  | 72°C                | 5 min  | 5 min  |

High quality products





## Total PCR Time



4 kb: Genomic DNA; 7 kb and 10 kb: Plasmid DNA

## TransStart<sup>®</sup> KD Plus DNA Polymerase

|            |          |             |
|------------|----------|-------------|
|            | AP301-01 | 100 units   |
| dNTPs-free | AP301-02 | 200 units   |
|            | AP301-03 | 6x200 units |
| dNTPs      | AP301-11 | 100 units   |
| (2.5 mM)   | AP301-12 | 200 units   |
|            | AP301-13 | 6x200 units |

### Concentration

1 unit/μl

### Contents

- TransStart<sup>®</sup> KD Plus DNA Polymerase
- 5xTransStart<sup>®</sup> KD Plus Buffer (100 mM Tris-HCl pH 9.2; 50 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>; 200 mM KCl; 5 mM MgSO<sub>4</sub>; 10% Glycerol; others)
- 50 mM MgSO<sub>4</sub>
- 6xDNA Loading Buffer

### Storage

at -20°C for two years

### Description

TransStart<sup>®</sup> KD Plus DNA Polymerase is a genetically modified high fidelity DNA polymerase. This enzyme provides higher amplification capability than traditional *Pfu* DNA polymerase and fast amplification speed equal to *Taq* DNA polymerase (1 kb/min). Due to strong 3'-5' exonuclease activity, this enzyme offers 108-fold fidelity as compared to EasyTaq<sup>®</sup> DNA Polymerase.

- PCR products can be directly cloned into pEASY<sup>®</sup>-Blunt vectors.
- Amplification of genomic DNA fragment up to 15 kb.
- Amplification of plasmid DNA fragment up to 20 kb.

### Applications

- Fast, high specificity amplification
- High fidelity, high yield amplification

### Unit Definition

One unit of TransStart<sup>®</sup> KD Plus DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

TransStart<sup>®</sup> KD Plus DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity, >99% homogeneous measured by SDS-PAGE. Each batch of TransStart<sup>®</sup> KD Plus DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.




## PROTOCOL

### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| Template   | Variable | as required         |
| Forward Primer (10 µM)                                       | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                                       | 1 µl     | 0.2 µM              |
| 5x <i>TransStart</i> <sup>®</sup> <i>KD</i> Plus Buffer      | 10 µl    | 1x                  |
| 2.5 mM dNTPs   | 4 µl     | 0.2 mM              |
| <i>TransStart</i> <sup>®</sup> <i>KD</i> Plus DNA Polymerase | 1 µl     | 1 unit              |
| ddH <sub>2</sub> O   | Variable | -                   |
| Total volume   | 50 µl    | -                   |

### Thermal cycling conditions

|         |          |                |
|---------|----------|----------------|
| 94°C    | 2-5 min  | } 30-35 cycles |
| 94°C    | 30 sec   |                |
| 50-60°C | 30 sec   |                |
| 68°C    | 1 kb/min |                |
| 68°C    | 5-10 min |                |



M: 1Kb Plus DNA Ladder  
 Human cDNA as templates      Human Genomic DNA as templates      Plasmid DNA as templates  
 1: GAPDH 0.9 kb;                      5: Rhod 4.17 kb;                      8: Fang 14.7 kb  
 2: REPA 1.8 kb;                      6: β-globin 6.1 kb;  
 3: NCBP 2.5 kb;                      7: Factor IX 7.5 kb;  
 4: ACTR 3 kb;



## GC Enhancer

AG101-01

200  $\mu$ l

### Storage

at -20°C for two years

### Description

GC Enhancer can be used to increase sensitivity and specificity for GC/AT-rich template or complex template. The stock concentration is 10 $\times$ , and the working concentration can be varied between 0.5 $\times$  to 5 $\times$ .

### Applications

- Complex templates
- GC/AT-rich templates  
(50  $\mu$ l reaction volume)

| Volume of GC Enhancer ( $\mu$ l) | Final Concentration |
|----------------------------------|---------------------|
| 2.5                              | 0.5 $\times$        |
| 5                                | 1 $\times$          |
| 10                               | 2 $\times$          |
| 15                               | 3 $\times$          |
| 20                               | 4 $\times$          |
| 25                               | 5 $\times$          |

## PCR Stimulant

AG111-01

200  $\mu$ l

### Storage

at -20°C for two years

### Description

PCR Stimulant can be used to increase sensitivity and specificity for GC/AT-rich template or complex template. It is especially suitable for *Pfu* enzymes. The stock concentration is 5 $\times$ , and the working concentration can be varied between 0.5 $\times$  to 2.5 $\times$ .

### Applications

- Complex templates
- GC/AT-rich templates  
(50  $\mu$ l reaction volume)

| Volume of PCR Stimulant ( $\mu$ l) | Final Concentration |
|------------------------------------|---------------------|
| 5                                  | 0.5 $\times$        |
| 10                                 | 1 $\times$          |
| 20                                 | 2 $\times$          |
| 25                                 | 2.5 $\times$        |




# PCR SuperMix

PCR SuperMix is a ready-to-use mixture of DNA polymerase, salt, magnesium, dNTPs and other components for efficient PCR amplification. Only add template, primers and ddH<sub>2</sub>O to the SuperMix for PCR. If PCR SuperMix with dye is used, dye needs to be removed before cloning or sequencing.


## Applications

| PCR SuperMix   | Application  |
|--|--|
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix                          | routine PCR  |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix for PAGE                 | short fragment PCR   |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix                      | complex templates, TA cloning  |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix I  | GC/AT-rich templates, complex templates, long PCR, genomic DNA amplification (<15 kb), TA cloning    |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix II | GC/AT-rich templates, complex templates, long PCR, λDNA, cDNA, plasmid DNA amplification, TA cloning |
| 2× <i>EasyPfu</i> PCR SuperMix                                       | high fidelity PCR, blunt cloning, site-directed mutagenesis  |
| 2× <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix        | high fidelity PCR, fast PCR, blunt cloning, site-directed mutagenesis                                |

PCR Enzyme vs. PCR SuperMix



PCR SuperMix Selection Chart



High quality products





## 2×EasyTaq<sup>®</sup> PCR SuperMix

|            |          |         |
|------------|----------|---------|
|            | AS111-01 | 1 ml    |
| Mix (-dye) | AS111-02 | 5×1 ml  |
|            | AS111-03 | 15×1 ml |
|            | AS111-11 | 1 ml    |
| Mix (+dye) | AS111-12 | 5×1 ml  |
|            | AS111-13 | 15×1 ml |
|            | AS111-14 | 6×80 ml |

### Storage

at -20°C for two years

### Description

*EasyTaq*<sup>®</sup> PCR SuperMix is a ready-to-use mixture of *EasyTaq*<sup>®</sup> DNA Polymerase, dNTPs and optimized buffer. The SuperMix is provided at 2× concentration and used at 1× concentration by adding template, primers and H<sub>2</sub>O. PCR products are not suitable for PAGE.

- Extension rate is about 1-2 kb/min.
- Template-independent “A” can be generated at the 3' end of the PCR product. PCR products can be cloned into *pEASY*<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 4 kb.

### Application

Routine PCR

## PROTOCOL


### Reaction Components

| Component                                   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template                                    | Variable | as required         |
| Forward Primer (10 μM)                      | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                      | 1 μl     | 0.2 μM              |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix | 25 μl    | 1×                  |
| ddH <sub>2</sub> O                          | Variable | -                   |
| Total volume                                | 50 μl    | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |

M 1 2 3 4 5 M 1 2 3 4 5 M



M: 1Kb Plus DNA Ladder  
 Lane 1: CCRD 0.5 kb  
 Lane 2: BDNF 0.8 kb  
 Lane 3: Rhod 1.2 kb  
 Lane 4: Rhod 2 kb  
 Lane 5: Rhod 4.17 kb  
 50 ng of Human Genomic DNA as templates

2×*EasyTaq*<sup>®</sup> PCR SuperMix (+dye)

2×*EasyTaq*<sup>®</sup> PCR SuperMix (-dye)



## 2×EasyTaq<sup>®</sup> PCR SuperMix for PAGE

|            |          |         |
|------------|----------|---------|
|            | AS112-11 | 1 ml    |
| Mix (+dye) | AS112-12 | 5×1 ml  |
|            | AS112-13 | 15×1 ml |

### Storage

at -20°C for two years

### Description

*EasyTaq*<sup>®</sup> PCR SuperMix for PAGE is a ready-to-use mixture of *EasyTaq*<sup>®</sup> DNA Polymerase for PAGE, dNTPs and optimized buffer. The SuperMix for PAGE is provided at 2× concentration and used at 1× concentration by adding template, primers and H<sub>2</sub>O.

- Extension rate is about 1-2 kb/min.
- Unique buffer system compatible with PAGE.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be cloned into *pEASY*<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 3 kb.

### Application

Short fragment PCR

## PROTOCOL

### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| Template   | Variable | as required         |
| Forward Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix for PAGE | 25 μl    | 1×                  |
| ddH <sub>2</sub> O                                   | Variable | -                   |
| Total volume   | 50 μl    | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



## 2×*TransTaq*<sup>®</sup>-T PCR SuperMix

|            |          |        |
|------------|----------|--------|
| Mix (-dye) | AS122-01 | 1 ml   |
|            | AS122-02 | 5×1 ml |
| Mix (+dye) | AS122-11 | 1 ml   |
|            | AS122-12 | 5×1 ml |

### Storage

at -20°C for two years

### Description

*TransTaq*<sup>®</sup>-T PCR SuperMix is a ready-to-use mixture of *TransTaq*<sup>®</sup>-T DNA Polymerase, dNTPs and optimized buffer. The SuperMix is provided at 2× concentration and used at 1× concentration by adding template, primers and H<sub>2</sub>O. Efficiency of PCR products with “A” is equal to *EasyTaq*<sup>®</sup> DNA polymerase. It is more suitable for high fidelity TA cloning.

- *TransTaq*<sup>®</sup>-T PCR SuperMix offers 18-fold fidelity as compared to *EasyTaq*<sup>®</sup> DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent “A” can be generated at the 3' end of the PCR product. PCR products can be cloned into *pEASY*<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 8 kb.

### Applications

- Complex templates
- TA cloning

## PROTOCOL


### Reaction Components

| Component                                       | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 μM)                          | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                          | 1 μl     | 0.2 μM              |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix | 25 μl    | 1×                  |
| ddH <sub>2</sub> O                              | Variable | -                   |
| Total volume                                    | 50 μl    | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |

M 1 2 3 4 5 6 7 M 1 2 3 4 5 6 7 M



M: 1Kb Plus DNA Ladder  
 Lane 1: BDNF 0.8 kb; Lane 2: Rhod 1.2 kb;  
 Lane 3: β-globin 1.3 kb; Lane 4: Rhod 2.0 kb;  
 Lane 5: β-globin 3.0 kb; Lane 6: Rhod 4.17 kb;  
 Lane 7: Factor IX 7.5 kb  
 50 ng of Human Genomic DNA as templates

*TransTaq*<sup>®</sup>-T DNA Polymerase      2×*TransTaq*<sup>®</sup>-T PCR SuperMix



## 2×*TransTaq*<sup>®</sup> High Fidelity (HiFi) PCR SuperMix

|               |          |        |
|---------------|----------|--------|
| Mix I (-dye)  | AS131-01 | 1 ml   |
|               | AS131-02 | 5×1 ml |
| Mix II (-dye) | AS131-21 | 1 ml   |
|               | AS131-22 | 5×1 ml |

### Storage

at -20°C for two years

### Description

*TransTaq*<sup>®</sup> High Fidelity (HiFi) PCR SuperMix I or II is a ready-to-use mixture of *TransTaq*<sup>®</sup> High Fidelity (HiFi) DNA polymerase, dNTPs and optimized buffer. *TransTaq*<sup>®</sup> High Fidelity (HiFi) PCR SuperMix I is optimized for the amplification of genomic DNA and PCR SuperMix II is optimized for the amplification of λDNA, cDNA or plasmid DNA. The SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primers and H<sub>2</sub>O.

- *TransTaq*<sup>®</sup> High Fidelity (HiFi) PCR SuperMix offers 18-fold fidelity as compared to *EasyTaq*<sup>®</sup> DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be directly cloned into *pEASY*<sup>®</sup>-T vectors.
- Amplification of genomic DNA fragment up to 15 kb.

### Applications

- Complex templates
- GC/AT-rich templates
- Long PCR
- High yield PCR

## PROTOCOL


### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 μM)                            | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                            | 1 μl     | 0.2 μM              |
| 2× <i>TransTaq</i> <sup>®</sup> HiFi PCR SuperMix | 25 μl    | 1×                  |
| ddH <sub>2</sub> O                                | Variable | -                   |
| Total volume                                      | 50 μl    | -                   |

### Thermal cycling conditions


|         |            |                |
|---------|------------|----------------|
| 94°C    | 2-5 min    | } 30-35 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |





TransTaq<sup>®</sup> HIFI DNA Polymerase      2xTransTaq<sup>®</sup> HIFI PCR SuperMix I

M1: 1Kb Plus DNA Ladder  
 M2: Trans15K DNA Marker  
 Lane 1: BDNF 0.8 kb;      Lane 5: Rhod 4.17 kb;  
 Lane 2: Rhod 1.2 kb;      Lane 6: Factor IX 7.5 kb;  
 Lane 3: Rhod 2.0 kb;      Lane 7: Serum albumin 12.4 kb  
 Lane 4:  $\beta$ -globin 3.0 kb;  
 50 ng of Human Genomic DNA as templates



TransTaq<sup>®</sup> HIFI DNA Polymerase      2xTransTaq<sup>®</sup> HIFI PCR SuperMix II

M1: 1Kb Plus DNA Ladder  
 M2: Trans15K DNA Marker  
 Lane 1: REPA 1.8 kb;      Lane 5: Pol 6.8 kb;  
 Lane 2: NCBP 2.5 kb;      Lane 6: APC 8.5 kb;  
 Lane 3: HDP 3.5 kb;      Lane 7: Dynein 12.3 kb  
 Lane 4: VIN 4.6 kb;  
 Human cDNA as templates



# 2x*EasyPfu* PCR SuperMix

|            |          |        |
|------------|----------|--------|
| Mix (-dye) | AS211-01 | 1 ml   |
|            | AS211-02 | 5x1 ml |

### Storage

at -20°C for two years

### Description

*EasyPfu* PCR SuperMix is a ready-to-use mixture of *EasyPfu* DNA Polymerase, dNTPs and optimized buffer. The SuperMix is provided at 2x concentration and used at 1x concentration by adding template, primers and H<sub>2</sub>O.

- *EasyPfu* PCR SuperMix offers 18-fold fidelity as compared to *EasyTaq*<sup>®</sup> DNA Polymerase.
- Extension rate is about 0.5 kb/min.
- PCR products can be directly cloned into *pEASY*<sup>®</sup>-Blunt vectors.
- Amplification of genomic DNA fragment up to 6 kb.
- Amplification of plasmid DNA fragment up to 10 kb.

### Applications

- High fidelity PCR
- Blunt end cloning
- Site-directed mutagenesis

## PROTOCOL


### Reaction Components

| Component                      | Volume   | Final Concentration |
|--------------------------------|----------|---------------------|
| Template                       | Variable | as required         |
| Forward Primer (10 μM)         | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)         | 1 μl     | 0.2 μM              |
| 2x <i>EasyPfu</i> PCR SuperMix | 25 μl    | 1x                  |
| ddH <sub>2</sub> O             | Variable | -                   |
| Total volume                   | 50 μl    | -                   |

### Thermal cycling conditions

94°C 2-5 min  
 94°C 30 sec  
 50-60°C 30 sec  
 72°C 0.5 kb/min  
 72°C 5-10 min

30-35 cycles



*EasyPfu* DNA Polymerase                      2x*EasyPfu* PCR SuperMix

M: 1Kb Plus DNA Ladder  
 Lane 1: CCRD 0.5 kb;                      Lane 5: Rhod 2 kb;  
 Lane 2: BDNF 0.8 kb;                      Lane 6: β-globin 3 kb;  
 Lane 3: Rhod 1.2 kb;                      Lane 7: Rhod 4.17 kb;  
 Lane 4: β-globin 1.3 kb;                      Lane 8: β-globin 4.1 kb  
 50 ng of Human Genomic DNA as templates

High quality products



## 2×*TransStart*<sup>®</sup> *FastPfu* PCR SuperMix

|            |          |        |
|------------|----------|--------|
| Mix (-dye) | AS221-01 | 1 ml   |
|            | AS221-02 | 5×1 ml |

### Storage

at -20°C for two years

### Description

*TransStart*<sup>®</sup> *FastPfu* PCR SuperMix is a ready-to-use mixture of *TransStart*<sup>®</sup> *FastPfu* DNA polymerase, dNTPs, and optimized buffer. The SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primers and H<sub>2</sub>O.

- *TransStart*<sup>®</sup> *FastPfu* PCR SuperMix offers 54-fold fidelity as compared to *EasyTaq*<sup>®</sup> DNA Polymerase.
- Extension rate is about 2-4 kb/min.
- PCR products can be directly cloned into *pEASY*<sup>®</sup>-Blunt vectors.
- Amplification of genomic DNA fragment up to 15 kb.
- Amplification of plasmid DNA fragment up to 20 kb.

### Applications

- High fidelity PCR
- High yield PCR
- Fast PCR
- Blunt end cloning
- Site-directed mutagenesis
- Complex templates

## PROTOCOL


### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 μM)  | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)  | 1 μl     | 0.2 μM              |
| 2× <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix | 25 μl    | 1×                  |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total volume  | 50 μl    | -                   |

### Thermal cycling conditions

94°C 2-5 min  
 94°C 20 sec  
 50-60°C 20 sec  
 72°C 2-4 kb/min  
 72°C 5-10 min

30-35 cycles



*TransStart*<sup>®</sup> *FastPfu* DNA Polymerase      2×*TransStart*<sup>®</sup> *FastPfu* PCR SuperMix

M: 1Kb Plus DNA Ladder  
 Lane 1: β-globin 1.3 kb Human genomic DNA  
 Lane 2: Rhod 2.0 kb Human genomic DNA  
 Lane 3: NCBP 2.5 kb Human cDNA  
 Lane 4: VIN 4.6 kb Human cDNA  
 Lane 5: Pol 6.8 kb Human cDNA  
 Lane 6: LN 10.0 kb Plasmid DNA



# Direct PCR

*TransDirect*<sup>®</sup> PCR uses a proprietary formulated lysis reagent to release nucleic acids from a variety of fresh or frozen animal cells/tissues and plant tissues. Unpurified DNA is used as template for PCR using *2xTransDirect*<sup>®</sup> PCR SuperMix which has extremely high resistance to PCR inhibitors found in animal tissues, plant tissues and blood.

| Direct PCR Kit  | Application  |
|---|--|
| <i>TransDirect</i> <sup>®</sup> Animal Tissue PCR Kit | Mammalian cell cultures, saliva, hair shaft, animal tissues, blood   |
| <i>TransDirect</i> <sup>®</sup> Plant Tissue PCR Kit  | Low polysaccharides, low polyphenols plant tissues   |
| <i>TransDirect</i> <sup>®</sup> Blood PCR Kit         | Fresh or frozen blood stored in EDTA, heparin, or citric acid<br>Fresh or dried blood without anticoagulant<br>Human oral epithelial cells |

## *TransDirect*<sup>®</sup> Animal Tissue PCR Kit

|          |                               |
|----------|-------------------------------|
| AD201-01 | 100 rxns (20 µl per reaction) |
| AD201-02 | 500 rxns (20 µl per reaction) |

### Storage

at -20°C for two years

### Description

*TransDirect*<sup>®</sup> Animal Tissue PCR Kit uses a unique lysis buffer to lyse animal tissues (fresh or frozen) and blood. The resulting lysate without purification can be directly used as PCR template. *2xTransDirect*<sup>®</sup> PCR SuperMix (+dye) is highly resistant to various PCR inhibitors present in animal tissues. PCR product can be directly used for gel electrophoresis.

### Applications

- Direct amplification from unpurified lysate. Suitable for high throughput applications.
- Suitable for mammalian cells, saliva, hair shaft, animal tissues and blood.
- Amplification of genomic DNA fragment up to 3 kb.

### Kit Contents

| Component   | AD201-01 | AD201-02 |
|---|----------|----------|
| AD1 Buffer  | 4 ml     | 20 ml    |
| AD2 Buffer  | 1 ml     | 5 ml     |
| AD3 Buffer  | 4 ml     | 2×10 ml  |
| <i>2xTransDirect</i> <sup>®</sup> PCR SuperMix (+dye) | 1 ml     | 5×1 ml   |
| ddH <sub>2</sub> O                                    | 5 ml     | 25 ml    |

### Materials

| Material        | Amount                |
|-----------------|-----------------------|
| Mammalian Cells | ≤10 <sup>6</sup> cell |
| Hair shaft      | ≤10 mg                |
| Mouse Tail      | ≤0.5 cm               |
| Mouse Ear       | ≤0.5 cm <sup>2</sup>  |
| Saliva          | ≤10 µl                |
| Animal Tissues  | ≤10 mg                |
| Blood           | ≤10 µl                |

High quality products



## PROTOCOL

### Genomic DNA extraction


- Mix 40  $\mu$ l of AD1 buffer with 10  $\mu$ l of AD2 buffer. For more samples, premix AD1 buffer with AD2 buffer at a ratio of 4:1. The mixture can be stored up to 2 hours at room temperature.
- Sample treatment
  - Mammalian Cells  
Pellet the cells by centrifugation and remove the supernatant. Add the mixture of AD1/AD2, mix thoroughly by pipetting up and down.
  - Saliva  
Directly add saliva into the mixture of AD1/AD2, mix thoroughly by pipetting up and down.
  - Hair Shafts  
Cut hair into pieces, add the mixture of AD1/AD2, mix thoroughly by pipetting up and down.
  - Animal Tissues  
Cut up tissues with sterile scissors or blade, add the mixture of AD1/AD2, mix thoroughly by pipetting up and down.
  - Blood  
Directly add blood into the mixture of AD1/AD2, mix thoroughly by pipetting up and down.
- Incubate at room temperature for 10 minutes, followed by at 95°C for 3 minutes (for hard-to-lyse tissues, like hair, we suggest incubating at 55°C for 10 minutes, followed by at 95°C for 3 minutes).
- Add 40  $\mu$ l of AD3 buffer, mix well. The lysate can be used as PCR template or stored at 4°C for three months or at -20°C for six months.

### PCR

| Component   | Volume                     | Final Concentration |
|---|----------------------------|---------------------|
| Unpurified Lysate                                       | Variable ( $\leq 4 \mu$ l) | as required         |
| Forward Primer (10 $\mu$ M)                             | 0.4 $\mu$ l                | 0.2 $\mu$ M         |
| Reverse Primer (10 $\mu$ M)                             | 0.4 $\mu$ l                | 0.2 $\mu$ M         |
| 2 $\times$ TransDirect <sup>®</sup> PCR SuperMix (+dye) | 10 $\mu$ l                 | 1 $\times$          |
| ddH <sub>2</sub> O                                      | Variable                   | -                   |
| Total volume  | 20 $\mu$ l                 | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 5-10 min   | } 35-40 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



M: Trans2K<sup>®</sup> Plus II DNA Marker

Lane 1: Hair 0.8 kb  
 Lane 2: Saliva 0.8 kb  
 Lane 3: HeLa cell 0.8 kb  
 Lane 4: HeLa cell 2 kb  
 Lane 5: HeLa cell 3 kb  
 Lane 6: Mouse ear 0.86 kb  
 Lane 7: Mouse ear 1.8 kb  
 Lane 8: Mouse ear 3 kb  
 Lane 9: Drosophila 0.42 kb  
 Lane 10: Drosophila 2 kb  
 Lane 11: Nematode 0.9 kb  
 Lane 12: Shrimp 1.1 kb  
 Lane 13: Crab 1 kb  
 Lane 14: Razor clam 0.56 kb





# TransDirect<sup>®</sup> Plant Tissue PCR Kit

|          |                               |
|----------|-------------------------------|
| AD301-01 | 100 rxns (20 µl per reaction) |
| AD301-02 | 500 rxns (20 µl per reaction) |

## Storage

at -20°C for two years

## Description

*TransDirect<sup>®</sup>* Plant Tissue PCR Kit uses a unique lysis buffer to lyse plant tissues (fresh or frozen). The resulting lysate without purification can be directly used as PCR template. *2×TransDirect<sup>®</sup>* PCR SuperMix (+dye) is highly resistant to various PCR inhibitors present in plant tissues. PCR product can be directly used for gel electrophoresis.

## Applications

- Direct amplification from unpurified lysate. Suitable for high throughput applications.
- Amplification of genomic DNA fragment up to 2 kb.

## Kit Contents

| Component  | AD301-01 | AD301-02 |
|--|----------|----------|
| PD1 Buffer   | 4 ml     | 20 ml    |
| PD2 Buffer   | 4 ml     | 20 ml    |
| <i>2×TransDirect<sup>®</sup></i> PCR SuperMix (+dye) | 1 ml     | 5×1 ml   |
| ddH <sub>2</sub> O                                   | 5 ml     | 25 ml    |

## PROTOCOL

### Genomic DNA Extraction

1. Cut 5 mg or 0.5 cm<sup>2</sup> plant tissues and add it to a tube containing 40 µl of PD1 buffer, vortex.
2. Incubate at 95°C for 10 minutes (for hard-to-lyse tissues, we suggest incubating at 95°C for 30 minutes).
3. Add 40 µl of PD2 buffer and vortex to mix. The lysate can be used as PCR template or stored at 4°C for three months or at -20°C for six months.

### PCR

| Component  | Volume           | Final Concentration |
|--|------------------|---------------------|
| Unpurified Lysate                                    | Variable (≤4 µl) | as required         |
| Forward Primer (10 µM)                               | 0.4 µl           | 0.2 µM              |
| Reverse Primer (10 µM)                               | 0.4 µl           | 0.2 µM              |
| <i>2×TransDirect<sup>®</sup></i> PCR SuperMix (+dye) | 10 µl            | 1×                  |
| ddH <sub>2</sub> O                                   | Variable         | -                   |
| Total volume   | 20 µl            | -                   |



### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 5-10 min   |                |
| 94°C    | 30 sec     | } 35-40 cycles |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



M: *Trans2K<sup>®</sup>* Plus II DNA Marker

|                        |                             |                        |
|------------------------|-----------------------------|------------------------|
| Lane 1: Corn 0.4 kb    | Lane 8: Soybean 1.5 kb      | Lane 15: Cotton 0.3 kb |
| Lane 2: Corn 0.8 kb    | Lane 9: Arabidopsis 0.5 kb  | Lane 16: Cotton 1 kb   |
| Lane 3: Corn 1.5 kb    | Lane 10: Arabidopsis 0.9 kb | Lane 17: Cotton 1.6 kb |
| Lane 4: Wheat 0.4 kb   | Lane 11: Arabidopsis 1.5 kb | Lane 18: Rice 0.3 kb   |
| Lane 5: Wheat 0.9 kb   | Lane 12: Tobacco 0.5 kb     | Lane 19: Rice 0.9 kb   |
| Lane 6: Wheat 1.5 kb   | Lane 13: Tobacco 0.9 kb     | Lane 20: Rice 1.9 kb   |
| Lane 7: Soybean 0.9 kb | Lane 14: Tobacco 1.5 kb     |                        |

## *TransDirect<sup>®</sup>* Blood PCR Kit

|          |                               |
|----------|-------------------------------|
| AD401-01 | 100 rxns (20 µl per reaction) |
| AD401-02 | 500 rxns (20 µl per reaction) |

### Storage

at -20°C for two years

### Description

*TransDirect<sup>®</sup>* Blood PCR Kit is designed for DNA amplification from whole blood without DNA extraction. *2×TransDirect<sup>®</sup>* PCR SuperMix (+dye) is highly resistant to various PCR inhibitors present in blood.

### Applications

- Fresh or frozen blood stored in EDTA, heparin or citric acid
- Fresh or dried blood without anticoagulant
- Human oral epithelial cells
- Amplification of genomic DNA fragment up to 4 kb

### Kit Contents

| Component  | AD401-01 | AD401-02 |
|--|----------|----------|
| <i>2×TransDirect<sup>®</sup></i> PCR SuperMix (+dye) | 1 ml     | 5×1 ml   |
| ddH <sub>2</sub> O                                   | 5 ml     | 25 ml    |



## PROTOCOL

### PCR

| Component                          | Volume                            | Final Concentration |
|------------------------------------|-----------------------------------|---------------------|
| Blood                              | Variable ( $\leq 1 \mu\text{l}$ ) | as required         |
| Forward Primer (10 $\mu\text{M}$ ) | 0.4 $\mu\text{l}$                 | 0.2 $\mu\text{M}$   |
| Reverse Primer (10 $\mu\text{M}$ ) | 0.4 $\mu\text{l}$                 | 0.2 $\mu\text{M}$   |
| 2xTransDirect® PCR SuperMix (+dye) | 10 $\mu\text{l}$                  | 1x                  |
| ddH <sub>2</sub> O                 | Variable                          | -                   |
| Total volume                       | 20 $\mu\text{l}$                  | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 94°C    | 5 min      | } 30-40 cycles |
| 94°C    | 30 sec     |                |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |

human frozen blood (EDTA anticoagulated) as templates

human frozen blood (heparin anticoagulated) as templates

human fresh blood (without anticoagulant) as templates

0.5  $\mu\text{l}$  of blood in 20  $\mu\text{l}$  reaction.

M: 1Kb Plus DNA Ladder  
 1: Hdt gene 0.32 kb  
 2: Hmt gene 0.5 kb  
 3: BDNF 0.8 kb  
 4: Rhod 1.2 kb  
 5:  $\beta$ -globin 1.3 kb  
 6: Rhod 2.0 kb  
 7:  $\beta$ -globin 3.0 kb  
 8: Rhod 4.17 kb  
 9:  $\beta$ -globin 4.1 kb

0.5  $\mu\text{l}$  of diluted chicken blood (sodium citrate anticoagulated) as templates to amplify a 0.25 kb fragment in 20  $\mu\text{l}$  reaction.  
 M: Trans2K® DNA Marker  
 Lane 1: Blood  
 Lane 2: 1:10 dilution  
 Lane 3: 1:100 dilution


0.5  $\mu\text{l}$  of 1:40 diluted mouse blood (heparin anticoagulated) as templates to amplify Neo and MAP genes in 20  $\mu\text{l}$  reaction.  
 M: Trans2K® DNA Marker  
 1: Neo 0.25 kb  
 2: MAP 0.6 kb

Human oral epithelial cells as templates to amplify Hdt gene (0.32 kb) in 20  $\mu\text{l}$  reaction.  
 M: Trans2K® DNA Marker




# RT-PCR

Trans reverse transcriptases




| Products   | Size of cDNA products | Temperature | Sensitivity | Fidelity | GC-rich or Complex template |
|--|-----------------------|-------------|-------------|----------|-----------------------------|
| EasyScript® RT   | ≤8 kb                 | 42°C        | •           | •        | •                           |
| TransScript® RT  | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript® II RT   | ≤15 kb                | 42°C-55°C   | •••         | ••       | •••                         |
| EasyScript® First-Strand cDNA Synthesis SuperMix   | ≤8 kb                 | 42°C        | •           | •        | •                           |
| EasyScript® One-Step gDNA Removal and cDNA Synthesis SuperMix                                    | ≤8 kb                 | 42°C        | •           | •        | •                           |
| TransScript® First-Strand cDNA Synthesis SuperMix  | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript® One-Step gDNA Removal and cDNA Synthesis SuperMix                                   | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript® Fly First-Strand cDNA Synthesis SuperMix  | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript®-Uni One-Step gDNA Removal and cDNA Synthesis SuperMix                               | ≤20 kb                | 42°C-65°C   | •••         | ••       | •••                         |
| TransScript®-Uni Cell to cDNA Synthesis SuperMix for qPCR  | ≤ 250 bp              | 42°C        | •••         | ••       | •••                         |
| TransScript® miRNA First-Strand cDNA Synthesis SuperMix  | ≤ 250 bp              | 42°C        | ••          | ••       | ••                          |
| TransScript® II First-Strand cDNA Synthesis SuperMix   | ≤15 kb                | 42°C-55°C   | •••         | ••       | •••                         |
| TransScript® II One-Step gDNA Removal and cDNA Synthesis SuperMix                                | ≤15 kb                | 42°C-55°C   | •••         | ••       | •••                         |
| TransScript® All-in-One First-Strand cDNA Synthesis SuperMix for PCR                             | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript® All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)    | ≤ 250 bp              | 42°C        | ••          | ••       | ••                          |
| TransScript® II All-in-One First-Strand cDNA Synthesis SuperMix for PCR                          | ≤15 kb                | 42°C-55°C   | •••         | ••       | •••                         |
| TransScript® II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) | ≤ 250 bp              | 42°C-55°C   | •••         | ••       | •••                         |
| TransScript® Two-Step RT-PCR SuperMix  | ≤12 kb                | 42°C        | ••          | ••       | ••                          |
| TransScript® II Two-Step RT-PCR SuperMix   | ≤15 kb                | 42°C-55°C   | •••         | ••       | •••                         |
| EasyScript® One-Step RT-PCR SuperMix   | ≤4 kb                 | 45°C        | •           | •        | •                           |
| TransScript® One-Step RT-PCR SuperMix  | ≤8 kb                 | 45°C        | ••          | ••       | ••                          |
| TransScript® II One-Step RT-PCR SuperMix   | ≤8 kb                 | 50°C        | •••         | ••       | •••                         |




• *Trans* RT (except one-step kits) use Anchored Oligo(dT) to increase cDNA yield and full length cDNA products.

Traditional Oligo(dT)<sub>12-18</sub> Primer



Anchored Oligo(dT) Primer



- Poly(A) tail can be a few hundreds nt long and Oligo(dT) Primer binds randomly.
- Lower efficiency because of long poly(A) tail.
- Less full length cDNA products.

- Anchored Oligo(dT) Primer only anneals at 5' end of the Poly(A) tail of mRNA.
- Higher efficiency because of Anchored Oligo(dT) Primer.
- More full length cDNA products.





# EasyScript<sup>®</sup> Reverse Transcriptase[M-MLV, RNase H<sup>-</sup>]

|          |                |
|----------|----------------|
| AE101-02 | 10,000 units   |
| AE101-03 | 5×10,000 units |

## Concentration

200 units/μl

## Contents

- EasyScript<sup>®</sup> RT
- 5×ES RT Buffer  
(375 mM KCl; 15 mM MgCl<sub>2</sub>;  
100 mM Tris-HCl pH 8.4)
- Anchored Oligo(dT)<sub>18</sub> Primer

## Storage

at -20°C for one year

## Description

EasyScript<sup>®</sup> Reverse Transcriptase is an engineered version of M-MLV reverse transcriptase with deficient RNase H activity. The enzyme is purified to near homogeneity from *E. coli* containing the modified M-MLV RT gene.

- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- Anchored Oligo(dT)<sub>18</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 8 kb.

## Applications

- First-strand cDNA synthesis
- Multiple copy gene detection

## Unit Definition

One unit of EasyScript<sup>®</sup> RT incorporates 1 nmol of deoxyribonucleotide into acid precipitable material in 10 minutes at 37°C using Poly(A)/Oligo(dT) as template/primer.

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume              |
|---|---------------------|
| Total RNA/mRNA  | 50 ng-5 μg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 μg/μl)<br>or Random Primer(N9) (0.1 μg/μl)<br>or GSP | 1 μl                |
| 10 mM dNTPs   | 1 μl                |
| 5×ES RT Buffer  | 4 μl                |
| Ribonuclease Inhibitor (50 units/μl)  | 0.5 μl              |
| EasyScript <sup>®</sup> RT  | 1 μl                |
| RNase-free Water  | to 20 μl            |

#### 2. Incubation

- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 30 minutes.
- For random primer, incubate at 25°C for 10 minutes, then at 42°C for 30 minutes.

#### 3. Incubate at 85°C for 5 seconds to inactivate enzymes.



### RT-PCR


#### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| Template   | Variable | as required         |
| Forward Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| Reverse Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| 2x <i>TransTaq</i> <sup>®</sup> HiFi PCR SuperMix II | 25 μl    | 1x                  |
| ddH <sub>2</sub> O                                   | Variable | -                   |
| Total volume   | 50 μl    | -                   |


#### Thermal cycling conditions

94°C            2-5 min  
 94°C            30 sec  
 50-60°C       30 sec  
 72°C           1-2 kb/min  
 72°C           5-10 min

} 35-40 cycles



M: 1Kb Plus DNA Ladder  
 Lanes 1-7: 0, 0.01, 0.1, 1, 10, 100,  
 1,000 pg human total RNA as templates



M: 1Kb Plus DNA Ladder  
 1: GAPDH 0.5 kb; 2: GAPDH 0.9 kb;  
 3: REPA 1.8 kb; 4: ACTR 3 kb;  
 5: ACTR 3.5 kb; 6: VIN 4.6 kb;  
 7: TSC 5.3 kb; 8: Pol 6.8 kb  
 Human cDNA as templates



# TransScript<sup>®</sup> Reverse Transcriptase[M-MLV, RNase H<sup>-</sup>]

|          |                |
|----------|----------------|
| AT101-02 | 10,000 units   |
| AT101-03 | 5×10,000 units |

## Concentration

200 units/μl

## Contents

- TransScript<sup>®</sup> RT
- 5×TS RT Buffer  
(250 mM KCl; 15 mM MgCl<sub>2</sub>;  
100 mM Tris-HCl pH 8.4)
- Anchored Oligo(dT)<sub>18</sub> Primer

## Storage

at -20°C for one year

## Description

TransScript<sup>®</sup> Reverse Transcriptase is a recombinant M-MLV reverse transcriptase with deficient RNase H activity.

- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- Anchored Oligo(dT)<sub>18</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 12 kb.

## Applications


- First-strand cDNA synthesis
- Multiple copy and low copy gene detection

## Unit Definition

One unit of TransScript<sup>®</sup> RT incorporates 1 nmol of deoxyribonucleotide into acid-precipitable material in 10 minutes at 37°C using Poly(A)/Oligo(dT) as template/primer.

## PROTOCOL

The suggested condition for the first-strand cDNA synthesis and RT-PCR are the same as described on pages 41-42.





# *TransScript*<sup>®</sup> II Reverse Transcriptase [M-MLV, RNase H<sup>-</sup>] (High Temperature RT)

AH101-02

10,000 units

## Concentration

200 units/μl

## Contents

- *TransScript*<sup>®</sup> II RT
- 10xTS II RT Buffer  
(500 mM KCl; 30 mM MgCl<sub>2</sub>;  
200 mM Tris-HCl pH 8.4)
- Anchored Oligo(dT)<sub>20</sub> Primer

## Storage

at -20°C for one year

## Description

*TransScript*<sup>®</sup> II Reverse Transcriptase is a recombinant M-MLV reverse transcriptase with deficient RNase H activity and increased thermostability. The enzyme is active at up to 55°C. It provides higher specificity, higher yield and more full-length cDNA products.

- Increased thermostability for more full-length cDNA products.
- Reaction temperature at 42°C-55°C.
- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- Anchored Oligo(dT)<sub>20</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 15 kb.

## Applications

- First-strand cDNA synthesis, cDNA library construction, 3' and 5' RACE
- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Unit Definition

One unit of *TransScript*<sup>®</sup> II RT incorporates 1 nmol of deoxyribonucleotide into acid-precipitable material in 10 minutes at 37°C using Poly(A)/Oligo(dT) as template/primer.

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume              |
|---|---------------------|
| Total RNA/mRNA  | 50 ng-5 μg/5-500 ng |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 μg/μl)<br>or Random Primer(N9) (0.1 μg/μl) | 1 μl                |
| or GSP  | 2 pmol              |
| 10 mM dNTPs   | 1 μl                |
| 10xTS II RT Buffer  | 2 μl                |
| Ribonuclease Inhibitor (50 units/μl)  | 0.5 μl              |
| <i>TransScript</i> <sup>®</sup> II RT   | 1 μl                |
| RNase-free Water  | to 20 μl            |

#### 2. Incubation


- For anchored oligo(dT)<sub>20</sub> primer or GSP, incubate at 50°C for 30 minutes.
- For random primer, incubate at 25°C for 10 minutes, then at 50°C for 30 minutes.
- For GC-rich or complex secondary structure RNA template, incubate at 55°C for 30 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.



## RT-PCR

The suggested reaction condition is the same as described on page 42.



M1: 1Kb Plus DNA Ladder  
 M2: *Trans15K* DNA Marker  
 1: REPA 1.8 kb; 2: NCBP 2.5 kb;  
 3: HDP 3.5 kb; 4: VIN 4.6 kb;  
 5: Pol 6.8 kb; 6: APC 8.5 kb;  
 7: Dynein 12.3 kb; 8: FAL 15.1 kb  
 Human cDNA as templates

## EasyScript<sup>®</sup> First-Strand cDNA Synthesis SuperMix

|          |                                    |
|----------|------------------------------------|
| AE301-02 | 50 rxns (20 $\mu$ l per reaction)  |
| AE301-03 | 100 rxns (20 $\mu$ l per reaction) |

### Storage

at -20°C for one year

### Description

*EasyScript<sup>®</sup>* First-Strand cDNA Synthesis SuperMix provides all the necessary components for cDNA synthesis from total RNA or mRNA. The cDNA is efficiently synthesized by *EasyScript<sup>®</sup>* RT/RI Enzyme Mix and 2 $\times$ ES Reaction Mix.

- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- Anchored Oligo(dT)<sub>18</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 8 kb.

### Application

Multiple copy gene detection

### Kit Contents

| Component  | AE301-02    | AE301-03    |
|--|-------------|-------------|
| <i>EasyScript<sup>®</sup></i> RT/RI Enzyme Mix                 | 50 $\mu$ l  | 100 $\mu$ l |
| 2 $\times$ ES Reaction Mix                                     | 500 $\mu$ l | 1 ml        |
| Random Primer(N9) (0.1 $\mu$ g/ $\mu$ l)                       | 50 $\mu$ l  | 100 $\mu$ l |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 $\mu$ g/ $\mu$ l) | 50 $\mu$ l  | 100 $\mu$ l |
| RNase-free Water   | 500 $\mu$ l | 1 ml        |





## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA  | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl)<br>or Random Primer(N9) (0.1 µg/µl) | 1 µl                 |
| or GSP  | 2 pmol               |
| 2×ES Reaction Mix   | 10 µl                |
| <i>EasyScript</i> <sup>®</sup> RT/RI Enzyme Mix   | 1 µl                 |
| RNase-free Water  | to 20 µl             |

#### 2. Incubation

- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

**The suggested reaction condition is the same as described on page 42.**



# EasyScript® One-Step gDNA Removal and cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AE311-02 | 50 rxns (20 µl per reaction)  |
| AE311-03 | 100 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

Unique genomic DNA remover is combined with *EasyScript*® First-Strand cDNA Synthesis SuperMix to achieve simultaneous genomic DNA removal and cDNA synthesis. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds.

- Simultaneous genomic DNA removal and cDNA synthesis in one tube to minimize RNA contamination.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- cDNA up to 8 kb.

## Application

Multiple copy gene detection

## Kit Contents

| Component   | AE311-02 | AE311-03 |
|---|----------|----------|
| <i>EasyScript</i> ® RT/RI Enzyme Mix                | 50 µl    | 100 µl   |
| gDNA Remover  | 50 µl    | 100 µl   |
| 2×ES Reaction Mix                                   | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA                                      | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 1 µl                 |
| or Random Primer(N9) (0.1 µg/µl)                    | 1 µl                 |
| or GSP  | 2 pmol               |
| 2×ES Reaction Mix                                   | 10 µl                |
| <i>EasyScript</i> ® RT/RI Enzyme Mix                | 1 µl                 |
| gDNA Remover  | 1 µl                 |
| RNase-free Water                                    | to 20 µl             |

#### 2. Incubation

- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

**The suggested reaction condition is the same as described on page 42.**



# TransScript® First-Strand cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AT301-02 | 50 rxns (20 µl per reaction)  |
| AT301-03 | 100 rxns (20 µl per reaction) |

### Storage

at -20°C for one year

### Description

TransScript® First-Strand cDNA Synthesis SuperMix provides all the necessary components for cDNA synthesis from total RNA or mRNA. The cDNA is efficiently synthesized by TransScript® RT/RI Enzyme Mix and 2×TS Reaction Mix.

- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- Anchored Oligo(dT)<sub>18</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 12 kb.

### Application

Multiple copy and low copy gene detection

### Kit Contents

| Component   | AT301-02 | AT301-03 |
|---|----------|----------|
| TransScript® RT/RI Enzyme Mix                       | 50 µl    | 100 µl   |
| 2×TS Reaction Mix                                   | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA  | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl)<br>or Random Primer(N9) (0.1 µg/µl) | 1 µl                 |
| or GSP  | 2 pmol               |
| 2×TS Reaction Mix   | 10 µl                |
| TransScript® RT/RI Enzyme Mix   | 1 µl                 |
| RNase-free Water  | to 20 µl             |

#### 2. Incubation

- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

**The suggested reaction condition is the same as described on page 42.**



# TransScript<sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AT311-02 | 50 rxns (20 µl per reaction)  |
| AT311-03 | 100 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

Unique genomic DNA remover is combined with *TransScript*<sup>®</sup> First-Strand cDNA Synthesis SuperMix to achieve simultaneous genomic DNA removal and cDNA synthesis. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds.

- Simultaneous genomic DNA removal and cDNA synthesis in one tube to minimize RNA contamination.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- cDNA up to 12 kb.

## Application

Multiple copy and low copy gene detection

## Kit Contents

| Component   | AT311-02 | AT311-03 |
|---|----------|----------|
| <i>TransScript</i> <sup>®</sup> RT/RI Enzyme Mix    | 50 µl    | 100 µl   |
| gDNA Remover  | 50 µl    | 100 µl   |
| 2×TS Reaction Mix                                   | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA                                      | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 1 µl                 |
| or Random Primer(N9) (0.1 µg/µl)                    | 1 µl                 |
| or GSP  | 2 pmol               |
| 2×TS Reaction Mix                                   | 10 µl                |
| <i>TransScript</i> <sup>®</sup> RT/RI Enzyme Mix    | 1 µl                 |
| gDNA Remover  | 1 µl                 |
| RNase-free Water                                    | to 20 µl             |


#### 2. Incubation

- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 42°C for 15 minutes (for qPCR) or incubate at 42°C for 30 minutes (for PCR).

#### 3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

The suggested reaction condition is the same as described on page 42.



M: *Trans2K*<sup>®</sup> Plus DNA Marker  
 PCR +/-gDNA remover  
 Lane 1: 200 ng Human Genomic DNA (-gDNA remover);  
 Lane 2: 200 ng Human Genomic DNA (+gDNA remover);  
 RT-PCR +/- gDNA remover  
 Lane 3: 100 ng Human total RNA (-gDNA remover);  
 Lane 4: 100 ng Human total RNA (+gDNA remover);  
 cDNA as template, 1 kb  
 Genomic DNA as template, 2 kb



# TransScript<sup>®</sup> Fly First-Strand cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AF301-02 | 50 rxns (20 µl per reaction)  |
| AF301-03 | 100 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

TransScript<sup>®</sup> Fly RT is generated by site mutations. It provides high affinity to RNA template with fast extension rate. The cDNA is efficiently synthesized by TransScript<sup>®</sup> Fly RT/RI Enzyme Mix and 2×TS Fly Reaction Mix. The entire reverse transcription can be completed within 5 minutes.

- 5 minutes reverse transcription.
- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- Anchored Oligo(dT)<sub>18</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 12 kb.

## Application

Multiple copy and low copy gene fast detection

## Kit Contents

| Component   | AF301-02 | AF301-03 |
|---|----------|----------|
| TransScript <sup>®</sup> Fly RT/RI Enzyme Mix       | 50 µl    | 100 µl   |
| 2×TS Fly Reaction Mix                               | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA                                      | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>18</sub> Primer (0.5 µg/µl) | 1 µl                 |
| or Random Primer(N9) (0.1 µg/µl)                    | 1 µl                 |
| or GSP  | 2 pmol               |
| 2×TS Fly Reaction Mix                               | 10 µl                |
| TransScript <sup>®</sup> Fly RT/RI Enzyme Mix       | 1 µl                 |
| RNase-free Water                                    | to 20 µl             |

#### 2. Incubation


- For anchored oligo(dT)<sub>18</sub> primer or GSP, incubate at 42°C for 5 minutes.
- For random primer, incubate at 25°C for 10 minutes, then at 42°C for 5 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

**The suggested reaction condition is the same as described on page 42.**





## *TransScript*<sup>®</sup>-Uni One-Step gDNA Removal and cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AU311-02 | 50 rxns (20 µl per reaction)  |
| AU311-03 | 100 rxns (20 µl per reaction) |

### Storage

at -20°C for one year

### Description

*TransScript*<sup>®</sup>-Uni RT is an improved version of M-MLV reverse transcriptase with broad range of reaction temperature (42°C-65°C) and higher thermostability. The suggested reaction temperature is 50°C. The SuperMix contains reagents for simultaneous genomic DNA removal and cDNA synthesis. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds.

- Broad range reaction temperature (42°C-65°C) .
- Simultaneous genomic DNA removal and cDNA synthesis in one tube to minimize RNA contamination.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- cDNA up to 20 kb.

### Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template
- cDNA library construction, primer extension, 3' and 5' RACE

**Kit Contents**

| Component   | AU311-02 | AU311-03 |
|---|----------|----------|
| <i>TransScript</i> <sup>®</sup> -Uni RT/RI Enzyme Mix | 50 µl    | 100 µl   |
| gDNA Remover  | 50 µl    | 100 µl   |
| 2xTS-Uni Reaction Mix                                 | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                         | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg/µl)   | 50 µl    | 100 µl   |
| RNase-free Water                                      | 500 µl   | 1 ml     |

**PROTOCOL****First-Strand cDNA synthesis**

## 1. Reaction Components

| Component   | Volume               |
|---|----------------------|
| Total RNA/mRNA  | 50 ng -5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg/µl)<br>or Random Primer(N9) (0.1 µg/µl) | 1 µl                 |
| or GSP  | 2 pmol               |
| 2xTS-Uni Reaction Mix   | 10 µl                |
| gDNA Remover  | 1 µl                 |
| <i>TransScript</i> <sup>®</sup> -Uni RT/RI Enzyme Mix                                   | 1 µl                 |
| RNase-free Water  | to 20 µl             |


## 2. Incubation

- For anchored oligo(dT)<sub>20</sub> primer or GSP, incubate at 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, at incubate 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For GC-rich or complex secondary structure RNA template, better yield can be obtained by optimizing the reaction temperature.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

**RT-PCR**

**The suggested reaction condition is the same as described on page 42.**



M: *Trans2K*<sup>®</sup> Plus II DNA Marker  
 Lane 1: Company A kit, 42°C  
 Lane 2: Company B kit, 42°C  
 Lane 3: *TransScript*<sup>®</sup>-Uni, 42°C  
 Lane 4: *TransScript*<sup>®</sup>-Uni, 50°C

Lane 5: *TransScript*<sup>®</sup>-Uni, 55°C  
 Lane 6: *TransScript*<sup>®</sup>-Uni, 60°C  
 Lane 7: *TransScript*<sup>®</sup>-Uni, 65°C  
 Human cDNA as template, VIN 4.6 kb



# TransScript<sup>®</sup>-Uni Cell to cDNA Synthesis SuperMix for qPCR

AC301-01

25 rxns

## Storage

at -20°C for one year

## Description

TransScript<sup>®</sup>-Uni Cell to cDNA Synthesis SuperMix for qPCR uses a unique lysis buffer to lyse cells. The resulting lysate (without purification) can be directly used as template for reverse transcription. Unique genomic DNA remover is combined with TransScript<sup>®</sup>-Uni RT/RI Enzyme Mix to achieve simultaneous genomic DNA removal and cDNA synthesis in one tube. This kit is suitable to generate qPCR-ready cDNA directly from cells.


## Application

Multiple copy and low copy gene detection


## Kit Contents

| Component                                      | AC301-01  |
|--|-----------|
| C to C Lysis Buffer                            | 2×1.25 ml |
| TransScript <sup>®</sup> -Uni RT/RI Enzyme Mix | 12.5 μl   |
| gDNA Remover                                   | 12.5 μl   |
| 2×TS-Uni Reaction Mix                          | 250 μl    |
| Oligo(dT)/RP Mix                               | 25 μl     |
| RNase-free Water                               | 250 μl    |

Amplification Plots



Standard Curve



## Cell Types Tested with the Kit

|          |            |            |
|----------|------------|------------|
| A549     | Hep G2     | SGC-7901   |
| CHO-K1   | K-562      | Sp2/0-Ag14 |
| HEK-293  | MCF7       | Vero       |
| HEK-293T | MDA-MB-231 | WI-38      |
| HeLa     | P815       |            |



## PROTOCOL

### Cell Lysis

1. Add 50  $\mu$ l of C to C Lysis Buffer to each well ( $5 \times 10^2$ - $5 \times 10^4$  cells), incubate at room temperature (22°C-25°C) for 5 minutes.
2. Mix by pipetting up and down. Transfer the lysate into a microcentrifuge tube. Incubate at 75°C for 5 minutes, then place the tube on ice.

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume        |
|---|---------------|
| Cell Lysate   | 2 $\mu$ l     |
| Oligo(dT)/RP Mix                                      | 1 $\mu$ l     |
| 2 $\times$ TS-Uni Reaction Mix                        | 10 $\mu$ l    |
| gDNA Remover  | 0.5 $\mu$ l   |
| <i>TransScript</i> <sup>®</sup> -Uni RT/RI Enzyme Mix | 0.5 $\mu$ l   |
| RNase-free Water                                      | to 20 $\mu$ l |

2. Gently mix and incubate at 42°C for 15 minutes.
3. Incubate at 85°C for 5 seconds to inactivate *TransScript*<sup>®</sup>-Uni RT/RI Enzyme Mix and gDNA Remover.

#### Suggested qPCR conditions (20 $\mu$ l reaction volume)

| Component   | Volume      | Final Concentration |
|---|-------------|---------------------|
| Template  | Variable    | as required         |
| Forward Primer (10 $\mu$ M)   | 0.4 $\mu$ l | 0.2 $\mu$ M         |
| Reverse Primer (10 $\mu$ M)   | 0.4 $\mu$ l | 0.2 $\mu$ M         |
| 2 $\times$ <i>TransStart</i> <sup>®</sup> Top/Tip Green qPCR SuperMix | 10 $\mu$ l  | 1 $\times$          |
| Passive Reference Dye (50 $\times$ ) (optional)                       | 0.4 $\mu$ l | 1 $\times$          |
| ddH <sub>2</sub> O  | Variable    | -                   |
| Total Volume  | 20 $\mu$ l  | -                   |

#### Thermal cycling conditions (three-step)

|         |         |                |
|---------|---------|----------------|
| 94°C    | 30 sec  | } 40-45 cycles |
| 94°C    | 5 sec   |                |
| 50-60°C | 15 sec* |                |
| 72°C    | 10 sec* |                |

Dissociation Stage

#### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 94°C | 30 sec  | } 40-45 cycles |
| 94°C | 5 sec   |                |
| 60°C | 30 sec* |                |

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

\* For ABI Prism7700/7900, the time is 30 seconds.

\* For ABI Prism7000/7300, the time is 31 seconds.

\* For ABI Prism7500, the time is 34 seconds.

\* For ABI Viia 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**



# TransScript® miRNA First-Strand cDNA Synthesis SuperMix

AT351-01

20 rxns ( 20 µl per reaction)

## Storage

at -20°C for one year

## Description

*TransScript*® miRNA First-Strand cDNA Synthesis SuperMix provides all the necessary components for cDNA synthesis from miRNA template. High efficient poly(A) tail addition and first-strand cDNA synthesis are performed by *TransScript*® miRNA RT Enzyme Mix (containing tailing enzyme and RT enzyme) and 2×TS miRNA Reaction Mix.

- Optimized enzyme and buffer system for high efficient cDNA synthesis.
- One-step Poly(A) tailing and cDNA synthesis.

## Application

miRNA synthesis

## Kit Contents

| Component                                | AT351-01 |
|--|----------|
| <i>TransScript</i> ® miRNA RT Enzyme Mix | 20 µl    |
| 2×TS miRNA Reaction Mix                  | 200 µl   |
| Universal miRNA qPCR Primer (10 µM)      | 200 µl   |
| RNase-free Water                         | 1 ml     |

## PROTOCOL

### Tail addition and First-Strand cDNA synthesis

#### 1. Reaction Components

| Component                                | Volume   |
|--|----------|
| Total RNA/miRNA*                         | x µl     |
| <i>TransScript</i> ® miRNA RT Enzyme Mix | 1 µl     |
| 2×TS miRNA Reaction Mix                  | 10 µl    |
| RNase-free Water                         | to 20 µl |

\* Total RNA ≤5 µg. Since miRNA cannot be directly quantified by spectrophotometer, we suggest using 1-9 µl for 20 µl reaction.

- Mix gently, and incubate at 37°C for 1 hour.
- Incubate at 85°C for 5 seconds to inactivate RT Enzyme Mix.

### Suggested qPCR conditions (20 µl reaction volume)

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| cDNA* <sup>1</sup>                                 | Variable | as required         |
| Forward Primer (10 µM)* <sup>2</sup>               | 0.4 µl   | 0.2 µM              |
| Universal miRNA qPCR Primer (10 µM)                | 0.4 µl   | 0.2 µM              |
| 2× <i>TransStart</i> ® Tip/Top Green qPCR SuperMix | 10 µl    | 1×                  |
| Passive Reference Dye (50×) (optional)             | 0.4 µl   | 1×                  |
| ddH <sub>2</sub> O                                 | Variable | -                   |
| Total volume                                       | 20 µl    | -                   |

\*1. We suggest diluting the synthesized cDNA 5-10 folds.

\*2. Upstream primer is target miRNA specific primer, which will be designed by customers according to target miRNA.



**Thermal cycling conditions (three-step)**

|                    |         |                |
|--------------------|---------|----------------|
| 94°C               | 30 sec  |                |
| 94°C               | 5 sec   | } 40-45 cycles |
| 50-60°C            | 15 sec* |                |
| 72°C               | 10 sec* |                |
| Dissociation Stage |         |                |


**Thermal cycling conditions (two-step)**

|                    |         |                |
|--------------------|---------|----------------|
| 94°C               | 30 sec  |                |
| 94°C               | 5 sec   | } 40-45 cycles |
| 60°C               | 30 sec* |                |
| Dissociation Stage |         |                |

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**  
**Three-step qPCR is more suitable for higher sensitivity assay.**



**Principle of miRNA Detection**



# TransScript<sup>®</sup> II First-Strand cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AH301-02 | 50 rxns (20 µl per reaction)  |
| AH301-03 | 100 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

*TransScript<sup>®</sup> II First-Strand cDNA Synthesis SuperMix* provides all the necessary components for cDNA synthesis from total RNA or mRNA. The cDNA is efficiently synthesized by *TransScript<sup>®</sup> II RT/RI Enzyme Mix* and *2×TS II Reaction Mix*.

- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.
- Anchored Oligo(dT)<sub>20</sub> Primer for higher yield and more full length cDNA.
- cDNA up to 15 kb.

## Applications

- cDNA library construction, 3' and 5' RACE
- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Kit Contents

| Component   | AH301-02 | AH301-03 |
|---|----------|----------|
| <i>TransScript<sup>®</sup> II RT/RI Enzyme Mix</i>  | 50 µl    | 100 µl   |
| <i>2×TS II Reaction Mix</i>                         | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component  | Volume                 |
|--|------------------------|
| Total RNA/mRNA   | 50 ng-5 µg/5-500 ng    |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg /µl)<br>or Random Primer(N9) (0.1 µg/µl)<br>or GSP | 1 µl<br>1 µl<br>2 pmol |
| <i>2×TS II Reaction Mix</i>  | 10 µl                  |
| <i>TransScript<sup>®</sup> II RT/RI Enzyme Mix</i>   | 1 µl                   |
| RNase-free Water   | to 20 µl               |

#### 2. Incubation

- For anchored oligo(dT)<sub>20</sub> primer or GSP, incubate at 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For GC-rich or complex secondary structure RNA template, incubate at 55°C for 30 minutes.

#### 3. Incubate at 85°C for 5 seconds to inactivate enzymes.

## RT-PCR

**The suggested reaction condition is the same as described on page 42.**



# TransScript® II One-Step gDNA Removal and cDNA Synthesis SuperMix

|          |                               |
|----------|-------------------------------|
| AH311-02 | 50 rxns (20 µl per reaction)  |
| AH311-03 | 100 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

Unique genomic DNA remover is combined with *TransScript*® II First-Strand cDNA Synthesis SuperMix to achieve simultaneous genomic DNA removal and cDNA synthesis. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds.

- Simultaneous genomic DNA removal and cDNA synthesis in one tube to minimize RNA contamination. The product obtained from 15 minutes reaction is used for qPCR; the product obtained from 30 minutes reaction is used for PCR.

- cDNA up to 15 kb.

## Applications

- cDNA library construction, 3' and 5' RACE
- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Kit Contents

| Component   | AH311-02 | AH311-03 |
|---|----------|----------|
| <i>TransScript</i> ® II RT/RI Enzyme Mix            | 50 µl    | 100 µl   |
| gDNA Remover  | 50 µl    | 100 µl   |
| 2xTS II Reaction Mix                                | 500 µl   | 1 ml     |
| Random Primer(N9) (0.1 µg/µl)                       | 50 µl    | 100 µl   |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg/µl) | 50 µl    | 100 µl   |
| RNase-free Water                                    | 500 µl   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume              |
|---|---------------------|
| Total RNA/mRNA                                      | 50 ng-5 µg/5-500 ng |
| Anchored Oligo(dT) <sub>20</sub> Primer (0.5 µg/µl) | 1 µl                |
| or Random Primer(N9) (0.1 µg/µl)                    | 1 µl                |
| or GSP  | 2 pmol              |
| 2xTS II Reaction Mix                                | 10 µl               |
| <i>TransScript</i> ® II RT/RI Enzyme Mix            | 1 µl                |
| gDNA Remover  | 1 µl                |
| RNase-free Water                                    | to 20 µl            |

#### 2. Incubation

- For anchored oligo(dT)<sub>20</sub> primer or GSP, incubate at 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For random primer, incubate at 25°C for 10 minutes. After that, incubate at 50°C for 15 minutes (for qPCR) or incubate at 50°C for 30 minutes (for PCR).
- For GC-rich or complex secondary structure RNA template, incubate at 55°C for 30 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

The suggested reaction condition is the same as described on page 42.



# TransScript<sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for PCR

AT321-01

50 rxns (20 µl per reaction)

## Storage

at -20°C for one year

## Description

TransScript<sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for PCR provides all the necessary components for cDNA synthesis from total RNA or mRNA. The SuperMix is provided at 5× concentration and used at 1× concentration by adding RNA and H<sub>2</sub>O. The resulting cDNA is suitable for regular PCR, not for qPCR.

- One-tube format for simple and fast setup and reducing pipetting variability.
- The optimal ratio of oligo(dT)<sub>18</sub> primer to random primer(N9) for PCR ready cDNA.
- PCR ready cDNA in 30 minutes (unsuitable for qPCR).
- cDNA up to 12 kb.

## Application

Multiple copy and low copy gene detection

## Kit Contents

| Component  | AT321-01 |
|--|----------|
| 5×TransScript <sup>®</sup> All-in-One SuperMix for PCR | 200 µl   |
| RNase-free Water                                       | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component  | Volume              |
|--|---------------------|
| Total RNA/mRNA   | 50 ng-5 µg/5-500 ng |
| 5×TransScript <sup>®</sup> All-in-One SuperMix for PCR | 4 µl                |
| RNase-free Water                                       | to 20 µl            |

#### 2. Incubation

- For RNA template with poly(A)<sup>+</sup>, incubate at 42°C for 30 minutes.
- For RNA template without poly(A)<sup>+</sup>, incubate at 25°C for 10 minutes, then at 42°C for 30 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

The suggested reaction condition is the same as described on page 42.



# TransScript<sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)

|          |                               |
|----------|-------------------------------|
| AT341-01 | 50 rxns (20 µl per reaction)  |
| AT341-02 | 100 rxns (20 µl per reaction) |

**Storage**  
at -20°C for one year

**Description**

The kit provides all the necessary components for cDNA synthesis from total RNA or mRNA. It is provided at 5× concentration and used at 1× concentration by adding gDNA remover, RNA and H<sub>2</sub>O. Simultaneous genomic DNA removal and cDNA synthesis are performed. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds. The resulting cDNA is suitable for qPCR, not for regular PCR.

- Simultaneous genomic DNA removal and cDNA synthesis.
- The optimal ratio of oligo(dT)<sub>18</sub> primer to random primer(N9) for qPCR ready cDNA.
- qPCR ready cDNA in 15 minutes.
- cDNA up to 250 bp.

**Application**

Multiple copy and low copy gene detection

**Kit Contents**

| Component   | AT341-01 | AT341-02 |
|---|----------|----------|
| 5×TransScript <sup>®</sup> All-in-One SuperMix for qPCR               | 200 µl   | 400 µl   |
| gDNA Remover  | 50 µl    | 100 µl   |
| 5×TransScript <sup>®</sup> All-in-One No-RT Control SuperMix for qPCR | 20 µl    | 40 µl    |
| RNase-free Water  | 1 ml     | 2×1 ml   |

## PROTOCOL

**First-Strand cDNA synthesis**

1. Reaction Components

| Component   | Volume        |
|---|---------------|
| Total RNA/mRNA  | ≤1 µg/≤100 ng |
| 5×TransScript <sup>®</sup> All-in-One SuperMix for qPCR | 4 µl          |
| gDNA Remover  | 1 µl          |
| RNase-free Water  | to 20 µl      |

2. Incubate at 42°C for 15 minutes.
3. Incubate at 85°C for 5 seconds to inactivate enzymes.





## qPCR

### Reaction Components

| Component   | Volume      | Final Concentration |
|---|-------------|---------------------|
| Template  | Variable    | as required         |
| Forward Primer (10 $\mu$ M)   | 0.4 $\mu$ l | 0.2 $\mu$ M         |
| Reverse Primer (10 $\mu$ M)   | 0.4 $\mu$ l | 0.2 $\mu$ M         |
| 2 $\times$ <i>TransStart</i> <sup>®</sup> Top/Tip Green qPCR SuperMix | 10 $\mu$ l  | 1 $\times$          |
| Passive Reference Dye (50 $\times$ ) (optional)                       | 0.4 $\mu$ l | 1 $\times$          |
| ddH <sub>2</sub> O  | Variable    | -                   |
| Total Volume  | 20 $\mu$ l  | -                   |

### Thermal cycling conditions (three-step)

|         |         |                |
|---------|---------|----------------|
| 94°C    | 30 sec  | } 40-45 cycles |
| 94°C    | 5 sec   |                |
| 50-60°C | 15 sec* |                |
| 72°C    | 10 sec* |                |

Dissociation Stage

### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 94°C | 30 sec  | } 40-45 cycles |
| 94°C | 5 sec   |                |
| 60°C | 30 sec* |                |

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**



# TransScript® II All-in-One First-Strand cDNA Synthesis SuperMix for PCR

AH321-01

50 rxns (20 µl per reaction)

## Storage

at -20°C for one year

## Description

TransScript® II All-in-One First-Strand cDNA Synthesis SuperMix for PCR provides all the necessary components for cDNA synthesis from total RNA or mRNA. The SuperMix is provided at 5× concentration and used at 1× concentration by adding RNA and H<sub>2</sub>O. The resulting cDNA is suitable for regular PCR, not for qPCR.

- One-tube format for simple and fast setup and reduced pipetting variability.
- The optimal ratio of oligo(dT)<sub>20</sub> primer to random primer(N9) for PCR ready cDNA.
- PCR ready cDNA in 30 minutes (unsuitable for qPCR).
- cDNA up to 15 kb.

## Applications

- cDNA library construction, 3' and 5' RACE
- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Kit Contents

| Component                                     | AH321-01 |
|---|----------|
| 5×TransScript® II All-in-One SuperMix for PCR | 200 µl   |
| RNase-free Water                              | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component                                     | Volume              |
|---|---------------------|
| Total RNA/mRNA                                | 50 ng-5 µg/5-500 ng |
| 5×TransScript® II All-in-One SuperMix for PCR | 4 µl                |
| RNase-free Water                              | to 20 µl            |

#### 2. Incubation

- For RNA template with poly(A)<sup>+</sup>, incubate at 50°C for 30 minutes.
- For RNA template without poly(A)<sup>+</sup>, incubate at 25°C for 10 minutes, then at 50°C for 30 minutes.
- For GC-rich or complex secondary structure RNA template, incubate at 55°C for 30 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### RT-PCR

**The suggested reaction condition is the same as described on page 42.**



# TransScript<sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)

AH341-01

50 rxns (20 µl per reaction)

## Storage

at -20°C for one year

## Description

The kit provides all the necessary components for cDNA synthesis from total RNA or mRNA. It is provided at 5× concentration and used at 1× concentration by adding gDNA remover, RNA and H<sub>2</sub>O. Simultaneous genomic DNA removal and cDNA synthesis are performed. After cDNA synthesis, gDNA remover and reverse transcriptase are inactivated by heating at 85°C for 5 seconds. The resulting cDNA is suitable for qPCR, not for regular PCR.

- Simultaneous genomic DNA removal and cDNA synthesis.
- The optimal ratio of Oligo(dT)<sub>20</sub> Primer to random primer(N9) for qPCR ready cDNA.
- qPCR ready cDNA in 15 minutes.
- cDNA up to 250 bp.

## Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Kit Contents

| Component  | AH341-01 |
|--|----------|
| 5× <i>TransScript</i> <sup>®</sup> II All-in-One SuperMix for qPCR               | 200 µl   |
| gDNA Remover   | 50 µl    |
| 5× <i>TransScript</i> <sup>®</sup> II All-in-One No-RT Control SuperMix for qPCR | 20 µl    |
| RNase-free Water   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component  | Volume        |
|--|---------------|
| Total RNA/mRNA   | ≤1 µg/≤100 ng |
| 5× <i>TransScript</i> <sup>®</sup> II All-in-One SuperMix for qPCR | 4 µl          |
| gDNA Remover   | 1 µl          |
| RNase-free Water   | to 20 µl      |

#### 2. Incubate at 50°C for 15 minutes.

For GC-rich or complex secondary structure RNA template, incubate at 55°C for 15 minutes.

#### 3. Incubate at 85°C for 5 seconds to inactivate enzymes.

### qPCR

**The suggested reaction condition is the same as described on page 61.**



# TransScript® Two-Step RT-PCR SuperMix

AT401-01

50 rxns (20 µl per RT reaction)

80 rxns (50 µl per PCR)

## Storage

at -20°C for one year

## Description

*TransScript*® Two-Step RT-PCR SuperMix performs first-strand cDNA synthesis and PCR in two steps. 5×*TransScript*® All-in-One SuperMix for PCR is used for reverse transcription and 2×*TransTaq*® HiFi PCR SuperMix II is used for PCR.

- Amplification of fragment up to 12 kb.

## Application

Multiple copy and low copy gene detection

## Kit Contents

| Component   | AT401-01 |
|---|----------|
| 5× <i>TransScript</i> ® All-in-One SuperMix for PCR | 200 µl   |
| 2× <i>TransTaq</i> ® HiFi PCR SuperMix II           | 2×1 ml   |
| RNase-free Water                                    | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

The suggested reaction condition is the same as described on page 59.

### RT-PCR

Reaction Components

| Component                                 | Volume   | Final Concentration |
|---|----------|---------------------|
| cDNA                                      | 2 µl     | as required         |
| Forward Primer (10 µM)                    | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                    | 1 µl     | 0.2 µM              |
| 2× <i>TransTaq</i> ® HiFi PCR SuperMix II | 25 µl    | 1×                  |
| ddH <sub>2</sub> O                        | to 50 µl | Not applicable      |

The suggested reaction condition is the same as described on page 42.



# TransScript<sup>®</sup> II Two-Step RT-PCR SuperMix

AH401-01

50 rxns (20 µl per RT reaction)

80 rxns (50 µl per PCR)

## Storage

at -20°C for one year

## Description

*TransScript*<sup>®</sup> II Two-Step RT-PCR SuperMix performs first-strand cDNA synthesis and PCR in two steps. *5×TransScript*<sup>®</sup> II All-in-One SuperMix for PCR is used for reverse transcription and *2×TransTaq*<sup>®</sup> HiFi PCR SuperMix II is used for PCR.

- Amplification of fragment up to 15 kb.

## Applications

- cDNA library construction, 3' and 5' RACE
- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

## Kit Contents

| Component  | AH401-01 |
|--|----------|
| <i>5×TransScript</i> <sup>®</sup> II All-in-One SuperMix for PCR | 200 µl   |
| <i>2×TransTaq</i> <sup>®</sup> HiFi PCR SuperMix II              | 2×1 ml   |
| RNase-free Water   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

The suggested reaction condition is the same as described on page 62.

### RT-PCR

Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| cDNA  | 2 µl     | as required         |
| Forward Primer (10 µM)                              | 1 µl     | 0.2 µM              |
| Reverse Primer (10 µM)                              | 1 µl     | 0.2 µM              |
| <i>2×TransTaq</i> <sup>®</sup> HiFi PCR SuperMix II | 25 µl    | 1×                  |
| ddH <sub>2</sub> O                                  | to 50 µl | Not applicable      |

The suggested reaction condition is the same as described on page 42.



# EasyScript® One-Step RT-PCR SuperMix

Mix (+dye) AE411-02 200 rxns (20 µl per reaction)

### Storage

at -20°C for one year

### Description

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. *EasyScript*® RT and *TransTaq*® HiFi DNA Polymerase are used in the kit.

- Amplification of fragment up to 4 kb.

### Application

Multiple copy gene detection

### Kit Contents

| Component                               | AE411-02 |
|---|----------|
| <i>EasyScript</i> ® One-Step Enzyme Mix | 80 µl    |
| 2xOne-Step Reaction Mix                 | 2x1 ml   |
| RNase-free Water                        | 2x1 ml   |

## PROTOCOL


### Reaction Components

| Component                               | Volume    | Final Concentration |
|---|-----------|---------------------|
| RNA Template                            | 1 pg~1 µg | as required         |
| Forward GSP (10 µM)                     | 0.4 µl    | 0.2 µM              |
| Reverse GSP (10 µM)                     | 0.4 µl    | 0.2 µM              |
| 2xOne-Step Reaction Mix                 | 10 µl     | 1x                  |
| <i>EasyScript</i> ® One-Step Enzyme Mix | 0.4 µl    | -                   |
| RNase-free Water                        | to 20 µl  | -                   |


### Thermal cycling conditions

45°C 15-30 min  
 94°C 2-5 min  
 94°C 30 sec  
 50-60°C 30 sec  
 72°C 1-2 kb/min  
 72°C 5-10 min

35-40 cycles



RT-PCR with *EasyScript*® One-Step RT-PCR SuperMix  
 M: 1Kb Plus DNA Ladder  
 1: β-actin 0.5 kb 2: BACH1 1.0 kb  
 3: REPA 1.8 kb 4: ACTR 3.0 kb  
 5: VIN 4.6 kb



RT-PCR with *EasyScript*® One-Step RT-PCR SuperMix to amplify β-actin using human total RNA as templates  
 M: *Trans2K*® DNA Marker  
 Lanes 1-9: 0 pg, 0.1 pg, 1 pg, 10 pg, 100 pg, 1,000 pg, 10 ng, 100 ng, 1,000 ng





# TransScript<sup>®</sup> One-Step RT-PCR SuperMix

|            |          |                               |
|------------|----------|-------------------------------|
| Mix (+dye) | AT411-02 | 200 rxns (20 µl per reaction) |
|------------|----------|-------------------------------|

## Storage

at -20°C for one year

## Description

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. *TransScript<sup>®</sup>* RT and *TransTaq<sup>®</sup>* HiFi DNA Polymerase are used in the kit.

- Amplification of fragment up to 8 kb.

## Application

Multiple copy and low copy gene detection

## Kit Contents

| Component  | AT411-02 |
|--|----------|
| <i>TransScript<sup>®</sup></i> One-Step Enzyme Mix | 80 µl    |
| 2×One-Step Reaction Mix                            | 2×1 ml   |
| RNase-free Water                                   | 2×1 ml   |


## PROTOCOL

### Reaction Components

| Component  | Volume    | Final Concentration |
|--|-----------|---------------------|
| RNA Template                                       | 1 pg~1 µg | as required         |
| Forward GSP (10 µM)                                | 0.4 µl    | 0.2 µM              |
| Reverse GSP (10 µM)                                | 0.4 µl    | 0.2 µM              |
| 2×One-Step Reaction Mix                            | 10 µl     | 1×                  |
| <i>TransScript<sup>®</sup></i> One-Step Enzyme Mix | 0.4 µl    | -                   |
| RNase-free Water                                   | to 20 µl  | -                   |

### Thermal cycling conditions

|         |            |                |
|---------|------------|----------------|
| 45°C    | 15-30 min  |                |
| 94°C    | 2-5 min    |                |
| 94°C    | 30 sec     | } 35-40 cycles |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



# *TransScript*<sup>®</sup> II One-Step RT-PCR SuperMix

Mix (+dye)      AH411-02      200 rxns (20  $\mu$ l per reaction)

### Storage

at -20°C for one year

### Description

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. *TransScript*<sup>®</sup> II RT and *TransTaq*<sup>®</sup> HiFi DNA Polymerase are used in the kit.

- Amplification of fragment up to 8 kb.

### Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

### Kit Contents

| Component  | AH411-02   |
|--|------------|
| <i>TransScript</i> <sup>®</sup> II One-Step Enzyme Mix | 80 $\mu$ l |
| 2xOne-Step Reaction Mix                                | 2x1 ml     |
| RNase-free Water                                       | 2x1 ml     |

High quality products




## PROTOCOL

### Reaction Components

| Component  | Volume    | Final Concentration |
|--|-----------|---------------------|
| RNA Template   | 1 pg~1 µg | as required         |
| Forward GSP (10 µM)                                    | 0.4 µl    | 0.2 µM              |
| Reverse GSP (10 µM)                                    | 0.4 µl    | 0.2 µM              |
| 2×One-Step Reaction Mix                                | 10 µl     | 1×                  |
| <i>TransScript</i> <sup>®</sup> II One-Step Enzyme Mix | 0.4 µl    | -                   |
| RNase-free Water                                       | to 20 µl  | -                   |

### Thermal cycling conditions


|         |            |                |
|---------|------------|----------------|
| 50°C    | 15-30 min  |                |
| 94°C    | 2-5 min    |                |
| 94°C    | 30 sec     | } 35-40 cycles |
| 50-60°C | 30 sec     |                |
| 72°C    | 1-2 kb/min |                |
| 72°C    | 5-10 min   |                |



RT-PCR with *TransScript*<sup>®</sup> II One-Step RT-PCR SuperMix

M: 1Kb Plus DNA Ladder

1: GAPDH 0.9 kb                      5: VIN 4.6 kb  
 2: REPA 1.8 kb                      6: Pol 6.8 kb  
 3: NCBP 2.5 kb                      7: APC 8.5 kb  
 4: HDP 3.5 kb



RT-PCR with *TransScript*<sup>®</sup> II One-Step RT-PCR SuperMix to amplify  $\beta$ -actin using human total RNA as templates

M: *Trans2K*<sup>®</sup> DNA Marker

Lanes 1-9: 0 pg, 0.1 pg, 1 pg, 10 pg, 100 pg, 1,000 pg, 10 ng, 100 ng, 1,000 ng



# Ribonuclease Inhibitor

|          |               |
|----------|---------------|
| AI101-01 | 2,000 units   |
| AI101-02 | 5×2,000 units |

**Concentration**

50 units/μl

**Storage**

at -20°C for one year

**Description**

Ribonuclease Inhibitor is a recombinant protein purified from *E. coli* strain carrying human placenta ribonuclease inhibitor gene. Ribonuclease Inhibitor specifically inhibits RNase A, RNase B, and RNase C. It is not effective against RNase 1, RNase T1, S1 nuclease, RNase H and aspergillus-originated RNase. It has no inhibition effect on DNA Polymerase, AMV, M-MLV, SP6, T7 and T3 RNA Polymerases.

**Unit Definition**

One unit is defined as the amount of enzyme required to inhibit 5 ng RNase A by 50%.

**Applications**

*In vitro* inhibition of ribonuclease, cDNA synthesis and *in vitro* transcription and translation.




# qPCR and qRT-PCR SuperMix


## Basic principle of real-time quantitative PCR

Real-time qPCR is a PCR method used to amplify and simultaneously quantify target DNA molecules. Two methods are frequently used for qPCR: double-strand DNA-binding dyes (e.g. SYBR Green I) or fluorescent reporter probes (e.g. TaqMan). In both cases, fluorescence signals are detected during the exponential phase.

Amplification Plots



Dissociation Curve



## TransStart<sup>®</sup> Green qPCR SuperMix

|          |         |
|----------|---------|
| AQ101-01 | 1 ml    |
| AQ101-02 | 5×1 ml  |
| AQ101-03 | 15×1 ml |

### Contents

- 2×TransStart<sup>®</sup> Green qPCR SuperMix
- Passive Reference Dye (50×)

### Storage

at -20°C in dark for one year

### Description

TransStart<sup>®</sup> Green qPCR SuperMix is a ready-to-use qPCR cocktail containing all components, except primer and template. It contains TransStart<sup>®</sup> Taq DNA Polymerase, SYBR Green I, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primer, passive reference dye (optional) and ddH<sub>2</sub>O.

- TransStart<sup>®</sup> Taq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data.
- Double cation (K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.

### Passive Reference Dye

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50×)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000



- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## PROTOCOL

### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 µM)                                | 0.4 µl   | 0.2 µM              |
| Reverse Primer (10 µM)                                | 0.4 µl   | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix | 10 µl    | 1×                  |
| Passive Reference Dye (50×) (optional)                | 0.4 µl   | 1×                  |
| ddH <sub>2</sub> O                                    | Variable | -                   |
| Total Volume  | 20 µl    | -                   |

### Thermal cycling conditions (three-step)

|         |         |                |
|---------|---------|----------------|
| 94°C    | 30 sec  |                |
| 94°C    | 5 sec   | } 40-45 cycles |
| 50-60°C | 15 sec* |                |
| 72°C    | 10 sec* |                |

Dissociation Stage

### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 94°C | 30 sec  |                |
| 94°C | 5 sec   | } 40-45 cycles |
| 60°C | 30 sec* |                |

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**





# TransStart<sup>®</sup> Green qPCR SuperMix UDG

|          |         |
|----------|---------|
| AQ111-01 | 1 ml    |
| AQ111-02 | 5×1 ml  |
| AQ111-03 | 15×1 ml |

## Contents

- 2×TransStart<sup>®</sup> Green qPCR SuperMix UDG
- Passive Reference Dye (50×)

## Storage

at -20°C in dark for one year

## Description

TransStart<sup>®</sup> Green qPCR SuperMix UDG is a ready-to-use qPCR cocktail containing all components, except primer and template. It contains TransStart<sup>®</sup> Taq DNA Polymerase, UDG, SYBR Green I, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primer, passive reference dye (optional) and ddH<sub>2</sub>O.

- TransStart<sup>®</sup> Taq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data.
- Double cation (K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.
- UDG and dUTP avoid cross contamination.

## Passive Reference Dye

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50×)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## PROTOCOL

### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 μM)                            | 0.4 μl   | 0.2 μM              |
| Reverse Primer (10 μM)                            | 0.4 μl   | 0.2 μM              |
| 2×TransStart <sup>®</sup> Green qPCR SuperMix UDG | 10 μl    | 1×                  |
| Passive Reference Dye (50×) (optional)            | 0.4 μl   | 1×                  |
| ddH <sub>2</sub> O                                | Variable | -                   |
| Total Volume                                      | 20 μl    | -                   |



**Thermal cycling conditions (three-step)**

|         |                           |                |
|---------|---------------------------|----------------|
| 50°C    | 2 min (UDG Incubation)    |                |
| 94°C    | 10 min (UDG Inactivation) |                |
| 94°C    | 5 sec                     | } 40-45 cycles |
| 50-60°C | 15 sec*                   |                |
| 72°C    | 10 sec*                   |                |

Dissociation Stage

**Thermal cycling conditions (two-step)**

|      |                           |                |
|------|---------------------------|----------------|
| 50°C | 2 min (UDG Incubation)    |                |
| 94°C | 10 min (UDG Inactivation) |                |
| 94°C | 5 sec                     | } 40-45 cycles |
| 60°C | 30 sec*                   |                |

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**

# TransStart® Top Green qPCR SuperMix

|          |         |
|----------|---------|
| AQ131-01 | 1 ml    |
| AQ131-02 | 5×1 ml  |
| AQ131-03 | 15×1 ml |
| AQ131-04 | 25×1 ml |

**Contents**

- 2×TransStart® Top Green qPCR SuperMix
- Passive Reference Dye (50×)

**Storage**

at -20°C in dark for one year

**Description**

TransStart® Top Green qPCR SuperMix is a ready-to-use qPCR cocktail containing all components, except primer and template. It contains TransStart® TopTaq DNA Polymerase, SYBR Green I, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primer, passive reference dye (optional) and ddH<sub>2</sub>O.

- TransStart® TopTaq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data.
- Double cation (K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.

**Passive Reference Dye**

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus



- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Plkoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## PROTOCOL

### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 µM)                                    | 0.4 µl   | 0.2 µM              |
| Reverse Primer (10 µM)                                    | 0.4 µl   | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Top Green qPCR SuperMix | 10 µl    | 1×                  |
| Passive Reference Dye (50x) (optional)                    | 0.4 µl   | 1×                  |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total Volume  | 20 µl    | -                   |

The suggested reaction condition is the same as described on page 72.

## *TransStart*<sup>®</sup> Tip Green qPCR SuperMix

|          |         |
|----------|---------|
| AQ141-01 | 1 ml    |
| AQ141-02 | 5×1 ml  |
| AQ141-03 | 15×1 ml |
| AQ141-04 | 25×1 ml |

### Contents

- 2× *TransStart*<sup>®</sup> Tip Green qPCR SuperMix
- Passive Reference Dye (50x)

### Storage

at -20°C in dark for one year

### Description

*TransStart*<sup>®</sup> Tip Green qPCR SuperMix is a ready-to-use qPCR cocktail. It contains a novel *TransStart*<sup>®</sup> *TipTaq* DNA Polymerase, unique hot start reagents (DNA binding proteins combined with unique chemical), optimized double cation buffer, SYBR Green I, dNTPs, PCR Enhancer and PCR stabilizer. qPCR SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primer, passive reference dye (optional) and ddH<sub>2</sub>O.

- A combination of chemical blocking technique with *TransStart*<sup>®</sup> hot start technique to achieve complete blocking. Compared with double blocking *TransStart*<sup>®</sup> *TopTaq*, this method provides higher sensitivity, higher specificity, better amplification.
- Double cation (K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.



### Passive Reference Dye

- Passive Reference Dye I (50x)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## PROTOCOL

### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| Template  | Variable | as required         |
| Forward Primer (10 µM)                                    | 0.4 µl   | 0.2 µM              |
| Reverse Primer (10 µM)                                    | 0.4 µl   | 0.2 µM              |
| 2x <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix | 10 µl    | 1x                  |
| Passive Reference Dye (50x) (optional)                    | 0.4 µl   | 1x                  |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total Volume  | 20 µl    | -                   |

The suggested reaction condition is the same as described on page 72.

## *TransScript*<sup>®</sup> Green Two-Step qRT-PCR SuperMix

|          |                                 |
|----------|---------------------------------|
| AQ201-01 | 50 rxns (20 µl per RT reaction) |
|          | 300 rxns (20 µl per qPCR)       |

### Storage

at -20°C in dark for one year

### Description

*TransScript*<sup>®</sup> Green Two-Step qRT-PCR SuperMix contains all the necessary reagents for gDNA removal, cDNA synthesis and qPCR.

- gDNA remover and 5x *TransScript*<sup>®</sup> All-in-One SuperMix for qPCR are provided for simultaneous gDNA removal and cDNA synthesis.
- *TransStart*<sup>®</sup> Tip Green qPCR SuperMix is provided for qPCR.
- 5x *TransScript*<sup>®</sup> All-in-One No-RT Control SuperMix for qPCR is provided for experimental control.
- Passive reference dyes are provided for different qPCR instruments.

### Application

Multiple copy and low copy gene detection

### Passive Reference Dye

- Passive Reference Dye I (50x)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus



- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

### Kit Contents

| Component   | AQ201-01 |
|---|----------|
| 5× <i>TransScript</i> <sup>®</sup> All-in-One SuperMix for qPCR               | 200 µl   |
| gDNA Remover  | 50 µl    |
| 5× <i>TransScript</i> <sup>®</sup> All-in-One No-RT Control SuperMix for qPCR | 20 µl    |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix                     | 3×1 ml   |
| Passive Reference Dye (50x)   | 120 µl   |
| RNase-free Water  | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume          |
|---|-----------------|
| Total RNA/mRNA  | ≤1 µg / ≤100 ng |
| 5× <i>TransScript</i> <sup>®</sup> All-in-one SuperMix for qPCR | 4 µl            |
| gDNA Remover  | 1 µl            |
| RNase-free Water  | to 20 µl        |

2. Incubate at 42°C for 15 minutes.

3. Incubate at 85°C for 5 seconds to inactivate enzymes.

**The suggested reaction components and condition for qPCR are the same as described on page 72.**



# TransScript® Green miRNA Two-Step qRT-PCR SuperMix

AQ202-01

20 rxns (20 µl per RT reaction)

500 rxns (20 µl per qPCR)

## Storage

at -20°C in dark for one year

## Description

*TransScript®* Green miRNA Two-Step qRT-PCR SuperMix provides all the necessary components for miRNA detection. High efficient poly(A) tail addition and first-strand cDNA synthesis are performed by *TransScript®* miRNA RT Enzyme Mix (containing tailing enzyme and RT enzyme) and 2×TS miRNA Reaction Mix. *TransStart®* Tip Green qPCR SuperMix is provided for miRNA detection.

- One-Step poly(A) tailing and cDNA synthesis.
- Passive reference dyes are provided for different qPCR instruments.

## Application

Multiple copy and low copy gene detection

## Passive Reference Dye

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50×)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## Kit Contents

| Component                                     | AQ202-01 |
|---|----------|
| <i>TransScript®</i> miRNA RT Enzyme Mix       | 20 µl    |
| 2×TS miRNA Reaction Mix                       | 200 µl   |
| Universal miRNA qPCR Primer (10 µM)           | 200 µl   |
| 2× <i>TransStart®</i> Tip Green qPCR SuperMix | 5×1 ml   |
| Passive Reference Dye (50×)                   | 200 µl   |
| RNase-free Water                              | 1 ml     |





## PROTOCOL

### Tail addition and First-Strand cDNA synthesis

#### 1. Reaction Components

| Component   | Volume        |
|---|---------------|
| Total RNA/miRNA*                                    | x $\mu$ l     |
| <i>TransScript</i> <sup>®</sup> miRNA RT Enzyme Mix | 1 $\mu$ l     |
| 2 $\times$ TS miRNA Reaction Mix                    | 10 $\mu$ l    |
| RNase-free Water                                    | to 20 $\mu$ l |

\* Total RNA  $\leq$  5  $\mu$ g. Since miRNA cannot be directly quantified by spectrophotometer, we suggest using 1-9  $\mu$ l for 20  $\mu$ l reaction.

2. Mix gently, and incubate at 37°C for 1 hour.

3. Incubate at 85°C for 5 seconds to inactivate RT Enzyme Mix.

**The suggested reaction components and condition for qPCR are the same as described on page 55-56.**

## *TransScript*<sup>®</sup> II Green Two-Step qRT-PCR SuperMix

|          |                                      |
|----------|--------------------------------------|
| AQ301-01 | 50 rxns (20 $\mu$ l per RT reaction) |
|          | 300 rxns (20 $\mu$ l per qPCR)       |

### Storage

at -20°C in dark for one year

### Description

*TransScript*<sup>®</sup> II Green Two-Step qRT-PCR SuperMix contains all the necessary reagents for gDNA removal, cDNA synthesis and qPCR.

- gDNA remover and 5 $\times$ *TransScript*<sup>®</sup> II All-in-One SuperMix for qPCR are provided for simultaneous gDNA removal and cDNA synthesis.
- *TransStart*<sup>®</sup> Tip Green qPCR SuperMix is provided for qPCR.
- 5 $\times$ *TransScript*<sup>®</sup> II All-in-One No-RT Control SuperMix for qPCR is provided for experimental control.
- Passive reference dyes are provided for different qPCR instruments.

### Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template

### Passive Reference Dye

- Passive Reference Dye I (50 $\times$ )  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50 $\times$ )  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q



### Kit Contents

| Component  | AQ301-01 |
|--|----------|
| 5× <i>TransScript</i> <sup>®</sup> II All-in-One SuperMix for qPCR               | 200 µl   |
| gDNA Remover   | 50 µl    |
| 5× <i>TransScript</i> <sup>®</sup> II All-in-One No-RT Control SuperMix for qPCR | 20 µl    |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix                        | 3×1 ml   |
| Passive Reference Dye (50×)  | 120 µl   |
| RNase-free Water   | 1 ml     |

## PROTOCOL

### First-Strand cDNA synthesis

#### 1. Reaction Components

| Component  | Volume          |
|--|-----------------|
| Total RNA/mRNA   | ≤1 µg / ≤100 ng |
| 5× <i>TransScript</i> <sup>®</sup> II All-in-one SuperMix for qPCR | 4 µl            |
| gDNA Remover   | 1 µl            |
| RNase-free Water   | to 20 µl        |

#### 2. Incubation at 50°C for 15 minutes.

For GC-rich or complex secondary structure RNA template, incubate at 55°C for 15 minutes.

#### 3. Incubate at 85°C for 5 seconds to inactivate enzymes.

**The suggested reaction components and condition for qPCR are the same as described on page 72.**

## *TransScript*<sup>®</sup> Green One-Step qRT-PCR SuperMix

|          |                               |
|----------|-------------------------------|
| AQ211-01 | 100 rxns (20 µl per reaction) |
| AQ211-02 | 400 rxns (20 µl per reaction) |

### Storage

at -20°C in dark for one year

### Description

*TransScript*<sup>®</sup> Green One-Step qRT-PCR SuperMix combines the first-strand cDNA synthesis and qPCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene specific primers can be used for this kit. *TransScript*<sup>®</sup> Green One-Step qRT-PCR SuperMix contains all the necessary reagents for cDNA synthesis and qPCR except total RNA/mRNA template and gene specific primers.

- 5 minutes cDNA synthesis
- cDNA synthesis and qPCR are performed in a single tube using gene specific primers with total RNA or mRNA as templates.
- Passive reference dyes are provided for different qPCR instruments.

### Applications

- Multiple copy and low copy gene detection
- Viral RNA and trace RNA detection



### Passive Reference Dye

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50×)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

### Kit Contents

| Component   | AQ211-01 | AQ211-02 |
|---|----------|----------|
| <i>TransScript</i> <sup>®</sup> One-Step RT/RI Enzyme Mix | 40 µl    | 160 µl   |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix | 1 ml     | 4x1 ml   |
| Passive Reference Dye (50×)                               | 40 µl    | 160 µl   |
| RNase-free Water  | 1 ml     | 4x1 ml   |

## PROTOCOL

### Reaction Components

| Component   | Volume    | Final Concentration |
|---|-----------|---------------------|
| RNA Template  | 1 pg~1 µg | as required         |
| Forward GSP (10 µM)                                       | 0.4 µl    | 0.2 µM              |
| Reverse GSP (10 µM)                                       | 0.4 µl    | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix | 10 µl     | 1×                  |
| <i>TransScript</i> <sup>®</sup> One-Step RT/RI Enzyme Mix | 0.4 µl    | -                   |
| Passive Reference Dye (50×) (optional)                    | 0.4 µl    | 1×                  |
| RNase-free Water  | Variable  | -                   |
| Total volume  | 20 µl     | -                   |

### Thermal cycling conditions (three-step)

|                    |         |                |
|--------------------|---------|----------------|
| 45°C               | 5 min   |                |
| 94°C               | 30 sec  |                |
| 94°C               | 5 sec   | } 40-45 cycles |
| 50-60°C            | 15 sec* |                |
| 72°C               | 10 sec* |                |
| Dissociation Stage |         |                |

### Thermal cycling conditions (two-step)

|                    |         |                |
|--------------------|---------|----------------|
| 45°C               | 5 min   |                |
| 94°C               | 30 sec  |                |
| 94°C               | 5 sec   | } 40-45 cycles |
| 60°C               | 30 sec* |                |
| Dissociation Stage |         |                |



Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**

## TransScript<sup>®</sup> II Green One-Step qRT-PCR SuperMix

|          |                               |
|----------|-------------------------------|
| AQ311-01 | 100 rxns (20 µl per reaction) |
| AQ311-02 | 400 rxns (20 µl per reaction) |

### Storage

at -20°C in dark for one year

### Description

*TransScript<sup>®</sup> II Green One-Step qRT-PCR SuperMix* combines the first-strand cDNA synthesis and qPCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene specific primers can be used for this kit. *TransScript<sup>®</sup> II Green One-Step qRT-PCR SuperMix* contains all the necessary reagents for cDNA synthesis and qPCR except total RNA/mRNA template and gene specific primers.

- 5 minutes cDNA synthesis
- cDNA synthesis and qPCR are performed in a single tube using gene specific primers with total RNA or mRNA as templates.
- Passive reference dyes are provided for different qPCR instruments.

### Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template
- Viral RNA and trace RNA detection

### Passive Reference Dye

- Passive Reference Dye I (50x)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q



### Kit Contents

| Component  | AQ311-01 | AQ311-02 |
|--|----------|----------|
| <i>TransScript</i> <sup>®</sup> II One-Step RT/RI Enzyme Mix | 40 µl    | 160 µl   |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix    | 1 ml     | 4×1 ml   |
| Passive Reference Dye (50×)                                  | 40 µl    | 160 µl   |
| RNase-free Water   | 1 ml     | 4×1 ml   |

## PROTOCOL

### Reaction Components

| Component  | Volume    | Final Concentration |
|--|-----------|---------------------|
| RNA Template   | 1 pg~1 µg | as required         |
| Forward GSP (10 µM)  | 0.4 µl    | 0.2 µM              |
| Reverse GSP (10 µM)  | 0.4 µl    | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix    | 10 µl     | 1×                  |
| <i>TransScript</i> <sup>®</sup> II One-Step RT/RI Enzyme Mix | 0.4 µl    | -                   |
| Passive Reference Dye (50×) (optional)                       | 0.4 µl    | 1×                  |
| RNase-free Water   | Variable  | -                   |
| Total volume   | 20 µl     | -                   |

### Thermal cycling conditions (three-step)

|         |         |                |
|---------|---------|----------------|
| 50°C    | 5 min   |                |
| 94°C    | 30 sec  |                |
| 94°C    | 5 sec   | } 40-45 cycles |
| 50-60°C | 15 sec* |                |
| 72°C    | 10 sec* |                |

Dissociation Stage

### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 50°C | 5 min   |                |
| 94°C | 30 sec  |                |
| 94°C | 5 sec   | } 40-45 cycles |
| 60°C | 30 sec* |                |

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

**Two-step qPCR is more suitable for higher specificity assay.**

**Three-step qPCR is more suitable for higher sensitivity assay.**



# TransStart<sup>®</sup> Probe qPCR SuperMix

|          |         |
|----------|---------|
| AQ401-01 | 1 ml    |
| AQ401-02 | 5x1 ml  |
| AQ401-03 | 15x1 ml |

## Contents

- 2x TransStart<sup>®</sup> Probe qPCR SuperMix
- Passive Reference Dye (50x)

## Storage

at -20°C for one year

## Description

TransStart<sup>®</sup> Probe qPCR SuperMix is a ready-to-use qPCR cocktail containing all components, except probe, primer and template. It contains TransStart<sup>®</sup> Taq DNA Polymerase, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2x concentration and can be used at 1x concentration by adding template, primer, probe, passive reference dye (optional) and ddH<sub>2</sub>O.

- TransStart<sup>®</sup> Taq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data.
- Double cation (K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.

## Passive Reference Dye

- Passive Reference Dye I (50x)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## PROTOCOL

### Reaction Components

| Component                                      | Volume   | Final Concentration |
|--|----------|---------------------|
| Template                                       | Variable | as required         |
| Forward Primer (10 μM)                         | 0.4 μl   | 0.2 μM              |
| Reverse Primer (10 μM)                         | 0.4 μl   | 0.2 μM              |
| Probe (10 μM)                                  | 0.4 μl   | 0.2 μM              |
| 2x TransStart <sup>®</sup> Probe qPCR SuperMix | 10 μl    | 1x                  |
| Passive Reference Dye (50x) (optional)         | 0.4 μl   | 1x                  |
| ddH <sub>2</sub> O                             | Variable | -                   |
| Total volume                                   | 20 μl    | -                   |



**Thermal cycling conditions (two-step)**

|      |         |                |
|------|---------|----------------|
| 94°C | 30 sec  | } 40-45 cycles |
| 94°C | 5 sec   |                |
| 60°C | 30 sec* |                |

For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.

## TransScript<sup>®</sup> Probe One-Step qRT-PCR SuperMix

|          |                               |
|----------|-------------------------------|
| AQ221-01 | 100 rxns (20 µl per reaction) |
| AQ221-02 | 400 rxns (20 µl per reaction) |

**Storage**

at -20°C for one year

**Description**

*TransScript<sup>®</sup> Probe One-Step qRT-PCR SuperMix* combines the first-strand cDNA synthesis and qPCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene specific primers can be used for this kit. *TransScript<sup>®</sup> Probe One-Step qRT-PCR SuperMix* contains all the necessary reagents for cDNA synthesis and qPCR except probe, total RNA/mRNA template and gene specific primers.

- 5 minutes cDNA synthesis
- cDNA synthesis and qPCR are performed in a single tube using gene specific primers with total RNA or mRNA as templates.
- Passive reference dyes are provided for different qPCR instruments.

**Applications**

- Multiple copy and low copy gene detection
- Viral RNA and trace RNA detection

**Passive Reference Dye**

- Passive Reference Dye I (50x)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50x)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q



### Kit Contents

| Component   | AQ221-01 | AQ221-02 |
|---|----------|----------|
| <i>TransScript</i> <sup>®</sup> One-Step RT/RI Enzyme Mix | 40 µl    | 160 µl   |
| 2× <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix     | 1 ml     | 4×1 ml   |
| Passive Reference Dye (50×)                               | 40 µl    | 160 µl   |
| RNase-free Water  | 1 ml     | 4×1 ml   |

## PROTOCOL

### Reaction Components

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| RNA Template  | Variable | as required         |
| Forward GSP (10 µM)                                       | 0.4 µl   | 0.2 µM              |
| Reverse GSP (10 µM)                                       | 0.4 µl   | 0.2 µM              |
| Probe (10 µM)   | 0.4 µl   | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix     | 10 µl    | 1×                  |
| <i>TransScript</i> <sup>®</sup> One-Step RT/RI Enzyme Mix | 0.4 µl   | -                   |
| Passive Reference Dye (50×) (optional)                    | 0.4 µl   | 1×                  |
| RNase-free Water  | Variable | -                   |
| Total volume  | 20 µl    | -                   |

### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 45°C | 5 min   |                |
| 94°C | 30 sec  |                |
| 94°C | 5 sec   | } 40-45 cycles |
| 60°C | 30 sec* |                |

For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.



# TransScript<sup>®</sup> II Probe One-Step qRT-PCR SuperMix

|          |                               |
|----------|-------------------------------|
| AQ321-01 | 100 rxns (20 µl per reaction) |
| AQ321-02 | 400 rxns (20 µl per reaction) |

## Storage

at -20°C for one year

## Description

*TransScript*<sup>®</sup> II Probe One-Step qRT-PCR SuperMix combines the first-strand cDNA synthesis and qPCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene specific primers can be used for this kit. *TransScript*<sup>®</sup> II Probe One-Step qRT-PCR SuperMix contains all the necessary reagents for cDNA synthesis and qPCR except probe, total RNA/mRNA template and gene specific primers.

- 5 minutes cDNA synthesis
- cDNA synthesis and qPCR are performed in a single tube using gene specific primers with total RNA or mRNA as templates.
- Passive reference dyes are provided for different qPCR instruments.

## Applications

- Multiple copy and low copy gene detection
- GC-rich or complex secondary structure RNA template
- Viral RNA and trace RNA detection

## Passive Reference Dye

- Passive Reference Dye I (50×)  
ABI Prism7000/7300/7700/7900, Eppendorf, ABI Step One, ABI Step One Plus
- Passive Reference Dye II (50×)  
ABI Prism7500, ABI Prism7500 Fast, ABI Q6, ABI Quant Studio 6/7 Flex, ABI ViiA 7, Stratagene Mx3000/Mx3005P, Corbett Rotor Gene 3000
- No Passive Reference Dye  
Roche LightCycler480, Roche Light Cycler96, MJ Research Chromo4, Opticon (II), Takara TP800, Bio-Rad iCycler iQ, iQ5, Bio-Rad CFX96, Bio-Rad C1000 Thermal Cycler, Thermo Pikoreal 96, Corbett Rotor Gene 6000, Corbett Rotor Gene G, Corbett Rotor Gene Q

## Kit Contents

| Component  | AQ321-01 | AQ321-02 |
|--|----------|----------|
| <i>TransScript</i> <sup>®</sup> II One-Step RT/RI Enzyme Mix | 40 µl    | 160 µl   |
| 2× <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix        | 1 ml     | 4×1 ml   |
| Passive Reference Dye (50×)                                  | 40 µl    | 160 µl   |
| RNase-free Water   | 1 ml     | 4×1 ml   |



## PROTOCOL

### Reaction Components

| Component  | Volume   | Final Concentration |
|--|----------|---------------------|
| RNA Template   | Variable | as required         |
| Forward GSP (10 µM)  | 0.4 µl   | 0.2 µM              |
| Reverse GSP (10 µM)  | 0.4 µl   | 0.2 µM              |
| Probe (10 µM)  | 0.4 µl   | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix        | 10 µl    | 1×                  |
| <i>TransScript</i> <sup>®</sup> II One-Step RT/RI Enzyme Mix | 0.4 µl   | -                   |
| Passive Reference Dye (50×) (optional)                       | 0.4 µl   | 1×                  |
| RNase-free Water   | Variable | -                   |
| Total volume   | 20 µl    | -                   |

### Thermal cycling conditions (two-step)

|      |         |                |
|------|---------|----------------|
| 50°C | 5 min   |                |
| 94°C | 30 sec  |                |
| 94°C | 5 sec   | } 40-45 cycles |
| 60°C | 30 sec* |                |

For ABI qPCR instrument, we suggest using the following signal collecting time:

- \* For ABI Prism7700/7900, the time is 30 seconds.
- \* For ABI Prism7000/7300, the time is 31 seconds.
- \* For ABI Prism7500, the time is 34 seconds.
- \* For ABI ViiA 7, the time is at least 19 seconds.



## High Pure dNTPs

|        |          |        |
|--------|----------|--------|
| 2.5 mM | AD101-01 | 1 ml   |
|        | AD101-02 | 5x1 ml |
| 10 mM  | AD101-11 | 1 ml   |
|        | AD101-12 | 5x1 ml |

**Storage**

at -20°C for two years

**Description**

High Pure dNTPs is an equal molar solution of high quality dATP, dCTP, dGTP, and dTTP with purity up to 99%. It is suitable for PCR, qPCR, DNA sequencing and cDNA synthesis.

---

## Chapter 2 DNA Molecular Weight Standards

---

|  |     |
|--|-----|
| <i>Trans2K</i> <sup>®</sup> DNA Marker .....         | 092 |
| <i>Trans2K</i> <sup>®</sup> Plus DNA Marker .....    | 092 |
| <i>Trans2K</i> <sup>®</sup> Plus II DNA Marker ..... | 092 |
| <i>Trans5K</i> DNA Marker .....                      | 093 |
| <i>Trans15K</i> DNA Marker .....                     | 093 |
| 1Kb DNA Ladder .....                                 | 093 |
| 1Kb Plus DNA Ladder.....                             | 094 |
| 100bp DNA Ladder .....                               | 094 |
| 100bp Plus DNA Ladder .....                          | 094 |
| 100bp Plus II DNA Ladder .....                       | 095 |
| GelStain .....                                       | 095 |
| Agarose .....  | 095 |















TransGen provides a broad range of double-strand DNA molecular weight markers for conventional electrophoresis. All DNA markers are generated from restriction enzymes digested plasmids. All DNA markers are in ready-to-use format.

DNA Marker Selection Guide

| DNA Marker                                     | Agarose | DNA Marker               | Agarose |
|--|---------|--------------------------|---------|
| <i>Trans2K</i> <sup>®</sup> DNA Marker         | 1.5%    | 1Kb DNA Ladder           | 1.0%    |
| <i>Trans2K</i> <sup>®</sup> Plus DNA Marker    | 1.5%    | 1Kb Plus DNA Ladder      | 1.0%    |
| <i>Trans2K</i> <sup>®</sup> Plus II DNA Marker | 1.5%    | 100bp DNA Ladder         | 2.0%    |
| <i>Trans5K</i> DNA Marker                      | 1.5%    | 100bp Plus DNA Ladder    | 2.0%    |
| <i>Trans15K</i> DNA Marker                     | 0.8%    | 100bp Plus II DNA Ladder | 2.0%    |



## Trans2K<sup>®</sup> DNA Marker

|          |          |
|----------|----------|
| BM101-01 | 500 µl   |
| BM101-02 | 5×500 µl |

### Concentration

0.07 mg/ml


### Band Size

100 bp, 250 bp, 500 bp, 750 bp (100 ng/5 µl, the double intensity band), 1,000 bp, 2,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

Trans2K<sup>®</sup> DNA Marker



1.5% TAE agarose gel stained with GelStain (loading volume: 5 µl)

## Trans2K<sup>®</sup> Plus DNA Marker

|          |          |
|----------|----------|
| BM111-01 | 500 µl   |
| BM111-02 | 5×500 µl |

### Concentration

0.09 mg/ml


### Band Size

100 bp, 250 bp, 500 bp, 750 bp (100 ng/5 µl, the double intensity band), 1,000 bp, 2,000 bp, 3,000 bp, 5,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

Trans2K<sup>®</sup> Plus DNA Marker



1.5% TAE agarose gel stained with GelStain (loading volume: 5 µl)

## Trans2K<sup>®</sup> Plus II DNA Marker

|          |          |
|----------|----------|
| BM121-01 | 500 µl   |
| BM121-02 | 5×500 µl |

### Concentration

0.10 mg/ml


### Band Size

100 bp, 250 bp, 500 bp, 750 bp (100 ng/5 µl, the double intensity band), 1,000 bp, 2,000 bp, 3,000 bp, 5,000 bp, 8,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

Trans2K<sup>®</sup> Plus II DNA Marker



1.5% TAE agarose gel stained with GelStain (loading volume: 5 µl)

High quality products



## Trans5K DNA Marker

|          |                        |
|----------|------------------------|
| BM141-01 | 500 $\mu$ l            |
| BM141-02 | 5 $\times$ 500 $\mu$ l |

### Concentration

0.095 mg/ml


### Band Size

300 bp, 500 bp, 800 bp, 1,000 bp, 1,500 bp (125 ng/5  $\mu$ l, the double intensity band), 2,000 bp, 3,000 bp, 5,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

Trans5K DNA Marker



1.5% TAE agarose gel stained with GelStain (loading volume: 5  $\mu$ l)

## Trans15K DNA Marker

|          |                        |
|----------|------------------------|
| BM161-01 | 500 $\mu$ l            |
| BM161-02 | 5 $\times$ 500 $\mu$ l |

### Concentration

0.09 mg/ml


### Band Size

500 bp, 1,000 bp, 1,500 bp, 3,000 bp, 5,000 bp (100 ng/5  $\mu$ l, the double intensity band), 7,500 bp, 10,000 bp, 15,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

Trans15K DNA Marker



0.8% TAE agarose gel stained with GelStain (loading volume: 5  $\mu$ l)

## 1Kb DNA Ladder

|          |                        |
|----------|------------------------|
| BM201-01 | 500 $\mu$ l            |
| BM201-02 | 5 $\times$ 500 $\mu$ l |

### Concentration

0.09 mg/ml


### Band Size

1,000 bp, 2,000 bp, 3,000 bp, 4,000 bp (100 ng/5  $\mu$ l, the double intensity band), 5,000 bp, 6,000 bp, 8,000 bp, 10,000 bp.

### Storage

at 4°C for six months; at -20°C for two years

1Kb DNA Ladder



1.0% TAE agarose gel stained with GelStain (loading volume: 5  $\mu$ l)

# 1Kb Plus DNA Ladder

|          |          |
|----------|----------|
| BM211-01 | 500 µl   |
| BM211-02 | 5x500 µl |

**Concentration**

0.13 mg/ml


**Band Size**

300 bp, 500 bp, 800 bp, 1,000 bp,  
1,500 bp, 2,000 bp (100 ng/5 µl, the double intensity band), 3,000 bp, 4,000 bp,  
5,000 bp, 6,000 bp, 8,000 bp, 10,000 bp.

**Storage**

at 4°C for six months; at -20°C for two years

1Kb Plus DNA Ladder



1.0% TAE agarose gel stained with GelStain (loading volume: 5 µl)

# 100bp DNA Ladder

|          |          |
|----------|----------|
| BM301-01 | 500 µl   |
| BM301-02 | 5x500 µl |

**Concentration**

0.12 mg/ml


**Band Size**

100 bp, 200 bp, 300 bp, 400 bp,  
500 bp (100 ng/5 µl, the double intensity band),  
600 bp, 700 bp, 800 bp,  
900 bp, 1,000 bp, 1,500 bp.

**Storage**

at 4°C for six months; at -20°C for two years

100bp DNA Ladder



2.0% TAE agarose gel stained with GelStain (loading volume: 5 µl)

# 100bp Plus DNA Ladder

|          |          |
|----------|----------|
| BM311-01 | 500 µl   |
| BM311-02 | 5x500 µl |

**Concentration**

0.13 mg/ml


**Band Size**

100 bp, 200 bp, 300 bp, 400 bp,  
500 bp (100 ng/5 µl, the double intensity band),  
600 bp, 700 bp, 800 bp,  
900 bp, 1,000 bp, 1,200 bp, 1,500 bp.

**Storage**

at 4°C for six months; at -20°C for two years

100bp Plus DNA Ladder



2.0% TAE agarose gel stained with GelStain (loading volume: 5 µl)

High quality products



## 100bp Plus II DNA Ladder

|          |                        |
|----------|------------------------|
| BM321-01 | 500 $\mu$ l            |
| BM321-02 | 5 $\times$ 500 $\mu$ l |

### Concentration


0.15 mg/ml

### Band Size

100 bp, 200 bp, 300 bp, 400 bp,  
500 bp (100 ng/5  $\mu$ l, the double intensity band),  
600 bp, 700 bp, 800 bp,  
900 bp, 1,000 bp, 1,500 bp, 2,000 bp,  
3,000 bp, 5,000 bp.

### Storage

at 4°C for six months; at -20°C for two years



## GelStain

|          |             |
|----------|-------------|
| GS101-01 | 500 $\mu$ l |
| GS101-02 | 1 ml        |

### Concentration

10,000 $\times$

### Storage

at 4°C in dark for one year

### Description

GelStain is a sensitive, stable and safe staining reagent for DNA/RNA. GelStain uses the same wavelength as ethidium bromide (EB), and it is more sensitive than EB.

### Highlights

- No toxicity: GelStain is a specific form of oily macromolecules, which is incapable of entering cells via the cell membrane.
- High sensitivity: GelStain provides high sensitivity, which can detect low amount of DNA even at 10-20 ng.
- Exceptional stability: GelStain can be heated or microwaved.
- Signal to noise ratio: strong fluorescent signal from samples, weak from background.
- Like EB, GelStain can be used before electrophoresis or after electrophoresis. No destaining is needed.
- No optical setting change: standard EB filter and SYBR filter can be used.

## Agarose

|          |       |
|----------|-------|
| GS201-01 | 100 g |
|----------|-------|

### Storage

at room temperature for two years

### Description

Extremely pure, molecular biology grade Agarose from TransGen is free of DNase, RNase and protease. This product is suitable for routine analysis of nucleic acids by gel electrophoresis and blotting.

| % of Agarose | Resolution (bp) |
|--------------|-----------------|
| 0.5%         | 1,000 ~ 30,000  |
| 0.7%         | 800 ~ 12,000    |
| 1.0%         | 500 ~ 10,000    |
| 1.2%         | 400 ~ 7,000     |
| 1.5%         | 200 ~ 3,000     |
| 2.0%         | 50 ~ 2,000      |

---

## Chapter 3 Cloning and Mutagenesis System

---

### Cloning Vectors

|  |     |
|--|-----|
| <i>pEASY</i> <sup>®</sup> -T1 Cloning Kit .....                        | 99  |
| <i>pEASY</i> <sup>®</sup> -Blunt Cloning Kit .....                     | 102 |
| <i>pEASY</i> <sup>®</sup> -T1 Simple Cloning Kit .....                 | 103 |
| <i>pEASY</i> <sup>®</sup> -Blunt Simple Cloning Kit .....              | 104 |
| <i>pEASY</i> <sup>®</sup> -T3 Cloning Kit .....                        | 105 |
| <i>pEASY</i> <sup>®</sup> -Blunt3 Cloning Kit .....                    | 106 |
| <i>pEASY</i> <sup>®</sup> -T5 Zero Cloning Kit .....                   | 107 |
| <i>pEASY</i> <sup>®</sup> -Blunt Zero Cloning Kit .....                | 108 |
| <i>pEASY</i> <sup>®</sup> -Uni Seamless Cloning and Assembly Kit ..... | 109 |









### Cloning Competent Cells

|  |     |
|--|-----|
| <i>Trans10</i> Chemically Competent Cell .....                   | 112 |
| <i>Trans5<math>\alpha</math></i> Chemically Competent Cell ..... | 112 |
| <i>Trans109</i> Chemically Competent Cell .....                  | 113 |
| <i>Trans110</i> Chemically Competent Cell .....                  | 113 |
| <i>Trans1-Blue</i> Chemically Competent Cell .....               | 113 |
| <i>Trans2-Blue</i> Chemically Competent Cell .....               | 114 |
| <i>Trans1-T1</i> Phage Resistant Chemically Competent Cell ..... | 114 |
| DMT Chemically Competent Cell .....                              | 114 |
| <i>TransStbl3</i> Chemically Competent Cell .....                | 115 |
| <i>TransDB3.1</i> Chemically Competent Cell .....                | 115 |

### Mutagenesis System

|  |     |
|--|-----|
| <i>Fast</i> Mutagenesis System .....           | 116 |
| <i>Fast MultiSite</i> Mutagenesis System ..... | 117 |





## Advantage of *pEASY*<sup>®</sup> vectors

- Fast: 5 minutes at room temperature.
- Simple: only add PCR products.
- High efficient: up to 90% clones with correct insert.

## Feature and application of *pEASY*<sup>®</sup> cloning vectors (MCS=multi-cloning site)

| Name                                    | Amp <sup>+</sup> | Kan <sup>+</sup> | <i>In vitro</i> transcription | Sequencing primer  | Characteristics  | Application   |
|---|------------------|------------------|-------------------------------|--|--|---------------|
| <i>pEASY</i> <sup>®</sup> -T1           | +                | +                | T7 Promoter                   | M13 Forward Primer;<br>M13 Reverse Primer;<br>T7 Promoter                  | Dual resistance, MCS   | TA cloning    |
| <i>pEASY</i> <sup>®</sup> -Blunt        | +                | +                | T7 Promoter                   | M13 Forward Primer;<br>M13 Reverse Primer;<br>T7 Promoter                  | Dual resistance, MCS   | Blunt cloning |
| <i>pEASY</i> <sup>®</sup> -T1 Simple    | +                | +                | T7 Promoter                   | M13 Forward Primer;<br>SR Primer   | Dual resistance, No MCS  | TA cloning    |
| <i>pEASY</i> <sup>®</sup> -Blunt Simple | +                | +                | T7 Promoter                   | M13 Forward Primer;<br>SR Primer   | Dual resistance, No MCS  | Blunt cloning |
| <i>pEASY</i> <sup>®</sup> -T3           | +                | -                | T7/SP6 Promoter               | M13 Forward Primer;<br>M13 Reverse Primer;<br>T7 Promoter;<br>SP6 Promoter | Dual <i>EcoR</i> I, Dual <i>Not</i> I restriction enzyme cut sites | TA cloning    |
| <i>pEASY</i> <sup>®</sup> -Blunt3       | +                | -                | T7/SP6 Promoter               | M13 Forward Primer;<br>M13 Reverse Primer;<br>T7 Promoter;<br>SP6 Promoter | Dual <i>EcoR</i> I, Dual <i>Not</i> I restriction enzyme cut sites | Blunt cloning |
| <i>pEASY</i> <sup>®</sup> -T5 Zero      | +                | +                | T3/T7 Promoter                | M13 Forward Primer;<br>M13 Reverse Primer                                  | Dual resistance, Zero background                                   | TA cloning    |
| <i>pEASY</i> <sup>®</sup> -Blunt Zero   | +                | +                | T3/T7 Promoter                | M13 Forward Primer;<br>M13 Reverse Primer                                  | Dual resistance, Zero background                                   | Blunt cloning |

## General notes for cloning using *pEASY*<sup>®</sup> vectors

- Do not add 5' phosphates to the PCR primers. PCR products with 5' phosphates will not be cloned into *pEASY*<sup>®</sup> vector.
- Choose the right PCR enzymes for TA cloning or blunt cloning.
- To clone diluted PCR products, increase the amount of PCR products or concentrate the PCR products.
- To clone PCR products with multi-bands, gel purify the products before cloning.
- Cloning reaction time cannot be more than 30 minutes.



# pEASY<sup>®</sup>-T1 Cloning Kit

|          |         |
|----------|---------|
| CT101-01 | 20 rxns |
| CT101-02 | 60 rxns |

## Storage

*Trans1*-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description


pEASY<sup>®</sup>-T1 Cloning Kit is designed for cloning and sequencing *Taq*-amplified PCR products.

- 5 minutes fast ligation of *Taq*-amplified PCR products.
- Kanamycin and Ampicillin resistance genes for selection.
- Easy blue/white selection.
- T7 promoter, M13 forward and M13 reverse primers for sequencing.
- T7 promoter for *in vitro* transcription.
- *Trans1*-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.

## Kit Contents

| Component   | CT101-01  | CT101-02  |
|---|-----------|-----------|
| pEASY <sup>®</sup> -T1 Cloning Vector (10 ng/μl)            | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                  | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                     | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                  | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                  | 50 μl     | 150 μl    |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

## pEASY<sup>®</sup>-T1 Cloning Vector Map



LacZα fragment: bases 1-544

M13 reverse priming site: bases 205-221

Multiple cloning site: bases 234-354

T7 promoter priming site: bases 361-380


M13 forward priming site: bases 387-403

f1 origin: bases 545-982

Kanamycin resistance ORF: bases 1,316-2,110

Ampicillin resistance ORF: bases 2,128-2,988

pUC origin: bases 3,133-3,806



## PROTOCOL

### Suggested cloning reaction condition

- Optimal amount of insert  
Molar ratio of vector to insert = 1:7 (1 kb, ~20 ng; 2 kb, ~40 ng)
- Optimal volume of vector: 1  $\mu$ l
- Optimal reaction volume: 3~5  $\mu$ l
- Optimal incubation time
  - 0.1~1 kb (including 1 kb): 5~10 minutes
  - 1~2 kb (including 2 kb): 10~15 minutes
  - 2~3 kb (including 3 kb): 15~20 minutes
  - $\geq$ 3 kb: 20~30 minutes  
Use the maximum incubation time if the insert is gel purified PCR product.
- Optimal incubation temperature: for most PCR inserts, the optimal temperature is about 25°C; for some PCR inserts, optimal results can be achieved with higher temperature (up to 37°C).

### Transformation

- Add the ligated products to 50  $\mu$ l of *Trans1*-T1 Phage Resistant Chemically Competent Cell and mix gently (do not mix by pipetting up and down).
- Incubate on ice for 20~30 minutes.
- Heat-shock the cells at 42°C for 30 seconds.
- Immediately place the tube on ice for 2 minutes.
- Add 250  $\mu$ l of room temperature SOC or LB medium. Shake the tube at 37°C (200 rpm) for 1 hour.
- In the meantime, mix 8  $\mu$ l of 500 mM IPTG with 40  $\mu$ l of 20 mg/ml X-gal. Spread them evenly onto a selective LB plate. Place the plate at 37°C for 30 minutes.
- Spread 200  $\mu$ l or all transformants on the pre-warmed plate. Incubate at 37°C overnight.

### Analysis of positive clones

- Transfer 5~10 white or light blue colonies into 10  $\mu$ l ddH<sub>2</sub>O and vortex.
- Use 1  $\mu$ l of the mixture as template for 25  $\mu$ l PCR using M13 forward and M13 reverse primers.
- PCR
 

|      |          |   |           |
|------|----------|---|-----------|
| 94°C | 10 min   | } | 30 cycles |
| 94°C | 30 sec   |   |           |
| 55°C | 30 sec   |   |           |
| 72°C | x min*   |   |           |
| 72°C | 5-10 min |   |           |

\* (depends on the insert size and PCR enzymes)

the PCR product size from vector self-ligation is 199 bp.

- Analyze positive clones by restriction enzyme digestion and DNA sequencing.

**PCR for control insert (700 bp)**

| Component                  | Volume   | Final Concentration |
|----------------------------|----------|---------------------|
| Control Template (5 ng/μl) | 1 μl     | 0.1 ng/μl           |
| Control Primers (10 μM)    | 1 μl     | 0.2 μM              |
| 2×EasyTaq® PCR SuperMix    | 25 μl    | 1×                  |
| ddH <sub>2</sub> O         | Variable | -                   |
| Total volume               | 50 μl    | -                   |

**Thermal cycling conditions**

|      |         |             |
|------|---------|-------------|
| 94°C | 2-5 min | } 30 cycles |
| 94°C | 30 sec  |             |
| 55°C | 30 sec  |             |
| 72°C | 1 min   |             |
| 72°C | 10 min  |             |

Ligate 1 μl of control PCR insert with 1 μl vector. Hundreds of colonies should be produced with cloning efficiency over 90%.

**General notes for cloning are the same as described on page 98.**

# pEASY<sup>®</sup>-Blunt Cloning Kit

|          |         |
|----------|---------|
| CB101-01 | 20 rxns |
| CB101-02 | 60 rxns |

### Storage

*Trans1*-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

### Description


pEASY<sup>®</sup>-Blunt Cloning Kit is designed for cloning and sequencing *Pfu*-amplified PCR products.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- Kanamycin and Ampicillin resistance genes for selection.
- Easy blue/white selection.
- T7 promoter, M13 forward and M13 reverse primers for sequencing.
- T7 promoter for *in vitro* transcription.
- *Trans1*-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.


### Kit Contents

| Component   | CB101-01  | CB101-02  |
|---|-----------|-----------|
| pEASY <sup>®</sup> -Blunt Cloning Vector (10 ng/μl)         | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                  | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                     | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                  | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                  | 50 μl     | 150 μl    |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

### pEASY<sup>®</sup>-Blunt Cloning Vector Map



LacZα fragment: bases 1-545  
 Multiple cloning site: bases 234-355  
 M13 reverse priming site: bases 205-221  
 T7 promoter priming site: bases 362-381  
 M13 forward priming site: bases 388-404  
 f1 origin: bases 546-983  
 Kanamycin resistance ORF: bases 1,317-2,111  
 Ampicillin resistance ORF: bases 2,129-2,989  
 pUC origin: bases 3,134-3,807



## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100, except the PCR product size from vector self-ligation is 200 bp. General notes for cloning are the same as described on page 98.



# pEASY<sup>®</sup>-T1 Simple Cloning Kit

|          |         |
|----------|---------|
| CT111-01 | 20 rxns |
| CT111-02 | 60 rxns |

## Storage

*Trans1-T1* Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description


pEASY<sup>®</sup>-T1 Simple Cloning Vector eliminates the multi-cloning sites of pEASY<sup>®</sup>-T1 Cloning Vector. It is designed for cloning and sequencing *Taq*-amplified PCR products.

- 5 minutes fast ligation of *Taq*-amplified PCR products.
- Kanamycin and Ampicillin resistance genes for selection.
- Easy blue/white selection.
- SR primer and M13 forward primer for sequencing.
- T7 promoter for *in vitro* transcription.
- *Trans1-T1* Phage Resistant Chemically Competent Cell, high transformation efficiency ( $>10^9$  cfu/μg pUC19 DNA) and fast growing.


## Kit Contents

| Component  | CT111-01  | CT111-02  |
|--|-----------|-----------|
| pEASY <sup>®</sup> -T1 Simple Cloning Vector (10 ng/μl)    | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                 | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                    | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                 | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                 | 50 μl     | 150 μl    |
| SR Primer (10 μM)  | 50 μl     | 150 μl    |
| <i>Trans1-T1</i> Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

## pEASY<sup>®</sup>-T1 Simple Cloning Vector Map



*LacZ* fragment: bases 1-445  
M13 reverse priming site: bases 205-221  
T7 promoter priming site: bases 262-281  
M13 forward priming site: bases 288-304  
f1 origin: bases 446-883  
Kanamycin resistance ORF: bases 1,217-2,011  
Ampicillin resistance ORF: bases 2,029-2,889  
pUC origin: bases 3,034-3,707



## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100, except the PCR product size from vector self-ligation is 100 bp. General notes for cloning are the same as described on page 98.



# pEASY<sup>®</sup>-Blunt Simple Cloning Kit

|          |         |
|----------|---------|
| CB111-01 | 20 rxns |
| CB111-02 | 60 rxns |

### Storage

*Trans1*-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

### Description


pEASY<sup>®</sup>-Blunt Simple Cloning Vector eliminates the multi-cloning sites of pEASY<sup>®</sup>-Blunt Cloning Vector. It is designed for cloning and sequencing *Pfu*-amplified PCR products.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- Kanamycin and Ampicillin resistance genes for selection.
- Easy blue/white selection.
- SR primer and M13 forward primer for sequencing.
- T7 promoter for *in vitro* transcription.
- *Trans1*-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.

### Kit Contents

| Component   | CB111-01  | CB111-02  |
|---|-----------|-----------|
| pEASY <sup>®</sup> -Blunt Simple Cloning Vector (10 ng/μl)  | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                  | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                     | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                  | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                  | 50 μl     | 150 μl    |
| SR Primer (10 μM)   | 50 μl     | 150 μl    |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

### pEASY<sup>®</sup>-Blunt Simple Cloning Vector Map



*LacZ*α fragment: bases 1-446  
M13 reverse priming site: bases 205-221  
T7 promoter priming site: bases 263-282  
M13 forward priming site: bases 289-305  
f1 origin: bases 447-884  
Kanamycin resistance ORF: bases 1,218-2,012  
Ampicillin resistance ORF: bases 2,030-2,890  
pUC origin: bases 3,035-3,708

SR Primer: CAG GCT TTA CAC TTT ATG CTT  
M13 Reverse Primer: CA GGA AAC AGC TAT GAC

T7 Promoter: AA GGG CAG CTT CAA TTC  
M13 Forward Primer: ACT GGC CGT CGT TTT ACA

PCR Product: CC CTA TAG TGA GTC GTA TTA  
GG GAT ATC ACT CAG CAT AAT

## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100, except the PCR product size from vector self-ligation is 101 bp. General notes for cloning are the same as described on page 98.



# pEASY<sup>®</sup>-T3 Cloning Kit

|          |         |
|----------|---------|
| CT301-01 | 20 rxns |
| CT301-02 | 60 rxns |

## Storage

*Trans1-T1* Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description

pEASY<sup>®</sup>-T3 Cloning Vector provides dual *EcoR* I and dual *Not* I restriction sites. It is designed for cloning and sequencing *Taq*-amplified PCR products. The cloned insert can be released from a single enzyme digestion.

- 5 minutes fast ligation of *Taq*-amplified PCR products.
- Ampicillin resistance gene for selection.
- Easy blue/white selection.
- T7 promoter, SP6 promoter, M13 forward and M13 reverse primers for sequencing.
- T7 promoter and SP6 promoter for *in vitro* transcription.
- *Trans1-T1* Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.


## Kit Contents

| Component  | CT301-01  | CT301-02  |
|--|-----------|-----------|
| pEASY <sup>®</sup> -T3 Cloning Vector (10 ng/μl)           | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                 | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                    | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                 | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                 | 50 μl     | 150 μl    |
| <i>Trans1-T1</i> Phage Resistant Chemically Competent Cell | 10x100 μl | 30x100 μl |

## pEASY<sup>®</sup>-T3 Cloning Vector Map



*Lac* operon sequence: bases 2,860-3,020,190-419  
 Multiple cloning site: bases 10-152  
 SP6 priming site: bases 163-182  
 M13 reverse priming site: bases 200-216  
*LacZ* start codon: base 204  
*Lac* operator: bases 224-240  
 pUC origin: bases 543-1,216  
 Ampicillin resistance ORF (c): bases 1,361-2,221  
 f1 origin: bases 2,421-2,858  
 M13 forward priming site: bases 3,000-3,016  
 T7 promoter priming site: bases 3,023-3  
 (c) = complementary strand



## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100, except the PCR product size from vector self-ligation is 253 bp. General notes for cloning are the same as described on page 98.

# pEASY<sup>®</sup>-Blunt3 Cloning Kit

|          |         |
|----------|---------|
| CB301-01 | 20 rxns |
| CB301-02 | 60 rxns |

### Storage

Trans1-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

### Description


pEASY<sup>®</sup>-Blunt3 Cloning Vector provides dual *EcoR* I and dual *Not* I enzyme digestion sites. It is designed for cloning and sequencing *Pfu*-amplified PCR products. The cloned insert can be released from a single enzyme digestion.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- Ampicillin resistance gene for selection.
- Easy blue/white selection.
- T7 promoter, SP6 promoter, M13 forward and M13 reverse primers for sequencing.
- T7 promoter and SP6 promoter for *in vitro* transcription.
- Trans1-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.


### Kit Contents

| Component  | CB301-01  | CB301-02  |
|--|-----------|-----------|
| pEASY <sup>®</sup> -Blunt3 Cloning Vector (10 ng/μl) | 20 μl     | 3x20 μl   |
| Control Template (5 ng/μl)                           | 5 μl      | 5 μl      |
| Control Primers (10 μM)                              | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                           | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                           | 50 μl     | 150 μl    |
| Trans1-T1 Phage Resistant Chemically Competent Cell  | 10x100 μl | 30x100 μl |

### pEASY<sup>®</sup>-Blunt3 Cloning Vector Map



*Lac* operon sequence: bases 2,861-3,021,191-420  
 Multiple cloning site: bases 10-153  
 SP6 priming site: bases 164-183  
 M13 reverse priming site: bases 201-217  
*LacZ* start codon: base 205  
*Lac* operator: bases 225-241  
 pUC origin: bases 544-1,217  
 Ampicillin resistance ORF (c): bases 1,362-2,222  
 f1 origin: bases 2,422-2,859  
 M13 forward priming site: bases 3,001-3,017  
 T7 promoter priming site: bases 3,024-3  
 (c) = complementary strand



## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100, except the PCR product size from vector self-ligation is 254 bp. General notes for cloning are the same as described on page 98.



# pEASY<sup>®</sup>-T5 Zero Cloning Kit

|          |         |
|----------|---------|
| CT501-01 | 20 rxns |
| CT501-02 | 60 rxns |

## Storage

*Trans1-T1* Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description


pEASY<sup>®</sup>-T5 Zero Cloning Vector contains a suicide gene. Ligation of PCR fragment disrupts the expression of the gene. Cells that contain non-recombinant vector are killed upon plating. Therefore, blue/white selection is not required.

- 5 minutes fast ligation of *Taq*-amplified PCR products.
- High cloning efficiency. Positive clones up to 100%.
- No blue/white selection needed.
- Suitable for short and large fragment cloning.
- Kanamycin and Ampicillin resistance genes for selection.
- M13 forward primer and M13 reverse primer for sequencing.
- T3 promoter and T7 promoter for *in vitro* transcription.
- *Trans1-T1* Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.


## Kit Contents

| Component  | CT501-01  | CT501-02  |
|--|-----------|-----------|
| pEASY <sup>®</sup> -T5 Zero Cloning Vector (10 ng/μl)      | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                 | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                    | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                 | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                 | 50 μl     | 150 μl    |
| <i>Trans1-T1</i> Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

## pEASY<sup>®</sup>-T5 Zero Cloning Vector Map



*LacZα* fragment: bases 217-809  
 M13 reverse priming site: bases 207-223  
 T7 promoter priming site: bases 327-346  
 M13 Forward priming site: bases 353-369  
 Kanamycin resistance ORF: bases 1,158-1,952  
 Ampicillin resistance ORF (c): bases 2,202-3,062  
 pUC origin: bases 3,160-3,833  
 (c) = complementary strand



## PROTOCOL

**Protocols for cloning, transformation and analysis are the same as described on page 100. General notes for cloning are the same as described on page 98.**

# pEASY<sup>®</sup>-Blunt Zero Cloning Kit

|          |         |
|----------|---------|
| CB501-01 | 20 rxns |
| CB501-02 | 60 rxns |

## Storage

*Trans1*-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description


pEASY<sup>®</sup>-Blunt Zero Cloning Vector contains a suicide gene. Ligation of PCR fragment disrupts the expression of the gene. Cells that contain non-recombinant vector are killed upon plating. Therefore, blue/white selection is not required.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- High cloning efficiency. Positive clones up to 100%.
- No blue/white selection needed.
- Suitable for short and large fragment cloning.
- Kanamycin and Ampicillin resistance genes for selection.
- M13 forward primer and M13 reverse primer for sequencing.
- T3 promoter and T7 promoter for *in vitro* transcription.
- *Trans1*-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.


## Kit Contents

| Component   | CB501-01  | CB501-02  |
|---|-----------|-----------|
| pEASY <sup>®</sup> -Blunt Zero Cloning Vector (10 ng/μl)    | 20 μl     | 3×20 μl   |
| Control Template (5 ng/μl)                                  | 5 μl      | 5 μl      |
| Control Primers (10 μM)                                     | 5 μl      | 5 μl      |
| M13 Forward Primer (10 μM)                                  | 50 μl     | 150 μl    |
| M13 Reverse Primer (10 μM)                                  | 50 μl     | 150 μl    |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | 10×100 μl | 30×100 μl |

## pEASY<sup>®</sup>-Blunt Zero Cloning Vector Map



LacZα fragment: bases 217-810  
 M13 reverse priming site: bases 205-221  
 T7 promoter priming site: bases 328-347  
 M13 Forward priming site: bases 354-370  
 Kanamycin resistance ORF: bases 1,159-1,953  
 Ampicillin resistance ORF (c): bases 2,203-3,063  
 pUC origin: bases 3,161-3,834  
 (c) = complementary strand



## PROTOCOL

Protocols for cloning, transformation and analysis are the same as described on page 100. General notes for cloning are the same as described on page 98.



# *pEASY*<sup>®</sup>-Uni Seamless Cloning and Assembly Kit

CU101-01

10 rxns

## Storage

*Trans1*-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for one year

## Description

This kit takes advantage of proprietary assembly mix and homologous recombination. This kit can achieve directional cloning of PCR fragments that share 15-25 bp overlapping sequences into any linearized vector.


- Fast: 15 minutes.
- Broad: no restriction enzyme digestions. Can be cloned into any sites.
- High efficiency: up to 95% cloning efficiency.
- Seamless: no extra sequences introduced; up to 5 fragments assembly.

## Kit Contents

| Component   | CU101-01 |
|---|----------|
| 2xAssembly Mix  | 50 µl    |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | 5×100 µl |
| Linearized pUC19 Control Vector (10 ng/µl)                  | 3 µl     |
| Control Insert (1 kb, 20 ng/µl)                             | 3 µl     |

## Principle

1. Prepare linearized vector by PCR/Enzyme digestion
2. PCR amplify inserts with 15-25 bp overlapping sequences
3. Mix vector, DNA fragments and Assembly Mix together, incubate at 50°C for 15 minutes



4. Transformation

## PROTOCOL

### Cloning

#### Preparation of Vector and Inserts

A: Preparation of Vector

- (1) Enzyme digestion: digest plasmid vector with restriction enzyme(s) to generate the linearized vector. Purify the digested vector using Gel Extraction Kit (Cat. No. EG101).
- (2) PCR amplification: prepare the linearized vector by high-fidelity DNA polymerase. If a single expected band is generated, use PCR Purification Kit (Cat. No. EP101) to purify the product. Otherwise, use Gel Extraction Kit to recover the product.


In order to increase the positive cloning efficiency, we suggest using DMT enzyme to digest plasmid template before PCR purification or gel extraction. Add DMT enzyme (Cat. No. GD111) after PCR amplification (1 µl of DMT enzyme for a 50 µl PCR system), and incubate at 37°C for 30 minutes.

B: Preparation of Inserts


(1) Forward primer (5'-3'): 15-25 nt homology of linearized vector + 20-25 nt target specific sequence.

Reverse primer (5'-3'): 15-25 nt homology of linearized vector + 20-25 nt target specific sequence.

Example:



(2) Primers for multiple fragments



(3) We suggest using high-fidelity DNA polymerases to generate both the linear vector and fragments.

(4) Reaction conditions

- Use 0.2-0.4 µM ( final concentration) primers for PCR.
- Use 60-68°C as annealing temperature.

(5) Purification of target DNA fragments

- To increase the cloning efficiency, if the recombinant vector has the same selection marker as the parental plasmid for PCR fragments, pretreat the PCR fragments with DMT enzyme before purification.
- If product is single band, we recommend using PCR Purification Kit (Cat. No. EP101) to purify your fragments.
- If products are multibands, we recommend using Gel Extraction Kit (Cat. No. EG101) to recover your fragments.





### Setting up the cloning reaction

| Component                   | Volume   |
|-----------------------------|----------|
| 2×Assembly Mix              | 5 µl     |
| Linearized vector(5-100 ng) | x µl*    |
| Inserts                     | y µl*    |
| ddH <sub>2</sub> O          | to 10 µl |

\* In a 10 µl system, we recommend using 0.01-0.025 pmols of vector and insert respectively, for optimal cloning efficiency, use 1:2 (vector: insert) molar ratio. pmols= (weight in ng)/(base pairs×0.65 kDa)

For example

100 ng of 2,000 bp insert is equal to  $100/(2,000 \times 0.65)$ , which is about 0.08 pmols. 100 ng of 5,000 bp insert is equal to  $100/(5,000 \times 0.65)$ , which is about 0.03 pmols. Gently mix and incubate at 50°C for 15 minutes. Place it on ice for a few seconds. The reaction mixture can be directly used for transformation or stored at -20°C.

### Transformation

- (1) Thaw a vial of *Trans1*-T1 Phage Resistant Chemically Competent Cell on ice.
- (2) Transfer 2 µl of reaction mixture into 50 µl of *Trans1*-T1 Phage Resistant Chemically Competent Cell and mix gently by flicking the tube (do not vortex). Incubate on ice for 30 minutes.
- (3) Heat-shock at 42°C for 30 seconds, and immediately place on ice for 2 minutes.
- (4) Add 450 µl of room temperature SOC/LB medium. Incubate at 37°C for 1 hour at 250 rpm.
- (5) Pre-warm LB plate containing the appropriate selection antibiotic at 37°C.
- (6) Spread 100 µl of cells on the selection plate and incubate overnight at 37°C.

### Analysis of Positive Clones

#### Analyzing positive clones by PCR

- (1) Pick single colony into 10 µl of sterile water. Mix by vortexing or pipetting up and down.
- (2) Add 1 µl of mixture into 25 µl of PCR system. Identify the positive clones by appropriate forward and reverse primer.

#### Analyzing positive clones by restriction enzyme digestion

Pick several single colony and culture them overnight in LB medium containing the appropriate selection antibiotic. Isolate plasmid DNA by *EasyPure*® Plasmid MiniPrep Kit. Analyze the plasmids by restriction enzyme digestion.

### Sequencing

Perform sequence analysis using vector universal primers

#### Cloning reaction for control insert

| Component                       | Volume |
|---------------------------------|--------|
| 2×Assembly Mix                  | 5 µl   |
| Linearized pUC19 Control Vector | 1 µl   |
| Control Insert                  | 1 µl   |
| ddH <sub>2</sub> O              | 3 µl   |

Reaction conditions, transformation and analysis of positive clones are the same as described above.

# Cloning Competent Cells

## Selection Guide

| Name               | Cat. No. | Transformation Efficiency  | Blue/White Selection/<br>( <i>lacZ</i> ΔM15) | Low Recombination Rate<br>( <i>recA</i> ) | High Quality Plasmid DNA Prepared (endA1) | Cloning of Toxic Gene | Phage Resistance |
|--------------------|----------|----------------------------|--|---|---|-----------------------|------------------|
| <i>Trans10</i>     | CD101    | 10 <sup>8</sup> cfu/μg DNA | •  | •   | •   | •                     | —                |
| <i>Trans5α</i>     | CD201    | 10 <sup>8</sup> cfu/μg DNA | •  | •   | •   | —                     | —                |
| <i>Trans109</i>    | CD301    | 10 <sup>8</sup> cfu/μg DNA | •  | ••  | ••  | —                     | —                |
| <i>Trans110</i>    | CD311    | 10 <sup>8</sup> cfu/μg DNA | •  | ••  | ••  | —                     | —                |
| <i>Trans1-Blue</i> | CD401    | 10 <sup>8</sup> cfu/μg DNA | •  | •   | •   | —                     | —                |
| <i>Trans2-Blue</i> | CD411    | 10 <sup>9</sup> cfu/μg DNA | •  | •   | •   | —                     | —                |
| <i>Trans1-T1</i>   | CD501    | 10 <sup>9</sup> cfu/μg DNA | •  | •   | •   | —                     | •                |
| DMT                | CD511    | 10 <sup>8</sup> cfu/μg DNA | •  | •   | •   | —                     | •                |
| <i>TransStbl3</i>  | CD521    | 10 <sup>8</sup> cfu/μg DNA | —  | •   | —   | —                     | —                |
| <i>TransDB3.1</i>  | CD531    | 10 <sup>8</sup> cfu/μg DNA | —  | •   | •   | •                     | —                |

## *Trans10* Chemically Competent Cell

|          |           |
|----------|-----------|
| CD101-01 | 10×100 μl |
| CD101-02 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency: >10<sup>8</sup> cfu/μg (pUC19 DNA).
- Str<sup>R</sup>.
- Blue/white selection.
- Toxic gene cloning and stable replication of plasmid DNA.

### Genotype

F<sup>-</sup> *mcrA* Δ(*mrr-hsdRMS-mcrBC*) φ80*lacZ*ΔM15 Δ*lacX74 recA1 araD139* Δ(*ara-leu*)7697 *galU galK rpsL* (Str<sup>R</sup>) *endA1 nupG*

## *Trans5α* Chemically Competent Cell

|          |           |
|----------|-----------|
| CD201-01 | 10×100 μl |
| CD201-02 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency: >10<sup>8</sup> cfu/μg (pUC19 DNA).
- Reduced recombination of cloned DNA.
- Blue/white selection.

### Genotype

F<sup>-</sup> φ80 *lacZ*ΔM15 Δ(*lacZYA-argF*) U169 *endA1 recA1 hsdR17* (*r<sub>k</sub><sup>-</sup>, m<sub>k</sub><sup>+</sup>*) *supE44λ- thi-1 gyrA96 relA1 phoA*

High quality products



## Trans109 Chemically Competent Cell

|          |           |
|----------|-----------|
| CD301-02 | 10×100 µl |
| CD301-03 | 20×100 µl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^8$  cfu/µg (pUC19 DNA).
- The lowest homologous recombination is favorable for plasmid DNA preparation.
- Routine cloning.
- Blue/white selection.

### Genotype

*endA1 recA1 gyrA96 thi-1 hsdR17 (r<sub>k</sub><sup>-</sup>, m<sub>k</sub><sup>+</sup>) relA1 supE44 Δ(lac-proAB) [F'*traD36 proAB lacI*<sup>q</sup>ZΔM15]*

## Trans110 Chemically Competent Cell

|          |           |
|----------|-----------|
| CD311-02 | 10×100 µl |
|----------|-----------|

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency:  $>10^8$  cfu/µg (pUC19 DNA).
- Unmethylated DNA due to *dam*<sup>-</sup>/*dcm*<sup>-</sup>.
- Str<sup>R</sup>.

### Genotype

*rpsL (Str<sup>R</sup>) thr leu thi-1 lacY galK galT ara tonA tsx dam dcm supE44 Δ(lac-proAB) /F' [traD36 proAB lacI<sup>q</sup> lacZΔM15]*

## Trans1-Blue Chemically Competent Cell

|          |           |
|----------|-----------|
| CD401-02 | 10×100 µl |
| CD401-03 | 20×100 µl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^8$  cfu/µg (pUC19 DNA).
- Tet<sup>R</sup>.
- Blue/white selection.

### Genotype

*recA1 endA1 gyrA96 thi-1 hsdR17 supE44 relA1 lac [F' proAB lacI<sup>q</sup>ZΔM15: Tn10 (Tet<sup>R</sup>)]*

## Trans2-Blue Chemically Competent Cell

|          |           |
|----------|-----------|
| CD411-02 | 10×100 µl |
| CD411-03 | 20×100 µl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^9$  cfu/µg (pUC19 DNA).
- Suitable for larger plasmid transformation.
- Reduced preference for plasmid size, suitable for library construction.
- Tet<sup>R</sup> and Cam<sup>R</sup>.
- Blue/white selection.

### Genotype

Tet<sup>R</sup>Δ(*mcrA*)183 Hte[F' {*proAB lacI<sup>q</sup> lacZ*ΔM15 *Tn10*(Tet<sup>R</sup>) *Amy Cam<sup>R</sup>*}]  
Δ(*mcrCB-hsdSMR-mrr*)173 *endA1 supE44 thi-1 recA1 gyrA96 relA1*

## Trans1-T1 Phage Resistant Chemically Competent Cell

|          |           |
|----------|-----------|
| CD501-01 | 5×100 µl  |
| CD501-02 | 10×100 µl |
| CD501-03 | 20×100 µl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^9$  cfu/µg (pUC19 DNA).
- Fast-growing, colonies are visible in 8–9 hours.
- Resistance to T1 and T5 phage.
- Blue/white selection.

### Genotype

F<sup>-</sup> φ80(*lacZ*)ΔM15 Δ*lacX74 hsdR*(r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>) Δ*recA1398 endA1 tonA*

## DMT Chemically Competent Cell

|          |          |
|----------|----------|
| CD511-01 | 10×50 µl |
| CD511-02 | 20×50 µl |

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^8$  cfu/µg (pUC19 DNA).
- Resistance to T1 and T5 phage.
- *In vivo* digestion of methylated DNA, suitable for site-directed mutagenesis.

### Genotype

F<sup>-</sup> φ80 *lacZ*ΔM15 Δ(*lacZYA-argF*)U169 *recA1 endA1 hsdR17*(r<sub>K</sub><sup>-</sup>, m<sub>K</sub><sup>+</sup>)  
*phoA supE44 thi-1 gyrA96 relA1 tonA*

High quality products



## TransStbl3 Chemically Competent Cell

CD521-01

10×100 μl

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^8$  cfu/μg (pUC19 DNA).
- Suitable for lentivirus and retrovirus plasmid vectors transformation.
- Str<sup>R</sup>
- Reduced the frequency of homologous recombination of long terminal repeats.

### Genotype

F<sup>-</sup> *mcrB mrr hsdS20*(r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>-</sup>) *recA13 supE44 ara-14 galK2 lacY1 proA2 rpsL20* (Str<sup>R</sup>) *xyl-5 λ-leu mtl-1*

## TransDB3.1 Chemically Competent Cell

CD531-01

10×100 μl

### Storage

at -70°C for six months

### Characteristics

- High transformation efficiency:  $>10^8$  cfu/μg (pUC19 DNA).
- Transformation and propagation of plasmids containing the *ccdB* gene.
- Str<sup>R</sup>.

### Genotype

F<sup>-</sup> *gyrA462 endA1 Δ(sr1-recA) mcrB mrr hsdS20*(r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>-</sup>) *supE44ara-14 galK2 lacY1 proA2 rpsL20*(Sm<sup>R</sup>) *xyl-5 λ-leu mtl1*





# Fast MultiSite Mutagenesis System

FM201-01

10 rxns

## Storage

DMT Chemically Competent Cell at  $-70^{\circ}\text{C}$  for six months; others at  $-20^{\circ}\text{C}$  for two years

## Description


*Fast MultiSite Mutagenesis System* is used for generating mutated PCR fragments by introducing mutation sites on overlapping regions. High fidelity *TransStart<sup>®</sup> FastPfu* PCR SuperMix is included for amplification. This kit uses proprietary assembly mix and homologous recombination to seamlessly assemble up to six mutagenesis fragments.

- **Fast:** Amplified with fast & high-fidelity *2xTransStart<sup>®</sup> FastPfu* PCR SuperMix; only 15 minutes for recombination.
- **Flexible:** Able to be cloned into any site to realize single-site/ multi-site, continuous/non-continuous mutagenesis.
- **Efficient:** >90% mutagenesis efficiency.

## Kit Contents

| Component  | FM201-01            |
|--|---------------------|
| <i>2xTransStart<sup>®</sup> FastPfu</i> PCR SuperMix | 1 ml                |
| DMT Enzyme (10 units/ $\mu\text{l}$ )                | 30 $\mu\text{l}$    |
| <i>2xAssembly Mix</i>                                | 50 $\mu\text{l}$    |
| DMT Chemically Competent Cell                        | 10x50 $\mu\text{l}$ |
| ddH <sub>2</sub> O                                   | 1 ml                |

## Cloning principle






PROTOCOL

**Preparation of multisite mutagenic fragment**


**(1) Primer Design**

- Both primers contain overlapping region at the 5' ends and extension region at the 3' ends, with mutation site on overlapping region, as shown in figure 1.
- Primer length: Both primers (forward and reverse) should be approximately at 25-40 nucleotides in length, excluding the mutation site. Primers should have an overlapping region of 15-25 nucleotides and have an extension region of at least 10 nucleotides.



**(2) Preparation of mutated fragment**

- The mutation sites are located on one pair of primers, as shown in figure 2.
- The mutation sites are located on multiple pairs of primers, with F1/R1, F2/R2, ..., Fn/Rn for amplification, as shown in figure 3.





### PCR System

| Component   | Volume        | Final Concentration |
|---|---------------|---------------------|
| Plasmid   | 1-10 ng       | as required         |
| Forward Primer (10 $\mu$ M)                             | 1 $\mu$ l     | 0.2 $\mu$ M         |
| Reverse Primer (10 $\mu$ M)                             | 1 $\mu$ l     | 0.2 $\mu$ M         |
| 2 $\times$ TransStart <sup>®</sup> FastPfu PCR SuperMix | 25 $\mu$ l    | 1 $\times$          |
| ddH <sub>2</sub> O                                      | to 50 $\mu$ l | Not applicable      |

#### PCR

|                         |            |                           |
|-------------------------|------------|---------------------------|
| 95°C                    | 3 min      | } 25 cycles <sup>*2</sup> |
| 95°C                    | 20 sec     |                           |
| 55°C-65°C <sup>*1</sup> | 20 sec     |                           |
| 72°C                    | 2-4 kb/min |                           |
| 72°C                    | 5-10 min   |                           |

#### Notes

\*1. Annealing temperature depends on primers.

\*2. We suggest performing 25 cycles for PCR. For low yield PCR products, we suggest using 30 cycles.

#### Electrophoresis Analysis

Amplified PCR products can be checked by electrophoresis with 10  $\mu$ l of PCR product on a 1% agarose gel.

#### (3) Digestion of PCR Product with DMT

Add 1  $\mu$ l of DMT enzyme into PCR product, mix thoroughly and incubate at 37°C for 1 hour.

#### (4) Purification of PCR products

For PCR product with the single expected band, we suggest using PCR Purification Kit to purify PCR products; for PCR product with multibands, we suggest using Quick Gel Extraction Kit to purify PCR products.

#### Assembly of Mutated Fragments

| Component               | Volume        |
|-------------------------|---------------|
| 2 $\times$ Assembly Mix | 5 $\mu$ l     |
| Amplified fragment A    | x $\mu$ l*    |
| Amplified fragment B    | y $\mu$ l*    |
| .....                   | .....         |
| Amplified fragment N    | z $\mu$ l*    |
| ddH <sub>2</sub> O      | to 10 $\mu$ l |

\*Suggested amount is 20-150 ng

Gently mix and perform reaction at 50°C for 15 minutes. After reaction, transfer the reaction tube on ice for a few seconds.

#### Transformation

- (1) Add 2  $\mu$ l of assembly products into 50  $\mu$ l of DMT Chemically Competent Cell (DNA should be added immediately after thawing the cells on ice) and mix by tapping gently. Incubate on ice for 20-30 minutes.
- (2) Heat-shock at 42°C for exactly 45 seconds, quickly remove from 42°C water bath and place on ice for 2 minutes.
- (3) Add 250  $\mu$ l of SOC or LB medium (pre-warm to room temperature), and incubate at 37°C for 1 hour with shaking at 200 rpm.
- (4) Pre-warm a selective plate at 37°C for 30 minutes.
- (5) Spread 100-200  $\mu$ l of transformants on the plate and incubate at 37°C overnight.

#### Positive Clone Analysis

Analyze the clones by sequencing.

---

## Chapter 4 Nucleic Acid Purification

---

### Genomic DNA Purification

|   |     |
|---|-----|
| <i>BloodZol</i> .....   | 121 |
| <i>PlantZol</i> .....   | 122 |
| <i>EasyPure</i> <sup>®</sup> Genomic DNA Kit .....                          | 123 |
| <i>EasyPure</i> <sup>®</sup> Plant Genomic DNA Kit .....                    | 125 |
| <i>EasyPure</i> <sup>®</sup> Blood Genomic DNA Kit .....                    | 126 |
| <i>EasyPure</i> <sup>®</sup> Marine Animal Genomic DNA Kit .....            | 127 |
| <i>EasyPure</i> <sup>®</sup> Bacteria Genomic DNA Kit .....                 | 128 |
| <i>EasyPure</i> <sup>®</sup> Food and Fodder Security Genomic DNA Kit ..... | 129 |
| <i>EasyPure</i> <sup>®</sup> Micro Genomic DNA Kit .....                    | 131 |

### Plasmid DNA Purification and *E.coli* Medium

|  |     |
|--|-----|
| <i>EasyPure</i> <sup>®</sup> Plasmid MiniPrep Kit .....        | 132 |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MiniPrep Kit ..... | 133 |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MaxiPrep Kit ..... | 134 |
| <i>ArtMedia</i> <sup>®</sup> Plasmid Culture .....             | 135 |

### DNA Purification

|   |     |
|---|-----|
| <i>EasyPure</i> <sup>®</sup> PCR Purification Kit .....     | 136 |
| <i>EasyPure</i> <sup>®</sup> Quick Gel Extraction Kit ..... | 137 |

### RNA Purification

|   |     |
|---|-----|
| <i>TransZol</i> .....                                   | 138 |
| <i>TransZol</i> Up .....                                | 139 |
| <i>TransZol</i> Plant .....                             | 140 |
| <i>EasyPure</i> <sup>®</sup> RNA Kit .....              | 141 |
| <i>EasyPure</i> <sup>®</sup> Viral DNA/RNA Kit .....    | 142 |
| <i>EasyPure</i> <sup>®</sup> Plant RNA Kit .....        | 143 |
| <i>EasyPure</i> <sup>®</sup> Blood RNA Kit .....        | 144 |
| <i>TransZol</i> Up Plus RNA Kit .....                   | 145 |
| <i>EasyPure</i> <sup>®</sup> miRNA Kit .....            | 146 |
| <i>EasyPure</i> <sup>®</sup> RNA Purification Kit ..... | 147 |
| <i>RNAhold</i> <sup>®</sup> .....                       | 148 |




# BloodZol

|          |                  |
|----------|------------------|
| EE131-01 | For 50 ml blood  |
| EE131-02 | For 200 ml blood |

## Storage

Proteinase K solution at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures




## Description

*BloodZol* provides an easy and fast method to isolate high quality genomic DNA from 0.1-20 ml of fresh or frozen blood. Isolated DNA is free of contaminants and enzyme inhibitors. Red Cell Lysis Buffer is provided to remove non-nucleated red cells and reduce hemoglobin contamination. Genomic DNA is precipitated with isopropanol.

- High quality, free of contaminants and inhibitors.
- Suitable for EDTA, sodium citrate and heparin-anticoagulated fresh and frozen blood.
- No organic solvents.
- Isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

## Kit Contents

| Component                   | EE131-01 | EE131-02 |
|-----------------------------|----------|----------|
| Red Cell Lysis Buffer (RCL) | 125 ml   | 2×250 ml |
| Lysis Buffer 3 (LB3)        | 30 ml    | 120 ml   |
| Elution Buffer (EB)         | 25 ml    | 80 ml    |
| Proteinase K (20 mg/ml)     | 250 µl   | 1 ml     |



## DNA yield from different samples

| Blood             | Amount | Yield  |
|-------------------|--------|--------|
| Human whole blood | 400 µl | ~30 µg |
| Mouse whole blood | 400 µl | ~20 µg |

# PlantZol

EE141-01

100 ml

### Storage



at room temperature (15-25°C) for one year

### Description

*PlantZol* provides an easy and fast method to isolate high quality plant genomic DNA. Plant tissue is disrupted by grinding in liquid nitrogen. DNA is released with detergent. DNA is separated from other components by centrifugation and precipitated with isopropanol. *PlantZol* is suitable to isolate DNA from plants rich in polysaccharide and polyphenol.

- Isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

### Procedures



### DNA yield from different fresh plant leaves (100 mg)

| Material           | Yield  |
|--------------------|--------|
| Tobacco leaves     | ~20 µg |
| Wheat leaves       | ~35 µg |
| Rape leaves        | ~9 µg  |
| Rice leaves        | ~29 µg |
| Soybean leaves     | ~16 µg |
| Arabidopsis leaves | ~28 µg |
| Maize leaves       | ~22 µg |
| Tomato leaves      | ~7 µg  |



# EasyPure<sup>®</sup> Genomic DNA Kit

|              |          |          |
|--------------|----------|----------|
| RNase A      | EE101-01 | 50 rxns  |
|              | EE101-02 | 200 rxns |
| RNase A-free | EE101-11 | 50 rxns  |
|              | EE101-12 | 200 rxns |

## Storage

RNase A and Proteinase K solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Description

EasyPure<sup>®</sup> Genomic DNA Kit provides a simple and convenient way to isolate high quality genomic DNA from a variety of mammalian cells, tissues, *E.coli* and yeast. Cells and tissues are enzymatically lysed. DNA is bound to silica-based column. The isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

- DNA yield up to 15 µg.
- Complete removal of contaminants and inhibitors.
- Column based purification, no organic extraction or ethanol precipitation.

## Kit Contents


| Component                                  | EE101-01        | EE101-02          |
|--|-----------------|-------------------|
|  | EE101-11        | EE101-12          |
| Lysis Buffer 2 (LB2)                       | 6 ml            | 24 ml             |
| Binding Buffer 2 (BB2)                     | 28 ml           | 110 ml            |
| Clean Buffer 2 (CB2)                       | 55 ml           | 2×110 ml          |
| Wash Buffer 2 (WB2)                        | 12 ml           | 2×22 ml           |
| Elution Buffer (EB)                        | 25 ml           | 80 ml             |
| RNase A (20 mg/ml)                         | 1 ml (EE101-01) | 4×1 ml (EE101-02) |
|  | 0 (EE101-11)    | 0 (EE101-12)      |
| Proteinase K (20 mg/ml)                    | 1 ml            | 4×1 ml            |
| Genomic Spin Columns with Collection Tubes | 50 each         | 200 each          |

## Sample Requirement

| Material            | Amount                    |
|---------------------|---------------------------|
| Mammalian Cells     | 1-5×10 <sup>6</sup> cells |
| Mammalian Tissues   | ≤25 mg                    |
| Mouse Tail          | 0.5 cm sections           |
| <i>E.coli</i> Cells | ≤2×10 <sup>9</sup> cells  |
| Yeast Cells         | ≤5×10 <sup>7</sup> cells  |

## Procedures





**DNA yield from different mouse tissues**

| Tissue | Amount | Yield   |
|--------|--------|---------|
| Heart  | 25 mg  | ~5 µg   |
| Liver  | 25 mg  | ~10 µg  |
| Spleen | 25 mg  | ~12 µg  |
| Lung   | 25 mg  | ~5 µg   |
| Kidney | 25 mg  | ~10 µg  |
| Muscle | 25 mg  | ~2.5 µg |





# EasyPure<sup>®</sup> Plant Genomic DNA Kit

|              |          |          |
|--------------|----------|----------|
| RNase A      | EE111-01 | 50 rxns  |
|              | EE111-02 | 200 rxns |
| RNase A-free | EE111-11 | 50 rxns  |
|              | EE111-12 | 200 rxns |

## Storage

RNase A at -20°C for one year; others at room temperature (15-25°C) for one year

## Description

EasyPure<sup>®</sup> Plant Genomic DNA Kit provides a simple and convenient way to isolate high quality genomic DNA from plant tissues (up to 100 mg). The isolated genomic DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

- DNA yield up to 15 µg.
- Complete removal of pigment, polysaccharides and other impurities.
- Column based purification, no organic extraction or ethanol precipitation.

## Kit Contents

| Component                                  | EE111-01        | EE111-02          |
|--|-----------------|-------------------|
|  | EE111-11        | EE111-12          |
| Resuspension Buffer 1 (RB1)                | 25 ml           | 100 ml            |
| Precipitation Buffer 1 (PB1)               | 6 ml            | 24 ml             |
| Binding Buffer 1 (BB1)                     | 10 ml           | 2×20 ml           |
| Wash Buffer 1 (WB1)                        | 12 ml           | 2×24 ml           |
| Elution Buffer (EB)                        | 25 ml           | 80 ml             |
| RNase A (10 mg/ml)                         | 1 ml (EE111-01) | 4×1 ml (EE111-02) |
|  | 0 (EE111-11)    | 0 (EE111-12)      |
| Genomic Spin Columns with Collection Tubes | 50 each         | 200 each          |

## Procedures

Prepare plant lysate and resuspend in Resuspension Buffer 1 (RB1)

Precipitate pigment, phenol and polysaccharide with Precipitation Buffer 1 (PB1)



Add Binding Buffer 1 (BB1) and ethanol to the plant lysate



Apply lysate to a Spin Column




Wash the column once with Clean Buffer 1 (CB1)



Wash the column twice with Wash Buffer 1 (WB1)



Elute DNA with Elution Buffer (EB) or ddH<sub>2</sub>O



## DNA yield from different fresh plant leaves (100 mg)

| Material            | Yield  |
|---------------------|--------|
| Tobacco leaves      | ~10 µg |
| Wheat leaves        | ~7 µg  |
| Sweet potato leaves | ~9 µg  |
| Pepper leaves       | ~6 µg  |
| Rape leaves         | ~5 µg  |
| Rice leaves         | ~8 µg  |
| Soybean leaves      | ~7 µg  |

# EasyPure<sup>®</sup> Blood Genomic DNA Kit


|              |          |          |
|--------------|----------|----------|
| RNase A      | EE121-01 | 50 rxns  |
|              | EE121-02 | 200 rxns |
| RNase A-free | EE121-11 | 50 rxns  |
|              | EE121-12 | 200 rxns |

### Storage

RNase A and Proteinase K solutions at -20°C for one year; others at room temperature (15-25°C) for one year

### Procedures

Release DNA by Proteinase K digestion



### Description

EasyPure<sup>®</sup> Blood Genomic DNA Kit provides a simple and convenient way to isolate high quality genomic DNA from 5-250 µl of fresh or frozen blood. Whole blood is incubated with binding/lysis buffer to release DNA. DNA is bound to silica-based column. The isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

- Simple and fast, red cell lysis buffer is no longer needed.
- Complete removal of contaminants and inhibitors.
- DNA yield up to 40 µg.
- Column based purification, no organic extraction or ethanol precipitation.
- Suitable for EDTA, sodium citrate and heparin-anticoagulated fresh or frozen blood in a volume of 5 to 250 µl.

### Kit Contents

| Component                                  | EE121-01          | EE121-02          |
|--|-------------------|-------------------|
|  | EE121-11          | EE121-12          |
| Binding Buffer 3 (BB3)                     | 30 ml             | 110 ml            |
| Clean Buffer 3 (CB3)                       | 6 ml              | 24 ml             |
| Wash Buffer 3 (WB3)                        | 12 ml             | 2×22 ml           |
| Elution Buffer (EB)                        | 25 ml             | 80 ml             |
| RNase A (20 mg/ml)                         | 500 µl (EE121-01) | 2×1 ml (EE121-02) |
|  | 0 (EE121-11)      | 0 (EE121-12)      |
| Proteinase K (20 mg/ml)                    | 1 ml              | 4×1 ml            |
| Genomic Spin Columns with Collection Tubes | 50 each           | 200 each          |



### DNA yield from different samples

| Material             | Volume | DNA yield |
|----------------------|--------|-----------|
| Human whole blood    | 100 µl | ~6 µg     |
| Mouse whole blood    | 100 µl | ~6 µg     |
| Bullfrog whole blood | 20 µl  | ~20 µg    |
| Chicken whole blood  | 20 µl  | ~29 µg    |



# EasyPure<sup>®</sup> Marine Animal Genomic DNA Kit

|              |          |         |
|--------------|----------|---------|
| RNase A      | EE151-01 | 50 rxns |
| RNase A-free | EE151-11 | 50 rxns |

## Storage

RNase A and Proteinase K solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures

Prepare lysate using Lysis Buffer 8 (LB8) and Proteinase K, RNase A

Add Binding Buffer 8 (BB8) and ethanol to the lysate



Apply lysate to a Spin Column



Wash the column twice with Clean Buffer 8 (CB8)



Wash the column twice with Wash Buffer 8 (WB8)



Elute DNA with Elution Buffer (EB) or ddH<sub>2</sub>O


## Description

EasyPure<sup>®</sup> Marine Animal Genomic DNA Kit provides a simple and convenient way to isolate high quality genomic DNA from up to 30 mg marine animals. DNA is bound to silica-based column. The isolated DNA is suitable for PCR, restriction enzyme digestion and Southern blot.

- DNA yield up to 40 µg.
- Complete removal of contaminants and inhibitors.
- Column based purification, no organic extraction or ethanol precipitation.

## Kit Contents

| Component                                  | EE151-01        | EE151-11     |
|--|-----------------|--------------|
| Lysis Buffer 8 (LB8)                       | 12 ml           |              |
| Binding Buffer 8 (BB8)                     | 9 ml            |              |
| Clean Buffer 8 (CB8)                       | 12 ml           |              |
| Wash Buffer 8 (WB8)                        | 12 ml           |              |
| Elution Buffer (EB)                        | 25 ml           |              |
| RNase A (10 mg/ml)                         | 1 ml (EE151-01) | 0 (EE151-11) |
| Proteinase K (20 mg/ml)                    | 1 ml            |              |
| Genomic Spin Columns with Collection Tubes | 50 each         |              |



M: Trans15K<sup>®</sup> DNA Marker  
 Lane 1: white clam  
 Lane 2: razor clam  
 Lane 3: oyster  
 Lane 4: crab  
 Lane 5: king prawn  
 Lane 6: scallop

## DNA yield from different animal tissues

| Material           | Amount | DNA yield |
|--------------------|--------|-----------|
| Scallop            | 30 mg  | ~30 µg    |
| Razor clam         | 30 mg  | ~26 µg    |
| Small-sized shrimp | 30 mg  | ~10 µg    |
| Mantis shrimp      | 30 mg  | ~16 µg    |
| Crab               | 30 mg  | ~2.5 µg   |
| Oyster             | 30 mg  | ~38 µg    |
| White clam         | 30 mg  | ~50 µg    |

# EasyPure<sup>®</sup> Bacteria Genomic DNA Kit

|              |          |         |
|--------------|----------|---------|
| RNase A      | EE161-01 | 50 rxns |
| RNase A-free | EE161-11 | 50 rxns |

### Storage

RNase A and Proteinase K solutions at -20°C for one year; others at room temperature (15-25°C) for one year

### Description

EasyPure<sup>®</sup> Bacteria Genomic DNA Kit uses lysozyme and moderate lysis buffer to lyse cells. Proteinase K is used for protein digestion and RNase A used for RNA digestion. DNA is specifically bound to silica-based column in hypersaline condition, and DNA is eluted by low salt and high pH solution. This kit is suitable for isolating high quality genomic DNA from Gram-positive and Gram-negative bacteria. The isolated DNA is suitable for PCR, restriction enzyme digestion, and Southern blot.


- Fast: the whole process can be completed in 50 minutes
- High yield: DNA yield up to 20 µg

### Procedures

Resuspend Gram-Positive Bacteria in Resuspension Buffer 11 (RB11) and Lysozyme


Prepare Lysate in Lysis Buffer 11 (LB11) and Proteinase K, RNase A

Prepare Gram-Negative Bacteria Lysate in Lysis Buffer 11 (LB11) and Proteinase K, RNase A



### Kit Contents

| Component                                  | EE161-01        |
|--|-----------------|
|  | EE161-11        |
| Resuspension Buffer 11 (RB11)              | 12 ml           |
| Lysis Buffer 11 (LB11)                     | 6 ml            |
| Binding Buffer 11 (BB11)                   | 10 ml           |
| Clean Buffer 11 (CB11)                     | 55 ml           |
| Wash Buffer 11 (WB11)                      | 12 ml           |
| Elution Buffer (EB)                        | 25 ml           |
| RNase A (10 mg/ml)                         | 1 ml (EE161-01) |
|  | 0 (EE161-11)    |
| Proteinase K (20 mg/ml)                    | 1 ml            |
| Genomic Spin Columns With Collection Tubes | 50 each         |



Extraction from Gram-positive Bacteria

- 1: *Streptomyces coelicolor*
  - 2: *Staphylococcus aureus*
  - 3: *Lactobacillus acidophilus*
  - 4: *Bacillus subtilis*
- M: Trans2K<sup>®</sup> Plus II DNA Marker

Extraction from Gram-negative Bacteria

- 1: *Escherichia coli*
  - 2: *Citrobacter freundii*
  - 3: *Pseudomonas fluorescens*
- M: Trans2K<sup>®</sup> Plus II DNA Marker



# EasyPure<sup>®</sup> Food and Fodder Security Genomic DNA Kit


EE171-01

50 rxns

## Storage

Proteinase K and Carrier RNA solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures



## Description

This kit uses modified cetyltrimethylammonium bromide (CTAB) lysis method to lysis cells. DNA is bound to high-adsorption silica-based column and eluted with elution buffer without phenol/chloroform. This kit is designed for total DNA extraction from highly processed food material due to high temperature, or/and extreme pH. It is also suitable to isolate trace amount of animal DNA from fodder. The purified DNA can be used for the detection of genetically modified organisms, animal species in food and fodder.


- Strong lysis, fast extraction
- High purity, high efficiency DNA isolation

## Kit Contents


| Component                                  | EE171-01 |
|--|----------|
| Resuspension Buffer 13 (RB13)              | 180 ml   |
| Lysis Buffer 13 (LB13)                     | 30 ml    |
| Precipitation Buffer 13 (PB13)             | 12 ml    |
| Precipitation Buffer 14 (PB14)             | 18 ml    |
| Dissolving Buffer (DB)                     | 6 ml     |
| Binding Buffer 13 (BB13)                   | 10 ml    |
| Carrier RNA (1 µg/µl)                      | 50 µl    |
| Clean Buffer 13 (CB13)                     | 55 ml    |
| Wash Buffer 13 (WB13)                      | 12 ml    |
| Elution Buffer (EB)                        | 10 ml    |
| Proteinase K (20 mg/ml)                    | 1 ml     |
| Genomic Spin Columns With Collection Tubes | 50 each  |

## Sample Requirement

| Material   | Amount |
|--|--------|
| Seeds and flour  | 200 mg |
| Liquid processed food (e.g. soybean sauce, soybean milk) | 20 ml  |
| Oil (e.g. soy oil, rapeseed oil)                         | 20 ml  |
| Processed food (e.g. instant noodle, chips, ketchup)     | 200 mg |
| Cocoa nuts, chocolate                                    | 200 mg |
| Raw meat (e.g. beef, lamb, pork)                         | 200 mg |
| Meat-derived processed food                              | 200 mg |
| Fodder for cattle and sheep                              | 200 mg |



1: genomic DNA from soybean sauce  
M: *Trans2K*<sup>®</sup> Plus II DNA Marker



Amplify plant 18S rDNA from isolated soybean sauce genomic DNA  
1: positive control (soybean genomic DNA)  
2: genomic DNA from soybean sauce  
3: negative control  
M: *Trans2K*<sup>®</sup> Plus II DNA Marker



# EasyPure<sup>®</sup> Micro Genomic DNA Kit

EE181-01

50 rxns

## Storage

Proteinase K and Carrier RNA solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures

Enzymatic digestion

(Lysis Buffer 14+Proteinase K)

Add Binding Buffer 14 (BB14)  
and Carrier RNA



Apply lysate to a Spin Column



Wash the column twice with  
Clean Buffer 14 (CB14)



Wash the column twice with  
Wash Buffer 14 (WB14)



Elute DNA with Elution  
Buffer (EB) or ddH<sub>2</sub>O



## Description

EasyPure<sup>®</sup> Micro Genomic DNA Kit uses enzyme digestion method to lyse samples. The unique lysis buffer in this kit can efficiently lyse small volume of cells from a variety of materials including blood, dried blood spots, serum/plasma, mouthwash, hair follicles, tissues, microdissected tissues. DNA from the lysate will bind to silica-based column and elute with elution buffer. The isolated DNA is suitable for PCR, restriction enzyme digestion, and other downstream applications.

## Kit Contents

| Component                                  | EE181-01 |
|--|----------|
| Lysis Buffer 14 (LB14)                     | 6 ml     |
| Binding Buffer 14 (BB14)                   | 28 ml    |
| Clean Buffer 14 (CB14)                     | 28 ml    |
| Wash Buffer 14 (WB14)                      | 12 ml    |
| Elution Buffer (EB)                        | 5 ml     |
| Carrier RNA (1 µg/µl)                      | 55 µl    |
| Proteinase K (20 mg/ml)                    | 1 ml     |
| Genomic Spin Columns with Collection Tubes | 50 each  |

## Sample Requirement

| Material               | Amount                                   |
|------------------------|--|
| Cultured cells         | 1×10 <sup>4</sup> -10 <sup>8</sup> cells |
| Tissues                | ≤10 mg                                   |
| Microdissected tissues | ≤10 mg                                   |
| Formalin fixed tissues | ≤10 mg                                   |
| <i>E. coli</i>         | ≤1×10 <sup>9</sup> cells                 |
| Anti-coagulant blood   | 1-50 µl                                  |
| Serum/plasma           | 50-250 µl                                |
| Mouthwash              | 2-20 ml                                  |
| Dried blood spots      | 5 mm <sup>2</sup> -100 mm <sup>2</sup>   |
| Hair follicles         | 1-20 pieces                              |



# EasyPure<sup>®</sup> Plasmid MiniPrep Kit

|          |          |
|----------|----------|
| EM101-01 | 50 rxns  |
| EM101-02 | 200 rxns |

### Storage

RNase A at -20°C for one year; others at room temperature (15-25°C) for one year

### Description


EasyPure<sup>®</sup> Plasmid MiniPrep Kit uses a modified alkaline lysis method to isolate high-quality plasmid DNA from ≤20 ml (LB) or ≤4 ml (*ArtMedia*<sup>®</sup> Plasmid Culture) of bacterial culture. Unique formulated lysis buffer and neutralization buffer permit error-free visual identification of complete bacterial cell lysis and neutralization. The purified plasmid DNA is suitable for a variety of molecular biology applications, including restriction enzyme digestion, ligation, transformation and DNA sequencing.

- Simple and fast: the whole procedure can be performed in 20 minutes.
- High yield: DNA yield up to 40 µg.
- Error-free visualization: colored buffers to visualize lysis and neutralization.

### Kit Contents

| Component                                       | EM101-01 | EM101-02   |
|---|----------|------------|
| Resuspension Buffer (RB)                        | 15 ml    | 60 ml      |
| Lysis Buffer (LB, Blue)                         | 15 ml    | 60 ml      |
| Neutralization Buffer (NB, Yellow)              | 20 ml    | 80 ml      |
| Wash Buffer (WB)                                | 10 ml    | 2×20 ml    |
| Elution Buffer (EB)                             | 5 ml     | 10 ml      |
| RNase A (10 mg/ml)                              | 150 µl   | 600 µl     |
| Mini-Plasmid Spin Columns with Collection Tubes | 50 each  | 2×100 each |

### Procedures





# EasyPure<sup>®</sup> HiPure Plasmid MiniPrep Kit

EM111-01

50 rxns

## Storage

RNase A at -20°C for one year; others at room temperature (15-25°C) for one year

## Description


EasyPure<sup>®</sup> HiPure Plasmid MiniPrep Kit provides an efficient way to isolate high yield (up to 40 µg) and high quality plasmid DNA from ≤20 ml (LB) or ≤4 ml (*ArtMedia*<sup>®</sup> Plasmid Culture) of bacterial culture. Unique formulated lysis buffer and neutralization buffer permit error-free visual identification of complete bacterial cell lysis and neutralization. Endotoxin is removed by a simple incubation on column with a novel buffer. The purified plasmid DNA is suitable for a variety of molecular biology applications, including restriction enzyme digestion, ligation, transformation, DNA sequencing, and transfection.

- Fast: the whole procedure can be performed in 20 minutes.
- Simple: endotoxin is removed on column.
- High yield: DNA yield up to 40 µg.
- Error-free visualization: colored buffers to visualize lysis and neutralization.

## Kit Contents

| Component                                       | EM111-01 |
|---|----------|
| Resuspension Buffer (RB)                        | 15 ml    |
| Lysis Buffer (LB, Blue)                         | 15 ml    |
| Neutralization Buffer (NB, Yellow)              | 20 ml    |
| ToxinOut Buffer (TB)                            | 15 ml    |
| Wash Buffer (WB)                                | 10 ml    |
| Elution Buffer (EB)                             | 5 ml     |
| RNase A (10 mg/ml)                              | 150 µl   |
| Mini-Plasmid Spin Columns with Collection Tubes | 50 each  |

## Procedures



# EasyPure<sup>®</sup> HiPure Plasmid MaxiPrep Kit

EM121-01

10 rxns

## Storage

RNase A at -20°C for one year; others at room temperature (15-25°C) for one year

## Description


EasyPure<sup>®</sup> HiPure Plasmid MaxiPrep Kit uses a modified alkaline lysis method to isolate high quality plasmid DNA from ≤ 500 ml (LB) or ≤100 ml (*ArtMedia*<sup>®</sup> Plasmid Culture) of bacterial culture. Unique formulated lysis buffer and neutralization buffer permit error-free visual identification of complete bacterial cell lysis and neutralization. Endotoxin is removed by a simple incubation on column with a novel buffer. The purified DNA is suitable for a variety of molecular biology applications including restriction enzyme digestion, ligation, transformation, DNA sequencing, and transfection.

- Fast: the whole procedure can be performed in one hour.
- Simple: endotoxin is removed on column.
- High yield: DNA yield up to 1 mg.
- Error-free visualization: colored buffers to visualize lysis and neutralization.

## Kit Contents

| Component                                       | EM121-01 |
|---|----------|
| Resuspension Buffer (RB)                        | 120 ml   |
| Lysis Buffer (LB, Blue)                         | 120 ml   |
| Neutralization Buffer (NB, Yellow)              | 160 ml   |
| ToxinOut Buffer (TB)                            | 60 ml    |
| Wash Buffer (WB)                                | 25 ml    |
| Elution Buffer (EB)                             | 30 ml    |
| RNase A (10 mg/ml)                              | 1.2 ml   |
| Maxi-Plasmid Spin Columns with Collection Tubes | 10 each  |

## Procedures





# ArtMedia<sup>®</sup> Plasmid Culture

EM201-01

95 ml+5 ml

## Storage

at 2-8°C for six months

## Description

*ArtMedia<sup>®</sup>* Plasmid Culture is an enriched bacteria growth medium, which is suitable for growing various *E.coli* strains. It improves bacterial growth rate, increases cell density and obtains high yields of plasmid DNA. Under the same culture condition, *ArtMedia<sup>®</sup>* Plasmid Culture produces 3-7 folds as much of plasmid DNA as compared with traditional LB medium.

## Kit Contents

| Component | EM201-01 |
|-----------|----------|
| AM1       | 95 ml    |
| AM2       | 5 ml     |

## Suitable bacterial strains

*Trans1-T1, Trans5α, Trans10, Trans109, Trans110, Trans1-Blue, Trans2-Blue, etc.*

## DNA yield

| Medium                                      | Volume | DNA yield |
|---|--------|-----------|
| LB  | 1 ml   | ~2 µg     |
| <i>ArtMedia<sup>®</sup></i> Plasmid Culture | 0.5 ml | ~7 µg     |
| <i>ArtMedia<sup>®</sup></i> Plasmid Culture | 1 ml   | ~14 µg    |
| <i>ArtMedia<sup>®</sup></i> Plasmid Culture | 2 ml   | ~28 µg    |


# EasyPure<sup>®</sup> PCR Purification Kit

|          |          |
|----------|----------|
| EP101-01 | 50 rxns  |
| EP101-02 | 200 rxns |

### Storage

at room temperature (15-25°C) for one year

### Procedures




### Description

EasyPure<sup>®</sup> PCR Purification Kit provides a simple and fast method to purify PCR product and enzyme-digested DNA. DNA is specifically bound to silica-based column. This kit can effectively remove impurities, including proteins, organic compounds, inorganic salt ion and primers. The purified DNA is suitable for restriction enzyme digestion, ligation, transformation and sequencing.

- Effective removal of primers, dNTPs, enzymes and inorganic salt ion.
- 95%-100% recoveries for PCR fragments of 100 bp to 10 kb.
- 5 minutes procedure.
- Purified DNA ideal for using in all molecular biology experiments, including restriction enzyme digestion, ligation and sequencing.

### Kit Contents

| Component                              | EP101-01 | EP101-02   |
|--|----------|------------|
| Binding Buffer (BB)                    | 30 ml    | 120 ml     |
| Wash Buffer (WB)                       | 10 ml    | 2x20 ml    |
| Elution Buffer (EB)                    | 5 ml     | 10 ml      |
| PCR Spin Columns with Collection Tubes | 50 each  | 2x100 each |






# EasyPure<sup>®</sup> Quick Gel Extraction Kit

|          |          |
|----------|----------|
| EG101-01 | 50 rxns  |
| EG101-02 | 200 rxns |

## Storage

at room temperature (15-25°C) for one year

## Procedures




## Description

EasyPure<sup>®</sup> Quick Gel Extraction Kit is designed for rapid purification and recovery of DNA from TAE or TBE agarose gel. DNA is specifically bound to a silica-based column. The purified DNA is suitable for a variety of molecular biology applications, including restriction enzyme digestion, ligation, cloning, and DNA sequencing.

- DNA fragments size of 100 bp to 10 kb.
- Colored GSB solution (yellow) to monitor gel dissolving efficiency.
- Less than 20 minutes procedures.

## Kit Contents

| Component                               | EG101-01 | EG101-02   |
|---|----------|------------|
| Gel Solubilization Buffer (GSB, Yellow) | 30 ml    | 120 ml     |
| Wash Buffer (WB)                        | 10 ml    | 2×20 ml    |
| Elution Buffer (EB)                     | 5 ml     | 10 ml      |
| Gel Spin Columns with Collection Tubes  | 50 each  | 2×100 each |



# TransZol

ET101-01

100 ml

### Storage


at 4°C in dark for one year

### Description

*TransZol* is a ready-to-use reagent for the isolation of total RNA from cells and tissues. *TransZol* combines phenol and guanidine thiocyanate in a mono-phase solution to inhibit RNase. After lysis and centrifugation, RNA remains in the aqueous phase and others in the interphase or organic phase. RNA is precipitated by addition of isopropanol.


- Isolate RNA from a variety of species: animal, plant, yeast, bacteria and virus.
- The whole procedure can be completed in one hour.
- Simultaneous isolation of RNA, DNA and protein from the same sample.
- Pink solution for easy visualizing different phases.
- Unique dissolving solution for long-term RNA storage.

### Procedures




### Kit Contents

| Component               | ET101-01 |
|-------------------------|----------|
| <i>TransZol</i>         | 100 ml   |
| RNA Dissolving Solution | 15 ml    |



total RNA from mouse liver



total RNA from tobacco leaves





# TransZol Up


ET111-01

100 ml

## Storage

at 4°C in dark for one year

## Procedures




## Description

*TransZol Up* is a ready-to-use reagent for the isolation of total RNA from cells and tissues. Unique lysis buffer is used to disrupt cells. After centrifugation, the solution is separated into an upper colorless aqueous phase containing RNA and a lower pink organic phase. RNA is precipitated and recovered with isopropanol. Proteins can be recovered from organic phase with isopropanol. Compared with other total RNA extraction reagents, *TransZol Up* provides a powerful lysis buffer to extract RNA from a variety of species.


- Suitable for isolating RNA from a variety of species including animal, plant and bacteria.
- Superior lysis capability and higher RNA yield.
- The whole procedure can be completed in one hour.
- Pink solution for easy visualizing different phases.
- Unique dissolving solution for long-term RNA storage.

## Kit Contents

| Component               | ET111-01 |
|-------------------------|----------|
| <i>TransZol Up</i>      | 100 ml   |
| RNA Dissolving Solution | 15 ml    |



*TransZol Up* isolates RNA from rat liver



*TransZol Up* isolates RNA from HeLa cells

## RNA yield from different samples

| Material    | Amount                  | RNA yield |
|-------------|-------------------------|-----------|
| Tobacco     | 100 mg                  | ~10 µg    |
| Human blood | 200 µl                  | ~2 µg     |
| HeLa cell   | 2×10 <sup>6</sup> cells | ~10 µg    |
| Mouse liver | 100 mg                  | ~16 µg    |
| Rat liver   | 100 mg                  | ~28 µg    |

# TransZol Plant


ET121-01

100 ml

### Storage

at room temperature (15-25°C) for one year

### Procedures




### Description

*TransZol Plant* is a ready-to-use reagent for the isolation of total RNA from polysaccharide-rich and/or polyphenol plant tissues, such as champignon, banana fruit, mango fruit, potato, carrot, sansevieria. It uses a modified CTAB method to lyse samples and phenol/chloroform to remove proteins and others impurities. It is also suitable for the isolation of total RNA from animal tissues like fat, connective tissues etc.

- Superior lysis capability and higher RNA yield.
- The whole procedure can be completed in one hour.
- Pink solution for easy visualizing different phases.
- Unique dissolving solution for long-term RNA storage.

### Kit Contents

| Component               | ET121-01 |
|-------------------------|----------|
| TP I Buffer             | 100 ml   |
| TP II Buffer            | 100 ml   |
| RNA Dissolving Solution | 15 ml    |



### RNA yield from different samples

| Material                    | Amount | RNA yield |
|-----------------------------|--------|-----------|
| Papaya                      | 100 mg | ~7 µg     |
| Banana                      | 100 mg | ~8.5 µg   |
| Apple                       | 100 mg | ~4 µg     |
| Chinese yam                 | 100 mg | ~9 µg     |
| Pear                        | 100 mg | ~1.5 µg   |
| Chlorophytum comosum leaves | 100 mg | ~7.5 µg   |
| Potato tuber                | 100 mg | ~7 µg     |
| Pine needle                 | 100 mg | ~2.5 µg   |



# EasyPure<sup>®</sup> RNA Kit


ER101-01

50 rxns

## Storage

Proteinase K and DNase I solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures



## Description


EasyPure<sup>®</sup> RNA Kit provides a simple and fast column based method to isolate total RNA from animal cells, animal tissues, bacteria and yeast. Cells and tissues are enzymatically lysed. DNA is digested with DNase I. RNA is bound to silica membrane. After washing, high quality RNA is eluted from the column. RNA is free of protein contamination, and is suitable for RT-PCR, qRT-PCR and Northern blot.

## Kit Contents

| Component                              | ER101-01   |
|--|------------|
| Binding Buffer 4 (BB4)                 | 40 ml      |
| Clean Buffer 4 (CB4)                   | 60 ml      |
| Wash Buffer 4 (WB4)                    | 12 ml      |
| Proteinase K (20 mg/ml)                | 1 ml       |
| DNase I (3 units/μl)                   | 1500 units |
| DNase I Reaction Buffer                | 4×1 ml     |
| RNase-free Water                       | 10 ml      |
| RNase-free Tube (1.5 ml)               | 50 each    |
| RNA Spin Columns with Collection Tubes | 50 each    |

## Sample Requirement

| Material       | Amount             | volume of BB4/β-ME |
|----------------|--------------------|--------------------|
| Animal cell    | ≤5×10 <sup>6</sup> | 0.3-0.6 ml         |
| Animal tissue  | ≤20 mg             | 0.3-0.6 ml         |
| Bacterial cell | ≤1×10 <sup>9</sup> | 0.35 ml            |



# EasyPure<sup>®</sup> Viral DNA/RNA Kit

ER201-01

50 rxns


## Storage

Carrier RNA and Proteinase K solutions at -20°C for one year; others at room temperature (15-25°C) for one year

## Description

EasyPure<sup>®</sup> Viral DNA/RNA Kit provides a simple and fast column based method to isolate viral DNA/RNA from up to 200 µl of plasma, serum, body fluid and mammalian cell supernatant. Samples are lysed with unique lysis buffer and DNA/RNA is enriched by carrier RNA. DNA/RNA is bound to silica membrane. After washing, high quality DNA/RNA is eluted from the column. DNA/RNA is free of protein contamination, and is suitable for PCR, RT-PCR, qPCR and qRT-PCR.

## Procedures



## Kit Contents

| Component                              | ER201-01 |
|--|----------|
| Binding Buffer 5 (BB5)                 | 15 ml    |
| Wash Buffer 5 (WB5)                    | 12 ml    |
| Proteinase K (20 mg/ml)                | 1 ml     |
| Carrier RNA (1 µg/µl)                  | 310 µl   |
| RNase-free Water                       | 10 ml    |
| RNase-free Tube (1.5 ml)               | 50 each  |
| RNA Spin Columns with Collection Tubes | 50 each  |



# EasyPure<sup>®</sup> Plant RNA Kit


ER301-01

50 rxns

## Storage

DNase I at -20°C for one year; others at room temperature (15-25°C) for one year

## Procedures



## Description


*EasyPure<sup>®</sup>* Plant RNA Kit provides a simple and fast column based method to isolate RNA from plant tissue. Samples are lysed with detergent to inactivate RNase. DNA is digested with DNase I. RNA is bound to silica membrane. After washing, high quality RNA is eluted from the column. RNA is free of protein contamination, and is suitable for RT-PCR, qRT-PCR, Microarray analysis and Northern blot.

## Kit Contents

| Component                              | ER301-01   |
|--|------------|
| Binding Buffer 6 (BB6)                 | 60 ml      |
| Wash Buffer 6 (WB6)                    | 12 ml      |
| Clean Buffer 6 (CB6)                   | 60 ml      |
| DNase I (3 units/ $\mu$ l)             | 1500 units |
| DNase I Reaction Buffer                | 4x1 ml     |
| RNase-free Water                       | 10 ml      |
| RNase-free Tube (1.5 ml)               | 50 each    |
| RNA Spin Columns with Collection Tubes | 50 each    |

## Sample Requirement

| Material      | Volume of BB6/ $\beta$ -ME |
|---------------|----------------------------|
| $\leq 100$ mg | 0.5 ml                     |
| 100-200 mg    | 1 ml                       |

M: *Trans2k<sup>®</sup>* Plus II DNA Marker

1: corn leaves

2: wheat leaves

3: soybean leaves

# EasyPure<sup>®</sup> Blood RNA Kit

ER401-01

50 rxns


## Storage

DNase I at -20°C for one year; others at room temperature (15-25°C) for one year

## Description

EasyPure<sup>®</sup> Blood RNA Kit provides a simple and fast column based method to isolate total RNA from 50 µl-1.5 ml of fresh or anticoagulated blood. Blood is lysed and DNA is digested with DNase I. RNA is bound to silica membrane. After washing, high quality RNA is eluted. Purified RNA is suitable for RT-PCR, qRT-PCR and Northern blot.

## Procedures



## Kit Contents

| Component                              | ER401-01   |
|--|------------|
| Red Cell Lysis Buffer 2 (RCL2)         | 125 ml     |
| Binding Buffer 7 (BB7)                 | 40 ml      |
| Clean Buffer 7 (CB7)                   | 60 ml      |
| Wash Buffer 7 (WB7)                    | 12 ml      |
| DNase I (3 units/µl)                   | 1500 units |
| DNase I Reaction Buffer                | 4×1 ml     |
| RNase-free Water                       | 10 ml      |
| RNA Spin Columns with Collection Tubes | 50 each    |
| RNase-free Tube (1.5 ml)               | 50 each    |

## Sample Requirement

Fresh or anticoagulated blood can be kept at 4°C for one week. Do not freeze blood sample. Blood sample should be extracted as soon as possible and mixed well before use.

| Amount of Blood | Volume of BB7 |
|-----------------|---------------|
| <500 µl         | 300 µl        |
| 500 µl-1.5 ml   | 600 µl        |



# TransZol Up Plus RNA Kit

ER501-01

100 rxns

## Storage

*TransZol Up* at 4°C in dark for one year, others at room temperature (15°C-25°C) for one year

## Description


*TransZol Up Plus* RNA Kit is a ready-to-use reagent for the isolation of total RNA from cells and tissues. After lysis and centrifugation, the solution separates into an upper colorless aqueous phase (containing RNA), intermediate phase and a lower pink organic phase. RNA is specifically bound to silica-based spin column. It is a new modified version of *TransZol Up*. Compared with other total RNA extraction methods, *TransZol Up Plus* RNA Kit provides powerful lysis and easy column based purification.

- Wide application: suitable for isolating RNA from a variety of species including human, animal, plant and bacteria.
- Powerful lysis capability: complete lysis, higher RNA yield and higher purity.
- Rapid extraction: the whole procedure can be completed in one hour.
- Visible operation: pink solution for easy visualizing different phases.

## Kit Contents

| Component                              | ER501-01 |
|--|----------|
| <i>TransZol Up</i>                     | 100 ml   |
| Clean Buffer 9 (CB9)                   | 110 ml   |
| Wash Buffer 9 (WB9)                    | 24 ml    |
| RNase-free Water                       | 40 ml    |
| RNase-free Tube (1.5 ml)               | 100 each |
| RNA Spin Columns with Collection Tubes | 100 each |

## Procedures





# EasyPure<sup>®</sup> miRNA Kit

ER601-01

50 rxns

### Storage

LB10 at 4°C in dark for one year; others at room temperature (15-25°C) for one year

### Description

EasyPure<sup>®</sup> miRNA Kit provides a simple and fast column based method to isolate small RNA ( $\leq 200$  nt) from cells, tissues, fresh blood and virus. Samples are lysed with lysis buffer. The addition of chloroform to the sample separates the solution into an upper colorless aqueous phase containing RNA, an interphase and a lower organic phase. High molecular RNA (28S rRNA, 18S rRNA, mRNA) is bound to a silica membrane. Small RNA in the flow-through can be bound to a miRNA spin column.


### Kit Contents

| Component                                | ER601-01 |
|--|----------|
| Lysis Buffer 10 (LB10)                   | 55 ml    |
| Wash Buffer 10 (WB10)                    | 12 ml    |
| RNA Spin Columns with Collection Tubes   | 50 each  |
| miRNA Spin Columns with Collection Tubes | 50 each  |
| RNase-free Tube (1.5 ml)                 | 50 each  |
| RNase-free Water                         | 10 ml    |

### Sample Requirement

| Material    | Amount                |
|-------------|-----------------------|
| Tissue      | 50-100 mg             |
| Cell        | $1 \times 10^7$ cells |
| Fresh Blood | 50-200 $\mu$ l        |

### Procedures





# EasyPure<sup>®</sup> RNA Purification Kit

ER701-01

25 rxns

## Storage

at room temperature (15°C- 25°C) for one year


## Description

EasyPure<sup>®</sup> RNA Purification Kit uses silica-based spin column for specific RNA binding. The kit can be used for RNA purification from DNase I-treated total RNA, *in vitro* transcription product, RNA-labelled product, synthetic RNA. This kit permits effective removal of proteins, organic chemicals, inorganic salt ion and other impurities. Purified RNA is suitable for RT-PCR, qRT-PCR, Northern blot and other applications.

## Kit Contents

| Component                              | ER701-01 |
|--|----------|
| Binding Buffer 12 (BB12)               | 10 ml    |
| Wash Buffer 12 (WB12)                  | 8 ml     |
| RNase-free Water                       | 1.5 ml   |
| RNase-free Tube (1.5 ml)               | 25 each  |
| RNA Spin Columns with Collection Tubes | 25 each  |

## Procedures



# RNAhold<sup>®</sup>

EH101-01

100 ml

## Storage

at room temperature for one year


## Description

*RNAhold<sup>®</sup>* is an aqueous, nontoxic tissue preservation solution. It can inactivate RNase and keep RNA intact by permeating cells and tissues. Cells and tissues can be stored at this solution for one week at room temperature without RNA degradation. It can be used for RNA preservation with bacteria, cells and most fresh animal tissues.

- Immediate RNase inactivation.
- Sample can be stored at room temperature for 1 week, 2-8°C for 1 month, -20°C or -80°C for long term storage.
- Ideal for field sample collection.

## Note

Tissues stored in *RNAhold<sup>®</sup>* solution can freeze and thaw at least 20 times without significantly affecting the yield or the integrity of the recoverable RNA.



M: *Trans2K<sup>®</sup>* Plus II DNA Marker

1: HeLa cells stored at 37°C for 1 day with *RNAhold<sup>®</sup>*

2: HeLa cells stored at room temperature for 1 week with *RNAhold<sup>®</sup>*

3: HeLa cells stored at room temperature for 1 week without *RNAhold<sup>®</sup>*

---

## Chapter 5 Gene Expression

---

### Prokaryotic Expression Vectors

*pEASY*<sup>®</sup>-Blunt E1 Expression Kit .....151

*pEASY*<sup>®</sup>-Blunt E2 Expression Kit .....154

### Expression Medium

*ArtMedia*<sup>®</sup> Protein Expression .....155

### Expression Competent Cells

BL21(DE3) Chemically Competent Cell .....156

BL21(DE3) pLysS Chemically Competent Cell .....156

*Transetta*(DE3) Chemically Competent Cell .....157





*TransB*(DE3) Chemically Competent Cell .....157

BL21 Chemically Competent Cell .....157

### Mammalian Expression Vectors

*pEASY*<sup>®</sup>-Blunt M2 Expression Kit .....158

*pEASY*<sup>®</sup>-Blunt M3 Expression Kit .....161



### Feature and application of *pEASY*<sup>®</sup> expression vectors

| Name                                | Amp <sup>+</sup> | Promoter     | Sequencing primer                              | Characteristics                             | Application            |
|-------------------------------------|------------------|--------------|--|---|------------------------|
| <i>pEASY</i> <sup>®</sup> -Blunt E1 | +                | <i>T7lac</i> | T7 Promoter Primer;<br>T7 Terminator Primer    | N-terminal 6xHis tag                        | Prokaryotic Expression |
| <i>pEASY</i> <sup>®</sup> -Blunt E2 | +                | <i>T7lac</i> | T7 Promoter Primer;<br>T7 Terminator Primer    | C-terminal 6xHis tag                        | Prokaryotic Expression |
| <i>pEASY</i> <sup>®</sup> -Blunt M2 | +                | Enhanced CMV | CMV Forward Primer;<br>TK PolyA Reverse Primer | C-terminal Myc tag;<br>Neomycin resistance  | Mammalian Expression   |
| <i>pEASY</i> <sup>®</sup> -Blunt M3 | +                | Enhanced CMV | CMV Forward Primer;<br>TK PolyA Reverse Primer | C-terminal Flag tag;<br>Neomycin resistance | Mammalian Expression   |



# pEASY<sup>®</sup>-Blunt E1 Expression Kit

CE111-01

10 rxns

## Storage

*Trans*1-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description



pEASY<sup>®</sup>-Blunt E1 Expression Vector is constructed from pET vector, it utilizes a highly efficient, five-minute blunt cloning strategy to clone PCR product into high-efficient expression vector. The size of control insert is 750 bp, and expressed target protein is about 27 kDa.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- Ampicillin resistance.
- T7 promoter primer and T7 terminator primer for sequencing.
- Bacteriophage T7 *lac* promoter for high level expression.
- N-terminal 6xHis tag for easy purification.
- *Trans*1-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency ( $>10^9$  cfu/ $\mu$ g pUC19 DNA) and fast growing.
- E1 Expression Plasmid included as negative control.

## Kit Contents

| Component   | CE111-01               |
|---|------------------------|
| pEASY <sup>®</sup> -Blunt E1 Expression Vector (15 ng/ $\mu$ l) | 10 $\mu$ l             |
| EControl Template (5 ng/ $\mu$ l)                               | 10 $\mu$ l             |
| EControl Forward Primer (10 $\mu$ M)                            | 10 $\mu$ l             |
| EControl Reverse Primer (10 $\mu$ M)                            | 10 $\mu$ l             |
| E1 Expression Plasmid (Negative Control) (15 ng/ $\mu$ l)       | 10 $\mu$ l             |
| T7 Promoter Primer (10 $\mu$ M)                                 | 50 $\mu$ l             |
| T7 Terminator Primer (10 $\mu$ M)                               | 50 $\mu$ l             |
| <i>Trans</i> 1-T1 Phage Resistant Chemically Competent Cell     | 5 $\times$ 100 $\mu$ l |

## pEASY<sup>®</sup>-Blunt E1 Prokaryotic Expression Vector Map



T7 promoter: bases 209-225  
 T7 transcription start: base 226  
 Lac operator(lacO): bases 228-252  
 RBS: bases 282-288  
 His-Tag coding sequence: bases 309-326  
 T7 terminator: bases 436-482  
 Ampicillin resistance ORF: bases 907-1,767  
 pBR322 origin: bases 1,922-2,541  
 ROP ORF: bases 2,953-3,144  
 LacI ORF: bases 4,459-5,547



## PROTOCOL

### Cloning reaction

- (1) Primer requirement: primer cannot be phosphorylated
- (2) PCR Enzyme: high fidelity *Pfu* DNA polymerase
- (3) Reaction conditions: for higher cloning efficiency, we recommend 5-10 minutes post PCR 72°C extension. After PCR, use agarose gel electrophoresis to verify the quality and quantity of PCR product.

### Suggested cloning reaction conditions

1. Optimal amount of insert  
Molar ratio of vector and insert = 1:7 (1 kb, ~20 ng; 2 kb, ~40 ng)
2. Optimal volume of vector: 1 µl
3. Optimal reaction volume: 3~5 µl
4. Optimal incubation time
  - (1) 0.1~1 kb (including 1 kb): 5~10 minutes
  - (2) 1~2 kb (including 2 kb): 10~15 minutes
  - (3) 2~3 kb (including 3 kb): 15~20 minutes
  - (4) ≥ 3 kb: 20~30 minutes  
Use the maximum incubation time if the insert is gel purified PCR product.
5. Optimal incubation temperature: for most PCR inserts, the optimal temperature is about 25°C; for some PCR inserts, optimal results can be achieved with higher temperature (up to 37°C).

### Transformation

1. Add the ligated products to 50 µl of *Trans1-T1* Phage Resistant Chemically Competent Cell and mix gently (do not mix by pipetting up and down).
2. Incubate on ice for 20-30 minutes.
3. Heat-shock the cells at 42°C for 30 seconds.
4. Immediately place the tube on ice for 2 minutes.
5. Add 250 µl of room temperature SOC or LB medium. Shake the tube at 37°C (200 rpm) for 1 hour.
6. Pre-warm a selective LB plate at 37°C for 30 minutes.
7. Spread 200 µl or all transformants on the pre-warmed plate. Incubate at 37°C overnight.

### Analysis of positive clones

1. Transfer 5~10 colonies into 10 µl ddH<sub>2</sub>O.
2. Use 1 µl of the mixture as template for 25 µl PCR using T7 promoter primer and gene reverse primer, or gene forward primer and T7 terminator primer.
3. PCR
 

|               |        |   |           |
|---------------|--------|---|-----------|
| 94°C          | 10 min | } | 30 cycles |
| 94°C          | 30 sec |   |           |
| 55°C          | 30 sec |   |           |
| 72°C          | X min* |   |           |
| 72°C 5-10 min |        |   |           |
- \*(depends on the insert size and PCR enzymes)
4. Analyze positive clones by restriction enzyme digestion and DNA sequencing.





### Target gene expression

#### 1. Competent cell

BL21(DE3) competent cell series are suitable for prokaryotic protein expression.

#### 2. Protein expression

##### Method 1

- Pick single colony and transfer into 5 ml of LB/Amp<sup>+</sup> medium and shake at 37°C (250 rpm) until OD<sub>600</sub> close to 0.5.
- Add IPTG to a final concentration of 0.5-1 mM and shake at 37°C for 3-5 hours.
- Remove a 500 µl aliquot during different time course and centrifuge at the maximum speed.

##### Method 2

- Pick single colony and transfer into 5 ml of *ArtMedia*<sup>®</sup> Protein Expression/Amp<sup>+</sup> medium, incubate at 37°C overnight.

#### 3. Check expression

Aspirate the supernatant and use the pellets for SDS-PAGE.

### PCR for control insert (750 bp)

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| EControl Template (5 ng/µl)                                   | 1 µl     | 0.1 ng/µl           |
| EControl Forward Primer (10 µM)                               | 1 µl     | 0.2 µM              |
| EControl Reverse Primer (10 µM)                               | 1 µl     | 0.2 µM              |
| 2× <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix | 25 µl    | 1×                  |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total Volume  | 50 µl    | -                   |

### Thermal cycling conditions

|      |         |             |
|------|---------|-------------|
| 94°C | 2-5 min | } 30 cycles |
| 94°C | 20 sec  |             |
| 55°C | 20 sec  |             |
| 72°C | 30 sec  |             |
| 72°C | 10 min  |             |

Ligate 1 µl of control PCR insert with 1 µl vector. Hundreds of colonies should be produced with cloning efficiency over 90%.

# pEASY<sup>®</sup>-Blunt E2 Expression Kit

CE211-01

10 rxns

### Storage

*Trans*1-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

### Description


pEASY<sup>®</sup>-Blunt E2 Expression Vector is constructed from pET vector. It utilizes a highly efficient, five-minute blunt cloning strategy to clone PCR product into high-efficient expression vector. The size of control insert is 750 bp, and expressed target protein is about 27 kDa.

- 5 minutes fast ligation of *Pfu*-amplified PCR products.
- Ampicillin resistance.
- T7 promoter primer and T7 terminator primer for sequencing.
- Bacteriophage T7lac promoter for high level expression.
- C-terminal 6xHis tag for easy purification.
- *Trans*1-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.
- E2 Expression Plasmid included as negative control.


### Kit Contents

| Component   | CE211-01 |
|---|----------|
| pEASY <sup>®</sup> -Blunt E2 Expression Vector (15 ng/μl)   | 10 μl    |
| EControl Template (5 ng/μl)                                 | 10 μl    |
| EControl Forward Primer (10 μM)                             | 10 μl    |
| EControl Reverse Primer (10 μM)                             | 10 μl    |
| E2 Expression Plasmid (Negative Control) (15 ng/μl)         | 10 μl    |
| T7 Promoter Primer (10 μM)                                  | 50 μl    |
| T7 Terminator Primer (10 μM)                                | 50 μl    |
| <i>Trans</i> 1-T1 Phage Resistant Chemically Competent Cell | 5×100 μl |

pEASY<sup>®</sup>-Blunt E2 Prokaryotic Expression Vector Map



T7 promoter: bases 5,117-5,133  
 T7 transcription start: base 5,134  
 Lac operator(lacO): bases 5,136-5,160  
 RBS: bases 5,190-5,196  
 His-Tag coding sequence: bases 5,238-5,255  
 T7 terminator: bases 5,323-5,369  
 ROP ORF: bases 2,648-2,839  
 LacI ORF: bases 3,651-4,739  
 pBR origin: bases 1,614-2,233  
 Ampicillin resistance ORF: bases 599-1,459  
 f1 origin: bases 13-450



High quality products



## PROTOCOL

Protocols for cloning, transformation, analysis and expression are the same as described on page 152-153.

# ArtMedia<sup>®</sup> Protein Expression

|          |            |
|----------|------------|
| CP101-01 | 95 ml+5 ml |
|----------|------------|

### Storage

at 2-8°C for six months

### Description

ArtMedia<sup>®</sup> Protein Expression is designed for higher protein yield with much less hands-on time. Protein is induced automatically without time-consuming OD monitoring and IPTG induction steps. Simply inoculate prepared ArtMedia<sup>®</sup> with colonies, grow the culture overnight and harvest cells for protein purification.

### Kit Contents


| Component | CP101-01 |
|-----------|----------|
| AM3       | 95 ml    |
| AM4       | 5 ml     |

### Suitable expression vectors


Lactose operons expression vectors: *pEASY<sup>®</sup>-Blunt E1*, *pEASY<sup>®</sup>-Blunt E2*, *pET*, *pGEX*, *pMAL*

### Strains

BL21 competent cell series.



Target protein: 36 kDa  
 Vector: *pGEX-5X-3 (tac promoter)*  
 Strain: BL21  
 M: *ProteinRuler<sup>®</sup> II*  
 Lane 1: LB medium only  
 Lane 2: LB, OD<sub>600</sub>=0.5, induced with 1 mM IPTG, 37°C for 12 hours  
 Lane 3: 37°C for 12 hours with ArtMedia<sup>®</sup>  
 Protein Expression



Target protein: 70 kDa  
 Vector: *pEASY<sup>®</sup>-Blunt E1 (T7lac promoter)*  
 Strain: *Transtetta(DE3)*  
 M: *ProteinRuler<sup>®</sup> II*  
 Lane 1: LB medium only  
 Lane 2: LB, OD<sub>600</sub>=0.5, induced with 1 mM IPTG, 37°C for 12 hours  
 Lane 3: 37°C for 12 hours with ArtMedia<sup>®</sup>  
 Protein Expression

### Notes

- AM4 may have a slight precipitate, which will not affect performance. If precipitate is observed, warm the bottle in a 37°C water bath to dissolve the precipitate.
- Add the whole volume of AM4 to AM3 for complete medium. Store the complete medium at 4°C up to 1 month.

# Expression Competent Cells

## Selection Guide

| Name                   | Cat. No. | Transformation Efficiency  | Application  |
|------------------------|----------|----------------------------|--|
| BL21(DE3)              | CD601    | 10 <sup>7</sup> cfu/μg DNA | High expression of non-toxic protein   |
| BL21(DE3) pLysS        | CD701    | 10 <sup>7</sup> cfu/μg DNA | High expression of toxic protein and non-toxic protein, low background   |
| <i>Transetta</i> (DE3) | CD801    | 10 <sup>7</sup> cfu/μg DNA | Contains tRNAs corresponding to 6 rare codons, application to eukaryotic gene expression                         |
| <i>TransB</i> (DE3)    | CD811    | 10 <sup>7</sup> cfu/μg DNA | Conducive to the formation of the correctly folded protein with disulfide, enhance the solubility of the protein |
| BL21                   | CD901    | 10 <sup>7</sup> cfu/μg DNA | High expression of toxic protein   |

## BL21(DE3) Chemically Competent Cell

|          |           |
|----------|-----------|
| CD601-02 | 10×100 μl |
| CD601-03 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency: >10<sup>7</sup> cfu/μg (pUC19 DNA).
- DE3 strains contains the λDE3 lysogen that carries the gene for T7 RNA polymerase.
- Suitable for T7 and T7lac such as pET, pEASY®.
- Suitable for high expression of non-toxic protein.
- Control plasmid I (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 25 kDa.

### Genotype

F<sup>-</sup> ompT hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup>m<sub>B</sub><sup>-</sup>) gal dcm(DE3)

## BL21(DE3) pLysS Chemically Competent Cell

|          |           |
|----------|-----------|
| CD701-02 | 10×100 μl |
| CD701-03 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency: >10<sup>7</sup> cfu/μg (pUC19 DNA).
- Cam<sup>r</sup>.
- Contains pLysS plasmid that expresses the T7 lysozyme gene to reduce the background of the target gene's expression without disturbing IPTG functioning.
- Suitable for non-toxic and toxic protein expression.
- Control plasmid I (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 25 kDa.

### Genotype

F<sup>-</sup> ompT hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup>m<sub>B</sub><sup>-</sup>) dcm(DE3) gal pLysS(Cam<sup>r</sup>)

High quality products



## Transetta(DE3) Chemically Competent Cell

|          |           |
|----------|-----------|
| CD801-02 | 10×100 μl |
| CD801-03 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency:  $>10^7$  cfu/μg (pUC19 DNA).
- Cam<sup>R</sup>.
- tRNAs for 6 rare codons AUA, AGG, AGA, CUA, CCC, GGA. Enhance the expression level of proteins in the prokaryotic system.
- Control plasmid I (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 25 kDa.

### Genotype

F<sup>-</sup> *ompT hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup> m<sub>B</sub><sup>-</sup>) gal dcm* (DE3) pRARE (argU, argW, ileX, glyT, leuW, proL) (Cam<sup>R</sup>)

## TransB(DE3) Chemically Competent Cell

|          |           |
|----------|-----------|
| CD811-02 | 10×100 μl |
|----------|-----------|

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency:  $>10^7$  cfu/μg (pUC19 DNA).
- Kan<sup>R</sup> and Tet<sup>R</sup>.
- Thioredoxin reductase (*trxB*) and glutathione reductase (*gor*) mutation greatly facilitates cytoplasmic disulfide bond formation.
- Control plasmid I (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 25 kDa.

### Genotype

F<sup>-</sup> *ompT hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup> m<sub>B</sub><sup>-</sup>) gal dcm lacY1 ahpC* (DE3) *gor522::Tn10 trxB* (Kan<sup>R</sup>, Tet<sup>R</sup>)

## BL21 Chemically Competent Cell

|          |           |
|----------|-----------|
| CD901-02 | 10×100 μl |
| CD901-03 | 20×100 μl |

### Storage

at -70°C for six months

### Characteristics

- Transformation efficiency:  $>10^7$  cfu/μg (pUC19 DNA).
- Tet<sup>R</sup>.
- Tight expression control ideal for toxic protein expression.
- Control plasmid II (Amp<sup>r</sup>) is used for detection of expression function of cell. The protein size is about 26 kDa.

### Genotype

*E. coli* B F<sup>-</sup> *dcm ompT hsdS (r<sub>B</sub><sup>-</sup> m<sub>B</sub><sup>-</sup>) gal araB::T7RNAP-tetA*

# pEASY<sup>®</sup>-Blunt M2 Expression Kit

CM211-01

10 rxns

### Storage

*Trans*1-T1 Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

### Description


pEASY<sup>®</sup>-Blunt M2 Expression Vector utilizes a highly efficient, five-minute blunt cloning strategy to clone PCR product into high-efficient expression vector under regulation of enhanced CMV promoter. The size of control insert is 750 bp, and expressed target protein is about 27 kDa.

- 5 minutes fast ligation of *Pfu*-amplified products.
- Enhanced CMV promoter for higher protein yield.
- CMV forward primer and TK polyA reverse primer for sequencing.
- C-terminal Myc tag for protein detection.
- Neomycin resistance gene for stable cell line selection.
- *Trans*1-T1 Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.
- M2 Expression Plasmid included as negative control.

### Kit Contents

| Component   | CM211-01 |
|---|----------|
| pEASY <sup>®</sup> -Blunt M2 Expression Vector (15 ng/μl)   | 10 μl    |
| MControl Template (5 ng/μl)                                 | 10 μl    |
| MControl Forward Primer (10 μM)                             | 10 μl    |
| MControl Reverse Primer (10 μM)                             | 10 μl    |
| M2 Expression Plasmid (Negative Control) (15 ng/μl)         | 10 μl    |
| CMV Forward Primer (10 μM)                                  | 50 μl    |
| TK PolyA Reverse Primer (10 μM)                             | 50 μl    |
| <i>Trans</i> 1-T1 Phage Resistant Chemically Competent Cell | 5×100 μl |

### pEASY<sup>®</sup>-Blunt M2 Mammalian Expression Vector Map



CMV promoter: bases 4,740-5,327  
 CMV forward primer binding site: bases 5,277-5,297  
 Cloning site: bases 6-7  
 TK PolyA reverse primer binding site: bases 73-91  
 TK polyadenylation signal: bases 66-337  
 f1 replication origin: bases 373-801  
 SV40 early promoter: bases 828-1,136  
 Neomycin resistance gene: bases 1,211-2,005  
 SV40 polyadenylation signal: bases 2,181-2,311  
 pUC origin: bases 2,694-3,367  
 Ampicillin (*bla*) resistance gene(c): bases 3,512-4,372  
*bla* promoter(c): bases 4,373-4,471  
 Myc epitope: bases 19-48  
 (c) = complementary strand

CMV Forward Primer Binding Site  
 TGA CGC AAA TGG GCG GTA GGC GTG TAC GGT GGG AGG TCT ATA TAA GCA GAG CTC GTT  
 ACT GCG TTT ACC CGC CAT CCG CAC ATG CCA CCC TCC AGA TAT ATT CGT CTC GAG CAA

TAG TGA ACC GTC AGA TCG CCT GGA GAC GCC ATC CAC GCT GTT TTG ACC TCC ATA GAA GAC ACC GGG ACC GAT CCA GCC TCC GGA CTC TAG AGG ATC  
 ATC ACT TGG CAG TCT AGC GGA CCT CTG CGG TAG GTG CGA CAA AAC TGG AGG TAT CTT CTG TGG CCC TGG CTA GGT CGG AGG CCT GAG ATC TCC TAG

Myc epitope  
 GCC CTT AAG GGC GAT CCG GAA CAA AAA CTC ATC TCA GAA GAG GAT CTG  
 CGG GAA PCR Product TTC CCG CTA GGC CTT GTT TTT GAG TAG AGT CTT CTC CTA GAC

TK PolyA TK PolyA Reverse Primer Binding Site  
 TAG TAA TGA GTT TAA ACG GGG GAG GCT AAC TGA AAC ACG GAA GGA GAC AAT ACC GGA AGG AAC CCG CG.....TAG C  
 ATC ATT ACT CAA ATT TGC CCC CTC CGA TTG ACT TTG TGC CTT CCT CTG TTA TGG CCT TCC TTC CTT GG.....ATC G

High quality products



## PROTOCOL

### Cloning reaction

1. PCR primer design
  - (1) Do not add 5' phosphates to PCR primers.
  - (2) Forward primer with Kozak consensus sequence: (G/A)NNATGN.
2. Using high fidelity *Pfu* DNA polymerase
3. Add PCR products and vector in a tube, mix gently and incubate at room temperature (20-37°C) for 5 minutes.

### Suggested cloning reaction conditions

1. Optimal amount of insert
 

Molar ratio of vector to insert = 1:7 (1 kb, ~20 ng; 2 kb, 40 ng)
2. Optimal volume of vector: 1  $\mu$ l
3. Optimal reaction volume: about 3~5  $\mu$ l
4. Optimal incubation time
  - (1) 0.1~1 kb (including 1 kb): 5~10 minutes
  - (2) 1~2 kb (including 2 kb): 10~15 minutes
  - (3) 2~3 kb (including 3 kb): 15~20 minutes
  - (4)  $\geq$ 3 kb: 20~30 minutes

Use the maximum incubation time if the insert is gel purified PCR product.
5. Optimal incubation temperature: for most PCR inserts, the optimal temperature is about 25°C; for some PCR inserts, optimal results can be achieved with higher temperature (up to 37°C).

### Transformation

1. Add the ligated products to 50  $\mu$ l of *Trans1*-T1 Phage Resistant Chemically Competent Cell and mix gently (do not mix by pipetting up and down).
2. Incubate on ice for 20-30 minutes.
3. Heat-shock the cells at 42°C for 30 seconds.
4. Immediately place the tube on ice for 2 minutes.
5. Add 250  $\mu$ l of room temperature SOC or LB medium. Shake the tube at 37°C (200 rpm) for 1 hour.
6. Pre-warm a selective LB plate at 37°C for 30 minutes.
7. Spread 200  $\mu$ l or all transformants on the pre-warmed plate. Incubate at 37°C overnight.

### Analysis of positive clones

1. Transfer 5-10 colonies into 10  $\mu$ l ddH<sub>2</sub>O.
2. Use 1  $\mu$ l of the mixture as template for 25  $\mu$ l PCR using CMV forward primer and gene reverse primer, or gene forward primer and TK polyA reverse primer.

### 3. PCR

|      |          |             |
|------|----------|-------------|
| 94°C | 10 min   | } 30 cycles |
| 94°C | 30 sec   |             |
| 55°C | 30 sec   |             |
| 72°C | X min*   |             |
| 72°C | 5-10 min |             |

\*(depends on the insert size and PCR enzymes)

4. Analyze positive clones by restriction enzyme digestion and DNA sequencing.

**Transfection**

See *TransLipid*<sup>®</sup> HL Transfection Reagent for the detailed protocol.

**Detection of target protein**

Anti-Myc antibody for the detection of proteins with Myc tag.

**Selection of stable cell lines**

G418 to select stable cell lines.

**PCR for control insert (750 bp)**

| Component   | Volume   | Final Concentration |
|---|----------|---------------------|
| MControl Template (5 ng/μl)                                   | 1 μl     | 0.1 ng/μl           |
| MControl Forward Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| MControl Reverse Primer (10 μM)                               | 1 μl     | 0.2 μM              |
| 2x <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix | 25 μl    | 1x                  |
| ddH <sub>2</sub> O  | Variable | -                   |
| Total Volume  | 50 μl    | -                   |

**Thermal cycling conditions**

|      |         |             |
|------|---------|-------------|
| 94°C | 2-5 min | } 30 cycles |
| 94°C | 20 sec  |             |
| 55°C | 20 sec  |             |
| 72°C | 30 sec  |             |
| 72°C | 10 min  |             |

Ligate 1 μl of control PCR insert with 1 μl vector. Hundreds of colonies should be produced with cloning efficiency over 90%.





# pEASY<sup>®</sup>-Blunt M3 Expression Kit

CM311-01

10 rxns

## Storage

*Trans1-T1* Phage Resistant Chemically Competent Cell at -70°C for six months; others at -20°C for nine months

## Description


pEASY<sup>®</sup>-Blunt M3 Expression Vector utilizes a highly efficient, five-minute blunt cloning strategy to clone PCR product into high-efficient expression vector under regulation of enhanced CMV promoter. The size of control insert is 750 bp, and expressed target protein is about 27 kDa.

- 5 minutes fast ligation of *Pfu*-amplified products.
- Enhanced CMV promoter for higher protein yield.
- CMV forward primer and TK polyA reverse primer for sequencing.
- C-terminal Flag-tag for protein detection.
- Neomycin resistance gene for stable cell line selection.
- *Trans1-T1* Phage Resistant Chemically Competent Cell, high transformation efficiency (>10<sup>9</sup> cfu/μg pUC19 DNA) and fast growing.
- M3 Expression Plasmid included as negative control.

## Kit Contents

| Component  | CM311-01 |
|--|----------|
| pEASY <sup>®</sup> -Blunt M3 Expression Vector (15 ng/μl)  | 10 μl    |
| MControl Template (5 ng/μl)                                | 10 μl    |
| MControl Forward Primer (10 μM)                            | 10 μl    |
| MControl Reverse Primer (10 μM)                            | 10 μl    |
| M3 Expression Plasmid (Negative Control) (15 ng/μl)        | 10 μl    |
| CMV Forward Primer (10 μM)                                 | 50 μl    |
| TK PolyA Reverse Primer (10 μM)                            | 50 μl    |
| <i>Trans1-T1</i> Phage Resistant Chemically Competent Cell | 5x100 μl |

## pEASY<sup>®</sup>-Blunt M3 Mammalian Expression Vector Map



CMV promoter: bases 4,734-5,321  
 CMV forward primer binding site: bases 5,271-5,291  
 Cloning site: bases 6-7  
 TK PolyA reverse primer binding site: bases 67-85  
 TK polyadenylation signal: bases 60-331  
 SV40 early promoter: bases 822-1,130  
 Neomycin resistance gene: bases 1,205-1,999  
 SV40 polyadenylation signal: bases 2,175-2,305  
 pUC origin: bases 2,688-3,361  
 Ampicillin (bla) resistance gene(c): bases 3,506-4,366  
 bla promoter(c): bases 4,367-4,465  
 Flag epitope: bases 19-42  
 (c) = complementary strand

### CMV Forward Primer Binding Site

TGA CGC AAA TGG GCG GTA GGC GTG TAC GGT GGG AGG TCT ATA TAA GCA GAG CTC GTT  
 ACT GCG TTT ACC CGC CAT CCG CAC ATG CCA CCC TCC AGA TAT ATT CGT CTC GAG CAA

TAG TGA ACC GTC AGA TCG CCT GGA GAC GCC ATC CAC GCT GTT TTG ACC TCC ATA GAA GAC ACC GGG ACC GAT  
 ATC ACT TGG CAG TCT AGC GGA CCT CTG CGG TAG GTG CGA CAA AAC TGG AGG TAT CTT CTG TGG CCC TGG CTA

CCA GCC TCC GGA CTC TAG AGG ATC GCC CTT AAG GGC GAT CCG GAT TAC AAG GAC GAT GAC GAT AAG GAA TTC  
 GGT CGG AGG CCT GAG ATC TOC TAG CGG GAA TTC CCG CTA GGC CTA ATG TTT CTG CTA CTG CTA TTT CTT AAG

### TK PolyA TK PolyA Reverse Primer Binding Site

TAG TAA TGA GTT TAA ACS GGG GAG GCT AAC TGA AAC ACG GAA GGA GAC AAT ACC GGA AGG AAC CCG CG.....TAG C  
 ATC ATT ACT CAA ATT TGC CCC CTC CGA TTG ACT TTG TGG CTT CCT CTG TTA TGG CCT TCC TTC CTT GG.....ATC G

## PROTOCOL

**Protocols for cloning and transformation are the same as described on page 159.**

### Analysis of positive clones

1. Transfer 5-10 colonies into 10  $\mu$ l ddH<sub>2</sub>O.
2. Use 1  $\mu$ l of the mixture as template for 25  $\mu$ l PCR using CMV forward primer and gene reverse primer, or gene forward primer and TK polyA reverse primer.
3. PCR
 

|      |          |           |
|------|----------|-----------|
| 94°C | 2-5 min  |           |
| 94°C | 30 sec   | }         |
| 55°C | 30 sec   |           |
| 72°C | X min*   |           |
|      |          | 30 cycles |
| 72°C | 5-10 min |           |

\*(depends on the insert size and PCR enzymes)

4. Analyze positive clones by restriction enzyme digestion and DNA sequencing.

### Transfection

See *TransLipid*<sup>®</sup> HL Transfection Reagent for the detailed protocol.

### Detection of target protein

Anti-Flag antibody for the detection of proteins with Flag tag.

### Selection of stable cell lines

G418 to select stable cell lines.

---

## Chapter 6 Protein Extraction, Purification and Detection

---

### Protein Extraction

|  |     |
|--|-----|
| <i>ProteinExt</i> <sup>™</sup> Mammalian Total Protein Extraction Kit .....                      | 164 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Nuclear and Cytoplasmic Protein<br>Extraction Kit ..... | 165 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Membrane Protein<br>Extraction Kit .....                | 166 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit for<br>Cultured Cells .....  | 167 |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit<br>for Tissue .....          | 168 |

### Protein Purification

|  |     |
|--|-----|
| <i>ProteinIso</i> <sup>®</sup> Ni-NTA Resin .....    | 169 |
| <i>ProteinIso</i> <sup>®</sup> Ni-IDA Resin .....    | 171 |
| <i>ProteinIso</i> <sup>®</sup> GST Resin .....       | 173 |
| <i>ProteinIso</i> <sup>®</sup> Protein A Resin ..... | 175 |
| <i>ProteinIso</i> <sup>®</sup> Protein G Resin ..... | 177 |

### Unstained Protein Marker

|   |     |
|---|-----|
| <i>ProteinRuler</i> <sup>®</sup> I .....  | 180 |
| <i>ProteinRuler</i> <sup>®</sup> II ..... | 180 |
| <i>ProteinRuler</i> <sup>®</sup> IV ..... | 181 |

### Prestained Protein Marker

|  |     |
|--|-----|
| <i>Blue Plus</i> <sup>®</sup> Protein Marker .....     | 182 |
| <i>Blue Plus</i> <sup>®</sup> II Protein Marker .....  | 182 |
| <i>Blue Plus</i> <sup>®</sup> III Protein Marker ..... | 183 |
| <i>Blue Plus</i> <sup>®</sup> IV Protein Marker .....  | 183 |

### Western Protein Marker

|   |     |
|---|-----|
| <i>EasySee</i> <sup>®</sup> Western Marker .....    | 184 |
| <i>EasySee</i> <sup>®</sup> II Western Marker ..... | 185 |

### Related Products

|   |     |
|---|-----|
| <i>EasySee</i> <sup>®</sup> Western Blot Kit .....            | 186 |
| 6×Protein Loading Buffer .....                                | 186 |
| Easy Protein Quantitative Kit (Bradford) .....                | 187 |
| Easy II Protein Quantitative Kit (BCA) .....                  | 188 |
| <i>ProteinEle</i> <sup>™</sup> Precast Tris-Glycine Gel ..... | 190 |

# ProteinExt™ Mammalian Total Protein Extraction Kit

DE101-01

100 ml

### Storage

at -20°C for one year


### Description

*ProteinExt™* Mammalian Total Protein Extraction Kit provides a fast and efficient method to extract total proteins (cytoplasmic, membrane and nuclear proteins) from mammalian cells and tissues without ultracentrifugation. The extracted proteins are suitable for SDS-PAGE, Western blot, ELISA, and other functional assays.

### Kit Content

| Component                              | DE101-01 |
|--|----------|
| Total Protein Extraction Buffer (TPEB) | 100 ml   |

### Procedures





# ProteinExt™ Mammalian Nuclear and Cytoplasmic Protein Extraction Kit

DE201-01

50 rxns

## Storage

at 2-8°C for one year


## Description

ProteinExt™ Mammalian Nuclear and Cytoplasmic Protein Extraction Kit provides a fast and efficient method to extract nuclear and cytoplasmic proteins from mammalian cells and tissues. Native proteins can be obtained within 80 minutes without ultracentrifugation. The extracted proteins are suitable for a variety of downstream applications, including SDS-PAGE, Western blot, ELISA, enzyme-activity assays, immunoprecipitation and transcription factor activity analysis.

## Kit Contents

| Component  | DE201-01 |
|--|----------|
| Cytoplasmic Protein Extraction Buffer I (CPEB I)   | 75 ml    |
| Cytoplasmic Protein Extraction Buffer II (CPEB II) | 3 ml     |
| Nuclear Protein Extraction Buffer (NPEB)           | 25 ml    |

## Procedures



# ProteinExt™ Mammalian Membrane Protein Extraction Kit

DE301-01

50 rxns

## Storage

at 2-8°C for one year


## Description

ProteinExt™ Mammalian Membrane Extraction Kit provides a fast and efficient method to extract membrane proteins from mammalian cells and tissues. Native proteins can be obtained within 70 minutes without ultracentrifugation. Membrane proteins with at least 1-2 transmembrane domains are typically extracted with an efficiency of up to 90%. Cross-contamination of cytosolic proteins into the membrane fraction is usually less than 10%. The extracted proteins are suitable for a variety of downstream applications, including SDS-PAGE, Western blot, ELISA, and enzyme-activity assays.

## Kit Contents

| Component   | DE301-01 |
|---|----------|
| Membrane Protein Extraction Buffer I (MPEB I)     | 50 ml    |
| Membrane Protein Extraction Buffer II (MPEB II)   | 7.5 ml   |
| Membrane Protein Extraction Buffer III (MPEB III) | 15 ml    |

## Procedures





# *ProteinExt*<sup>TM</sup> Mammalian Mitochondria Isolation Kit for Cultured Cells

DE401-01

50 rxns

## Storage

MSB at -20°C for one year, others at 2-8°C for one year


## Description

*ProteinExt*<sup>TM</sup> Mammalian Mitochondria Isolation Kit for Cultured Cells provides a fast and efficient method to isolate mitochondria from cultured mammalian cells. This kit provides two options for the separation of mitochondria from cytosolic components: a reagent-based method or homogenization-based method. Reagent-based method uses a mild procedure to process single or multiple samples. The isolated mitochondria is suitable for a variety of downstream applications, including protein analysis, apoptosis, signal transduction and metabolic assays.

## Kit Contents

| Component                                   | DE401-01 |
|---|----------|
| Mitochondria Isolation Buffer I (MIB I)     | 50 ml    |
| Mitochondria Isolation Buffer II (MIB II)   | 500 µl   |
| Mitochondria Isolation Buffer III (MIB III) | 65 ml    |
| Mitochondria Storage Buffer (MSB)           | 4 ml     |

## Procedures



# ProteinExt™ Mammalian Mitochondria Isolation Kit for Tissue

DE501-01

50 rxns

### Storage

MSB at -20°C for one year, others at 2-8°C for one year


### Description

ProteinExt™ Mammalian Mitochondria Isolation Kit for Tissue provides a fast and efficient isolation of mitochondria from tissues with simple procedure. This kit provides two options for the separation of mitochondria from cytosolic components: a reagent-based method or homogenization-based method. Reagent-based method uses a mild procedure to process single or multiple samples. The isolated mitochondria is suitable for a variety of downstream applications, including protein analysis, apoptosis, signal transduction and metabolic assays.

### Kit Contents

| Component                                   | DE501-01 |
|---|----------|
| Mitochondria Isolation Buffer I (MIB I)     | 50 ml    |
| Mitochondria Isolation Buffer II (MIB II)   | 500 µl   |
| Mitochondria Isolation Buffer III (MIB III) | 65 ml    |
| Mitochondria Storage Buffer (MSB)           | 4 ml     |
| Bovine Serum Albumin (BSA)                  | 500 mg   |

### Procedures







# ProteinIso<sup>®</sup> Ni-NTA Resin

|          |       |
|----------|-------|
| DP101-01 | 5 ml  |
| DP101-02 | 25 ml |

## Storage

at 2~8°C (20% ethanol) for two years

## Description

*ProteinIso*<sup>®</sup> Ni-NTA Resin allows rapid affinity purification of His-tagged proteins. The His-tagged proteins bind to Ni<sup>2+</sup> cations, which are immobilized on the Ni-NTA resin by 4 metal-chelating sites. After wash, the target proteins are recovered by gradient elution. The resin can be used for both native and denatured protein purification.

## Resin Specifications

|  |                         |
|--|-------------------------|
| Resin                                      | Cross-linked 6% agarose |
| Ligand                                     | NTA                     |
| Shape                                      | sphere                  |
| Pore size                                  | 45~165 μm               |
| Binding capacity                           | 10~20 mg/ml wet gel     |
| Recommended flow rate                      | <300 cm/h               |
| Highest resistance of atmospheric pressure | 0.3 Mpa                 |
| pH stability                               | 3-13                    |

## PROTOCOL


### Notes

- Samples should be centrifuged and filtrated with 0.45 μm filter before loading.
- **Equilibration Buffer for soluble protein**  
300 mM NaCl, 50 mM sodium phosphate buffer, 10 mM imidazole, 10 mM Tris-HCl pH 8.0
- **Equilibration Buffer for inclusion body**  
6 M GuHCl, 100 mM sodium phosphate buffer, 10 mM Tris-HCl pH 8.0; or 8 M urea, 100 mM sodium phosphate buffer, 10 mM Tris-HCl pH 8.0

### Procedures

1. Prepare Ni-NTA purification column
  - (1) Thoroughly resuspend the Ni-NTA resin to achieve a homogeneous suspension of the resin in the 20% ethanol storage buffer.
  - (2) Immediately transfer the resin into a purification column. Ensure that the bottom of the column is plugged with a stopper. Close the valve of the column. Allow the resin to settle.
  - (3) Equilibrate the column with 5~10 bed volume of equilibration buffer.
2. Prepare samples  
To avoid blocking column, samples should be centrifuged and filtrated with 0.45 μm filter before loading.
3. Load samples and wash  
Load samples and wash with 5~10 bed volume of equilibration buffer and collect the flow-through in one tube.
4. Elute  
Elute target proteins with imidazole or low pH buffer.
5. Regeneration of Ni-NTA resin
  - (1) Wash the column/resin with 2 bed volume of 6 M GuHCl, 0.2 M acetic acid
  - (2) 5 bed volume of deionized water
  - (3) 3 bed volume of 2% SDS
  - (4) 1 bed volume of 25% ethanol
  - (5) 1 bed volume of 50% ethanol
  - (6) 1 bed volume of 75% ethanol

- (7) 5 bed volume of 100% ethanol
- (8) 1 bed volume of 75% ethanol
- (9) 1 bed volume of 50% ethanol
- (10) 1 bed volume of 25% ethanol
- (11) 1 bed volume of deionized water
- (12) 5 bed volume of 100 mM EDTA, pH 8.0
- (13) 10 bed volume of deionized water
- (14) 5 bed volume of 100 mM NiSO<sub>4</sub>
- (15) Store column/resin in 20% ethanol



M: *ProteinRuler*<sup>®</sup> II

Lane 1: Sample

Lane 2: Flow-through

Lane 3: Wash

Lane 4: Elution (imidazole 40 mM)

Lane 5: Elution (imidazole 80 mM)

Lane 6: Elution (imidazole 120 mM)

Lane 7: Elution (imidazole 160 mM)

Lane 8: Elution (imidazole 200 mM)



# ProteinIso<sup>®</sup> Ni-IDA Resin

|          |       |
|----------|-------|
| DP111-01 | 5 ml  |
| DP111-02 | 25 ml |

## Storage

at 2~8°C (20% ethanol) for two years

## Description

*ProteinIso*<sup>®</sup> Ni-IDA Resin allows rapid affinity purification of His-tagged proteins. The His-tagged proteins bind to Ni<sup>2+</sup> cations, which are immobilized on the Ni-IDA resin by 3 metal-chelating sites. After wash, the target proteins are recovered by gradient elution. The resin can be used for both native and denatured protein purification.

## Resin Specifications

|  |                         |
|--|-------------------------|
| Resin                                      | Cross-linked 6% agarose |
| Ligand                                     | IDA                     |
| Shape                                      | sphere                  |
| Pore size                                  | 90 μm                   |
| Binding capacity                           | 20~40 mg/ml wet gel     |
| Recommended flow rate                      | <300 cm/h               |
| Highest resistance of atmospheric pressure | 0.3 Mpa                 |
| pH stability                               | 2~14                    |

## PROTOCOL


### Notes

- Samples should be centrifuged and filtrated with 0.45 μm filter before loading.
- **Equilibration Buffer for soluble protein**  
300 mM NaCl, 50 mM sodium phosphate buffer, 10 mM imidazole, 10 mM Tris-HCl pH 8.0
- **Equilibration Buffer for inclusion body**  
6 M GuHCl, 100 mM sodium phosphate buffer, 10 mM Tris-HCl pH 8.0; or 8 M urea, 100 mM sodium phosphate buffer, 10 mM Tris-HCl pH 8.0

### Procedures

1. Prepare Ni-IDA purification column
  - (1) Thoroughly resuspend the Ni-IDA resin to achieve a homogeneous suspension of the resin in the 20% ethanol storage buffer.
  - (2) Immediately transfer the resin into a purification column. Ensure that the bottom of the column is plugged with a stopper. Close the valve of the column. Allow the resin to settle.
  - (3) Equilibrate the column with 5~10 bed volume of equilibration buffer.
2. Prepare samples  
To avoid blocking column, samples should be centrifuged and filtrated with 0.45 μm filter before loading.
3. Load samples and wash  
Load samples and wash with 5~10 bed volume of equilibration buffer and collect the flow-through in one tube.
4. Elute  
Elute target proteins with imidazole or low pH buffer.
5. Regeneration of Ni-IDA resin
  - (1) Wash the column/resin with 2 bed volume of 6 M GuHCl, 0.2 M acetic acid
  - (2) 5 bed volume of deionized water
  - (3) 3 bed volume of 2% SDS
  - (4) 1 bed volume of 25% ethanol
  - (5) 1 bed volume of 50% ethanol
  - (6) 1 bed volume of 75% ethanol
  - (7) 5 bed volume of 100% ethanol

- (8) 1 bed volume of 75% ethanol
- (9) 1 bed volume of 50% ethanol
- (10) 1 bed volume of 25% ethanol
- (11) 1 bed volume of deionized water
- (12) 5 bed volume of 100 mM EDTA, pH 8.0
- (13) 10 bed volume of deionized water
- (14) 5 bed volume of 100 mM NiSO<sub>4</sub>
- (15) Store column/resin in 20% ethanol





# ProteinIso<sup>®</sup> GST Resin

DP201-01

10 ml

## Storage

at 2-8°C (20% ethanol) for two years

## Description

ProteinIso<sup>®</sup> GST Resin allows rapid affinity purification of GST-tagged proteins. GST fusion proteins expressed in *E.coli*, insect cells and mammalian cells can be purified with ProteinIso<sup>®</sup> GST Resin. The GST Resin is only suitable for soluble protein purification.

## Resin Specifications

|  |  |
|--|--|
| Resin                                      | Cross-linked 4% agarose  |
| Ligand                                     | glutathione  |
| Shape                                      | sphere   |
| Pore size                                  | 90 μm  |
| Support density                            | 8 mg GSH/ml wet gel  |
| Binding capacity                           | 10~12 mg/ml wet gel (MW 42 kDa)<br>3.5 mg/ml wet gel (rat liver) |
| Maximum flow rate (25°C)                   | 450 cm/h   |
| Recommended flow rate                      | <150 cm/h  |
| Highest resistance of atmospheric pressure | 0.3 Mpa  |
| pH stability                               | 3~10   |

## PROTOCOL


### Notes

- Samples should be centrifuged and filtrated with 0.45 μm filter before loading.
- To avoid cross-contamination, do not use the same medium to purify different proteins.
- **Equilibration Buffer**  
140 mM NaCl, 2.7 mM KCl, 10 mM Na<sub>2</sub>HPO<sub>4</sub>, 1.8 mM KH<sub>2</sub>PO<sub>4</sub>, pH 7.3
- **Elution Buffer**  
50 mM Tris-HCl pH 8.0, 10 mM reduced glutathione.

### Procedures

1. Prepare GST purification column
  - (1) Thoroughly resuspend the GST resin to achieve a homogeneous suspension of the resin in the 20% ethanol storage buffer.
  - (2) Immediately transfer the resin into a purification column. Ensure that the bottom of the column is plugged with a stopper. Close the valve of the column. Allow the resin to settle.
  - (3) Equilibrate the column with 5~10 bed volume of equilibration buffer.
2. Prepare samples  
To avoid blocking column, samples should be centrifuged and filtrated with 0.45 μm filter before loading.
3. Load samples and wash  
Load samples and wash with 5~10 bed volume of equilibration buffer and collect the flow-through in one tube.
4. Elute  
Elute target protein with elution buffer.
5. Regeneration of GST resin
  - (1) Wash the column/resin with 2 bed volume of 6 M GuHCl, 0.2 M acetic acid and then 5 bed volume of deionized water or PBS buffer.  
Or
  - (2) 3-4 bed volume of 70% ethanol or 30% isopropanol and then 3-5 bed volume of deionized water.  
Or

- (3) 2 bed volume of 10-100 mM NaOH and then 10 bed volume of deionized water.
- (4) Store column/resin in 20% ethanol.



— Target protein  
(82 kDa)

Lane 1: Sample  
 Lane 2: Flow-through  
 Lane 3: Wash  
 M: *ProteinRuler*<sup>®</sup> II  
 Lane 4-7: Elution (10 mM GSH)



# ProteinIso<sup>®</sup> Protein A Resin

DP301-01

5 ml

## Storage

at 2-8°C (20% ethanol) for two years

## Description

*ProteinIso*<sup>®</sup> Protein A Resin is an affinity chromatography resin with high binding capacity for IgG. The recombinant protein A ligand is coupled to highly cross-linked agarose. *ProteinIso*<sup>®</sup> Protein A Resin is suitable for purification of monoclonal antibody, polyclonal antibody and immunology complex, such as IP, Co-IP.

## Resin Specifications

|  |                            |
|--|----------------------------|
| Resin                                      | Cross-linked 4% agarose    |
| Ligand                                     | r-Protein A                |
| Shape                                      | sphere                     |
| Pore size                                  | 90 μm (45~165)             |
| Support density                            | 6 mg Protein A/ml wet gel  |
| Binding capacity                           | 40~50 mg h-IgG /ml wet gel |
| Maximum flow rate (25°C)                   | 300 cm/h                   |
| Recommended flow rate                      | <150 cm/h                  |
| Highest resistance of atmospheric pressure | 0.3 Mpa                    |
| pH stability                               | 3~10                       |

## PROTOCOL


### Notes

- Samples should be centrifuged and filtrated with 0.45 μm filter before loading.
- **Equilibration Buffer**  
20 mM PB, 150 mM KCl pH 7.0
- **Elution Buffer**  
20 mM citric acid pH 3.0-4.0;  
or 100 mM glycine pH 3.0;  
or 20 mM sodium acetate pH 3.0-4.0.
- **Neutralization Buffer**  
1 M Tris-HCl pH 9.0.

### Procedures

1. Prepare protein A purification column
  - (1) Thoroughly resuspend the protein A resin to achieve a homogeneous suspension of the resin in the 20% ethanol storage buffer.
  - (2) Immediately transfer the resin into a purification column. Ensure that the bottom of the column is plugged with a stopper. Close the valve of the column. Allow the resin to settle.
  - (3) Equilibrate the column with 5~10 bed volume of equilibration buffer.
2. Prepare samples  
To avoid blocking column, samples should be centrifuged and filtrated with 0.45 μm filter before loading.
3. Load samples and wash  
Load samples and wash with 5~10 bed volume of equilibration buffer and collect the flow-through in one tube.
4. Elute  
Elute antibodies with elution buffer.  
Collect the elution containing the target immunoglobulin and immediately neutralized to pH>7.0 with neutralization buffer.  
The elution conditions are closely related with binding strength and stability of antibody. When necessary, optimize the elution buffer.
5. Regeneration of Protein A Resin
  - (1) Wash the column/resin with 3~5 bed volume of 0.1 M citric acid or 0.1 M citric acid /20% ethanol and then 5 bed volume of PBS buffer (pH=7.0).  
Or

- (2) 3~5 bed volume of 0.05 M NaOH/1 M NaCl or 6 M GuHCl, and then 5 bed volume of deionized water.
- (3) Store column/resin in 20% ethanol.







# ProteinIso<sup>®</sup> Protein G Resin

DP401-01

5 ml

## Storage

at 2-8°C (20% ethanol) for two years

## Description

*ProteinIso*<sup>®</sup> Protein G Resin is an affinity chromatography resin with high binding capacity for IgG. The recombinant protein G ligand is coupled to highly cross-linked agarose. *ProteinIso*<sup>®</sup> Protein G Resin is suitable for purification of monoclonal antibody, polyclonal antibody and immunology complex, such as IP, Co-IP.

## Resin Specifications

|  |                            |
|--|----------------------------|
| Resin                                      | Cross-linked 4% agarose    |
| Ligand                                     | r-Protein G                |
| Shape                                      | sphere                     |
| Pore size                                  | 90 µm (45~165)             |
| Support density                            | 3 mg Protein G/ml wet gel  |
| Binding capacity                           | 20~25 mg h-IgG/ ml wet gel |
| Maximum flow rate (25°C)                   | 300 cm/h                   |
| Recommended flow rate                      | <150 cm/h                  |
| Highest resistance of atmospheric pressure | 0.3 Mpa                    |
| pH stability                               | 3~10                       |

## PROTOCOL

### Notes

- Samples should be centrifuged and filtrated with 0.45 µm filter before loading.
- **Equilibration Buffer**  
20 mM PB, 150 mM KCl pH 7.0
- **Elution Buffer**  
20 mM citric acid pH 3.0-4.0;  
or 100 mM glycine pH 3.0;  
or 20 mM sodium acetate pH 3.0-4.0.
- **Neutralization Buffer**  
1 M Tris-HCl pH 9.0.

### Procedures

1. Prepare protein G purification column
  - (1) Thoroughly resuspend the protein G resin to achieve a homogeneous suspension of the resin in the 20% ethanol storage buffer.
  - (2) Immediately transfer the resin into a purification column. Ensure that the bottom of the column is plugged with a stopper. Close the valve of the column. Allow the resin to settle.
  - (3) Equilibrate the column with 5~10 bed volume of equilibration buffer.
2. Prepare samples  
To avoid blocking column, samples should be centrifuged and filtrated with 0.45 µm filter before loading.
3. Load samples and wash  
Load samples and wash with 5~10 bed volume of equilibration buffer and collect the flow-through in one tube.
4. Elute  
Elute antibodies with elution buffer.  
Collect the elution containing the target immunoglobulin and immediately neutralize to pH>7.0 with neutralization buffer.  
The elution conditions are closely related with binding strength and stability of antibody. When necessary, optimize the elution buffer.


5. Regeneration of Protein G Resin

(1) Wash the column/resin with 3~5 bed volume of 0.1 M citric acid or 0.1 M citric acid /20% ethanol and then 5 bed volume of PBS buffer (pH=7.0).

Or

(2) 3~5 bed volume of 0.05 M NaOH/1 M NaCl or 6 M GuHCl, and then 5 bed volume of deionized water.

(3) Store column/resin in 20% ethanol.












**Affinity of Protein A/G for IgG Types**

| Sources | IgG subtypes      | Protein A binding capacity | Protein G binding capacity |
|---------|-------------------|----------------------------|----------------------------|
| human   | IgG <sub>1</sub>  | ++++                       | ++++                       |
|         | IgG <sub>2</sub>  | ++++                       | ++++                       |
|         | IgG <sub>3</sub>  | -                          | ++++                       |
|         | IgG <sub>4</sub>  | ++++                       | ++++                       |
| mouse   | IgG <sub>1</sub>  | +                          | ++++                       |
|         | IgG <sub>2a</sub> | ++++                       | ++++                       |
|         | IgG <sub>2b</sub> | +++                        | +++                        |
|         | IgG <sub>3</sub>  | ++                         | +++                        |
| rabbit  | IgG               | ++++                       | +++                        |
| goat    | IgG               | -                          | ++                         |
| horse   | IgG               | ++                         | ++++                       |
| dog     | IgG               | ++                         | +                          |
| cattle  | IgG               | ++                         | ++++                       |
| pig     | IgG               | +++                        | +++                        |
| monkey  | IgG               | ++++                       | ++++                       |



## Protein Marker Selection Guide

| Type                      | Name   | Cat. No. | SDS-PAGE | Western Blot | Monitor migration in SDS-PAGE | MW Range   |
|---------------------------|--|----------|----------|--------------|-------------------------------|------------|
| Unstained Protein Marker  | <i>ProteinRuler</i> <sup>®</sup> I               | DR101    | √        | -            | √                             | 12-80 kDa  |
|                           | <i>ProteinRuler</i> <sup>®</sup> II              | DR201    | √        | -            | √                             | 12-120 kDa |
|                           | <i>ProteinRuler</i> <sup>®</sup> IV              | DR401    | √        | -            | √                             | 30-200 kDa |
| Prestained Protein Marker | <i>Blue Plus</i> <sup>®</sup> Protein Marker     | DM101    | √        | √            | √                             | 14-100 kDa |
|                           | <i>Blue Plus</i> <sup>®</sup> II Protein Marker  | DM111    | √        | √            | √                             | 14-120 kDa |
|                           | <i>Blue Plus</i> <sup>®</sup> III Protein Marker | DM121    | √        | √            | √                             | 14-160 kDa |
|                           | <i>Blue Plus</i> <sup>®</sup> IV Protein Marker  | DM131    | √        | √            | √                             | 10-180 kDa |
| Western Protein Marker    | <i>EasySee</i> <sup>®</sup> Western Marker       | DM201    | √        | √            | √                             | 25-90 kDa  |
|                           | <i>EasySee</i> <sup>®</sup> II Western Marker    | DM211    | √        | √            | √                             | 30-150 kDa |



## ProteinRuler® I (12-80 kDa)

|          |        |
|----------|--------|
| DR101-01 | 250 µl |
| DR101-02 | 500 µl |

### Concentration

2 µg/5 µl for 12 kDa band; 1 µg/5 µl for 40 kDa band; 0.5 µg/5 µl for each of other bands


### Storage

at -20°C for two years

### Description

*ProteinRuler® I* is composed of five unstained recombinant proteins (20 kDa, 30 kDa, 40 kDa, 60 kDa, 80 kDa) and one blue prestained recombinant protein (12 kDa). The prestained band allows monitoring electrophoresis and membrane transfer. The 40 kDa band has doubled intensity to serve as a reference band.

- MW range from 12 to 80 kDa.
- Ready-to-use format, direct load on gels without heating.



12% Tris-glycine SDS gel (5 µl/well)

## ProteinRuler® II (12-120 kDa)

|          |        |
|----------|--------|
| DR201-01 | 250 µl |
| DR201-02 | 500 µl |

### Concentration

2 µg/5 µl for 12 kDa band; 1 µg/5 µl for 50 kDa band; 0.5 µg/5 µl for each of other bands


### Storage

at -20°C for two years

### Description

*ProteinRuler® II* is composed of eight unstained recombinant proteins (20 kDa, 30 kDa, 40 kDa, 50 kDa, 60 kDa, 80 kDa, 100 kDa, 120 kDa) and one blue prestained recombinant protein (12 kDa). The prestained band allows monitoring electrophoresis and membrane transfer. The 50 kDa band has doubled intensity to serve as a reference band.

- MW range from 12 to 120 kDa.
- Ready-to-use format, direct load on gels without heating.



12% Tris-glycine SDS gel (5 µl/well)



# ProteinRuler® IV

(30-200 kDa)

|          |        |
|----------|--------|
| DR401-01 | 250 µl |
| DR401-02 | 500 µl |

## Concentration

2 µg/5 µl for 30 kDa band; 1 µg/5 µl for 100 kDa band; 0.5 µg/5 µl for each of other bands

## Storage


at -20°C for two years

## Description

*ProteinRuler® IV* is composed of six unstained recombinant proteins (40 kDa, 60 kDa, 80 kDa, 100 kDa, 150 kDa, 200 kDa) and one blue prestained recombinant protein (30 kDa). The prestained band allows monitoring electrophoresis and membrane transfer. The 100 kDa band has doubled intensity to serve as a reference band.

- MW range from 30 to 200 kDa.
- Ready-to-use format, direct load on gels without heating.

*ProteinRuler® IV*  
(30-200 kDa)



8% Tris-glycine SDS gel (5 µl/well)

# Blue Plus<sup>®</sup> Protein Marker

(14-100 kDa)

|          |        |
|----------|--------|
| DM101-01 | 250 µl |
| DM101-02 | 500 µl |

### Concentration

about 2 µg/5 µl for each band

### Storage


at -20°C for two years

### Description

*Blue Plus*<sup>®</sup> Protein Marker is composed of seven prestained proteins ranging from 14 to 100 kDa. The protein of 50 kDa band is covalently coupled to orange dye. The protein of 14 kDa band is covalently coupled to yellow dye. The other five bands are covalently coupled to blue dye. This prestained protein marker is designed for monitoring the electrophoresis and membrane transfer.

- Five blue bands, one orange band and one yellow band.
- MW range from 14 to 100 kDa.
- Ready-to-use format, direct load on gels without heating.

*Blue Plus*<sup>®</sup> Protein Marker  
(14-100 kDa)



12% Tris-glycine SDS gel (5 µl/well)

# Blue Plus<sup>®</sup> II Protein Marker

(14-120 kDa)

|          |        |
|----------|--------|
| DM111-01 | 250 µl |
| DM111-02 | 500 µl |

### Concentration

about 2 µg/5 µl for each band

### Storage


at -20°C for two years

### Description

*Blue Plus*<sup>®</sup> II Protein Marker is composed of eight prestained proteins ranging from 14 to 120 kDa. The proteins of 50 kDa and 120 kDa bands are covalently coupled to orange dye. The protein of 14 kDa band is covalently coupled to yellow dye. The other five bands are covalently coupled to blue dye. This prestained protein marker is designed for monitoring the electrophoresis and membrane transfer.

- Five blue bands, two orange bands and one yellow band.
- MW range from 14 to 120 kDa.
- Ready-to-use format, direct load on gels without heating.

*Blue Plus*<sup>®</sup> II Protein Marker  
(14-120 kDa)



12% Tris-glycine SDS gel (5 µl/well)



# Blue Plus® III Protein Marker

(14-160 kDa)

|          |        |
|----------|--------|
| DM121-01 | 250 µl |
| DM121-02 | 500 µl |

**Concentration**

about 2 µg/5 µl for each band

**Storage**


at -20°C for two years

**Description**

Blue Plus® III Protein Marker is composed of nine prestained proteins ranging from 14 to 160 kDa. The proteins of 50 kDa, 120 kDa and 160 kDa bands are covalently coupled to orange dye. The protein of 14 kDa band is covalently coupled to yellow dye. The other five bands are covalently coupled to blue dye. This prestained protein marker is designed for monitoring the electrophoresis and membrane transfer.

- Five blue bands, three orange bands and one yellow band.
- MW range from 14 to 160 kDa.
- Ready-to-use format, direct load on gels without heating.

Blue Plus® III Protein Marker (14-160 kDa)



12% Tris-glycine SDS gel (5 µl/well)

# Blue Plus® IV Protein Marker

(10-180 kDa)

|          |        |
|----------|--------|
| DM131-01 | 250 µl |
| DM131-02 | 500 µl |

**Concentration**

about 2 µg/5 µl for each band

**Storage**


at -20°C for two years

**Description**

Blue Plus® IV Protein Marker is composed of ten prestained proteins ranging from 10 to 180 kDa. The protein of 80 kDa band is covalently coupled to orange dye. The protein of 25 kDa band is covalently coupled to green dye. The other eight bands are covalently coupled to blue dye. This prestained protein marker is designed for monitoring the electrophoresis and membrane transfer.

- Eight blue bands, one orange band and one green band.
- MW range from 10 to 180 kDa.
- Ready-to-use format, direct load on gels without heating.

Blue Plus® IV Protein Marker (10-180 kDa)



12% Tris-glycine SDS gel (5 µl/well)

# EasySee<sup>®</sup> Western Marker (25-90 kDa)

|   |                      |                                |
|---|----------------------|--------------------------------|
| without EasySee <sup>®</sup> Western Blot Kit | DM201-01<br>DM201-02 | 250 µl<br>500 µl               |
| with EasySee <sup>®</sup> Western Blot Kit    | DM201-11<br>DM201-12 | 250 µl+100 ml<br>500 µl+200 ml |


### Storage

at -20°C for two years

### Description

EasySee<sup>®</sup> Western Marker is composed of eight proteins ranging from 25 to 90 kDa. The 25 kDa, 45 kDa and 85 kDa bands are prestained allowing easy identification and monitoring electrophoresis and membrane transfer. Other five bands contain several IgG binding sites, allowing marker visualization using the same reagents and protocol for your target proteins. These no-dye-attached proteins provide more accurate molecular weight estimation.

- Three prestained bands for monitoring electrophoresis and membrane transfer.
- No label, no dye attached to other five recombinant protein bands.
- Five recombinant protein bands contain IgG binding sites, which can be developed with the standard Western blot substrates.
- Ready-to-use format, direct load on gels without heating.







# EasySee<sup>®</sup> II Western Marker

(30-150 kDa)

|   |          |                    |
|---|----------|--------------------|
| without EasySee <sup>®</sup> Western Blot Kit | DM211-01 | 250 $\mu$ l        |
|   | DM211-02 | 500 $\mu$ l        |
| with EasySee <sup>®</sup> Western Blot Kit    | DM211-11 | 250 $\mu$ l+100 ml |
|   | DM211-12 | 500 $\mu$ l+200 ml |


## Storage

at -20°C for two years

## Description

EasySee<sup>®</sup> II Western Marker is composed of seven proteins ranging from 30 to 150 kDa. The 30 kDa and 120 kDa bands are prestained allowing easy identification and monitoring electrophoresis and membrane transfer. Other five bands contain several IgG binding sites, allowing marker visualization using the same reagents and protocol for your target proteins. These no-dye-attached proteins provide more accurate molecular weight estimation.

- Two prestained bands for monitoring electrophoresis and membrane transfer.
- No label, no dye attached to other five recombinant protein bands.
- Five recombinant protein bands contain IgG binding sites, which can be developed with the standard Western blot substrates.
- Ready-to-use format, direct load on gels without heating.



# EasySee<sup>®</sup> Western Blot Kit

|          |        |
|----------|--------|
| DW101-01 | 100 ml |
| DW101-02 | 200 ml |

### Storage

at 2-8°C in dark for two years

### Description

EasySee<sup>®</sup> Western Blot Kit is optimized for enhanced chemiluminescence detection of western blots. The kit can be used in detection of horseradish peroxidase (HRP) conjugated secondary antibodies and related antigen.

- High sensitivity.
- Extended exposure time.
- High stability.


### Kit Contents

| Component                         | DW101-01 | DW101-02 |
|-----------------------------------|----------|----------|
| WB Solution A (Luminol)           | 50 ml    | 100 ml   |
| WB Solution B (Oxidant)           | 50 ml    | 100 ml   |
| WB Solution C (Light intensifier) | 150 µl   | 300 µl   |

## PROTOCOL

### Procedures

1. After electrophoresis, transfer proteins onto PVDF or NC membrane. Probe membrane with primary antibody followed by HRP-conjugated secondary antibody. Wash membrane for three times.
2. Mix equal volume of Solution A with Solution B. Add 0.05%-0.1% (v/v) of Solution C to the mixture (for example, add 1-2 µl of Solution C to 2 ml of Solution A + Solution B).
3. Add the mixed solution to the membrane (0.125 ml of reagent per cm<sup>2</sup> membrane). Incubate at room temperature for 1 minute.
4. Drain off excess solution from membrane. Do not let the membrane dry. Wrap the membrane with plastic wrap.
5. Develop image with x-Ray film.



# 6×Protein Loading Buffer

|          |        |
|----------|--------|
| DL101-02 | 5×1 ml |
|----------|--------|

### Storage

at -20°C for one year

### Description

6×Protein Loading Buffer is especially formulated for protein sample preparation used in SDS-PAGE. Prior to loading, add appropriate volume of 6×Protein Loading Buffer to protein sample to make its working concentration at 1×.

High quality products



# Easy Protein Quantitative Kit (Bradford)

DQ101-01

100 ml

## Storage

BSA Standard Solution at  $-20^{\circ}\text{C}$  for two years;  
Coomassie Brilliant Blue Solution at  $2-8^{\circ}\text{C}$   
in dark for two years

## Description

Easy Protein Quantitative Kit is a ready-to-use modified Bradford Coomassie-binding, colorimetric method for protein quantification. Under acidic condition, Coomassie Brilliant Blue G-250 binds to proteins providing an immediate shift in absorption maximum from 465 nm to 596 nm and a color change from brown to blue.


## Kit Contents

| Component                          | DQ101-01 |
|------------------------------------|----------|
| Coomassie Brilliant Blue Solution  | 100 ml   |
| BSA Standard Solution (0.22 mg/ml) | 4×1 ml   |

## PROTOCOL

### Procedures

1. Prior to use, equilibrate Coomassie Brilliant Blue Solution to room temperature and gently invert to mix well.
2. Transfer 0, 10, 30, 50, 70, 90, 100  $\mu\text{l}$  of BSA Standard Solution (0.22 mg/ml) into seven of 1.5 ml microfuge tubes, and add  $\text{H}_2\text{O}$  to a final volume of 100  $\mu\text{l}$ .
3. Transfer protein sample into a new 1.5 ml microfuge tube, and add  $\text{H}_2\text{O}$  to a final volume of 100  $\mu\text{l}$ .
4. Pipette 1.0 ml Coomassie Brilliant Blue Solution into each tube, mix thoroughly and incubate at room temperature for 5-10 minutes.
5. Measure the absorbance at 595 nm by spectrophotometer and record the value. Use the absorbance of sample without BSA as a blank control.
6. Plot the standard curve and calculate protein concentration in sample. Dilute the sample and re-measure it if the protein concentration falls out of the range of the standard curve.
7. The above procedures can be performed with microtiter-plate with 1/10 of the original volume.



# Easy II Protein Quantitative Kit (BCA)

DQ111-01 100 ml

## Storage

BSA Standard Solution at -20°C for two years;  
others at room temperature for one year

## Description

The BCA protein assay is one of the most commonly used methods for protein quantification. Under alkaline condition, the reduction of  $\text{Cu}^{2+}$  to  $\text{Cu}^+$  is realized by peptide bonds in proteins (biuret reaction). The amount of reduced copper is directly proportional to the amount of total proteins.

## Kit Contents

| Component                       | DQ111-01 |
|---------------------------------|----------|
| BCA Solution A                  | 100 ml   |
| BCA Solution B                  | 3 ml     |
| BSA Standard Solution (2 mg/ml) | 2×1 ml   |

## Interfering Substances



Certain substances are known to interfere with the BCA assay including those substances with reducing potential, chelating agents, and strong acids/bases. The following table shows the highest concentration of these substances in the protein sample buffer without interfering the BCA assay.

| Interfering Substances | Tolerant Concentration | Interfering Substances  | Tolerant Concentration |
|------------------------|------------------------|-------------------------|------------------------|
| Salts/Buffer           |                        | Detergents              |                        |
| HEPES (pH 7.9)         | 100 mM                 | NP-40                   | 5%                     |
| PIPES (pH 6.8)         | 100 mM                 | Triton X-100            | 5%                     |
| Ammonium sulfate       | 1.5 M                  | CHAPS, CHAPSO           | 5%                     |
| Sodium chloride        | 1 M                    | SDS                     | 5%                     |
| Sodium bicarbonate     | 100 mM                 | Tween 20                | 5%                     |
| MOPS (pH 7.2)          | 100 mM                 | Tween 60                | 5%                     |
| Sodium citrate         | 200 mM                 | Tween 80                | 5%                     |
| Tricine (pH 8.0)       | 25 mM                  | Mixture/Polar compounds |                        |
| Sodium acetate         | 200 mM                 | PMSF                    | 1 mM                   |
| Guanidine-HCl          | 4 M                    | Acetone                 | 10%                    |
| Tris                   | 250 mM                 | Ethanol                 | 10%                    |
| Chelating Agents       |                        | Glycerol                | 10%                    |
| EDTA                   | 10 mM                  | Urea                    | 3 M                    |
| Reducing Agents        |                        | DMSO                    | 10%                    |
| DTT                    | 1 mM                   | Sucrose                 | 40%                    |
| 2-Mercaptoethanol      | 0.01%                  |                         |                        |



# PROTOCOL

## Procedures



# ProteinEle™ Precast Tris-Glycine Gel

|          |             |
|----------|-------------|
| DG101-01 | 8%, 10/Box  |
| DG101-02 | 10%, 10/Box |
| DG101-03 | 12%, 10/Box |

### Storage

at 2-8°C for one year

### Description

*ProteinEle™* Precast Tris-Glycine Gel is a polyacrylamide gel used for native and denatured protein electrophoresis. It provides shorter running time and higher transfer efficiency.

### Highlights

- High stability: gel shelf life up to one year at 2-8°C.
- High resolution: superior band resolution on a broad range of native and denatured proteins.
- High reproducibility: consistent performance from gel to gel.

### Kit Contents

| Component | Specification                |
|-----------|------------------------------|
| DG101-01  | 8%, 1.0 mm, 12 well, 10/box  |
| DG101-02  | 10%, 1.0 mm, 12 well, 10/box |
| DG101-03  | 12%, 1.0 mm, 12 well, 10/box |


## PROTOCOL

### Note

If using the mini electrophoresis tank (Bio-Rad) or other similar electrophoresis tank, remove the sealing tape from electrophoresis apparatus, invert it and install again, then place the precast gel in the gel running tank. Otherwise, precast gel cannot be completely sealed with the outer surface of the sealing tape, which may result in electrophoresis buffer leak and therefore affect the result.

### Instruction

1. Take out the precast gel, remove the sealing tape near the bottom of the gel cassette. Place the gel in the gel running tank. Fill the gel wells with the electrophoresis buffer to immerse the chamber, gently pull out the comb from the chamber. Load the sample on the gel and run electrophoresis.
2. After electrophoresis is complete. Remove the Gel cassette from the gel running tank. To open the Gel cassette, insert a screwdriver into the gap between the two plastic plates between the gel. Apply pressure to separate them.



12% ProteinEle™ Precast Tris-Glycine Gel (5 µl/well)

---

## Chapter 7 Cell Culture and Detection

---

### Cell Culture

|  |     |
|--|-----|
| <i>TransSerum</i> <sup>®</sup> HQ Fetal Bovine Serum .....   | 192 |
| <i>TransLipid</i> <sup>®</sup> HL Transfection Reagent ..... | 193 |
| <i>TransIn</i> <sup>™</sup> EL Transfection Reagent .....    | 195 |
| Penicillin-Streptomycin (100×) .....                         | 197 |
| L-Glutamine (100×) .....                                     | 197 |
| Trypsin .....  | 197 |
| G418 .....   | 197 |
| PBS (1×) .....   | 198 |

### Cell Detection

|  |     |
|--|-----|
| <i>TransDetect</i> <sup>®</sup> Double-Luciferase Reporter Assay Kit .....                             | 198 |
| <i>TransDetect</i> <sup>®</sup> Cell Counting Kit (CCK) .....  | 199 |
| <i>TransDetect</i> <sup>®</sup> Annexin V-FITC/PI Cell Apoptosis<br>Detection Kit .....                | 201 |
| <i>TransDetect</i> <sup>®</sup> Annexin V-EGFP/PI Cell Apoptosis<br>Detection Kit .....                | 202 |
| <i>TransDetect</i> <sup>®</sup> <i>In Situ</i> Fluorescein TUNEL Cell<br>Apoptosis Detection Kit ..... | 203 |

# TransSerum<sup>®</sup> HQ Fetal Bovine Serum

FS101-02

500 ml

### Storage

at -20°C for five years

### Description

Fetal Bovine Serum is the most widely used undefined supplement in mammalian cell culture. *TransSerum*<sup>®</sup> HQ Fetal Bovine Serum is collected from healthy fetal blood, followed by filtration through 0.1 µm filters for three times.

- Low toxicity
- Promoting better cell growth
- Suitable for a broad range of cells

### Physical Properties

|                             |                   |
|-----------------------------|-------------------|
| pH                          | 6.9-7.8           |
| Osmolality                  | 280-340 mOsmol/kg |
| Total proteins              | 3.0-3.5 g/dl      |
| Albumin                     | 1.8-2.45 g/dl     |
| Endotoxin                   | <10 EU/ml         |
| IgG                         | <100 µg/ml        |
| Hemoglobin                  | <20 mg/dl         |
| Test of germ, mycete, phage | Negative          |
| Test of mycoplasma          | Negative          |
| Test of virus               | Negative          |

### Note

This product is not heat inactivated. If needed, incubate thawed FBS at 56°C water bath for 30 minutes to inactivate complement.

### Successfully cultured cell types with *TransSerum*<sup>®</sup> HQ Fetal Bovine Serum

|          |            |          |
|----------|------------|----------|
| 5637     | HEK-293T   | OV45     |
| 7721     | HeLa       | P815     |
| A2780    | Hep G2     | PANC-1   |
| A549     | HL-60      | PC-12    |
| B16-F10  | HT-29      | PT67     |
| BEL-7402 | Huh7       | RAW264.7 |
| BHK21    | Jurkat     | Sf9      |
| BGC-823  | K-562      | SGC-7901 |
| CEF      | L929       | SK-OV-3  |
| CEK      | LS 174T    | Sp2/0    |
| CHO-K1   | MCF7       | STO      |
| COS-1    | MDA-MB-231 | SW480    |
| COS-7    | MEF        | T24      |
| DF-1     | MGC803     | THP-1    |
| DLD-1    | MKN-28     | U87      |
| EJ       | MKN-45     | U937     |
| GLC-82   | MRC-5      | Vero     |
| HCT 116  | NIH3T3     | WEHI-3B  |
| HEK-293  | NRK        |          |





# TransLipid<sup>®</sup> HL Transfection Reagent

|          |           |
|----------|-----------|
| FT111-01 | 0.75 ml   |
| FT111-02 | 2×0.75 ml |

## Storage

at 2-8°C for one year

## Description

*TransLipid<sup>®</sup>* HL Transfection Reagent is a proprietary formulated cationic lipid that offers superior transfection efficiency and low cytotoxicity across a broad range of mammalian cell lines.


- High efficiency
- Low cytotoxicity
- Transfect DNA, RNA, siRNA
- Adherent or suspension cells
- Can be used in the presence of serum and antibiotics

## PROTOCOL

We suggest using DNA( $\mu$ g)/*TransLipid<sup>®</sup>* HL ( $\mu$ l) with ratio at 1:2-1:3. To obtain better transfection efficiency, we suggest using high density cell (70%-90% confluency).

### Plasmid DNA Transfection

24-well format as an example



### siRNA Transfection

Cells should be 30-50% confluency at the time of transfection. For 24-well plate, use 20 pmol of siRNA and 1  $\mu$ l of *TransLipid*<sup>®</sup> HL. The experimental procedure is the same as DNA transfection described above.

### Optimization of plasmid DNA and siRNA transfection

In order to achieve optimal combination of high transfection efficiency and low cytotoxicity, the ratio of DNA or siRNA to *TransLipid*<sup>®</sup> HL as well as the initial cell density for transfection could be optimized. DNA transfection can be optimized within the range of 1:2-1:5, it is recommended to use a range of 10 to 50 pmol of siRNA and 0.5 to 1.5  $\mu$ l of *TransLipid*<sup>®</sup> HL.

### Amount of culture medium, nucleic acid and *TransLipid*<sup>®</sup> HL in transfection of different cell culture plates

| Culture Vessel | Surface Area per Well | Volume of Plating Medium | Dilution Volume        | DNA Transfection |                                   | siRNA Transfection |                                   |
|----------------|-----------------------|--------------------------|------------------------|------------------|-----------------------------------|--------------------|-----------------------------------|
|                |                       |                          |                        | DNA              | <i>TransLipid</i> <sup>®</sup> HL | siRNA              | <i>TransLipid</i> <sup>®</sup> HL |
| 96-well        | 0.3 cm <sup>2</sup>   | 100 $\mu$ l              | 2 $\times$ 10 $\mu$ l  | 0.2 $\mu$ g      | 0.4-1 $\mu$ l                     | 5 pmol             | 0.25 $\mu$ l                      |
| 48-well        | 1 cm <sup>2</sup>     | 250 $\mu$ l              | 2 $\times$ 25 $\mu$ l  | 0.4 $\mu$ g      | 0.8-2 $\mu$ l                     | 10 pmol            | 0.5 $\mu$ l                       |
| 24-well        | 2 cm <sup>2</sup>     | 500 $\mu$ l              | 2 $\times$ 50 $\mu$ l  | 0.8 $\mu$ g      | 1.6-4 $\mu$ l                     | 20 pmol            | 1 $\mu$ l                         |
| 12-well        | 4 cm <sup>2</sup>     | 1 ml                     | 2 $\times$ 100 $\mu$ l | 1.6 $\mu$ g      | 3.2-8 $\mu$ l                     | 40 pmol            | 2 $\mu$ l                         |
| 6-well         | 10 cm <sup>2</sup>    | 2 ml                     | 2 $\times$ 250 $\mu$ l | 4 $\mu$ g        | 8-20 $\mu$ l                      | 100 pmol           | 5 $\mu$ l                         |
| 35 mm          | 10 cm <sup>2</sup>    | 2 ml                     | 2 $\times$ 250 $\mu$ l | 4 $\mu$ g        | 8-20 $\mu$ l                      | 100 pmol           | 5 $\mu$ l                         |
| 60 mm          | 20 cm <sup>2</sup>    | 5 ml                     | 2 $\times$ 0.5 ml      | 8 $\mu$ g        | 16-40 $\mu$ l                     | 200 pmol           | 10 $\mu$ l                        |
| 10 cm          | 60 cm <sup>2</sup>    | 10 ml                    | 2 $\times$ 1.5 ml      | 24 $\mu$ g       | 48-120 $\mu$ l                    | 600 pmol           | 30 $\mu$ l                        |
| T 25           | 25 cm <sup>2</sup>    | 6 ml                     | 2 $\times$ 0.625 ml    | 10 $\mu$ g       | 20-50 $\mu$ l                     | 250 pmol           | 12.5 $\mu$ l                      |
| T 75           | 75 cm <sup>2</sup>    | 20 ml                    | 2 $\times$ 1.875 ml    | 30 $\mu$ g       | 60-150 $\mu$ l                    | 800 pmol           | 40 $\mu$ l                        |

### Successfully transfected cell types with *TransLipid*<sup>®</sup> HL Transfection Reagent

|          |          |        |           |
|----------|----------|--------|-----------|
| A549     | HCT 116  | Huh7   | PC3       |
| B16-F10  | HeLa     | K-562  | PANC-1    |
| BEL-7402 | Hep G2   | MCF7   | PT67      |
| BHK21    | HEK-293  | MDS-L  | RAW 264.7 |
| CHO-K1   | HEK-293T | MRC-5  | STO       |
| COS-1    | HMC-1    | NIH3T3 | Vero      |
| COS-7    | HT-29    | NRK    |           |



# TransIn™ EL Transfection Reagent

|          |           |
|----------|-----------|
| FT201-01 | 0.75 ml   |
| FT201-02 | 2×0.75 ml |

## Storage

at 2-8°C for one year

## Description


*TransIn*™ EL Transfection Reagent is a non-liposomal formulation designed to transfect DNA into a wide variety of eukaryotic cell lines with high efficiency and low toxicity. Primary cells and other difficult-to-transfect cells can also be effectively transfected by this reagent.

- Non-liposomal transfection reagent
- High efficiency
- Low cytotoxicity
- Adherent or suspension cells
- Can be used in the presence of serum and antibiotics

## PROTOCOL

### Plasmid DNA Transfection

24-well format as an example



### Optimization of plasmid DNA transfection

In order to achieve optimal combination of high transfection efficiency and low cytotoxicity, the ratio of DNA to *TransIn*<sup>TM</sup> EL as well as the initial cell density for transfection could be optimized. Ratio of DNA (µg) and *TransIn*<sup>TM</sup> EL (µl) can be optimized within the range of 1:1-1:3.

### Amount of culture medium, nucleic acid and *TransIn*<sup>TM</sup> EL in transfection of different cell culture plates

| Culture Vessel | Surface Area per Well | Volume of Plating Medium | Dilution Volume | DNA    | <i>TransIn</i> <sup>TM</sup> EL |
|----------------|-----------------------|--------------------------|-----------------|--------|---------------------------------|
| 96-well        | 0.3 cm <sup>2</sup>   | 100 µl                   | 10 µl           | 0.2 µg | 0.2-0.6 µl                      |
| 48-well        | 1 cm <sup>2</sup>     | 250 µl                   | 25 µl           | 0.4 µg | 0.4-1.2 µl                      |
| 24-well        | 2 cm <sup>2</sup>     | 500 µl                   | 50 µl           | 0.8 µg | 0.8-2.4 µl                      |
| 12-well        | 4 cm <sup>2</sup>     | 1 ml                     | 100 µl          | 1.6 µg | 1.6-4.8 µl                      |
| 6-well         | 10 cm <sup>2</sup>    | 2 ml                     | 200 µl          | 4 µg   | 4-12 µl                         |
| 35 mm          | 10 cm <sup>2</sup>    | 2 ml                     | 200 µl          | 4 µg   | 4-12 µl                         |
| 60 mm          | 20 cm <sup>2</sup>    | 5 ml                     | 0.5 ml          | 8 µg   | 8-24 µl                         |
| 10 cm          | 60 cm <sup>2</sup>    | 10 ml                    | 1 ml            | 24 µg  | 24-72 µl                        |
| T 25           | 25 cm <sup>2</sup>    | 6 ml                     | 0.5 ml          | 10 µg  | 10-30 µl                        |
| T 75           | 75 cm <sup>2</sup>    | 20 ml                    | 1 ml            | 30 µg  | 30-90 µl                        |

### Successfully transfected cell types with *TransIn*<sup>TM</sup> EL Transfection Reagent

|         |          |            |          |
|---------|----------|------------|----------|
| A549    | HEK-293  | MARC-145   |          |
| B16-F10 | HEK-293T | MCF-7      | PT67     |
| BHK-21  | HeLa     | MEF        | SGC-7901 |
| CEF     | Hep G2   | MIA PaCa-2 | SH-SY5Y  |
| CHO     | HL-60    | NIH/3T3    | STO      |
| COS-1   | K562     | NRK        | Vero     |
| HCT-116 | L929     | P815       |          |



## Penicillin-Streptomycin (100x)

FG101-01

100 ml

### Storage

-20°C for one year

### Description

Penicillin-Streptomycin (100x) contains 10 kU/ml of penicillin and 10 mg/ml of streptomycin. The solution has been filter-sterilized. It can be used for cell culture at a final concentration of 1x.

### Note

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks, or at -20°C for one year.

## L-Glutamine (100x)

FG201-01

100 ml

### Storage

-20°C for one year

### Description

L-Glutamine (100x) contains 200 mM of L-Glutamine. The solution has been filter-sterilized. It can be used for cell culture at a final concentration of 1x.

### Note

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks, or at -20°C for one year.

## Trypsin

FG301-01 (+EDTA)

100 ml

FG301-11 (-EDTA)

100 ml

### Storage

-20°C for 18 months

### Description

Trypsin contains porcine trypsin (0.25%), EDTA (+/-), and phenol red. It does not contain calcium and magnesium ion. Trypsin solution has been filter-sterilized and it can be used for cell dissociation.

### Note

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks, or at -20°C for one year.

## G418

FG401-01

5 ml

### Storage

2-8°C for two years

### Description

G418 is an aminoglycoside antibiotic, which blocks polypeptide synthesis by interfering with the function of 80S ribosome. Due to its toxicity on prokaryotic and eukaryotic cells (including bacteria, fungi, plants and mammalian cells), it is widely used as a selective reagent for stable cell line construction. The resistance mechanism is based on that resistance gene (Neomycin) specifically express aminoglycoside phosphotransferase, which confers resistance on cells. Thus cells are able to grow in selective culture medium containing G418.

## PBS (1x)

FG701-01

500 ml

### Storage

at room temperature for two years

### Description

PBS (phosphate buffered saline) contains 1.06 mM  $\text{KH}_2\text{PO}_4$ , 155.17 mM NaCl, 2.97 mM  $\text{Na}_2\text{HPO}_4$  (pH 7.4). PBS has been widely used for a variety of cell culture applications, such as washing, dissociation and dilution. PBS is formulated without calcium, magnesium and phenol red. This product has been filter-sterilized.

## TransDetect<sup>®</sup> Double-Luciferase Reporter Assay Kit

FR201-01

50 rxns

FR201-02

200 rxns

### Storage

at -20°C in dark for one year

### Description

*TransDetect*<sup>®</sup> Double-Luciferase Reporter Assay Kit provides an efficient method to perform dual-reporter assays. In the assay, the activities of firefly (*Photinus pyralis*) and *Renilla* (sea pansy) luciferases are measured sequentially from a single sample. The firefly luciferase reporter is measured firstly by adding Luciferase Reaction Reagent to generate a luminescent signal. After quantifying the firefly luminescence, this reaction is quenched, and the *Renilla* luciferase reaction is initiated by adding Luciferase Reaction Reagent II to the same tube.

### Kit Contents

| Component                                   | FR201-01          | FR201-02          |
|---|-------------------|-------------------|
| Luciferase Reaction Buffer                  | 5 ml              | 20 ml             |
| Luciferase Reaction Substrate (Lyophilized) | 1 vial            | 4 vials           |
| Luciferase Reaction Buffer II               | 5 ml              | 20 ml             |
| Luciferase Reaction Substrate II (50x)      | 100 $\mu\text{l}$ | 400 $\mu\text{l}$ |
| Cell Lysis Buffer (5x)                      | 5 ml              | 20 ml             |

### Procedures

1. Prepare Luciferase Reaction Reagent  
Dissolve the lyophilized Luciferase Reaction Substrate (whole vial) by adding 5 ml of Luciferase Reaction Buffer to the vial. The prepared Luciferase Reaction Reagent can be aliquoted and stored at -20°C for one month or -70°C for one year.
2. Prepare Luciferase Reaction Reagent II  
Mix Luciferase Reaction Substrate II with Luciferase Reaction Buffer II with the ratio of 1:50. The prepared Luciferase Reaction Reagent II can be aliquoted and stored at -70°C for one month. For best results, the Luciferase Reaction Reagent II should be prepared fresh before each use.
3. Prepare 1xCell Lysis Buffer  
Mix 5xCell Lysis Buffer with ddH<sub>2</sub>O with the ratio of 1:4. 1xCell Lysis Buffer can be stored at 2-8°C for one month.
4. Prepare Cell Lysate  
Wash cells twice with PBS and add appropriate volume of 1xCell Lysis Buffer. Incubate at room temperature for 10 minutes and then scrape

High quality products



the cells. Centrifuge at 12,000xg for 10 minutes at 2-8°C. Transfer the supernatant to a new tube.

#### 5. Assay

Mix 20 µl of lysate with 100 µl of Luciferase Reaction Reagent. Place the tube in luminometer to measure the luminescent signal of firefly luciferase. Then, add 100 µl of Luciferase Reaction Reagent II to the tube, place the tube in luminometer to measure the luminescent signal of *Renilla* luciferase. Record the data.

## TransDetect<sup>®</sup> Cell Counting Kit (CCK)

|          |       |
|----------|-------|
| FC101-01 | 1 ml  |
| FC101-02 | 5 ml  |
| FC101-03 | 10 ml |
| FC101-04 | 30 ml |

### Storage

at 2-8°C in dark for one year or at -20°C in dark for two years

### Description

*TransDetect<sup>®</sup>* Cell Counting Kit (CCK) is designed for cell proliferation assays as well as cytotoxicity assays by utilizing a water-soluble tetrazolium salt. The salt can be reduced to an orange water-soluble formazan by mitochondrial dehydrogenase in the presence of an electron coupling reagent 1-Methoxy PMS. The amount of the formazan dye generated by dehydrogenases in cells is directly proportional to the number of living cells. Faster cell proliferation, lower cytotoxicity, and more cell number produce deeper color. The depth of the dye is directly proportional to the number of living cells. The toxicity of CCK solution is so low that the same cells can be used for other assays after the CCK assay is completed. Compared with MTT, XTT, MST and WST-1, this method provides higher sensitivity and broader linear range. It is suitable for drug screening, cell proliferation test, cytotoxicity assay and drug sensitivity test.

### Kit Content


| Component    | FC101-01 | FC101-02 | FC101-03 | FC101-04 |
|--------------|----------|----------|----------|----------|
| CCK Solution | 1 ml     | 5 ml     | 10 ml    | 30 ml    |

### Procedures

1. Inoculate cell suspension in a 96-well plate (100 µl/well). Pre-incubate the plate in a cell incubator at 37°C for 12-24 hours according to experimental need. In general, use  $2 \times 10^3$  cells per well for cell proliferation assays, use  $5 \times 10^3$  cells per well for cytotoxicity assays.
2. Add appropriate volume (0-10 µl) of substance to be tested to the plate, incubate for an appropriate length of time in the incubator.
3. Add CCK solution (equal to 1/10 of the media volume) to each well of the 96-well plate (e.g. add 10 µl of CCK solution for 100 µl of the media).
4. Incubate the plate for 1-4 hour in the incubator.
5. Measure the absorbance at 450 nm using a microplate reader.

**Notes**

1. The presence of phenol red has no effect on the result, but can increase the background absorption. Thus the blank absorbance should be subtracted.
2. The incubation time after adding CCK solution varies by the cell type and density. Perform initial experiments to determine the appropriate incubation time. Generally, lymphocyte has lower sensitivity, which requires longer incubation time or higher cell density.
3. Assays by this kit depend on the catalyzation of dehydrogenase in cells. If the substance to be tested has strong oxidativity or reductivity, fresh media should be changed prior to adding CCK solution.
4. Be careful not to introduce bubbles to the wells since they will interfere with absorbance value.
5. If there is no 450 nm optical filter, the filter with absorbance between 430 nm and 490 nm can be used.
6. For highly turbid cell suspension, 600 nm (or above 600 nm) can be used as a reference to perform dual wavelength measurement.
7. To stop the reaction, add 10  $\mu$ l of 0.1 M HCl or 1% w/v SDS solution into each well of 96-well plate, and store in dark. The absorbance will not change within 24 hours.
8. The toxicity of CCK solution is so low that the same cells can be used for other assays after the CCK assay is completed.
9. This kit should be stored in dark. For long-term storage, aliquot CCK solution and store at  $-20^{\circ}\text{C}$ .







# TransDetect<sup>®</sup> Annexin V-FITC/PI Cell Apoptosis Detection Kit

|          |         |
|----------|---------|
| FA101-01 | 25 rxns |
| FA101-02 | 50 rxns |

## Storage


at 2-8°C in dark for one year

## Description

The Annexin V-FITC/PI Cell Apoptosis Kit provides a rapid and sensitive method for early apoptosis detection. In normal cells, the membrane phospholipid phosphatidylserine (PS) is located on the cytoplasmic surface of the membrane. In apoptotic cells, PS is translocated from the inner to the outer leaflet of the plasma membrane, thereby exposing PS to the external cellular environment. The FITC-conjugated Annexin V, a Ca<sup>2+</sup> dependent phospholipid-binding protein, can bind specifically to the exposed PS. Propidium iodide (PI) is a nucleic acid binding dye, which binds tightly to the nucleic acids in the cells and stains the cells with red fluorescence. PI is impermeant to live cells and early apoptotic cells, so the combination of Annexin V-FITC and PI staining allows the differentiation among different stage of apoptotic cells and necrosis cells.

## Kit Contents

| Component                   | FA101-01 | FA101-02  |
|-----------------------------|----------|-----------|
| Annexin V-FITC              | 125 µl   | 250 µl    |
| Propidium iodide (PI)       | 125 µl   | 250 µl    |
| 1× Annexin V Binding Buffer | 12.5 ml  | 2×12.5 ml |



Apoptosis detection of Camptothecin treated Jurkat cells by Flow Cytometry

# TransDetect<sup>®</sup> Annexin V-EGFP/PI Cell Apoptosis Detection Kit

|          |         |
|----------|---------|
| FA111-01 | 25 rxns |
| FA111-02 | 50 rxns |

### Storage


at 2-8°C in dark for one year

### Description

The Annexin V-EGFP/PI Cell Apoptosis Kit provides a rapid and sensitive method for early apoptosis detection. In normal cells, the membrane phospholipid phosphatidylserine (PS) is located on the cytoplasmic surface of the membrane. In apoptotic cells, PS is translocated from the inner to the outer leaflet of the plasma membrane, thereby exposing PS to the external cellular environment. The EGFP-conjugated Annexin V, a Ca<sup>2+</sup> dependent phospholipid-binding protein, can bind specifically to the exposed PS. Propidium iodide (PI) is a nucleic acid binding dye, which binds tightly to the nucleic acids in the cells and stains the cells with red fluorescence. PI is impermeant to live cells and early apoptotic cells, so the combination of Annexin V-EGFP and PI staining allows the differentiation among different stage of apoptotic cells and necrosis cells. Compared with FITC, EGFP is brighter and more photo-stable. Because Annexin V-EGFP is a fusion protein with a 1:1 binding ratio of EGFP to PS, this kit can also be used for quantitative detection.

### Kit Contents

| Component                   | FA111-01 | FA111-02  |
|-----------------------------|----------|-----------|
| Annexin V-EGFP              | 125 µl   | 250 µl    |
| Propidium Iodide (PI)       | 125 µl   | 250 µl    |
| 1× Annexin V Binding Buffer | 12.5 ml  | 2×12.5 ml |



Apoptosis detection of Camptothecin treated Jurkat cells by Flow Cytometry



# *TransDetect*<sup>®</sup> *In Situ* Fluorescein TUNEL Cell Apoptosis Detection Kit

|          |         |
|----------|---------|
| FA201-01 | 25 rxns |
| FA201-02 | 50 rxns |

## Storage

TdT at -20°C for one year, 1×Labeling Solution at -20°C in dark for one year


## Description

*TransDetect*<sup>®</sup> *In Situ* Fluorescein TUNEL Cell Apoptosis Detection Kit provides a precise, simple and low-toxicity way to detect and quantify apoptotic cell death at single cell level in cells and tissues. TdT-mediated dUTP Nick-End Labeling (TUNEL) reaction preferentially labels DNA strand breaks generated during apoptosis with fluorescein-labeled dUTP. The fluorescein labeled DNA can be detected and quantified by fluorescence microscopy or flow cytometry. This kit can be used to detect apoptosis in paraffin-embedded tissue sections, cryopreserved tissue sections, cells cultured on chamber slides, cell smear and cell suspensions.

## Kit Contents


| Component           | FA201-01 | FA201-02  |
|---------------------|----------|-----------|
| TdT                 | 50 µl    | 100 µl    |
| 1×Labeling Solution | 1.25 ml  | 2×1.25 ml |

DNase I-treated



TransGen Company R

Camptothecin-treated



TransGen Company R

Apoptosis detection of DNase I or Camptothecin treated HeLa cells by fluorescence microscopy

---

## Chapter 8 Antibodies

---

### Primary Antibodies

|   |          |
|---|----------|
| <i>ProteinFind</i> <sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody        | .....205 |
| <i>ProteinFind</i> <sup>®</sup> Anti-DYKDDDDK Tag Mouse Monoclonal Antibody | .....205 |
| <i>ProteinFind</i> <sup>®</sup> Anti-HA Mouse Monoclonal Antibody           | .....206 |
| <i>ProteinFind</i> <sup>®</sup> Anti-V5 Mouse Monoclonal Antibody           | .....207 |
| <i>ProteinFind</i> <sup>®</sup> Anti-His Mouse Monoclonal Antibody          | .....207 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GST Mouse Monoclonal Antibody          | .....208 |
| <i>ProteinFind</i> <sup>®</sup> Anti-MBP Mouse Monoclonal Antibody          | .....208 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GFP Mouse Monoclonal Antibody          | .....209 |

### Control Antibodies

|  |          |
|--|----------|
| <i>ProteinFind</i> <sup>®</sup> Anti- $\beta$ -Tubulin Mouse Monoclonal Antibody | .....209 |
| <i>ProteinFind</i> <sup>®</sup> Anti- $\beta$ -Actin Mouse Monoclonal Antibody   | .....210 |
| <i>ProteinFind</i> <sup>®</sup> Anti-GAPDH Mouse Monoclonal Antibody             | .....210 |

### Secondary Antibodies

|  |          |
|--|----------|
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), HRP Conjugate   | .....211 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), FITC Conjugate  | .....211 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), PE Conjugate    | .....212 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), AF488 Conjugate | .....212 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), HRP Conjugate    | .....213 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), FITC Conjugate   | .....214 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), PE Conjugate     | .....214 |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), AF488 Conjugate  | .....215 |

### Related Products

|                           |          |
|---------------------------|----------|
| TMB ELISA Substrate       | .....216 |
| Super TMB ELISA Substrate | .....216 |



## ProteinFind<sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody

|          |             |
|----------|-------------|
| HT101-01 | 50 $\mu$ l  |
| HT101-02 | 100 $\mu$ l |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody is a purified monoclonal antibody that detects recombinant proteins containing the c-Myc (EQKLISEEDL) epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-DYKDDDDK Tag Mouse Monoclonal Antibody

|          |             |
|----------|-------------|
| HT201-01 | 50 $\mu$ l  |
| HT201-02 | 100 $\mu$ l |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-DYKDDDDK Tag Mouse Monoclonal Antibody is the same as FLAG antibody from Sigma. It is a purified monoclonal antibody that detects recombinant proteins containing the DYKDDDDK epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind® Anti-HA Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HT301-01 | 50 µl  |
| HT301-02 | 100 µl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind® Anti-HA Mouse Monoclonal Antibody is a purified monoclonal antibody that detects recombinant proteins containing the HA (YPYDVPDYA) epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.





## ProteinFind<sup>®</sup> Anti-V5 Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HT401-01 | 50 µl  |
| HT401-02 | 100 µl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-V5 Mouse Monoclonal Antibody is a purified monoclonal antibody that detects recombinant proteins containing the V5 (CGKPIPPELLGLDST) epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-His Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HT501-01 | 50 µl  |
| HT501-02 | 100 µl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-His Mouse Monoclonal Antibody is a purified monoclonal antibody that detects recombinant proteins containing the 6xHis (HHHHHH) epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-GST Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HT601-01 | 50 µl  |
| HT601-02 | 100 µl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-GST Mouse Monoclonal Antibody is a purified monoclonal antibody against yeast Y258 GST recombinant proteins that detects GST fusion proteins.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-MBP Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HT701-01 | 50 µl  |
| HT701-02 | 100 µl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-MBP Mouse Monoclonal Antibody is a purified monoclonal antibody that detects MBP fusion proteins.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IP: 1:100-500 dilution.







## ProteinFind<sup>®</sup> Anti-GFP Mouse Monoclonal Antibody

|          |             |
|----------|-------------|
| HT801-01 | 50 $\mu$ l  |
| HT801-02 | 100 $\mu$ l |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

ProteinFind<sup>®</sup> Anti-GFP Mouse Monoclonal Antibody is a purified monoclonal antibody that detects GFP fusion proteins.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti- $\beta$ -Tubulin Mouse Monoclonal Antibody

|          |             |
|----------|-------------|
| HC101-01 | 50 $\mu$ l  |
| HC101-02 | 100 $\mu$ l |

### Concentration

1 mg/ml

### Storage

at 2-8°C for one month; at -20°C for one year


### Description

Tubulin is an important component of the cytoskeleton. It is widely present in various mammalian cells and mainly consists of  $\alpha$ -tubulin and  $\beta$ -tubulin. The expression level of  $\beta$ -tubulin is relatively stable. It is widely used as expression control.

ProteinFind<sup>®</sup> Anti- $\beta$ -Tubulin Mouse Monoclonal Antibody is a purified monoclonal antibody that detects  $\beta$ -Tubulin in human, rat, mouse, goat and other species.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-β-Actin Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HC201-01 | 50 μl  |
| HC201-02 | 100 μl |

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

Actin is an important component of the cytoskeleton. It is widely present in various mammalian cells and mainly consists of the β-Actin. The expression level of β-Actin is relatively stable. It is widely used as expression control. ProteinFind<sup>®</sup> Anti-β-Actin Mouse Monoclonal Antibody is a purified monoclonal antibody that detects β-Actin in human, mouse, rabbit and other species.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Anti-GAPDH Mouse Monoclonal Antibody

|          |        |
|----------|--------|
| HC301-01 | 50 μl  |
| HC301-02 | 100 μl |

### Concentration

1 mg/ml

### Storage

at 2-8°C for one month; at -20°C for one year


### Description

GAPDH (glyceraldehyde-3-phosphate dehydrogenase) is a key enzyme for the glycolysis process. It is widely present in various cells and has been used as expression control.

ProteinFind<sup>®</sup> Anti-GAPDH Mouse Monoclonal Antibody is a purified monoclonal antibody that detects GAPDH in human, mouse, rabbit and other species.

### Suggested Dilution

- Western Blot: 1:500-10,000 dilution.
- ELISA: 1:500-5,000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.



High quality products



## *ProteinFind*<sup>®</sup> Goat Anti-Rabbit IgG(H+L), HRP Conjugate

HS101-01 100  $\mu$ l**Concentration**

1 mg/ml

**Storage**


at 2-8°C for one month; at -20°C for one year

**Description**

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Rabbit IgG(H+L) Antibody is a horseradish peroxidase (HRP) conjugated secondary antibody for ELISA and Western blot detection.

**Suggested Dilution**

- Western: 1:1000-10,000 dilution.
- ELISA: 1:1,000-5,000 dilution.



## *ProteinFind*<sup>®</sup> Goat Anti-Rabbit IgG(H+L), FITC Conjugate

HS111-01 100  $\mu$ l**Concentration**

2 mg/ml

**Storage**


at 2-8°C in dark for one year

**Description**

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Rabbit IgG(H+L) Antibody is conjugated with Fluorescein Isothiocyanate (FITC) dye under optimal conditions. FITC dye is a bright, yellow green-fluorescence dye with a maximal absorption wavelength at 490~495 nm and a maximal emission wavelength at 520~530 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

**Suggested Dilution**

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



## ProteinFind<sup>®</sup> Goat Anti-Rabbit IgG(H+L), PE Conjugate

HS121-01 100 µl

### Concentration

0.4 mg/ml

### Storage


at 2-8°C in dark for one year

### Description

Affinity purified ProteinFind<sup>®</sup> Goat Anti-Rabbit IgG(H+L) Antibody is conjugated with phycoerythrin (PE) dye under optimal conditions. PE dye is a natural fluorescent dye extracted from red algae with a maximal absorption wavelength at 488 nm and a maximal emission wavelength at 575 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

### Suggested Dilution

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



Detection of fusion protein expression in HEK-293 cell by IF

## ProteinFind<sup>®</sup> Goat Anti-Rabbit IgG(H+L), AF488 Conjugate

HS131-01 100 µl

### Concentration

1 mg/ml

### Storage


at 2-8°C in dark for one year

### Description

Affinity purified ProteinFind<sup>®</sup> Goat Anti-Rabbit IgG(H+L) Antibody is conjugated with Alexa Fluor 488 (AF488) dye under optimal conditions. AF488 dye is a bright, green-fluorescence dye with a maximal absorption wavelength at 495 nm and a maximal emission wavelength at 519 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

### Suggested Dilution

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



Detection of fusion protein expression in HEK-293 cell by IF

## *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L), HRP Conjugate

HS201-01 100 µl

### Concentration

1 mg/ml

### Storage


at 2-8°C for one month; at -20°C for one year

### Description

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L) Antibody is a horseradish peroxidase (HRP) conjugated secondary antibody for ELISA and Western blot detection.

### Suggested Dilution

- Western Blot: 1:1,000-10,000 dilution.
- ELISA: 1:1,000-5,000 dilution.



EasySee<sup>®</sup> Western Marker (5 µl/well)

## *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L), FITC Conjugate

HS211-01 100 µl

### Concentration

2 mg/ml

### Storage


at 2-8°C in dark for one year

### Description

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L) Antibody is conjugated with Fluorescein Isothiocyanate (FITC) dye under optimal conditions. FITC dye is a bright, yellow green-fluorescence dye with a maximal absorption wavelength at 490~495 nm and a maximal emission wavelength at 520~530 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

### Suggested Dilution

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



Detection of fusion protein expression in HEK-293 cell by IF

## *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L), PE Conjugate

HS221-01 100 µl

### Concentration

0.4 mg/ml

### Storage


at 2-8°C in dark for one year

### Description

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L) Antibody is conjugated with phycoerythrin (PE) dye under optimal conditions. PE dye is a natural fluorescent dye extracted from red algae with a maximal absorption wavelength at 488 nm and a maximal emission wavelength at 575 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

### Suggested Dilution

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



Detection of fusion protein expression in HEK-293 cell by IF

High quality products



## *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L), AF488 Conjugate

HS231-01 100  $\mu$ l**Concentration**

1 mg/ml

**Storage**


at 2-8°C in dark for one year

**Description**

Affinity purified *ProteinFind*<sup>®</sup> Goat Anti-Mouse IgG(H+L) Antibody is conjugated with Alexa Fluor 488 (AF488) dye under optimal conditions. AF488 dye is a bright, green-fluorescence dye with a maximal absorption wavelength at 495 nm and a maximal emission wavelength at 519 nm. This product has been optimized for use as a secondary antibody in immunofluorescence and flow cytometry.

**Suggested Dilution**

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.



Detection of fusion protein expression in HEK-293 cell by IF

## TMB ELISA Substrate

HE101-01 100 ml

### Storage

at 2-8°C in dark for one year

### Description

TMB ELISA Substrate is a ready-to-use chromogenic substrate for detection of horseradish peroxidase (HRP) activity. HRP can catalyze 3,3',5,5'-tetramethylbenzidine (TMB) to yield a blue color, the maximal absorbance is at 370 nm or 620-652 nm; however, upon addition of the stop solution, the solution turns to yellow and can be measured at 450 nm.

## Super TMB ELISA Substrate


HE111-01 100 ml

### Storage

at 2-8°C in dark for one year

### Description

Super TMB ELISA Substrate is a ready-to-use chromogenic substrate for detection of horseradish peroxidase (HRP) activity. HRP can catalyze 3,3',5,5'-tetramethylbenzidine (TMB) to yield a blue color, the maximal absorbance is at 370 nm or 620-652 nm. Upon addition of the stop solution, the solution turns to yellow and can be measured at 450 nm. This one-component method is 40-50% more sensitive than the traditional TMB ELISA method.



**Comparison of sensitivity between Super TMB and TMB ELISA Substrates**



---

## Chapter 9 Other Products

---

|  |     |
|--|-----|
| T4 DNA Ligase .....                        | 218 |
| DMT Enzyme .....                           | 219 |
| DNase I (RNase-free) .....                 | 219 |
| RNase A .....                              | 220 |
| Proteinase K .....                         | 220 |
| IPTG .....                                 | 220 |
| X-gal .....                                | 221 |
| Ampicillin .....                           | 221 |
| Kanamycin .....                            | 221 |
| Chloramphenicol .....                      | 221 |
| 6xDNA Loading Buffer .....                 | 222 |
| 2xRNA Loading Buffer .....                 | 222 |
| ddH <sub>2</sub> O .....                   | 222 |
| RNase-free Water .....                     | 222 |
| T7 High Efficiency Transcription Kit ..... | 223 |
| shRNA Synthesis Kit .....                  | 225 |

# T4 DNA Ligase

|          |              |
|----------|--------------|
| FL101-01 | 10,000 units |
| FL101-02 | 20,000 units |

## Concentration

200 units/μl

## Contents

- T4 DNA Ligase
- 5×T4 DNA Ligase Buffer

## Storage

at -20°C for one year

## Description

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5'-phosphate and 3'-hydroxyl termini in duplex DNA or RNA with blunt or cohesive end. The enzyme repairs single-strand nicks in duplex DNA, RNA or DNA/RNA hybrids but has no activity on single-strand nucleic acids. T4 DNA Ligase requires ATP as a cofactor.

## Source

*E.coli* strain carrying T4 DNA ligase gene

## Unit Definition

One unit is the amount of enzyme required to give 50% ligation of *Hind* III fragments of λDNA (5' DNA termini concentration of 0.12 μM, 200 μg/ml) in a total reaction volume of 20 μl in 30 minutes at 16°C in 1×T4 DNA Ligase Buffer.

## Quality Control

Functional absence of endonucleases and exonucleases activities

## Applications

- Cloning blunt end or cohesive end fragments.
- Ligation of synthetic linkers or adaptors.

## PROTOCOL

### Note

It is recommended to use a molar ratio of insert to vector at 3:1-10:1.

### Reaction Components

| Component              | Volume   | Final concentration |
|------------------------|----------|---------------------|
| Vector                 | Variable | as required         |
| Insert                 | Variable | as required         |
| 5×T4 DNA Ligase Buffer | 2 μl     | 1×                  |
| T4 DNA Ligase          | 0.5-1 μl | 100-200 units       |
| ddH <sub>2</sub> O     | Variable | -                   |
| Total volume           | 10 μl    | -                   |

### Reaction Conditions

- Cohesive ends ligation: incubate at 25°C for 10 minutes.
- Blunt ends ligation: incubate at 25°C for 2 hours, or overnight at 16°C.
- Cohesive and blunt ends ligation: incubate at 25°C for 2 hours.



## DMT Enzyme

GD111-01

200 units

### Concentration

10 units/μl

### Storage

at -20°C for two years

### Description

DMT cuts the sequence GmATC (G is methylated) but does not cut the sequence GATC (G is not methylated). This enzyme cuts DNA prepared from most commonly used *E.coli* strains (dam<sup>+</sup> strain), but does not cut PCR products.

### Source

An *E.coli* strain that carries the cloned DMT enzyme gene from *Diplococcus pneumoniae*.

### Unit Definition

One unit is the amount of enzyme required to completely digest 1 μg of pBR322 DNA (prepared from dam<sup>+</sup> strain) in 50 μl of reaction mixture in 1 hour at 37°C.

### Quality Control

Functional absence of endonucleases and exonucleases activities

### Applications

*In vitro* site-directed mutagenesis, digestion of methylated DNA.

## DNase I (RNase-free)

GD201-01

1,500 units

### Concentration

3 units/μl

### Contents

- DNase I
- 10×DNase I Reaction Buffer
- 200 mM EDTA

### Storage

at -20°C for one year

### Description

Deoxyribonuclease I (DNase I) is an endonuclease that degrades double- and single-strand DNA and chromatin. It functions by hydrolyzing phosphodiester linkages, producing mono and oligonucleotides with a 5'-phosphate and a 3'-hydroxyl group. Ribonuclease has been reduced to non-detectable levels. Its activity depends on Mg<sup>2+</sup> or Mn<sup>2+</sup> ion. DNase I with Mg<sup>2+</sup> randomly cuts double strand DNA at any sites, DNase I with Mn<sup>2+</sup> cuts double strand DNA at the same site to form sticky end with 1-2 nucleotide or form blunt end.

### Source

Purified from bovine pancreas.

### Unit Definition

One unit is the amount of enzyme required to completely degrade 1 μg pBR322 plasmid DNA in 10 minutes at 37°C.

### Applications

- DNase I footprinting
- Nick translation
- Remove DNA from RNA preparation

## RNase A

GE101-01

1 ml

**Concentration**

20 mg/ml

**Storage**

at -20°C for one year

**Description**

RNase A is a ribonuclease that cleaves single-strand RNA. It has no DNase activity.

**Source**

Bovine pancreas

**Activity**

&gt;60 U/mg

**Applications**

- Remove RNA from DNA samples
- RNase protection assay

## Proteinase K

GE201-01

1 ml

**Concentration**

20 mg/ml

**Storage**

at -20°C for one year

**Description**

Proteinase K is a nonspecific serine protease that will hydrolyze a variety of peptide bands. Proteinase K is active in a broad range of temperature and buffers. It cannot be inactivated by metal ions, chelating agents (e.g., EDTA), or detergents such as SDS.

- Active in a wide range of buffers and pH value.
- Incubation temperature: 55-65°C; optimal temperature: 58°C
- Incubation time: 15 minutes to 48 hours; optimal incubation time: 2 hours

**Source**

 Purified from *Tritirachium album*
**Applications**

- Preparation of DNA and RNA.
- Inactivation of RNase, DNase and enzymes.
- Isolation of genomic DNA.
- Isolation of RNA.

## IPTG

GF101-01

1 ml

**Concentration**

500 mM

**Storage**

at -20°C for six months

**Description**

Isopropylthio- $\beta$ -galactoside (IPTG) is an effective inducer of  $\beta$ -galactosidase activity. It is commonly used with X-gal to detect *lac* gene activity in cloning based on blue/white selection. It is also used as an inducer of protein expression in *lac* or *tac* promoter-regulated expression vectors.

High quality products



## X-Gal

GF201-01

1 ml

**Concentration**

20 mg/ml

**Storage**

at -20°C in dark for six months

**Description**

X-gal is a substrate of  $\beta$ -galactosidase. It is commonly used with IPTG to detect *lac* gene activity in cloning based on blue/white selection.

## Ampicillin

GG101-01

1 ml

**Concentration**

100 mg/ml

**Storage**

at -20°C for one year

**Description**

Extremely pure, molecular biology grade Ampicillin from TransGen can be used as a selective antibiotic for resistant bacteria.

## Kanamycin

GG201-01

1 ml

**Concentration**

50 mg/ml

**Storage**

at -20°C for one year

**Description**

Extremely pure, molecular biology grade Kanamycin from TransGen can be used as a selective antibiotic for resistant bacteria.

## Chloramphenicol

GG301-01

1 ml

**Concentration**

34 mg/ml

**Storage**

at -20°C for one year

**Description**

Extremely pure, molecular biology grade Chloramphenicol from TransGen can be used as a selective antibiotic for resistant bacteria.

## 6×DNA Loading Buffer

GH101-01

5×1 ml

### Storage

at -20°C for two years

### Description

6×DNA Loading Buffer is used as loading buffer in nucleic acid electrophoresis. Prior to loading, add appropriate volume of 6×DNA Loading Buffer to DNA sample to make its working concentration at 1×, and then load the DNA samples on the gel for electrophoresis.

## 2×RNA Loading Buffer

GH201-01

1 ml

### Storage

at 4°C for one month, at -20°C for two years

### Description

This product is used as loading buffer in RNA electrophoresis. It is suitable for denatured or native agarose gel and polyacrylamide gel electrophoresis. Prior to loading, add equal volume of the Loading Buffer to RNA sample (or RNA marker), heat at 70°C for 10 minutes, immediately chill on ice, and then load on the gel.

## ddH<sub>2</sub>O

GI101-01

25 ml

### Storage

at room temperature for two years

### Description

ddH<sub>2</sub>O is purified by reverse osmosis method. It is suitable for most molecular and cell biology applications.

## RNase-free Water

GI201-01

25 ml

### Storage

at room temperature for two years

### Description

RNase-free Water is prepared from deionized water incubated with 0.01% DEPC, then autoclaved to remove residual DEPC. It is suitable for RNA-related molecular biology applications.

High quality products



# T7 High Efficiency Transcription Kit

JT101-01

20 µl×25 rxns

## Storage

at -20°C for one year

## Description

T7 High Efficiency Transcription Kit is designed for *in vitro* RNA synthesis by T7 RNA Polymerase with supercoiled or linearized DNA templates. Up to 150 µg of RNA can be produced from a 20 µl reaction. Synthesized RNA can be used for *in vitro* translation, RNase protection assays, RNA splicing, and hybridization assays.

## Kit Contents

| Component                                  | JT101-01 |
|--|----------|
| T7 Transcription Enzyme Mix                | 50 µl    |
| 5×T7 Transcription Reaction Buffer         | 100 µl   |
| 10 mM NTP Mix                              | 200 µl   |
| DNase I (1 unit/µl)                        | 25 µl    |
| 500 mM EDTA (pH 8.0)                       | 25 µl    |
| RNase-free Water                           | 500 µl   |
| Control Transcription Template (0.5 µg/µl) | 10 µl    |


## PROTOCOL

### Notes

- RNase contamination should be avoided.
- Transcript produced from the control template is 2 kb.

## RNA Synthesis

### Principle of *In Vitro* Transcription



### Template Preparation

- Supercoiled plasmid DNA

Supercoiled plasmid DNA should contain a T7 promoter and an effective terminator. Termination efficiency varies with terminators. The following sequence is recommended.

T7 Promoter ———— Transcription template ———— Terminator

T7 Promoter: 5'-TAATACGACTCACTATAGGG<sup>\*</sup>-3' #: G/A

Terminator: 5'-TTCCATCTGTTTTCTTATCTGTTCTTTCATCTGTTCTTTTATCTGTTTGT-3'

- Linearized DNA  
Linearized plasmid DNA or PCR product, with T7 promoter and terminal sequences, can be used as template for *in vitro* transcription. We suggest to use 5'-overhang or blunt end restriction enzymes to generate the linearized templates, and avoid to use 3'-overhang restriction enzymes to generate the template. Digested linearized DNA should be purified.

### Transcription

- Reaction Components

| Component                       | Volume   |
|---------------------------------|----------|
| Template                        | 1 µg     |
| 5×Transcription Reaction Buffer | 4 µl     |
| 10 mM NTP Mix                   | 8 µl     |
| T7 Transcription Enzyme Mix     | 2 µl     |
| RNase-free Water                | to 20 µl |

- Mix thoroughly and incubate at 37°C for 2 hours.
- Add 1 µl of DNase I, incubate at 37°C for 15 minutes. Then add 1 µl of 500 mM EDTA (pH 8.0) to terminate reaction (immediately proceed to the following purification step after termination).

### Purification of Synthesized RNA

Please refer to *EasyPure*<sup>®</sup> RNA Purification Kit.

### Quantification and Analysis of synthesized RNA

- RNA concentration can be determined by ultraviolet light spectrophotometer.
- Transcripts of 0.1-1 kb can be run on denatured gel (6% acrylamide, 7 M urea). Use 1×TBE Buffer as the running buffer. (10×TBE Buffer: 0.9 M Tris Base, 0.9 M Boric Acid, 20 mM EDTA.)
- Transcripts of 0.5-5 kb can also be run on 1% formaldehyde denatured gel. Use 1×MOPS Buffer as the running buffer. (10×MOPS Buffer: 0.4 M MOPS (pH 7.0), 0.1 M Sodium Acetate, 10 mM EDTA.)
- For electrophoresis analysis, dilute 0.2-1 µg RNA with RNase-free water to make the total volume to 5 µl, add 5 µl of 2×RNA Loading Buffer and mix thoroughly, incubate at 70°C for 10 minutes and followed by incubation on ice for 2 minutes, then load samples on the gel. After electrophoresis, stain the gel.





# shRNA Synthesis Kit

JT111-01

20 µl×25 rxns

## Storage

at -20°C for one year

## Description

shRNA Synthesis Kit is designed for RNA synthesis by T7 RNA Polymerase with annealed dsDNA as template. Template contains Loop structure, thus synthesized RNA can be formed as shRNA. After RNA synthesis, DNA template can be digested by DNase I, followed by purification. Purified product can be directly used for cell transfection and other shRNA assays. This kit is suitable for production of RNA ≤ 100 nt. Up to 40 µg shRNA can be produced from a 20 µl reaction.

## Kit Contents

| Component  | JT111-01 |
|--|----------|
| shRNA Enzyme Mix                                       | 50 µl    |
| 5×shRNA Synthesis Buffer                               | 100 µl   |
| 10 mM NTP Mix  | 200 µl   |
| Annealing Buffer                                       | 450 µl   |
| DNase I (1 unit/µl)                                    | 25 µl    |
| 500 mM EDTA (pH 8.0)                                   | 25 µl    |
| RNase-free Water                                       | 500 µl   |
| Universal T7 Promoter Oligonucleotide (20 µM)          | 30 µl    |
| Control Transcription Template Oligonucleotide (20 µM) | 5 µl     |


## PROTOCOL

### Note

RNase contamination should be avoided.

## shRNA Synthesis

### Principle



### Oligonucleotide Template Design

Two Oligos are used for shRNA synthesis: Universal T7 Promoter Oligonucleotide (supplied with the kit) and Transcription Template Oligonucleotide (not supplied with the kit). Transcription Template Oligonucleotide contains sequences complementary to Universal T7 Promoter Oligonucleotide and sequences for shRNA.

Universal T7 Promoter Oligonucleotide (supplied with this kit)

5' -GAGTCCTGCAATTAATACGACTCACTATAG- 3'

Transcription Template Oligonucleotide (provided by users)

3' -CTCAGGACGTTAATTATGCTGAGTGATATCTTCTTCAGCACGACGAAGTATTCTTCGTCGTGCTGAAGAAGAA- 5'

Universal T7 Promoter  
Complement Sequence

Encoding siRNA Sequence Loop

Encoding siRNA  
Complement Sequence

### Preparation of Hybrid DNA Template

Dissolve Transcription Template Oligonucleotide with RNase-free Water or TE to a final concentration at 20  $\mu$ M.

- Reaction Components

| Component   | Volume     |
|---|------------|
| Universal T7 Promoter Oligonucleotide (20 $\mu$ M)  | 1 $\mu$ l  |
| Transcription Template Oligonucleotide (20 $\mu$ M) | 1 $\mu$ l  |
| Annealing Buffer                                    | 18 $\mu$ l |
| Total Volume  | 20 $\mu$ l |

Incubate at 95°C for 5 minutes, then slowly cool down to room temperature. (the prepared Hybrid DNA Template can be stored at 4°C for one week)

### shRNA Synthesis

- Reaction Components

| Component                | Volume        |
|--------------------------|---------------|
| Hybrid DNA Template      | 2 $\mu$ l     |
| 5xshRNA Synthesis Buffer | 4 $\mu$ l     |
| 10 mM NTP mix            | 8 $\mu$ l     |
| shRNA Enzyme Mix         | 2 $\mu$ l     |
| RNase-free Water         | to 20 $\mu$ l |

- Mix thoroughly and incubate at 37°C for 2 hours for shRNA synthesis.
- After shRNA synthesis reaction, add 1  $\mu$ l of DNase I to the reaction mixture and incubate at 37°C for 15 minutes to digest the templates. Stop the reaction by adding 1  $\mu$ l of 500 mM EDTA (pH8.0).

### Purification of shRNA

Please refer to *EasyPure*<sup>®</sup> RNA Purification Kit.

### Quantification and Analysis of shRNA

- shRNA concentration can be determined by ultraviolet light spectrophotometer.
- shRNA can be analyzed with denatured gel (12% acrylamide, 7 M urea.)



## Services

TransGen provides the following services with fast turnaround time and very competitive price. Please contact Customer Service Department for further information.

- PCR and qPCR
- RT-PCR and qRT-PCR
- Cloning
- Vector construction
- Plasmid DNA Mini, Midi, and Maxi Prep
- Genomic DNA Mini, Midi, and Maxi Prep
- Total RNA or/and mRNA Isolation
- Mutagenesis
- Protein expression and purification
- Establishment of stable cell line
- Lentivirus package
- Cell transfection
- Luciferase assay
- Other molecular and cell biology related services

**TransGen will keep all information confidential. All data and materials are the property of the purchasers. TransGen will not use, disclose, or publish the materials without the written consent of the purchasers.**

**PCR, RT-PCR, qPCR and qRT-PCR**

| <b>Products Name</b>  | <b>Catalog Number</b> | <b>Quantity</b> | <b>Page</b> |
|---|-----------------------|-----------------|-------------|
| <i>TransFast</i> <sup>®</sup> Taq DNA Polymerase  | AP101-01              | 500 units       | 6           |
|   | AP101-02              | 6×500 units     |             |
| <i>TransFast</i> <sup>®</sup> Taq DNA Polymerase<br>(with 2.5 mM dNTPs)                 | AP101-11              | 500 units       |             |
|   | AP101-12              | 6×500 units     |             |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase  | AP111-01              | 500 units       | 7           |
|   | AP111-02              | 6×500 units     |             |
|   | AP111-03              | 4×2,500 units   |             |
|   | AP111-04              | 10×5,000 units  |             |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase<br>(with 2.5 mM dNTPs)                       | AP111-11              | 500 units       |             |
|   | AP111-12              | 6×500 units     |             |
|   | AP111-13              | 4×2,500 units   |             |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase for PAGE                                     | AP112-01              | 2,500 units     | 9           |
|   | AP112-02              | 4×2,500 units   |             |
| <i>EasyTaq</i> <sup>®</sup> DNA Polymerase for PAGE<br>(with 2.5 mM dNTPs)              | AP112-11              | 2,500 units     |             |
|   | AP112-12              | 4×2,500 units   |             |
| <i>TransTaq</i> <sup>®</sup> -T DNA Polymerase  | AP122-01              | 250 units       | 10          |
|   | AP122-02              | 500 units       |             |
|   | AP122-03              | 6×500 units     |             |
| <i>TransTaq</i> <sup>®</sup> -T DNA Polymerase<br>(with 2.5 mM dNTPs)                   | AP122-11              | 250 units       |             |
|   | AP122-12              | 500 units       |             |
|   | AP122-13              | 6×500 units     |             |
| <i>TransTaq</i> <sup>®</sup> DNA Polymerase High Fidelity (HiFi)                        | AP131-01              | 250 units       | 11          |
|   | AP131-02              | 500 units       |             |
|   | AP131-03              | 6×500 units     |             |
| <i>TransTaq</i> <sup>®</sup> DNA Polymerase High Fidelity (HiFi)<br>(with 2.5 mM dNTPs) | AP131-11              | 250 units       |             |
|   | AP131-12              | 500 units       |             |
|   | AP131-13              | 6×500 units     |             |
| <i>TransStart</i> <sup>®</sup> Taq DNA Polymerase                                       | AP141-01              | 250 units       | 14          |
|   | AP141-02              | 500 units       |             |
|   | AP141-03              | 6×500 units     |             |
| <i>TransStart</i> <sup>®</sup> Taq DNA Polymerase<br>(with 2.5 mM dNTPs)                | AP141-11              | 250 units       |             |
|   | AP141-12              | 500 units       |             |
|   | AP141-13              | 6×500 units     |             |
| <i>TransStart</i> <sup>®</sup> TopTaq DNA Polymerase                                    | AP151-01              | 250 units       | 16          |
|   | AP151-02              | 500 units       |             |
|   | AP151-03              | 6×500 units     |             |
| <i>TransStart</i> <sup>®</sup> TopTaq DNA Polymerase<br>(with 2.5 mM dNTPs)             | AP151-11              | 250 units       |             |
|   | AP151-12              | 500 units       |             |
|   | AP151-13              | 6×500 units     |             |
| <i>EasyPfu</i> DNA Polymerase   | AP211-01              | 250 units       | 18          |
|   | AP211-02              | 500 units       |             |
|   | AP211-03              | 6×500 units     |             |

High quality products



| Products Name   | Catalog Number | Quantity    | Page |
|---|----------------|-------------|------|
| <i>EasyPfu</i> DNA Polymerase<br>(with 2.5 mM dNTPs)                                    | AP211-11       | 250 units   | 18   |
|   | AP211-12       | 500 units   |      |
|   | AP211-13       | 6×500 units |      |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> DNA Polymerase                            | AP221-01       | 250 units   | 19   |
|   | AP221-02       | 500 units   |      |
|   | AP221-03       | 6×500 units |      |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> DNA Polymerase<br>(with 2.5 mM dNTPs)     | AP221-11       | 250 units   | 21   |
|   | AP221-12       | 500 units   |      |
|   | AP221-13       | 6×500 units |      |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> Fly DNA Polymerase                        | AP231-01       | 250 units   | 21   |
|   | AP231-02       | 500 units   |      |
|   | AP231-03       | 6×500 units |      |
| <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> Fly DNA Polymerase<br>(with 2.5 mM dNTPs) | AP231-11       | 250 units   | 23   |
|   | AP231-12       | 500 units   |      |
|   | AP231-13       | 6×500 units |      |
| <i>TransStart</i> <sup>®</sup> <i>KD</i> Plus DNA Polymerase                            | AP301-01       | 100 units   | 23   |
|   | AP301-02       | 200 units   |      |
|   | AP301-03       | 6×200 units |      |
| <i>TransStart</i> <sup>®</sup> <i>KD</i> Plus DNA Polymerase<br>(with 2.5 mM dNTPs)     | AP301-11       | 100 units   | 25   |
|   | AP301-12       | 200 units   |      |
|   | AP301-13       | 6×200 units |      |
| GC Enhancer   | AG101-01       | 200 µl      | 25   |
| PCR Stimulant   | AG111-01       | 200 µl      |      |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix (-dye)                                      | AS111-01       | 1 ml        | 27   |
|   | AS111-02       | 5×1 ml      |      |
|   | AS111-03       | 15×1 ml     |      |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix (+dye)                                      | AS111-11       | 1 ml        | 28   |
|   | AS111-12       | 5×1 ml      |      |
|   | AS111-13       | 15×1 ml     |      |
|   | AS111-14       | 6×80 ml     |      |
| 2× <i>EasyTaq</i> <sup>®</sup> PCR SuperMix for PAGE (+dye)                             | AS112-11       | 1 ml        | 29   |
|   | AS112-12       | 5×1 ml      |      |
|   | AS112-13       | 15×1 ml     |      |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix (-dye)                                  | AS122-01       | 1 ml        | 29   |
|   | AS122-02       | 5×1 ml      |      |
| 2× <i>TransTaq</i> <sup>®</sup> -T PCR SuperMix (+dye)                                  | AS122-11       | 1 ml        | 30   |
|   | AS122-12       | 5×1 ml      |      |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix I<br>(-dye)           | AS131-01       | 1 ml        | 30   |
|   | AS131-02       | 5×1 ml      |      |
| 2× <i>TransTaq</i> <sup>®</sup> High Fidelity (HiFi) PCR SuperMix II<br>(-dye)          | AS131-21       | 1 ml        | 32   |
|   | AS131-22       | 5×1 ml      |      |
| 2× <i>EasyPfu</i> PCR SuperMix (-dye)   | AS211-01       | 1 ml        | 32   |
|   | AS211-02       | 5×1 ml      |      |

| <b>Products Name</b>  | <b>Catalog Number</b> | <b>Quantity</b>                 | <b>Page</b> |
|---|-----------------------|---------------------------------|-------------|
| <i>2xTransStart</i> <sup>®</sup> <i>FastPfu</i> PCR SuperMix (-dye)   | AS221-01              | 1 ml                            | 33          |
|   | AS221-02              | 5×1 ml                          |             |
| <i>TransDirect</i> <sup>®</sup> Animal Tissue PCR Kit   | AD201-01              | 100 rxns (20 µl per reaction)   | 34          |
|   | AD201-02              | 500 rxns (20 µl per reaction)   |             |
| <i>TransDirect</i> <sup>®</sup> Plant Tissue PCR Kit  | AD301-01              | 100 rxns (20 µl per reaction)   | 36          |
|   | AD301-02              | 500 rxns (20 µl per reaction)   |             |
| <i>TransDirect</i> <sup>®</sup> Blood PCR Kit   | AD401-01              | 100 rxns (20 µl per reaction)   | 37          |
|   | AD401-02              | 500 rxns (20 µl per reaction)   |             |
| <i>EasyScript</i> <sup>®</sup> Reverse Transcriptase  | AE101-02              | 10,000 units                    | 41          |
|   | AE101-03              | 5×10,000 units                  |             |
| <i>TransScript</i> <sup>®</sup> Reverse Transcriptase   | AT101-02              | 10,000 units                    | 43          |
|   | AT101-03              | 5×10,000 units                  |             |
| <i>TransScript</i> <sup>®</sup> II Reverse Transcriptase  | AH101-02              | 10,000 units                    | 44          |
| <i>EasyScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix   | AE301-02              | 50 rxns (20 µl per reaction)    | 45          |
|   | AE301-03              | 100 rxns (20 µl per reaction)   |             |
| <i>EasyScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix                                    | AE311-02              | 50 rxns (20 µl per reaction)    | 47          |
|   | AE311-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> First-Strand cDNA Synthesis SuperMix  | AT301-02              | 50 rxns (20 µl per reaction)    | 48          |
|   | AT301-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> One-Step gDNA Removal and cDNA Synthesis SuperMix                                   | AT311-02              | 50 rxns (20 µl per reaction)    | 49          |
|   | AT311-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> Fly First-Strand cDNA Synthesis SuperMix  | AF301-02              | 50 rxns (20 µl per reaction)    | 50          |
|   | AF301-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> -Uni One-Step gDNA Removal and cDNA Synthesis SuperMix                              | AU311-02              | 50 rxns (20 µl per reaction)    | 51          |
|   | AU311-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> -Uni Cell to cDNA Synthesis SuperMix for qPCR                                       | AC301-01              | 25 rxns                         | 53          |
| <i>TransScript</i> <sup>®</sup> miRNA First-Strand cDNA Synthesis SuperMix  | AT351-01              | 20 rxns (20 µl per reaction)    | 55          |
| <i>TransScript</i> <sup>®</sup> II First-Strand cDNA Synthesis SuperMix   | AH301-02              | 50 rxns (20 µl per reaction)    | 57          |
|   | AH301-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> II One-Step gDNA Removal and cDNA Synthesis SuperMix                                | AH311-02              | 50 rxns (20 µl per reaction)    | 58          |
|   | AH311-03              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for PCR                             | AT321-01              | 50 rxns (20 µl per reaction)    | 59          |
| <i>TransScript</i> <sup>®</sup> All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)    | AT341-01              | 50 rxns (20 µl per reaction)    | 60          |
|   | AT341-02              | 100 rxns (20 µl per reaction)   |             |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for PCR                          | AH321-01              | 50 rxns (20 µl per reaction)    | 62          |
| <i>TransScript</i> <sup>®</sup> II All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal) | AH341-01              | 50 rxns (20 µl per reaction)    | 63          |
| <i>TransScript</i> <sup>®</sup> Two-Step RT-PCR SuperMix  | AT401-01              | 50 rxns (20 µl per RT reaction) | 64          |
|   |                       | 80 rxns (50 µl per PCR)         |             |
| <i>TransScript</i> <sup>®</sup> II Two-Step RT-PCR SuperMix   | AH401-01              | 50 rxns (20 µl per RT reaction) | 65          |
|   |                       | 80 rxns (50 µl per PCR)         |             |
| <i>EasyScript</i> <sup>®</sup> One-Step RT-PCR SuperMix (+dye)  | AE411-02              | 200 rxns (20 µl per reaction)   | 66          |
| <i>TransScript</i> <sup>®</sup> One-Step RT-PCR SuperMix (+dye)   | AT411-02              | 200 rxns (20 µl per reaction)   | 67          |
| <i>TransScript</i> <sup>®</sup> II One-Step RT-PCR SuperMix (+dye)  | AH411-02              | 200 rxns (20 µl per reaction)   | 68          |
| Ribonuclease Inhibitor  | AI101-01              | 2,000 units                     | 70          |
|   | AI101-02              | 5×2,000 units                   |             |

High quality products



| Products Name   | Catalog Number | Quantity                        | Page |
|---|----------------|---------------------------------|------|
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix                    | AQ101-01       | 1 ml                            | 71   |
|   | AQ101-02       | 5×1 ml                          |      |
|   | AQ101-03       | 15×1 ml                         |      |
| <i>TransStart</i> <sup>®</sup> Green qPCR SuperMix UDG                | AQ111-01       | 1 ml                            | 73   |
|   | AQ111-02       | 5×1 ml                          |      |
|   | AQ111-03       | 15×1 ml                         |      |
| <i>TransStart</i> <sup>®</sup> Top Green qPCR SuperMix                | AQ131-01       | 1 ml                            | 74   |
|   | AQ131-02       | 5×1 ml                          |      |
|   | AQ131-03       | 15×1 ml                         |      |
|   | AQ131-04       | 25×1 ml                         |      |
| <i>TransStart</i> <sup>®</sup> Tip Green qPCR SuperMix                | AQ141-01       | 1 ml                            | 75   |
|   | AQ141-02       | 5×1 ml                          |      |
|   | AQ141-03       | 15×1 ml                         |      |
|   | AQ141-04       | 25×1 ml                         |      |
| <i>TransScript</i> <sup>®</sup> Green Two-Step qRT-PCR SuperMix       | AQ201-01       | 50 rxns (20 µl per RT reaction) | 76   |
|   |                | 300 rxns (20 µl per qPCR)       |      |
| <i>TransScript</i> <sup>®</sup> Green miRNA Two-Step qRT-PCR SuperMix | AQ202-01       | 20 rxns (20 µl per RT reaction) | 78   |
|   |                | 500 rxns (20 µl per qPCR)       |      |
| <i>TransScript</i> <sup>®</sup> II Green Two-Step qRT-PCR SuperMix    | AQ301-01       | 50 rxns (20 µl per RT reaction) | 79   |
|   |                | 300 rxns (20 µl per qPCR)       |      |
| <i>TransScript</i> <sup>®</sup> Green One-Step qRT-PCR SuperMix       | AQ211-01       | 100 rxns (20 µl per reaction)   | 80   |
|   | AQ211-02       | 400 rxns (20 µl per reaction)   |      |
| <i>TransScript</i> <sup>®</sup> II Green One-Step qRT-PCR SuperMix    | AQ311-01       | 100 rxns (20 µl per reaction)   | 82   |
|   | AQ311-02       | 400 rxns (20 µl per reaction)   |      |
| <i>TransStart</i> <sup>®</sup> Probe qPCR SuperMix                    | AQ401-01       | 1 ml                            | 84   |
|   | AQ401-02       | 5×1 ml                          |      |
|   | AQ401-03       | 15×1 ml                         |      |
| <i>TransScript</i> <sup>®</sup> Probe One-Step qRT-PCR SuperMix       | AQ221-01       | 100 rxns (20 µl per reaction)   | 85   |
|   | AQ221-02       | 400 rxns (20 µl per reaction)   |      |
| <i>TransScript</i> <sup>®</sup> II Probe One-Step qRT-PCR SuperMix    | AQ321-01       | 100 rxns (20 µl per reaction)   | 87   |
|   | AQ321-02       | 400 rxns (20 µl per reaction)   |      |
| High Pure dNTPs (2.5 mM)  | AD101-01       | 1 ml                            | 89   |
|   | AD101-02       | 5×1 ml                          |      |
| High Pure dNTPs (10 mM)   | AD101-11       | 1 ml                            | 89   |
|   | AD101-12       | 5×1 ml                          |      |

## DNA Molecular Weight Standards

| Products Name                                  | Catalog Number | Quantity | Page |
|--|----------------|----------|------|
| <i>Trans2K</i> <sup>®</sup> DNA Marker         | BM101-01       | 500 µl   | 92   |
|  | BM101-02       | 5×500 µl |      |
| <i>Trans2K</i> <sup>®</sup> Plus DNA Marker    | BM111-01       | 500 µl   | 92   |
|  | BM111-02       | 5×500 µl |      |
| <i>Trans2K</i> <sup>®</sup> Plus II DNA Marker | BM121-01       | 500 µl   | 92   |
|  | BM121-02       | 5×500 µl |      |

| <b>Products Name</b>        | <b>Catalog Number</b> | <b>Quantity</b> | <b>Page</b> |    |
|-----------------------------|-----------------------|-----------------|-------------|----|
| <i>Trans</i> 5K DNA Marker  | BM141-01              | 500 µl          | 93          |    |
|                             | BM141-02              | 5×500 µl        |             |    |
| <i>Trans</i> 15K DNA Marker | BM161-01              | 500 µl          |             |    |
|                             | BM161-02              | 5×500 µl        |             |    |
| 1Kb DNA Ladder              | BM201-01              | 500 µl          |             |    |
|                             | BM201-02              | 5×500 µl        |             |    |
| 1Kb Plus DNA Ladder         | BM211-01              | 500 µl          |             |    |
|                             | BM211-02              | 5×500 µl        |             |    |
| 100bp DNA Ladder            | BM301-01              | 500 µl          |             | 94 |
|                             | BM301-02              | 5×500 µl        |             |    |
| 100bp Plus DNA Ladder       | BM311-01              | 500 µl          |             |    |
|                             | BM311-02              | 5×500 µl        |             |    |
| 100bp Plus II DNA Ladder    | BM321-01              | 500 µl          | 95          |    |
|                             | BM321-02              | 5×500 µl        |             |    |
| GelStain                    | GS101-01              | 500 µl          |             |    |
|                             | GS101-02              | 1 ml            |             |    |
| Agarose                     | GS201-01              | 100 g           |             |    |

## Cloning and Mutagenesis System

| <b>Products Name</b>   | <b>Catalog Number</b> | <b>Quantity</b> | <b>Page</b> |
|--|-----------------------|-----------------|-------------|
| <i>pEASY</i> <sup>®</sup> -T1 Cloning Kit                        | CT101-01              | 20 rxns         | 99          |
|  | CT101-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -Blunt Cloning Kit                     | CB101-01              | 20 rxns         | 102         |
|  | CB101-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -T1 Simple Cloning Kit                 | CT111-01              | 20 rxns         | 103         |
|  | CT111-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -Blunt Simple Cloning Kit              | CB111-01              | 20 rxns         | 104         |
|  | CB111-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -T3 Cloning Kit                        | CT301-01              | 20 rxns         | 105         |
|  | CT301-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -Blunt3 Cloning Kit                    | CB301-01              | 20 rxns         | 106         |
|  | CB301-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -T5 Zero Cloning Kit                   | CT501-01              | 20 rxns         | 107         |
|  | CT501-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -Blunt Zero Cloning Kit                | CB501-01              | 20 rxns         | 108         |
|  | CB501-02              | 60 rxns         |             |
| <i>pEASY</i> <sup>®</sup> -Uni Seamless Cloning and Assembly Kit | CU101-01              | 10 rxns         | 109         |
| <i>Trans</i> 10 Chemically Competent Cell                        | CD101-01              | 10×100 µl       | 112         |
|  | CD101-02              | 20×100 µl       |             |
| <i>Trans</i> 5α Chemically Competent Cell                        | CD201-01              | 10×100 µl       |             |
|  | CD201-02              | 20×100 µl       |             |
| <i>Trans</i> 109 Chemically Competent Cell                       | CD301-02              | 10×100 µl       | 113         |
|  | CD301-03              | 20×100 µl       |             |
| <i>Trans</i> 110 Chemically Competent Cell                       |                       | CD311-02        |             |

High quality products





| Products Name   | Catalog Number | Quantity  | Page |
|---|----------------|-----------|------|
| <i>Trans1</i> -Blue Chemically Competent Cell               | CD401-02       | 10×100 µl | 113  |
|   | CD401-03       | 20×100 µl |      |
| <i>Trans2</i> -Blue Chemically Competent Cell               | CD411-02       | 10×100 µl |      |
|   | CD411-03       | 20×100 µl |      |
| <i>Trans1</i> -T1 Phage Resistant Chemically Competent Cell | CD501-01       | 5×100 µl  | 114  |
|   | CD501-02       | 10×100 µl |      |
|   | CD501-03       | 20×100 µl |      |
| DMT Chemically Competent Cell                               | CD511-01       | 10×50 µl  |      |
|   | CD511-02       | 20×50 µl  |      |
| <i>TransStbl3</i> Chemically Competent Cell                 | CD521-01       | 10×100 µl | 115  |
| <i>TransDB3.1</i> Chemically Competent Cell                 | CD531-01       | 10×100 µl |      |
| Fast Mutagenesis System                                     | FM111-01       | 10 rxns   | 116  |
|   | FM111-02       | 20 rxns   |      |
| Fast MultiSite Mutagenesis System                           | FM201-01       | 10 rxns   | 117  |

## Nucleic Acid Purification

| Products Name   | Catalog Number | Quantity         | Page |
|---|----------------|------------------|------|
| <i>BloodZol</i>   | EE131-01       | for 50 ml blood  | 121  |
|   | EE131-02       | for 200 ml blood |      |
| <i>PlantZol</i>   | EE141-01       | 100 ml           | 122  |
| <i>EasyPure</i> <sup>®</sup> Genomic DNA Kit (with RNase A)               | EE101-01       | 50 rxns          | 123  |
|   | EE101-02       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Genomic DNA Kit                              | EE101-11       | 50 rxns          |      |
|   | EE101-12       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Plant Genomic DNA Kit (with RNase A)         | EE111-01       | 50 rxns          | 125  |
|   | EE111-02       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Plant Genomic DNA Kit                        | EE111-11       | 50 rxns          |      |
|   | EE111-12       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Blood Genomic DNA Kit (with RNase A)         | EE121-01       | 50 rxns          | 126  |
|   | EE121-02       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Blood Genomic DNA Kit                        | EE121-11       | 50 rxns          |      |
|   | EE121-12       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> Marine Animal Genomic DNA Kit (with RNase A) | EE151-01       | 50 rxns          | 127  |
|   | EE151-11       | 50 rxns          |      |
| <i>EasyPure</i> <sup>®</sup> Bacteria Genomic DNA Kit (with RNase A)      | EE161-01       | 50 rxns          | 128  |
|   | EE161-11       | 50 rxns          |      |
| <i>EasyPure</i> <sup>®</sup> Food and Fodder Security Genomic DNA Kit     | EE171-01       | 50 rxns          | 129  |
| <i>EasyPure</i> <sup>®</sup> Micro Genomic DNA Kit                        | EE181-01       | 50 rxns          | 131  |
| <i>EasyPure</i> <sup>®</sup> Plasmid MiniPrep Kit                         | EM101-01       | 50 rxns          | 132  |
|   | EM101-02       | 200 rxns         |      |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MiniPrep Kit                  | EM111-01       | 50 rxns          | 133  |
| <i>EasyPure</i> <sup>®</sup> HiPure Plasmid MaxiPrep Kit                  | EM121-01       | 10 rxns          | 134  |
| <i>ArtMedia</i> <sup>®</sup> Plasmid Culture                              | EM201-01       | 95 ml+5 ml       | 135  |
| <i>EasyPure</i> <sup>®</sup> PCR Purification Kit                         | EP101-01       | 50 rxns          | 136  |
|   | EP101-02       | 200 rxns         |      |

| Products Name   | Catalog Number | Quantity | Page |
|---|----------------|----------|------|
| <i>EasyPure</i> <sup>®</sup> Quick Gel Extraction Kit | EG101-01       | 50 rxns  | 137  |
|   | EG101-02       | 200 rxns |      |
| <i>TransZol</i>                                       | ET101-01       | 100 ml   | 138  |
| <i>TransZol</i> Up                                    | ET111-01       | 100 ml   | 139  |
| <i>TransZol</i> Plant                                 | ET121-01       | 100 ml   | 140  |
| <i>EasyPure</i> <sup>®</sup> RNA Kit                  | ER101-01       | 50 rxns  | 141  |
| <i>EasyPure</i> <sup>®</sup> Viral DNA/RNA Kit        | ER201-01       | 50 rxns  | 142  |
| <i>EasyPure</i> <sup>®</sup> Plant RNA Kit            | ER301-01       | 50 rxns  | 143  |
| <i>EasyPure</i> <sup>®</sup> Blood RNA Kit            | ER401-01       | 50 rxns  | 144  |
| <i>TransZol</i> Up Plus RNA Kit                       | ER501-01       | 100 rxns | 145  |
| <i>EasyPure</i> <sup>®</sup> miRNA Kit                | ER601-01       | 50 rxns  | 146  |
| <i>EasyPure</i> <sup>®</sup> RNA Purification Kit     | ER701-01       | 25 rxns  | 147  |
| <i>RNAhold</i> <sup>®</sup>                           | EH101-01       | 100 ml   | 148  |

## Gene Expression

| Products Name                                      | Catalog Number | Quantity   | Page |
|--|----------------|------------|------|
| <i>pEASY</i> <sup>®</sup> -Blunt E1 Expression Kit | CE111-01       | 10 rxns    | 151  |
| <i>pEASY</i> <sup>®</sup> -Blunt E2 Expression Kit | CE211-01       | 10 rxns    | 154  |
| <i>ArtMedia</i> <sup>®</sup> Protein Expression    | CP101-01       | 95 ml+5 ml | 155  |
| BL21(DE3) Chemically Competent Cell                | CD601-02       | 10×100 µl  | 156  |
|  | CD601-03       | 20×100 µl  |      |
| BL21(DE3)pLysS Chemically Competent Cell           | CD701-02       | 10×100 µl  | 157  |
|  | CD701-03       | 20×100 µl  |      |
| <i>Transetta</i> (DE3) Chemically Competent Cell   | CD801-02       | 10×100 µl  | 157  |
|  | CD801-03       | 20×100 µl  |      |
| <i>TransB</i> (DE3) Chemically Competent Cell      | CD811-02       | 10×100 µl  | 157  |
| BL21 Chemically Competent Cell                     | CD901-02       | 10×100 µl  | 158  |
|  | CD901-03       | 20×100 µl  |      |
| <i>pEASY</i> <sup>®</sup> -Blunt M2 Expression Kit | CM211-01       | 10 rxns    | 158  |
| <i>pEASY</i> <sup>®</sup> -Blunt M3 Expression Kit | CM311-01       | 10 rxns    | 161  |

## Protein Extraction, Purification and Detection

| Products Name   | Catalog Number | Quantity | Page |
|---|----------------|----------|------|
| <i>ProteinExt</i> <sup>™</sup> Mammalian Total Protein Extraction Kit                   | DE101-01       | 100 ml   | 164  |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Nuclear and Cytoplasmic Protein Extraction Kit | DE201-01       | 50 rxns  | 165  |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Membrane Protein Extraction Kit                | DE301-01       | 50 rxns  | 166  |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit for Cultured Cells  | DE401-01       | 50 rxns  | 167  |
| <i>ProteinExt</i> <sup>™</sup> Mammalian Mitochondria Isolation Kit for Tissue          | DE501-01       | 50 rxns  | 168  |
| <i>ProteinIso</i> <sup>®</sup> Ni-NTA Resin   | DP101-01       | 5 ml     | 169  |
|   | DP101-02       | 25 ml    |      |
| <i>ProteinIso</i> <sup>®</sup> Ni-IDA Resin   | DP111-01       | 5 ml     | 171  |
|   | DP111-02       | 25 ml    |      |

High quality products



| Products Name  | Catalog Number | Quantity      | Page |
|--|----------------|---------------|------|
| <i>ProteinIso</i> <sup>®</sup> GST Resin   | DP201-01       | 10 ml         | 173  |
| <i>ProteinIso</i> <sup>®</sup> Protein A Resin   | DP301-01       | 5 ml          | 175  |
| <i>ProteinIso</i> <sup>®</sup> Protein G Resin   | DP401-01       | 5 ml          | 177  |
| <i>ProteinRuler</i> <sup>®</sup> I (12-80 kDa)   | DR101-01       | 250 µl        | 180  |
|  | DR101-02       | 500 µl        |      |
| <i>ProteinRuler</i> <sup>®</sup> II (12-120 kDa)   | DR201-01       | 250 µl        | 181  |
|  | DR201-02       | 500 µl        |      |
| <i>ProteinRuler</i> <sup>®</sup> IV (30-200 kDa)   | DR401-01       | 250 µl        | 182  |
|  | DR401-02       | 500 µl        |      |
| <i>Blue Plus</i> <sup>®</sup> Protein Marker (14-100 kDa)  | DM101-01       | 250 µl        | 183  |
|  | DM101-02       | 500 µl        |      |
| <i>Blue Plus</i> <sup>®</sup> II Protein Marker (14-120 kDa)   | DM111-01       | 250 µl        | 184  |
|  | DM111-02       | 500 µl        |      |
| <i>Blue Plus</i> <sup>®</sup> III Protein Marker (14-160 kDa)  | DM121-01       | 250 µl        | 185  |
|  | DM121-02       | 500 µl        |      |
| <i>Blue Plus</i> <sup>®</sup> IV Protein Marker (10-180 kDa)   | DM131-01       | 250 µl        | 186  |
|  | DM131-02       | 500 µl        |      |
| <i>EasySee</i> <sup>®</sup> Western Marker (25-90 kDa)   | DM201-01       | 250 µl        | 187  |
|  | DM201-02       | 500 µl        |      |
| <i>EasySee</i> <sup>®</sup> Western Marker<br>(with <i>EasySee</i> <sup>®</sup> Western Blot Kit)    | DM201-11       | 250 µl+100 ml | 188  |
|  | DM201-12       | 500 µl+200 ml |      |
| <i>EasySee</i> <sup>®</sup> II Western Marker<br>(30-150 kDa)  | DM211-01       | 250 µl        | 189  |
|  | DM211-02       | 500 µl        |      |
| <i>EasySee</i> <sup>®</sup> II Western Marker<br>(with <i>EasySee</i> <sup>®</sup> Western Blot Kit) | DM211-11       | 250 µl+100 ml | 190  |
|  | DM211-12       | 500 µl+200 ml |      |
| <i>EasySee</i> <sup>®</sup> Western Blot Kit   | DW101-01       | 100 ml        | 191  |
|  | DW101-02       | 200 ml        |      |
| 6× Protein Loading Buffer  | DL101-02       | 5×1 ml        | 192  |
| <i>Easy</i> Protein Quantitative Kit (Bradford)  | DQ101-01       | 100 ml        | 193  |
| <i>Easy</i> II Protein Quantitative Kit (BCA)  | DQ111-01       | 100 ml        | 194  |
| <i>ProteinEle</i> <sup>™</sup> Precast Tris-Glycine Gel  | DG101-01       | 8%, 10/Box    | 195  |
|  | DG101-02       | 10%, 10/Box   |      |
|  | DG101-03       | 12%, 10/Box   |      |

## Cell Culture and Detection

| Products Name  | Catalog Number | Quantity  | Page |
|--|----------------|-----------|------|
| <i>TransSerum</i> <sup>®</sup> HQ Fetal Bovine Serum   | FS101-02       | 500 ml    | 192  |
| <i>TransLipid</i> <sup>®</sup> HL Transfection Reagent | FT111-01       | 0.75 ml   | 193  |
|  | FT111-02       | 2×0.75 ml |      |
| <i>TransIn</i> <sup>™</sup> EL Transfection Reagent    | FT201-01       | 0.75 ml   | 194  |
|  | FT201-02       | 2×0.75 ml |      |
| Penicillin-Streptomycin (100×)                         | FG101-01       | 100 ml    | 197  |
| L-Glutamine (100×)                                     | FG201-01       | 100 ml    |      |
| Trypsin (+EDTA)  | FG301-01       | 100 ml    |      |
| Trypsin (-EDTA)  | FG301-11       | 100 ml    |      |
| G418   | FG401-01       | 5 ml      |      |

| Products Name   | Catalog Number | Quantity | Page |
|---|----------------|----------|------|
| PBS (1x)  | FG701-01       | 500 ml   |      |
| <i>TransDetect</i> <sup>®</sup> Double-Luciferase Reporter Assay Kit                          | FR201-01       | 50 rxns  | 198  |
|   | FR201-02       | 200 rxns |      |
| <i>TransDetect</i> <sup>®</sup> Cell Counting Kit (CCK)                                       | FC101-01       | 1 ml     | 199  |
|   | FC101-02       | 5 ml     |      |
|   | FC101-03       | 10 ml    |      |
|   | FC101-04       | 30 ml    |      |
| <i>TransDetect</i> <sup>®</sup> Annexin V-FITC/PI Cell Apoptosis Detection Kit                | FA101-01       | 25 rxns  | 201  |
|   | FA101-02       | 50 rxns  |      |
| <i>TransDetect</i> <sup>®</sup> Annexin V-EGFP/PI Cell Apoptosis Detection Kit                | FA111-01       | 25 rxns  | 202  |
|   | FA111-02       | 50 rxns  |      |
| <i>TransDetect</i> <sup>®</sup> <i>In Situ</i> Fluorescein TUNEL Cell Apoptosis Detection Kit | FA201-01       | 25 rxns  | 203  |
|   | FA201-02       | 50 rxns  |      |

## Antibodies

| Products Name   | Catalog Number | Quantity | Page |
|---|----------------|----------|------|
| <i>ProteinFind</i> <sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody        | HT101-01       | 50 µl    | 205  |
|   | HT101-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-DYKDDDDK Tag Mouse Monoclonal Antibody | HT201-01       | 50 µl    | 206  |
|   | HT201-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-HA Mouse Monoclonal Antibody           | HT301-01       | 50 µl    | 207  |
|   | HT301-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-V5 Mouse Monoclonal Antibody           | HT401-01       | 50 µl    | 208  |
|   | HT401-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-His Mouse Monoclonal Antibody          | HT501-01       | 50 µl    | 209  |
|   | HT501-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-GST Mouse Monoclonal Antibody          | HT601-01       | 50 µl    | 210  |
|   | HT601-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-MBP Mouse Monoclonal Antibody          | HT701-01       | 50 µl    | 211  |
|   | HT701-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-GFP Mouse Monoclonal Antibody          | HT801-01       | 50 µl    | 212  |
|   | HT801-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-β-Tubulin Mouse Monoclonal Antibody    | HC101-01       | 50 µl    | 213  |
|   | HC101-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-β-Actin Mouse Monoclonal Antibody      | HC201-01       | 50 µl    | 214  |
|   | HC201-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Anti-GAPDH Mouse Monoclonal Antibody        | HC301-01       | 50 µl    | 215  |
|   | HC301-02       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), HRP Conjugate    | HS101-01       | 100 µl   | 216  |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), FITC Conjugate   | HS111-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), PE Conjugate     | HS121-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Rabbit IgG(H+L), AF488 Conjugate  | HS131-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), HRP Conjugate     | HS201-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), FITC Conjugate    | HS211-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), PE Conjugate      | HS221-01       | 100 µl   |      |
| <i>ProteinFind</i> <sup>®</sup> Goat Anti-Mouse IgG(H+L), AF488 Conjugate   | HS231-01       | 100 µl   |      |
| TMB ELISA Substrate   | HE101-01       | 100 ml   |      |
| Super TMB ELISA Substrate   | HE111-01       | 100 ml   |      |

High quality products



## Other Products

| Products Name                        | Catalog Number | Quantity      | Page |
|--------------------------------------|----------------|---------------|------|
| T4 DNA Ligase                        | FL101-01       | 10,000 units  | 218  |
|                                      | FL101-02       | 20,000 units  |      |
| DMT Enzyme                           | GD111-01       | 200 units     | 219  |
| DNase I (RNase-free)                 | GD201-01       | 1,500 units   |      |
| RNase A                              | GE101-01       | 1 ml          | 220  |
| Proteinase K                         | GE201-01       | 1 ml          |      |
| IPTG                                 | GF101-01       | 1 ml          |      |
| X-gal                                | GF201-01       | 1 ml          |      |
| Ampicillin                           | GG101-01       | 1 ml          |      |
| Kanamycin                            | GG201-01       | 1 ml          |      |
| Chloramphenicol                      | GG301-01       | 1 ml          | 222  |
| 6×DNA Loading Buffer                 | GH101-01       | 5×1 ml        |      |
| 2×RNA Loading Buffer                 | GH201-01       | 1 ml          |      |
| ddH <sub>2</sub> O                   | GI101-01       | 25 ml         |      |
| RNase-free Water                     | GI201-01       | 25 ml         |      |
| T7 High Efficiency Transcription Kit | JT101-01       | 20 μl×25 rxns |      |
| shRNA Synthesis Kit                  | JT111-01       | 20 μl×25 rxns | 225  |



Phone: +86-10-57815027  
+86-10-57815087  
Fax: +86-10-57815011  
Customer Service: +86-400-898-0321  
E-mail: [trade@transgen.com.cn](mailto:trade@transgen.com.cn)  
Technical Service: [trans@transgen.com.cn](mailto:trans@transgen.com.cn)  
Website: [www.transgenbiotech.com](http://www.transgenbiotech.com)



WeChat

---

Trans, The Highest Quality Products For Life Sciences